# Carbon Dioxide Test Kit

1.25 to 25 mg/L, 2 to 40 mg/L, 5 to 100 mg/L CO<sub>2</sub>

For test kit 143601 (CA-23)

DOC326.98.00004

Additional copies available on [www.hach.com](http://www.hach.com)

## Test preparation

- Rinse labware with deionized water between tests.
- When titrating, count each drop of titrant. Hold the dropper vertically. Swirl after each drop is added.

**CAUTION:** Handle chemical standards and reagents carefully. Review Material Safety Data Sheets for safe handling, storage and disposal information.

## Required items

Description	Unit	Catalog no.
Measuring Tube	each	43800
Mixing Bottle	6/pkg	232706
Phenolphthalein Indicator Solution	15 mL (½ oz) SCDB <sup>1</sup>	189736
Sodium Hydroxide Solution, 0.01 N	100 mL MDB <sup>2</sup>	67132


<sup>1</sup> Self-contained dropping bottle

<sup>2</sup> Marked dropping bottle

## Optional items

Description	Unit	Catalog no.
Deionized Water	500 mL	27249

## Low range (1.25 to 25 mg/L) and medium range (2 to 40 mg/L) test procedure


**1. Low Range:** Fill the bottle to the 23-mL mark with sample.


**Medium Range:** Fill the bottle to the 15-mL mark with sample.

**2.** Add one drop of Phenolphthalein Indicator Solution.

**3.** Add Sodium Hydroxide Solution by drops. Count the drops until the color changes to light pink and persists for 30 seconds. Swirl to mix after each drop.


**4. Low Range:** Calculate the result. Each drop of Sodium Hydroxide Solution used in step 3 equals 1.25 mg/L carbon dioxide (CO<sub>2</sub>).

**Medium Range:** Calculate the result. Each drop of Sodium Hydroxide Solution used in step 3 equals 2 mg/L carbon dioxide (CO<sub>2</sub>).


## High range (5 to 100 mg/L) test procedure


**1.** Fill the plastic tube to the top with sample. Pour the sample into the bottle.

**2.** Add one drop of Phenolphthalein Indicator Solution.

**3.** Add Sodium Hydroxide Solution by drops. Count the drops until the color changes to light pink and persists for 30 seconds. Swirl to mix after each drop.

**4.** Calculate the result. Each drop of Sodium Hydroxide Solution used in step 3 equals 5 mg/L carbon dioxide (CO<sub>2</sub>).

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