# BEMER Therapy System Evo

System for targeted and full-body therapy  
by means of magnetic fields and light

---

## Instructions for use

EN | FR | MY

## Instructions for use

Translation of the original instructions for use

### **BEMER Int. AG**

Austrasse 15  
 LI-9495 Triesen  
 www.bemergroup.com  
 Tel.: +423 399 39 99  
 Fax: +423 399 39 98  
 Collection ID 2590, Version 4  
 Date: 06/2024

### **Emergo Consulting (UK) Limited**

c/o Cr360 — UL International  
 Compass House, Vision Park Histon  
 Cambridge CB24 9BZ  
 United Kingdom  
 Phone: +44(0) 1223 772 671  
 Email: UKRPvigilance@ul.com

## Table of Contents

<b>1. General information</b>	<b>6</b>
1.1 Instructions for use	7
1.2 Symbols	8
1.3 Liability	9
1.4 Copyright notice	9
1.5 Warranty	10
1.6 Notification of incidents	10
<b>2. Safety</b>	<b>11</b>
2.1 Intended use	11
2.1.1 Intended purpose	11
2.1.2 Medical applications	11
2.1.3 Indications for use	11
2.1.4 Contraindications	12
2.1.5 Potential side effects	12
2.1.6 Intended use	13
2.1.7 Intended users	13
2.1.8 Intended patient group	14
2.1.9 Body regions	14
2.1.10 Intended environment and field of application	14
2.2 Clinical benefits	14
2.3 General safety instructions	15
2.4 Medical notes	20
2.4.1 Medical notes on PEMF	20
2.4.2 Medical notes on LLLT	21
2.5 Notes for professional users of medical products	22
2.6 Safety labels	22
2.7 Protective equipment	23
2.8 Environmental protection	23
<b>3. Technical data</b>	<b>24</b>
3.1 Storage, transport, and operating conditions	24
3.2 Product-specific technical data	25
3.3 EMC conformity	32

<b>4. BEMER Therapy System Evo</b>	<b>37</b>
4.1 System description	37
4.2 System overview	38
4.3 Scope of delivery for sets, packages, and individual products	40
4.4 Main components of the BEMER Therapy System Evo	44
4.5 Applicators/application modules and their use	46
4.5.1 B.Body Evo (PEMF) – full-body applicator for universal use	46
4.5.2 B.Bed Evo (PEMF) – full-body applicator for use in bed	46
4.5.3 B.Pad Evo (PEMF) – universal applicator for targeted treatment	47
4.5.4 B.Spot Evo (PEMF) – applicator module for targeted treatment	48
4.5.5 B.Sit Evo (PEMF) – applicator for targeted treatment while sitting	48
4.5.6 B.Light Clear Evo/B.Light Restore Evo (LLLT) – for targeted light therapy	48
4.5.7 B.Grip Evo interchangeable adapter	49
4.6 Accessories	50
4.6.1 B.Box Evo Battery Stand	50
4.6.2 B.Grip Evo Attachment Strap	50
4.6.3 B.Box Evo Power Supply	51
4.6.4 B.Box Evo Car Adapter	51
4.6.5 B.Light Evo Safety Glasses	52
4.6.6 B.Box Evo Wall Mount	52
<b>5. Transport and storage</b>	<b>53</b>
5.1 Safety	53
5.2 Symbols on the packaging	53
5.3 Storage of the packaging	53
<b>6. Recommended use</b>	<b>54</b>
6.1 Recommendations for application of magnetic field therapy (PEMF)	54
6.1.1 The plus Signal	54
6.1.2 Full-body therapy	54
6.1.2.1 Intensity modes and sleep program (basic plan)	54
6.1.2.2 Sleep program	56
6.1.3 Targeted treatment	56
6.2 Recommended use of light therapy (LLLT)	57
6.2.1 B.Light Clear Evo	58

6.2.2 B.Light Restore Evo	58
<b>7. Commissioning</b>	<b>59</b>
7.1 Connecting the B.Box Evo Stand, Wall Mount, and Battery Stand	59
7.2 Connecting the power supply	60
7.3 Connecting the applicators and application modules	60
7.4 Switching on the B.Box Evo	60
<b>8. Using the BEMER Therapy System Evo</b>	<b>62</b>
8.1 Switch the unit on (start screen)	62
8.2 Main menu	63
8.2.1 Status bar	63
8.2.2 Settings (start screen)	64
8.2.3 Adjust the signal volume	64
8.2.4 Adjust the display brightness	65
8.2.5 Set the current time	66
8.2.6 Expert mode editor	66
8.2.7 System	67
8.2.7.1 Switch the ambient light on/off	68
8.2.7.2 Enable / disable energy saving mode	68
8.2.7.3 Select system language	68
8.2.7.4 Set the time format	68
8.2.7.5 System information	69
8.3 Main menu	69
8.3.1 Intensity mode	70
8.3.2 Program mode	71
8.3.3 Sleep program	72
8.3.4 Light therapy	74
<b>9. Cleaning and care</b>	<b>75</b>
9.1 Cleaning	75
9.2 Disinfection	76
<b>10. Disposal</b>	<b>77</b>
<b>11. Error messages and remedies</b>	<b>78</b>
<b>12. Meaning of the symbols on the labels (devices and packaging)</b>	<b>80</b>

## 1. General information

Thank you for purchasing our BEMER Therapy System Evo and for the confidence you have placed in us. The BEMER Therapy System Evo is a versatile and flexible product thanks to the various applicators. The BEMER Therapy System Evo is your daily companion, regardless of whether your goal is to prevent illness and maintain an active lifestyle or to support a prescribed course of treatment. (Please also refer to the information in chapter 2.)

Please read these instructions for use carefully before using the system for the first time. By doing so, you will avoid damage to the device and ensure that the warranty is not invalidated.

When purchasing the BEMER Therapy System Evo, please ensure that you are given an introduction to the system by an official and certified BEMER partner.

If you have further questions or require training, our customer service team will be happy to assist you.

The original language of these instructions for use is German.

### 1.1 Instructions for use

These instructions for use are an integral part of the BEMER Therapy System Evo. They enable the user to use the BEMER Therapy System Evo safely and efficiently.

The user must have carefully read and understood these instructions for use before commissioning the system. The basic prerequisite for safe use is compliance with all specified safety instructions.

In addition to the information contained in these instructions for use, the local accident prevention regulations and occupational safety regulations apply.

These instructions for use must be kept in the immediate vicinity of the BEMER Therapy System Evo and within reach of the operator at all times.

#### Product images

Your BEMER Therapy System Evo may differ from the illustrations in this document. However, all descriptions are designed to apply analogously. If components are described that are not included in the scope of delivery, they are marked as optional.

#### Brand names and trademarks

Product and / or company names mentioned in these instructions for use may be registered trademarks of the respective companies.

#### Gender-specific wording

For easier readability, masculine nouns/pronouns are sometimes used in this document. However, its content applies equally to all readers regardless of their sex or gender.

## 1.2 Symbols

### DANGER

DANGER indicates an imminent hazard which, unless avoided, could result in severe injuries or even death.

### WARNING

WARNING indicates a potential hazard which, unless avoided, could result in severe injuries or even death.

### CAUTION

CAUTION indicates a potential hazard which, unless avoided, could result in minor injuries.

### NOTE

NOTE indicates a potential hazard which, unless avoided, could result in damage to the product and / or data loss.

#### Warning sign

Safety symbol that indicates a risk or hazard.

#### Command sign

Safety symbol that prescribes a specific action.

#### Prohibition sign

Safety symbol that indicates a prohibited action.

#### Information

Indicates user tips and generally useful information for optimal use of the product.

Use of symbols	Description	Example
●	A bullet point describes an action (activity)	● Switch on the device.
→	The result of an action (activity)	→ e.g., a new operating window opens
—	List - in no particular order	— B.Box Evo — B.Box Evo — ...
(Cross reference)	Reference to a chapter or page	(chapter 5.1)
<BUTTON>	Refers to the actuation of an operating element	e.g. <Save>
Menu - Submenu	Specifies the menu path	Settings - Time

## 1.3 Liability

In the case of damage and/or defects resulting from improper installation, assembly, or use of the product, or from failure to observe the instructions for use and/or the safety instructions, any liability on the part of BEMER Int. AG may be reduced or ruled out altogether and our warranty obligations may cease to apply. In the aforementioned cases, no warranty claims will be accepted.

## 1.4 Copyright notice

All contents of these instructions for use, in particular texts, photographs, and graphics, are protected by copyright. This form of legal protection also applies to databases and similar facilities. No part of these instructions for use may be reproduced in any form outside the narrow framework of copyright law without the written permission of BEMER Int. AG.

Anyone who infringes copyright (e.g., copies images or texts without permission) may be liable to prosecution, receive a warning combined with a financial penalty, or be required to pay damages. We reserve the right to enforce our rights in such cases.

## 1.5 Warranty

### Warranty conditions

When you purchase BEMER products, you are given the opportunity to find out about our warranty conditions. You can also view the current warranty conditions at any time in the relevant section of our websites.

Your statutory rights in the event of defects are not restricted by our warranty terms; you can therefore assert them free of charge. Consequently, your existing statutory warranty rights in relation to our products shall remain unaffected by our warranty commitment. The manufacturer's warranty terms therefore do not infringe your statutory rights, but rather strengthen your legal position.

## 1.6 Notification of incidents

If serious incidents occur while using this product, both the manufacturer (BEMER Int. AG) and the competent authorities in the regions in which the user is based must be informed.

## 2. Safety

BEMER products may only be used for the purposes described in this chapter. Any use of the product beyond the scope of the intended use is regarded as improper.

### 2.1 Intended use

#### 2.1.1 Intended purpose

The B.Box Evo product, together with the PEMF applicators, the B.Light Clear Evo, and the B.Light Restore Evo, is part of the "BEMER Therapy System Evo".

The BEMER B.Box Evo is only used to generate the electrical signal that is used for the pulsed electromagnetic field therapy (PEMF) and the low-level light therapy (LLLT), and to control the individual programs.

In the context of pulsed electromagnetic field therapy (PEMF), the PEMF applicators are used in conjunction with the B.Box Evo to stimulate blood flow in the small and very small blood vessels (microcirculation), and to alleviate certain clinical conditions.

The B.Light Clear Evo emits light with wavelengths of 465 nm and 645 nm ( $\pm 20$  nm) and is used in the context of low-level light therapy (LLLT) in conjunction with the B.Box Evo for skin treatment on or near the skin surface.

The B.Light Restore Evo emits light with wavelengths of 645 nm and 860 nm ( $\pm 20$  nm) and is used in the context of low-level light therapy (LLLT) in conjunction with the B.Box Evo for skin treatment on or near the skin surface.

#### 2.1.2 Medical applications

The B.Box Evo serves as the interface between the applicators and the users for LLLT and PEMF applications.

#### 2.1.3 Indications for use

The B.Box Evo has no indications for use per se.

The approved indications for use for PEMF and LLLT are defined by the applicators.

PEMF therapy is an adjuvant application and cannot replace medically prescribed therapy.

In the case of pre-existing conditions, the described user groups are limited to the following underlying conditions, their consequences and / or accompanying symptoms:

- Impaired wound healing
- Degenerative diseases of the musculoskeletal system
- Polyneuropathy due to diabetes mellitus or after cancer treatment

- Chronic fatigue, e.g., associated with chronic stress or multiple sclerosis
- Acute and chronic pain

LLLT is an adjuvant application and cannot replace medically prescribed therapy. It supports the treatment of skin diseases and is also intended as a supplement for cosmetic treatment.

#### **Example applications for the B.Light Clear Evo:**

- Treatment of mild to moderate acne vulgaris
- Improvement of the general appearance of the skin
- To help reduce inflammation (acne vulgaris)

#### **Example applications for the B.Light Restore Evo:**

- Cosmetic use: to help reduce the appearance of wrinkles and fine lines, to improve skin texture
- To support wound healing
- To help reduce inflammation
- To support muscles and joints

#### **2.1.4 Contraindications**

The B.Box Evo has no contraindications per se.

Contraindications for PEMF and LLLT are defined by the applicators.

PEMF therapy is contraindicated for the following user groups:

- Individuals fitted with active medical implants (e.g., medication pumps, pacemakers)
- Recipients of organ transplants, allogeneic cell transplants, bone marrow or stem cell transplants in combination with immunosuppressive therapy (= intentional suppression of the immune system)

LLLT therapy (B.Light Clear Evo, B.Light Restore Evo) is contraindicated for the following applications:

- Do not apply directly to mucous membranes or in the vicinity of the eye

#### **2.1.5 Potential side effects**

The B.Box Evo has no potential side effects per se.

The potential side effects of PEMF und LLLT are defined by the applicators.

When PEMF therapy is applied via the corresponding PEMF applicators, the following short-term side effects may occur in very rare cases:

- Change in pulse rate
- Change in blood pressure

When using LLL therapy by means of the associated light applicator modules (B.Light Clear Evo, B.Light Restore Evo), the following localized and short-term skin reactions may occur in very rare cases:

- Redness of the skin (erythema)
- Itching
- Burning / stinging
- Dry skin
- Hyperpigmentation

#### **2.1.6 Intended use**

The intended use is defined by the applicators.

The PEMF applicators are used together with the B.Box Evo for regular systemic application and additional targeted treatment. Depending on the applicator, different forms of application are possible.

Systemic (regular application / see chapter 6, basic plan):

B.Body Evo and B.Bed Evo are used for full-body therapy in the prone state.

Localized (optional):

B.Pad Evo, B.Sit Evo, and B.Spot Evo are additionally used for targeted treatment of individual body regions.

The B.Grip Evo holder module (Class I) is required for use of the B.Spot Evo.

The LLLT application modules B.Light Clear Evo and B.Light Restore Evo are used in conjunction with the B.Box Evo for targeted local application of polychromatic light on or near the skin surface.

The B.Grip Evo holder module (Class I) is required for use of the B.Light Clear Evo.

#### **2.1.7 Intended users**

The system is intended for use by end users aged 14 years and older, as well as by medically trained personnel.

### 2.1.8 Intended patient group

The device is intended for use on people aged 14 years or older in accordance with the respective indications for use and contraindications.

Children under 14 years of age and individuals with impaired physical, sensory or mental capabilities must be supervised and / or instructed by a suitable person who accepts responsibility for their safety.

### 2.1.9 Body regions

The body region to be treated is defined by the applicators.

The whole-body applicators (B.Body Evo and B.Bed Evo) are designed for systemic application of the PEMF therapy. Due to the distribution of the coils, all regions of the body are reached by the magnetic field.

The local applicators (B.Spot Evo, B.Pad Evo and B.Sit Evo) are optionally used for targeted application of the magnetic field in specific areas of the body.

The LLLT applicators are optionally used for targeted application of LLLT on specific areas of the body.

### 2.1.10 Intended environment and field of application

The product is intended exclusively for the application of PEMF and LLLT in conjunction with the applicators by laypersons in a home environment and by professional users in a clinical setting.

The PEMF applicators and LLT applicator modules, in combination with the B.Box Evo, are intended for use by laypersons in a home environment and by medically trained personnel in a clinical or outpatient setting.

### 2.2 Clinical benefits

The B.Box Evo itself has no clinical utility.

Users of PEMF therapy benefit from improved blood flow to the capillary network, in particular in the smaller and smallest blood vessels, and thus from improved tissue supply, which is desirable for various health conditions.

The application of LLLT using the B.Light Clear Evo leads to an improvement in the skin's appearance, especially in cases of mild to moderate acne vulgaris.

The application of LLLT using the B.Light Restore Evo leads to an improvement in the skin's appearance.

### 2.3 General safety instructions

#### WARNING

#### Interference with active implants due to electromagnetic forces (PEMF)

Active implants (e.g., pacemakers, insulin pumps, etc.) may be disrupted by electromagnetic forces.

- Do not use the BEMER Therapy System Evo (PEMF) under any circumstances if you are a patient with an active implant.

#### Risk of strangulation due to loose cables

Loose cables and leads pose a risk of injury, e.g., caused by tripping or strangulation.

- Lay the cables flat to ensure they do not create a tripping hazard.
- Use the supplied fastening aids to lay the cables.

#### Risk of burns due to damaged or worn out B.Box Evo Battery Stand

The failure of safety mechanisms may lead to spontaneous combustion or explosions.

- Never open or use lithium batteries or accumulators that are swollen, deformed, outgassed, have "leaked", or are covered with a "greasy film" or external deposits in the area of the poles are damaged. Damaged batteries are associated with an increased hazard risk.
- Therefore, dispose of these batteries and rechargeable batteries immediately, preferably at an electronics store or recycling center – and in a manner that ensures that the staff there can accept them.
- Speak to the qualified personnel there and point out the damage.

#### Life-threatening hazards due to non-compliance with the safety instructions

Failure to follow the instructions for use correctly may result in operating errors and life-threatening hazards.

- Always read and familiarize yourself with the instructions for use included in the scope of delivery.
- Observe the safety notes.

**WARNING****Risk of electric shock due to use of a damaged power cable or a power cable not approved by the manufacturer**

Contact with exposed electrical parts or power cables not approved by the manufacturer may result in electric shock.

- Disconnect the device from the power supply.
- Only use the power cables approved by the manufacturer.

**Electric shock due to use of the device in a humid environment**

Water and electricity are a dangerous combination that can cause an electric shock.

- Do not use this device in a humid environment (e.g., in the bathroom or near a shower or swimming pool).
- Do not allow water to run into the unit.

**Electric shock due to incorrect electrical voltage**

Incorrect electrical voltage from the local power grid may cause electric shock and permanently damage the device.

- Check that the voltage specified on the device matches the local mains voltage before connecting the device to avoid the risk of electric shock or permanent damage to the device.

**Risk of infection due to transmission of disease carriers**

Shared use of the application modules can lead to the transmission of diseases.

- Between each use, clean the applicators using the cleaning and disinfecting agent recommended by the manufacturer.

**WARNING****Risk of burns / fire hazard due to overheated equipment**

Overheated unattended equipment may result in an increased fire hazard and therefore burns.

- To avoid the risk of fire or burns, do not leave the unit unattended when it is switched on.
- The device must not be used by individuals with reduced physical, sensory or mental faculties and / or who lack the necessary experience and knowledge, unless they are supervised or given instructions on how to use the equipment to avoid the risk of fire or burns.
- The device is not intended for use by children. Ensure that children are supervised and do not play with the device to avoid the risk of fire and burns.

**Risk of infection due to application on injured skin**

The use of contaminated application modules on injured skin can result in the transmission of diseases.

- Do not apply the application modules to injured skin.
- Between each use, clean and disinfect the applicators using the cleaning and disinfecting agent recommended by the manufacturer.

**CAUTION****Risk of blinding due to optical radiation**

If the B.Light Clear Evo or B.Light Restore Evo is used in the vicinity of the eye, there is a risk of injury to the retina.

- Always wear the supplied safety glasses when using the B.Light Evo application modules.

**Toxic skin reactions during application with ointments and medications**

Use of the B.Light Clear Evo and / or B.Light Restore Evo light application module in conjunction with light-intensive or light-reactive ointments or medications may result in toxic skin reactions.

- Do not use light therapy in combination with ointments or medications.

**CAUTION****Risk of burns due to high leakage currents**

Simultaneous contact with two metallic parts can result in high leakage currents, which can cause skin burns.

- Do not touch any metallic parts during use.

**Electric shock due to water entering the device**

Water ingress may cause a short circuit in the device and endanger the user.

- Disconnect live parts from the power supply before cleaning them.
- The electrical contacts on the B.Box Evo Rechargeable Battery (B.Box Evo Battery Stand) must not come into contact with liquids.

**Crushing injuries due to components with magnets**

When the B.Grip Evo holder module is placed in close proximity to the B.Spot Evo, B.Light Clear Evo, and B.Light Restore Evo application modules, skin pinching may occur due to the resulting magnetic attraction. Skin bruising may also occur when attaching the B.Box Evo to the B.Box Evo Stand or the B.Box Evo Battery Stand.

- Follow the instructions on the corresponding components and do not reach between the holder module and the application modules or the B.Box Evo and the B.Box Evo Stand or the B.Box Evo Battery Stand.

**Increase in vital signs due to incorrect operation of the device**

User errors may result in elevated vital signs, e.g., higher blood pressure, especially among untrained users.

- When purchasing the BEMER Therapy System Evo, ensure that you are given an introduction to the system by an official and certified BEMER partner.
- Always read and familiarize yourself with the instructions for use included in the scope of delivery.
- Observe the safety notes.

**CAUTION****Allergic reactions due to material incompatibility**

The materials used in the applicators may cause skin reactions due to intolerances.

- If this occurs, discontinue use of the therapy system and consult your doctor.

**NOTE****Increased surface temperature of the LLLT application modules due to high ambient temperature**

At ambient temperatures above 35°C, the surface of the light therapy application modules can reach 44°C.

- If the ambient temperature exceeds 35°C, allow the light therapy application modules to cool for at least 10 minutes between treatments.

**Performance degradation of portable RF communication devices due to electromagnetic forces**

The performance of the respective devices may be affected by electromagnetic forces.

- Do not use portable RF communication devices (including accessories such as antenna cables and external antennas) within 12 inches (30 cm) of any part of the BEMER Therapy System Evo, including the cables specified by the manufacturer.

**Using the device in an unsuitable environment**

Use of this device in the immediate vicinity of other equipment and / or in a humid environment can lead to malfunctions.

- Do not use this device when it is directly adjacent to or stacked on/under other equipment.
- If such use is unavoidable, monitor the device as well as the other equipment to ensure normal operation.
- Only use the BEMER Therapy System Evo in dry rooms.

## NOTE

**Electronic storage media can be corrupted or erased**

The integral magnets in the B.Box Evo Stand and some of the cable connectors are very strong. Electromagnetic fields can disrupt the functions of storage media (e.g., credit and EC cards, data carriers) and can erase their contents.

- Keep storage media of this kind away from the magnets.

**Property damage due to maintenance and repairs by unauthorized personnel**

Repairs and maintenance activities performed by unauthorized and / or unqualified individuals may result in material damage to the device.

- Maintenance and repair work may only be carried out by authorized personnel.

- There are no user-replaceable parts and no maintenance is required during the service life of the system.

**2.4 Medical notes****2.4.1 Medical notes on PEMF**

In case of atypical reactions to the PEMF therapy, a doctor should be consulted.

For conditions requiring immunosuppression that are not related to transplantation, e.g., autoimmune diseases or dermatological diseases, there are no contraindications for use of the PEMF therapy.

At the end of the first cycle, new users who regularly take blood thinning/anticoagulant agents or antihypertensive medications are advised to consult their attending physician to check for any change in their effectiveness.

If the following circumstances or complaints exist, the prior approval of the attending physician/specialist must be obtained before application of the PEMF therapy:

- Unexplained fever
- Infectious diseases
- Severe cardiac rhythm disorder
- Severe psychosis
- Uncontrolled seizure disorders (e.g., epilepsy)
- Long-term use of  $\beta$ -receptor antagonists
- Use of high-dose corticoid agents

- Long-term use of anticoagulants (coumarin derivatives)
- Ongoing use of prescription medications
- Pregnancy
- Tumor diseases

**2.4.2 Medical notes on LLLT**

In case of atypical reactions to LLLT, a doctor should be consulted.

If the following circumstances or complaints exist, the prior approval of the attending physician/specialist must be obtained before application of LLLT:

- Light-induced seizures (light sensitivity)
- Light-induced migraine headaches
- Ingestion or application of ointments (cosmetics), medications, or supplements known to cause photosensitivity
- Ongoing use of prescription medications
- Allergic reaction caused by light
- Tumor diseases
- Cancerous lesions on the skin
- Skin lesions caused by bacteria, viruses or fungi

### 2.5 Notes for professional users of medical products

Professional users must ensure that the respective employees are aware of and implement the applicable occupational health and safety requirements. Furthermore, they must ensure that all employees have read and understood the instructions for use.

Professional users must train their employees at regular intervals, inform them about potential hazards, and provide them with the necessary protective equipment.

Individuals who are being trained, instructed, or are considered trainees in general may only work on the BEMER Therapy System Evo under the constant supervision of an experienced person.

Work on electrical components may only be carried out by qualified personnel with appropriate training and in compliance with all applicable provisions of the accident prevention regulations. A safety check must be carried out by the operator at regular intervals.

### 2.6 Safety labels

Symbol	Definition	City
	All users must read the instructions for use	This symbol is shown on each product label
	Patients with active implants must not use the BEMER Therapy System Evo (PEMF)	This symbol is shown on the back of the B.Box Evo
	Warning of high leakage currents	This symbol is shown on the inside of the B.Grip Evo
	Warning of crushing injuries due to components containing magnets	This symbol is shown on the inside of the B.Grip Evo

### 2.7 Protective equipment

The Beauty Pack Evo includes the following protective equipment

- B.Light Evo Safety Glasses

The B.Light Evo Safety Glasses must be worn when using the B.Light Clear Evo or the B.Light Restore Evo.

### 2.8 Environmental protection

BEMER Int. AG manufactures state-of-the-art therapy systems in terms of safety and environmental protection. These therapy systems do not pose any danger to human health or the environment, provided that they are operated properly.

#### CAUTION

#### Risk to human health and the environment due to toxic materials

Environmentally harmful materials, which the BEMER Therapy System Evo may contain, pose a risk to human health and the environment.

- The BEMER Therapy System Evo must not be disposed of with industrial or household waste, either in whole or in part.

Batteries contain toxic heavy metals. They are subject to hazardous waste treatment and must be handed in at municipal collection points or disposed of by a specialist company.

### 3. Technical data

#### 3.1 Storage, transport, and operating conditions

Temperature range (operating)	+5 to 40°C
Humidity (operating)	15 to 90% (non-condensing)
Ambient air pressure (operating)	700 to 1060 hPa
Temperature range (storage, transport)	-25 to +70°C
Humidity (storage, transport)	10 to 90% (non-condensing)
Ambient air pressure (storage, transport)	500 to 1060 hPa
Time until the operating temperature range is reached, starting from the minimum transport temperature	~ 30 minutes
Time until the operating temperature range is reached, starting from the maximum transport temperature	~ 30 minutes

#### 3.2 Product-specific technical data

Article number	424000
Product designation	B.Box Evo
Product type	Control unit
Dimensions (L x W x D) in mm	210 x 150 x 43
Weight (g)	926
Type	Portable device
Surface material	PC/ABS, aluminum, glass
IP degree of protection	22
Protection class (IEC 61140)	SK II
Protection against electric shock	2MOPP Class II
EMC class (CISPR 11:2009)	Class B
Input voltage	100 to 240 V AC / 50 to 60 Hz
Output voltage	15VDC / 2A
Operating voltage (V)	15
Max. power (watts)	30
Display dimensions (inches)	7"
Display resolution (px)	1024 x 600
Display brightness (cd/m <sup>2</sup> )	450
Display viewing angle (°)	80
Display contrast ratio	800:1
Location of type plate	Rear of device

<b>Designation</b>	<b>B.Body Evo (full-body applicator)</b>	<b>B.Bed Evo (full-body applicator)</b>	<b>B.Pad Evo (local applicator)</b>
Article number	434300	434400	434100
Dimensions L x W x D in cm	152 x 59 x 3	190 x 90 x 1	133 x 15 x 0.8
Weight in kg	1.9	1.97	0.33
Number of copper coils	16	16	4
Average flux density	≈ 35 μT (max. level)	≈ 35 μT (max. level)	≈ 100 μT (max. level)
Average flux density plus	≈ 50 μT (max. level)	≈ 50 μT (max. level)	≈ 150 μT (max. level)
Number of external connections	1 x magnetic connector with flexible cable with PVC insulation	1 x magnetic connector with flexible cable with PVC insulation	1 x magnetic connector with flexible cable with PVC insulation
Composition of the surface material in contact with the body	100% PES	100% PES	100% PES
Cable length	250 cm	250 cm	250 cm
Protection against moisture	IP22	IP22	IP22
Device class	Applied part, type BF	Applied part, type BF	Applied part, type BF
Location of type plate	Rear of applicator	Rear of applicator	Rear of applicator

<b>Designation</b>	<b>B.Spot Evo (local applicator module)</b>	<b>B.Sit Evo (local applicator)</b>	<b>B.Grip Evo (holder module)</b>
Article number	434000	434200	454000
Dimensions L x W x D in cm	12.12 x 12.12 x 2.46	44 x 36.7 x 5	12 x 12 x 4.8
Weight in kg	0.185	1.88	0.213
Number of copper coils	1	1	-
Average flux density	≈ 100 μT (max. level)	≈ 100 μT (max. level)	-
Average flux density plus	≈ 150 μT (max. level)	≈ 150 μT (max. level)	-
Number of external connections	Rotation-free 5-pole tracks	1 x magnetic connector with flexible cable with PVC insulation	Rotation-free 5-pole tracks
Material composition (surface)	100% PC	Surface material in contact with the body 66% PES   12% Rayon   2% Spandex 20% TPU film	PC/ABS
Cable length	Cable on B.Grip Evo	250 cm	250 cm
Protection against moisture	IP22	IP22	IP22
Device class	Applied part, type BF	Applied part, type BF	-
Location of type plate	Rear of applicator module	Rear of applicator	Inside of holder module

<b>Designation</b>	<b>B.Light Clear Evo (light applicator module)</b>	<b>B.Light Restore Evo (light applicator module)</b>	<b>B.Box Evo Battery Stand</b>
Article number	434500	434600	454100
Dimensions L x W x D in cm	12.12 x 12.12 x 2.46	12.12 x 12.12 x 2.46	18.3 x 10.8 x 9.8
Weight in kg	0.12	0.12	0.498
Supply voltage	-	-	7.2 VDC
Wavelength range	465 nm and 645 nm ( $\pm 20$ nm)	645 nm and 860 nm ( $\pm 20$ nm)	-
Number of LEDs	100	100	-
Maximum radiation intensity (mW/cm <sup>2</sup> ) at the skin surface	465 nm: $\approx 0.8$ 645 nm: $\approx 1.2$	645 nm: $\approx 0.56$ 860 nm: $\approx 1.4$	-
Treatment area	(J/cm <sup>2</sup> ) on the skin during a 480 sec. treatment	(J/cm <sup>2</sup> ) on the skin during a 480 sec. treatment	-
Treatment dose (J/cm <sup>2</sup> ) at the skin surface per 480 sec. treatment	$\approx 1$	$\approx 1$	-
Material composition (surface)	100% PC	100% PC	80 % Al, 20% PC/ABS
Protection against moisture	IP22	IP22	IP22
Risk group	1	0	-
Device class	Applied part, type BF	Applied part, type BF	-
Battery type	-	-	Li-Ion
Battery capacity	-	-	48 Wh
Number of treatments	-	-	$\approx 50$ @ 8 minutes each
Number of external connections	-	-	5-pin connector
Location of type plate	Rear of applicator module	Rear of applicator module	Underside of B.Box Evo Battery Stand

<b>Designation</b>	<b>B.Box Evo Car Adapter</b>	<b>B.Box Evo Power Supply</b>	<b>B.Light Evo Safety Glasses</b>
Article number	444100	444000	454900
Dimensions L x W x D in cm	10 x 4.5 x 2.25	10 x 4.5 x 2.25	14.1 x 6.2 x 4.5
Weight in kg	0.14	0.14	0.031
Primary voltage	10 to 32 VDC / 4A	100 to 240 VAC / 50-60Hz	-
Number of external connections	1	-	-
Cable length	360 cm	360 cm	-
Power supply	15 V DC / 2A	15 V DC	-
Protection class	-	SK II	-
Protection against moisture	IP21	IP21	IP22
Lenses	-	-	Shade 2, GA 166 CE (anti-fog, scratch-resistant, 100% UV protection)
Material composition	PC/ABS	PC/ABS	Polycarbonate (PC)
Location of type plate	Underside of Car Adapter	Underside of B.Box Evo Power Supply	no type plate on safety glasses

<b>Designation</b>	<b>Extension cord for Evo applicators</b>	<b>Velcro extension strap for B.Pad Evo</b>	<b>B.Box Evo Stand</b>
Article number	454500	454400	454800
Dimensions L x W x D in cm	-	30.0 x 8.0 x 0.02	18 x 13 x 8
Weight in kg	0.118	0.031	0.48
Material composition	Nylon / PC / ABS	85% PES, 15% Spandex	Aluminium
Cable length, Ø in cm	250, 0.55	-	-
Protection against moisture	IP22	IP21	IP22
Location of type plate	Product packaging	Product packaging	Underside of B.Box Evo Stand

<b>Designation</b>	<b>B.Grip Evo Fastening Strap</b>	<b>B.Box Evo Wall Mount</b>	<b>B.Bed Evo Fastening Strap</b>
Article number	454200	454600	454700
Dimensions L x W x D in cm	80.2 x 9.2 x 1.2	17.9 x 12.3 x 1.2	Flexible length (elastic material)
Weight in kg	0.085	0.12	0.09
Material composition	50% PU, 40% PA, 10% SP	PC/ABS, sheet steel	Spandex
Location of type plate	no type plate on Fastening Strap	Underside of B.Box Evo Wall Mount	Product packaging

<b>Designation</b>	<b>B.Body Evo foot protection</b>
Article number	450500
Dimensions L x W x D in cm	67.5 x 60.0 x 0.02
Weight in kg	0.23
Material composition	PES
Location of type plate	Product packaging

### 3.3 EMC conformity

The BEMER therapy system Evo is intended for operation in an electromagnetic environment in which immunity to RF disturbances is uncontrolled.

#### Electromagnetic Emission

<b>Emission Measurements</b>	<b>Conformity</b>	<b>Electromagnetic Environment – Guidelines</b>
RF_emissions in accordance with CISPR 11	Group 1	The BEMER therapy system uses RF energy only for its internal function. As a result, its RF emission is very low and it is unlikely to disturb any adjacent electronic devices.
RF_emissions in accordance with CISPR 11	Class B	The BEMER therapy system is intended for use in all institutions including residential areas and those that are directly connected to a public grid that also supplies buildings used for residential purposes.
Harmonics in accordance with IEC 61000-3-2	Class A	
Voltage fluctuations/flicker in accordance with IEC 61000-3-3	Compliant	

### Recommended protective distances between portable and mobile RF telecommunication devices and the BEMER therapy system

Customers or users can help to avoid electromagnetic disturbances by complying with the minimum distance between portable and mobile RF telecommunication devices (transmitters) and the BEMER therapy system, depending on the output of the communication device as indicated in the table on the next page.

<b>Rated output of transmitter [W]</b>	<b>Protective distance, depending on transmission frequency [m]</b>			
	150 kHz to 80 MHz outside ISM bands	150 kHz to 80 MHz inside ISM bands	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d = \frac{3.5}{3} \sqrt{P} = 1.17\sqrt{P}$	$d = \frac{12}{10} \sqrt{P} = 1.2\sqrt{P}$	$d = \frac{12}{10} \sqrt{P} = 1.2\sqrt{P}$	$d = \frac{23}{10} \sqrt{P} = 2.3\sqrt{P}$
0,01	0,12	0,12	0,12	0,23
0,1	0,37	0,38	0,38	0,73
1	1,17	1,2	1,2	2,3
10	3,69	3,79	3,79	7,27
100	11,67	12	12	23

For transmitters where the maximum rated output is not listed in the table above, the distance can be determined by using the equation that belongs to the relevant column, in which P is the maximum rated output of the transmitter in watts (W) according to information from the transmitter manufacturer.

## Electromagnetic immunity

Emission measurements	IEC-60601 test level	Compliance level	
Discharge of static electricity in accordance with IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be made of wood or concrete or covered with ceramic tiles.
Electrical fast transient/burst immunity in accordance with IEC 61000-4-4	± 2 kV for power cables ± 1 kV for input/output cables	± 2 kV for power cables ± 1 kV for input/output cables	The quality of the supply voltage should correspond to that of a typical store or clinic environment.
Impulse voltage/surges in accordance with IEC 61000-4-5	± 1 kV voltage phase conductor – phase conductor ± 2 kV voltage phase conductor – ground	± 1 kV voltage phase conductor – phase conductor ± 2 kV voltage phase conductor – ground	
Magnetic field at supply frequency (50/60 Hz) in accordance with IEC 61000-4-8	3 A/m	3 A/m	
Voltage dips, short-term interruptions and fluctuations in the supply voltage in accordance with IEC 61000-4-11	< 5 % $U_T$ for 1/2 period (> 95 % dip) < 40 % $U_T$ for 10 periods (60 % dip) < 70 % $U_T$ for 25 periods (30 % dip) < 5 % $U_T$ 5 s (> 95 % dip)	< 5 % $U_T$ for 1/2 period (> 95 % dip) < 40 % $U_T$ for 10 periods (60 % dip) < 70 % $U_T$ for 25 periods (30 % dip) < 5 % $U_T$ 5 s (> 95 % dip)	The quality of the supply voltage should correspond to that of a typical store or clinic environment. If the users of the device require continuous operation, even in case of power outages, it is recommended to feed power to the device via an uninterruptible power supply or a battery.
Radiated fields in close proximity in accordance with IEC 61000-4-39	8 A/m at 30kHz 65 A/m at 134.2kHz 7.5 A/m at 13.56kHz	8 A/m at 30kHz 65 A/m at 134.2kHz 7.5 A/m at 13.56kHz	Avoid exposure to known sources of EMI (electromagnetic interference) such as diathermy, lithotripsy, electrocautery, RFID (Radio Frequency Identification), and electromagnetic security systems such as anti-theft/electronic article surveillance systems, metal detectors. Note that the presence of RFID devices may not be obvious. If such interference is suspected, reposition the equipment if possible, to maximize distances
<b>Remark: <math>U_T</math> is the AC mains voltage prior to applying the test level</b>			

Immunity tests	ICE-60601 test level	Compliance level
Conducted RF disturbance in accordance with IEC 61000-4-6	3 V <sub>eff</sub> 150kHz to 80MHz inside the ISM/ amateur band	3V <sub>eff</sub>
	6 V <sub>eff</sub> 150kHz to 80MHz inside the ISM/ amateur band	6V <sub>eff</sub>
Radiated RF disturbance in accordance with IEC 61000-4-3	10V/m 80 Mhz to 2.7GHz	10V/m 80 Mhz to 2.7GHz
<b>Immunity to RF telecommunication equipment pursuant to limit values of IEC 61000-4-3 tested and passed</b>		

**Electromagnetic Environment – Guidelines**

Recommended protective distance:  $d = 1.2 \sqrt{P}$

$d = 1.2 \sqrt{P}$  for 80 MHz to 800 MHz

$d = 2,3 \sqrt{P}$  for 800 MHz to 2,5 GHz

With P as rated output of the transmitter in watts (W) pursuant to information from the transmitter manufacturer and d as the recommended protective distance in meters (m). For all frequencies pursuant to an examination on site, the field strength of stationary radio transmitters is lower than the compliance level. Disturbances can occur in proximity to devices that are marked with the following symbol:

**Remark 1:** The higher value applies for 80 MHz and 800 MHz.

**Remark 2:** These guidelines may not apply in all situations. The radiation of electromagnetic waves is influenced by absorption and reflections from buildings, objects and people.

(a) Theoretically, the field strength of stationary transmitters such as the base stations of mobile telephones and mobile terrestrial radio equipment, amateur radio stations, AM and FM radio and television transmitters cannot be determined in advance. To determine the electromagnetic environment with respect to the stationary transmitter, a study of the location should be considered.

If the measured field strength at the location at which the BEMER therapy system is used exceeds the compliance level indicated above, the BEMER therapy system should be observed in order to document its intended function. If atypical cable properties are observed, additional measures, such as a change in alignment or a different location for the BEMER therapy system, may be necessary.

(b) The field strength should be lower than 3 V/m across the 150 kHz to 80 MHz frequency range.

## 4. BEMER Therapy System Evo

The applicators (B.Pad Evo, B.Body Evo, B.Bed Evo, and B.Sit Evo) and application modules (B.Spot Evo, B.Light Restore Evo, B.Light Clear Evo) are operated via the controller (B.Box Evo). The application modules additionally require the interchangeable adapter (B.Grip Evo).

Power is supplied either via the mains or the rechargeable B.Box Evo Battery Stand (optional accessory). The medically approved B.Box Evo Car Adapter also allows use of the system via the on-board power supply of a motor vehicle or boat. Only use the system when the vehicle is stationary and ensure the B.Box Evo is appropriately secured.

### 4.1 System description

The BEMER Therapy System Evo is a medical product intended for daily use on humans, which aims to stimulate microcirculation through the use of pulsed electromagnetic fields (PEMF). It is also used for treatment of the skin through exposure to light of specific wavelengths (low-level light therapy or LLLT).

The BEMER Therapy System Evo consists of various applicators (B.Body Evo, B.Bed Evo, B.Pad Evo and B.Sit Evo) and applicator modules (B.Spot Evo, B.Light Clear Evo und B.Light Restore Evo). The applicator modules are connected to the B.Box Evo controller together with the holder module (B.Grip Evo). Without the B.Box Evo controller, neither the individual applicators nor the application modules can be used.

## 4.2 System overview

### Controller with accessories

CE 0483

B.Box Evo

B.Box Evo Stand

B.Box Evo Car Adapter

B.Box Evo Power Supply

B.Box Evo Battery Stand

### Applicators

B.Grip Evo

B.Spot Evo (PEMF)

CE 0483

B.Light Clear Evo (LLLT)

CE 0483

B.Light Restore Evo (LLLT)

CE 0483

B.Body Evo (PEMF)

CE 0483

B.Bed Evo (PEMF)

CE 0483

B.Pad Evo (PEMF)

CE 0483

B.Sit Evo (PEMF)

CE 0483

B.Body Evo foot protection

CE

B.Bed Evo Attachment Strap

CE

Velcro extension strap for B.Pad Evo

Extension cord for Evo applicators

CE

### Accessories

B.Grip Evo Attachment Strap

B.Light Evo Safety Glasses

CE

B.Box Evo Wall Mount

### 4.3 Scope of delivery for sets, packages, and individual products

Product designation	BEMER Basic-Set Evo	BEMER Classic-Set Evo (B.Body)	BEMER Classic-Set Evo (B.Bed)	BEMER Premium-Set Evo	Beauty Pack Evo	Description
B.Box Evo (REF 424000)	X	X	X	X		Controller, high-resolution touchscreen display, operation of individual applicators (medical product, Class IIa)
B.Body Evo (PEMF) (REF 434300)	X	X		X		Comfortable full-body applicator with connection cable and magnetic connector for connection to the B.Box Evo (medical product, Class IIa)
B.Bed Evo (PEMF) (REF 434400) B.Bed Evo Attachment Strap (REF 454700)			X	X		Flat, breathable full-body applicator including B.Bed Evo Attachment Strap, with connection cable and magnetic connector for connection to the B.Box Evo (medical product, Class IIa)
B.Pad Evo (PEMF) (REF 434100) Velcro extension strap for B.Pad Evo (REF 454400)		X	X	X		Local applicator with connection cable and magnetic connector for connection to the B.Box Evo controller (medical product, Class IIa) Extension belt for the B.Pad Evo with Velcro attachment
B.Spot Evo (PEMF) (REF 434000)		X	X	X		Local applicator module, can be connected to the B.Box Evo controller using the B.Grip Evo holder module (medical product, Class IIa)
B.Light Clear Evo (LLLT) (REF 434500) B.Light Restore Evo (LLLT) (REF 434600) B.Light Evo Safety Glasses (REF 454900)		X	X		X	Light applicator modules can be connected to the B.Box Evo controller using the B.Grip Evo holder module (medical product, Class IIa) B.Light Evo Safety Glasses for use of the light applicator modules in the facial area
B.Grip Evo (REF 454000) B.Grip Evo Attachment Strap (REF 454200)		X	X	X		Holder module: connection cable and magnetic connector (B.Spot Evo, B.Light Clear Evo, B.Light Restore Evo) (medical product, class I) The B.Grip Evo Attachment Strap allows targeted positioning of an application module (B.Spot Evo, B.Light Clear Evo or B.Light Restore Evo) on a specific area of the body
B.Box Evo Power Supply (REF 444000)	X	X	X	X		Power supply for connection of the B.Box Evo to the domestic power supply
B.Box Evo Stand (REF 454800)	X					Stand to position the B.Box Evo on a level surface

### Products sold separately

Product designation	Description
B.Bed Evo (REF 434400)	Flat, breathable full-body applicator including B.Bed Evo Attachment Strap, with connection cable and magnetic connector for connection to the B.Box Evo
B.Bed Evo Attachment Strap (REF 454700)	
B.Sit Evo (REF 434200)	Local applicator with connection cable and magnetic connector for connection to the B.Box Evo
B.Pad Evo (PEMF) (REF 434100)	Local applicator with connection cable and magnetic connector for connection to the B.Box Evo controller
B.Body Evo (PEMF) (REF 434300)	Comfortable full-body applicator with connection cable and magnetic connector for connection to the B.Box Evo
B.Spot Evo (PEMF) (REF 434000)	Local applicator module, can be connected to the B.Box Evo controller using the B.Grip Evo holder module

### Accessories sold separately

Product designation	Description
B.Grip Evo (REF 454000)	Interchangeable adapter with connection cable and magnetic connector to connect the B.Spot Evo, B.Light Clear Evo and B.Light Restore Evo application module to the B.Box Evo
B.Grip Evo Attachment Strap (REF 454200)	The B.Grip Evo Attachment Strap allows targeted positioning of an application module (B.Spot Evo, B.Light Clear Evo or B.Light Restore Evo) on a specific area of the body

Product designation	Description
B.Box Evo Battery Stand (REF 454100)	Portable power supply for connection to the B.Box Evo
B.Box Evo Power Supply (REF 444000)	Power supply for connection of the B.Box Evo to the domestic power supply
B.Box Evo Car Adapter (REF 444100)	Power supply to connect the B.Box Evo to a vehicle battery (10 – 32 VDC)
B.Box Evo Stand (REF 454800)	Allows the B.Box Evo to be positioned on a level surface
B.Body Evo foot protection (REF 450500)	Cover to protect against soiling
Extension cord for Evo applicators (REF 454500)	Extension cable to extend the cable to a maximum length of 5m
B.Bed Evo Attachment Strap (REF 454700)	Strap and tensioning system for installing the B.Bed Evo on a mattress
B.Light Evo Safety Glasses (REF 454900)	Safety glasses for use of the light applicator modules in the facial area
Velcro extension strap for B.Pad Evo (REF 454400)	Extension belt for the B.Pad Evo with Velcro attachment
B.Box Evo Wall Mount (REF 454600)	Mounting bracket that allows the B.Box Evo to be attached to a wall
B.Box Evo Travel Bag (REF 455000)	Travel bag for safe transport of medical devices and accessories

#### 4.4 Main components of the BEMER Therapy System Evo

**1** The LED indicator shows whether the B.Box Evo is switched on (white light) or the B.Box Evo Battery Stand is charging (green light). In the charging state (green light), the device is also ready for operation. If the LED indicator is white, the B.Box Evo battery holder is not being charged (e.g., it is already fully charged) or no B.Box Evo Battery Stand is connected.

**2** On the top of the B.Box Evo there is a push-button switch to turn the device on and off.

- Short press (<2 seconds) puts the device into stand-by mode.
- Long press (>2 seconds) turns the device off completely.

In all states, the applicators are disconnected from the mains to avoid e-smog.

B.Box Evo with all functional elements

B.Box Evo with ambient light ring

**3** The magnetic connection ports (Applicator 1 and Applicator 2) for the application modules are located on the right side. Two applicators can be operated in parallel.

**4** On the front there is a high-resolution control panel which serves as the user interface for the controller. All settings for the B.Box Evo and its applicators are configured via this user interface. Furthermore, there is a light sensor on the front for automatic adjustment of the display brightness.

**5** The removable B.Box Evo Stand is also located here. This can be removed to allow use of the rechargeable B.Box Evo Battery Stand. The rechargeable B.Box Evo Battery Stand consists of a battery integrated into a special stand.

**6** The ambient light ring is located on the back of the B.Box Evo and indicates the device status via different light colors.

**7** The loudspeaker for playback of acoustic signals and the relaxation melody that accompanies the therapy (this can be optionally switched off) is located on the left side of the B.Box Evo.

**8** Below the speaker output is the mains connection, via which the B.Box Evo is supplied with power. When the B.Box Evo Battery Stand is connected, the B.Box Evo Power Supply also serves as a charger.

- 1** LED display
  - Green = B.Box Evo Battery Stand is charging
  - White = device is switched on
- 6** Ambient light ring
  - White = device is ready for operation
  - Blue = active therapy
  - Orange = warning
  - Flashing red = error

This illustration shows the special B.Box Evo Battery Stand that can be used for mains-independent use of the B.Box Evo and its applicators.

No.	Description
<b>1</b>	LED display
<b>2</b>	Main switch
<b>3</b>	Magnetic field connections for applicators
<b>4</b>	Control panel
<b>5</b>	B.Box Evo Stand or Battery Stand

No.	Description
<b>6</b>	Ambient light ring
<b>7</b>	Loudspeaker
<b>8</b>	Power supply unit for B.Box Evo
<b>9</b>	Product label

## 4.5 Applicators/application modules and their use

### 4.5.1 B.Body Evo (PEMF) – full-body applicator for universal use

#### NOTE

Risk of damage to the device due to use of incompatible applicators.

The use of incompatible applicators could result in extensive damage to the BEMER Therapy System Evo, following which its continued (safe) use may no longer be possible.

- Only use compatible applicators supplied by BEMER Int. AG

The B.Body Evo is an applicator for full-body therapy. It can be placed above or underneath the user.


The user's body should be positioned in the center of the applicator. The applicator contains coils that are arranged according to the anatomical shape of the body. The logo is located on the side that should be facing up.

The electromagnetic field is emitted from both sides of the applicator. The B.Body Evo, like the B.Bed Evo, can be used with the sleep program. The B.Body Evo has a microfiber surface that is easy to keep clean.

Inside the B.Body Evo there are a total of 16 coils, which are adapted to the human anatomy and stimulate several regions of the body simultaneously.

Top side of the B.Body Evo full-body applicator

Example application

### 4.5.2 B.Bed Evo (PEMF) – full-body applicator for use in bed

The B.Bed Evo is an applicator for full-body therapy. The material is specially designed for this type of use and can be positioned on the mattress without slipping. The B.Bed Evo is larger than the B.Body Evo and is compatible with all standard mattress sizes.


The B.Bed Evo can be placed underneath the user. The user's body should be positioned in the center of the applicator. The logo shows the end that should be nearest to the user's head. The electromagnetic field is emitted from both sides of the applicator. The B.Bed Evo, like the B.Body Evo, can be used with the sleep program.

The B.Bed Evo features a highly breathable surface for additional comfort while sleeping.

On the underside there are Velcro fasteners, which can be attached to the mattress using the supplied tension straps (B.Bed Evo Attachment Strap).

Top side of the B.Bed Evo full-body applicator

Example application

### 4.5.3 B.Pad Evo (PEMF) – universal applicator for targeted treatment


The B.Pad Evo is a medium-sized, flexible applicator module for targeted treatment of small areas of the body. If the B.Body Evo is not available, e.g., when traveling, the B.Pad Evo can be used as a substitute.

The B.Pad Evo features comfortable and breathable fabric on the inside and a microfiber surface on the outside, which is easy to clean.

The length of the B.Pad Evo can be extended using the supplied accessories (B.Pad Evo Extension Belt).

Top side of the B.Pad Evo

Example applications

#### 4.5.4 B.Spot Evo (PEMF) – applicator module for targeted treatment

The B.Spot Evo is an application module with a very focused treatment area. The B.Spot Evo can only be used in conjunction with the interchangeable adapter (B.Grip Evo).

It is rigid and must be attached to the universal B.Grip Evo holder module, and held against the body part to be treated, either manually or with the supplied attachment strap.

B.Spot Evo

Example applications

#### 4.5.5 B.Sit Evo (PEMF) – applicator for targeted treatment while sitting

The B.Sit Evo is a medium-sized applicator for targeted treatment of a seated patient. The electromagnetic field is emitted from both sides of the applicator. The soft memory foam provides a high level of comfort and the abrasion-resistant surface ensures excellent durability.

B.Sit Evo

Example applications

#### 4.5.6 B.Light Clear Evo/B.Light Restore Evo (LLLT) – for targeted light therapy

The B.Light Clear Evo and B.Light Restore Evo are application modules for the emission of monochromatic, non-coherent light. The B.Light Evo must be attached to the B.Grip Evo universal holder and held against the part of the body to be treated.

B.Light Clear Evo (left)

B.Light Restore Evo (right)


Example applications

#### 4.5.7 B.Grip Evo interchangeable adapter

The B.Grip Evo is an interchangeable adapter for the B.Spot Evo, B.Light Clear Evo, and B.Light Restore Evo application modules. The applicators are magnetically connected to the B.Grip Evo and automatically recognized by the control unit.

B.Grip Evo interchangeable adapter

Example application

## 4.6 Accessories

### 4.6.1 B.Box Evo Battery Stand

The B.Box Evo Battery Stand supplies the B.Box Evo with power and makes it independent of the mains. The B.Box Evo Battery Stand is connected to the B.Box Evo in place of the regular stand. Before use, remove the protective cover using a commercially available Phillips screwdriver, size PH 1. Remember to reattach the protective cover to the device if it will be used or transported without the B.Box Evo Battery Stand.

The B.Box Evo Battery Stand is charged via the control unit and the current battery status is shown on the display.

B.Box Evo Battery Stand

Example applications

### 4.6.2 B.Grip Evo Attachment Strap

The B.Grip Evo Attachment Strap is a double-layer textile strap for easy attachment of the B.Grip Evo holding module to a specific body part. The strap is designed so that it can be attached with just one hand and is flexible enough to adapt to different parts of the body. The strap can be clipped into the B.Grip Evo before use.

B.Grip Evo Attachment Strap

Example application

### 4.6.3 B.Box Evo Power Supply

The B.Box Evo Power Supply is designed for connection to a domestic power supply. It also serves as a charger for the optionally available B.Box Evo Battery Stand.

- The B.Box Evo Power Supply is a special, medically approved power supply unit that may only be used in conjunction with the B.Box Evo. Commercially available power adapters must not be connected to the device.

External medical power supply (2MOPP) for the B.Box Evo

### 4.6.4 B.Box Evo Car Adapter

The B.Box Evo Car Adapter is designed for connection to a 12 V DC power supply. It can also be used as a charger for the optionally available B.Box Evo Battery Stand.

B.Box Evo Car Adapter

- The B.Box Evo Car Adapter is a special, medically approved power supply unit that may only be used in conjunction with the B.Box Evo. Commercially available power adapters must not be connected to the device.

#### 4.6.5 B.Light Evo Safety Glasses

The B.Light Evo Safety Glasses should be worn to protect the eyes during light therapy with the B.Light Clear Evo and B.Light Restore Evo.

B.Light Evo Safety Glasses for use during light therapy

#### 4.6.6 B.Box Evo Wall Mount

The B.Box Evo Wall Mount consists of a metal profile that is protected against corrosion and can be mounted on a wall by means of two universal screws (Ø 4.5 x 35 mm) and two dowels (S6). The B.Box Evo is securely attached to the B.Box Evo Wall Mount via the magnets located inside the housing. Battery operation is not possible when using the B.Box Evo Wall Mount.

B.Box Evo Wall Mount for the B.Box Evo

## 5. Transport and storage

### 5.1 Safety

#### NOTE

##### **Risk of damage to the device due to incorrect transport or storage**

The BEMER therapy system may be damaged as a result of incorrect transport or storage.

- Check the packaging for signs of damage.
- Check that all components of the BEMER Therapy System Evo are undamaged prior to commissioning.
- Always store the BEMER therapy system in a dust-free and dry environment.

##### **Risk of damage to the device due to incorrect storage**

Incorrect storage conditions could result in damage to the therapy system and impair its function.

- Follow the storage instructions described in chapter 3.1.

### 5.2 Symbols on the packaging

See chapter 12 "Meaning of the symbols on the labels (devices and packaging)".

### 5.3 Storage of the packaging

BEMER Int. AG recommends keeping the original packaging. This can be used to return the product within the framework of warranty or repair services.

In addition, the packaging is ideal for storing the individual components of the BEMER Therapy System Evo.

## 6. Recommended use

### 6.1 Recommendations for application of magnetic field therapy (PEMF)

- For full-body therapy, the B.Body Evo and B.Bed Evo applicators are connected to the B.Box Evo controller and used.
- For targeted treatment, the B.Spot Evo applicator module is connected to the B.Box Evo via the B.Grip Evo holder module before use.
- Additional local applicators (B.Pad Evo and B.Sit Evo) can be connected directly to the B.Box Evo and used.

#### 6.1.1 The plus Signal

The plus Signal is a short, cyclic pulse that is modulated every 20 seconds in addition to the base signal. The intensity is amplified up to 150% in five steps within 165 ms. This additional pulse results in more intensive cell excitation during this period. The plus Signal cannot be used with the sleep program.

#### 6.1.2 Full-body therapy

During daytime, the full-body therapy is carried out using the B.Body Evo or B.Bed Evo applicators according to the basic plan and with varying intensities. During sleep, the treatment is carried out using the B.Bed Evo applicator (optional accessory) with the sleep program. Alternatively, the B.Body Evo can also be used.

##### 6.1.2.1 Intensity modes and sleep program (basic plan)


The basic plan structures the generalized full-body therapy as follows: initial application via the B.Body Evo or B.Bed Evo for 2 x 8 minutes daily. The first cycle lasts for six weeks and the second cycle starts from the seventh week of use. From Cycle 2, the sleep program is used. For first-time users with sleep problems, during Cycle 1 we recommend additionally using "LOW" intensity before sleeping. The plus Signal can also be used regularly during the morning treatment.

		Week	Intensity (morning)	Intensity (evening)	Sleep program (nighttime)
Cycle 1	Level	1	low*	low	Number of sleep programs per week
		2	low*	low	
		3	low*	low	
		4	medium*	medium	
		5	medium*	medium	
		6	medium*	medium	
Cycle 2	Level	7	low*	low	1
		8	medium*	medium	2
		9	medium*	medium	3
		10	medium*	medium	4
		11	medium*	medium	5
		12	medium*	medium	6
		13	medium*	medium	7
(*) = treatment with plus Signal					


### 6.1.2.2 Sleep program

The sleep program helps to stimulate vasomotion during sleep.

The sleep program runs from the program start (A) until the wake-up time (B). For the first two hours it runs at an intensity below the LOW level (7  $\mu$ T). This is followed by a rest period; the B.Box Evo switches to stand-by mode (C). During the two hours before the wake-up time, the treatment is again applied at LOW intensity.

Example of sleep program  
 A & B = two hours of treatment each  
 C = rest phase (no treatment)

The B.Bed Evo should be placed under the sheet

Example: The sleep program will start at 9.45 pm. Beginning at this time, the user is treated with magnetic fields for two hours. After two hours, the magnetic field therapy ends and does not resume until two hours before the set wake-up time. During these two hours, the magnetic field therapy resumes and then ends with an acoustic alarm signal (if enabled).

- The sleep program can only be activated if there is a minimum of six and a maximum of ten hours between the program start and end. This is indicated by a green or red marker. The intensity cannot be adjusted for the sleep program.

### 6.1.3 Targeted treatment

Targeted treatments can be applied two or three times a day. Starting with P1, you can switch to the next higher program level every two to three days until P3 is reached. Subsequently, the treatment is continued using P3.

#### User-defined programs

Each user can adapt the program sequences to suit his or her own needs. This is done in the expert mode editor (see chapter 8.2.6).

Program	Total duration (in minutes)	Level (intensity)	Description
P1	8	low	Low intensity for superficial body regions
P2	16	low / medium	Medium intensity for slightly deeper lying tissue
P3	20	medium / high	High intensity for deep lying body tissue

Program structure and parameters (total duration and intensity)

### 6.2 Recommended use of light therapy (LLLT)

B.Light Clear Evo and B.Light Restore Evo can be used daily. The application modules should be placed as close as possible to the area to be treated. The skin must be clean and dry (e.g., without make-up, sunscreen, or creams).

The intensity and application time for the B.Light Clear Evo and B.Light Restore Evo light applicators are explained below.

### 6.2.1 B.Light Clear Evo

An application module for complementary treatment of skin diseases, impaired wound healing, and cosmetic use. The light-emitting diodes produce a red and blue light with a strong photobiological effect. This results in additional formation of energy-rich compounds in the skin and connective tissue cells. Complementary use of this module is recommended for acne.

Please continue treatment until the acne lesions subside and the skin appears healthy and supple. If symptoms persist and / or excessive redness or discomfort occurs, a physician must be consulted before continuing the therapy.

### 6.2.2 B.Light Restore Evo


The B.Light Restore Evo module is an application module for complementary treatment and temporary relief of mild muscle and joint pain, arthritis and muscle spasms, and relief from stiffness. It also promotes the relaxation of muscle tissue and temporarily boosts the local circulatory system. In addition, it helps to reduce the appearance of wrinkles and fine lines.

Please continue treatment until the fine lines and wrinkles appear less pronounced and the firmness and smoothness of the skin has improved. If symptoms persist and / or excessive redness or discomfort occurs, a physician must be consulted before continuing the therapy.

## 7. Commissioning

### 7.1 Connecting the B.Box Evo Stand, Wall Mount, and Battery Stand

Take care to open the packaging of the BEMER Therapy System Evo in the correct place (see the "Upside" symbol on the packaging) and ensure that the contents of the packaging cannot fall out after opening. Having opened the packaging, first remove the B.Box Evo Stand or B.Box Evo Battery Stand and connect it to the B.Box Evo controller. Connect the B.Box Evo Battery Stand (Premium Set) or the B.Box Evo Stand (Basic Set) to the B.Box Evo. To do this, bring the B.Box Evo Battery Stand close to the B.Box Evo as shown – the magnetic holder will then ensure a secure attachment.

Attaching the B.Box Evo Stand or the B.Box Evo Battery Stand

Before commissioning the B.Box Evo Battery Stand, the battery protection cover must be removed from the B.Box Evo. To do this, use a Phillips screwdriver (PH 1) and be sure to replace the protective cover after removal of the B.Box Evo Battery Stand / B.Box Evo Stand. Keep the protective cover in a safe place.

#### Note on battery operation

As long as the B.Box Evo is supplied with voltage, the B.Box Evo Battery Stand will be continuously charged. The current charging status can be seen on the display. In addition, the green LED on the top of the unit lights up.

#### Note on operation with the B.Box Evo Wall Mount

The B.Box Evo Wall Mount is attached to the desired location on the wall using the two included screws (4.5 x 35 mm) and dowels. Ideally, the chosen mounting position should take account of the length of the applicator cables (2.5 meters). If longer distances are required between the B.Box Evo and the application modules, the extension cord for Evo applicators (optional accessory) can be used.

## 7.2 Connecting the power supply

Set up the B.Box Evo so that the cables from the Power Supply (or Car Adapter), applicators, and application modules are not under tension and lie flat on the floor without creating tripping hazards. Do not position the B.Box Evo in a manner that makes it difficult to disconnect the unit from the mains. Organize a suitable place for the treatment.

- Connect the mains adapter to the B.Box Evo Power Supply.
- Connect the B.Box Evo Power Supply to the B.Box Evo.

## 7.3 Connecting the applicators and application modules

The applicators and application modules are connected to the B.Box Evo at ports (1) and (2) as follows:

- Guide the connector (3) into the desired port (1) or (2); the magnetic retainers ensure a secure connection and correct polarization.

## 7.4 Switching on the B.Box Evo

The main switch is located on the top-right of the B.Box Evo and is used to place the device in the operating state. After pressing the main switch, the device performs a self-test within 15 seconds and then displays the start screen. The entire start-up process can take up to 30 seconds.

- Press **<main switch>**.


→ LED indicator lights up white.

(It lights up green as soon as the B.Box Evo Battery Stand is connected and charging.)


→ An audible start signal sounds.

→ The B.Box Evo start screen is loaded.

→ After approx. 30 seconds, the language settings screen appears.

B.Box Evo start screen

B.Box Evo main menu

- Press **<◀>** or **<▶>** to set the user language.
- Select the user language **<English>**.
- Press **<Save>** to save your settings.

→ The B.Box Evo main menu is displayed in the selected language.

- In the operating state, the B.Box Evo control unit can be set to energy saving mode by briefly pressing the main switch. In energy saving mode, the device is not completely turned off and can be “woken up” again by touching the display. A long press turns the device off completely.

## 8. Using the BEMER Therapy System Evo

### 8.1 Switch the unit on (start screen)

- Switch on the device via the main switch.

→ The start screen appears after a few seconds.


The user interface offers three different application modes for magnetic field therapy and one application mode for light therapy.

#### Magnetic field therapy

- 1** Intensity: 8 minutes of therapy with constant magnetic field strength
- 2** Program: 8–20 minutes of therapy with increasing magnetic field strength
- 3** Sleep program: 6–10 hour duration with a total of 4 hours of therapy

#### Light therapy

- 4** Light therapy: 8 minutes of therapy at three different intensity levels


B.Box Evo main menu

### 8.2 Main menu

After you have switched on the B.Box Evo, the main menu is displayed via the user interface. You can choose between the four different therapy variants. Six additional device status icons and the current time are displayed in the upper part of the main menu.

(For a more detailed description see chapter 8.2.1 “Status bar”)

Main menu of the B.Box Evo showing all information and operating elements

No.	Description
1	Port A1 / A2, indication and description of the connected application module
2	● Active therapy
3	Power save mode
4	Volume
5	Battery operation, charging status: icon with lightning bolt – battery is charging
6	Mains operation
7	Current time (in 12 h or 24 h format)

No.	Description
8	Intensity menu
9	Program menu
10	Sleep program menu
11	Light therapy menu
12	Settings menu
13	Switch between ports A1 / A2

Table 1: Description of the information and controls in the B.Box Evo main menu

#### 8.2.1 Status bar

Symbol	Description
● A1 ● A2	No applicator connected to A1 / A2
● A1 B.Body Evo ● A2 B.Pad Evo	B.Body Evo applicator connected to A1 / B.Pad connected to A2
● A1 B.Body Evo ● A2 B.Pad Evo	Active therapy via Port 1


Description of the symbols on the status bar of the B.Box Evo

### 8.2.2 Settings (start screen)

Before starting an application, the basic settings should be configured.


- Press **12** <Settings>.  
 → The settings menu opens.

Settings menu

### 8.2.3 Adjust the signal volume

- Press **1** <Volume and voiceover>.  
 → Settings – Volume and voiceover opens.

Settings: Volume control and voiceover

- Adjust the volume for signals, voiceover and music via **1** <Volume control>.

To increase the accessibility of the BEMER therapy system for visually impaired users,


- Enable the **2** <Voiceover> function.
- Press <Save> to save your changes and return to the settings menu.
- Press <Back> to exit the menu without saving your changes.

- **2** Voiceover is an operating aid (for visually impaired users). If voiceover is enabled, a pre-recorded voice describing the controls is heard.

If the relaxation music is disabled during an ongoing treatment, the voiceover function remains active if it was previously enabled.

### 8.2.4 Adjust the display brightness

- Press **2** <Display and brightness>.  
 → Settings – Display and brightness opens.


Settings – Display and brightness

- Adjust the brightness of the display using the **1** <Brightness control>.
- Enable or disable **2** <Auto Brightness>.
- Press <Save> to save your changes and return to the settings menu.
- Press <Back> to exit the menu without saving your changes.

- If necessary, enable or disable **2** <Auto Brightness>. In automatic mode, the display brightness is adjusted in accordance with the ambient light levels. If the control is grayed out, auto brightness mode is enabled.

### 8.2.5 Set the current time

- Press **3** <Time>.
- Settings – Time opens.

Settings – Time

- Press <◀> to adjust the hours.
- Increase or decrease the displayed hours via <+> and <->.
- Press <▶> to adjust the minutes.
- Increase or decrease the displayed minutes via <+> and <->.
- Press <Save> to save your changes and return to the settings menu.
- Press <Back> to exit the menu without saving your changes.

### 8.2.6 Expert mode editor


Via the expert mode editor, three custom programs (E1 to E3) can be created by the user. The intensity and time parameters are set.

- Press **4** <Expert mode editor>.
- Settings – Expert mode editor opens.

The program editor is used to configure user-defined program sequences and enable/disable the plus Signal.

- Select the expert program to edit: <E1>, <E2> or <E3>.
- Press the corresponding button.
- Expert mode – <Program name> opens.

Expert mode editor – detailed settings


- Press <◀> or <▶>, to select the treatment time.
- Increase or decrease the treatment intensity via <+> and <->.
- Enable or disable the **1** <Plus> signal.
- When the <Plus> signal is enabled, the <Plus> symbol appears in the **2** status bar.
- Press <Save> to save your changes and return to the settings menu.
- Press <Back> to exit the menu without saving your changes.

- The individually predefined expert programs can be selected via the PROGRAMS menu.

### 8.2.7 System

Basic system settings are configured here (ambient lighting, energy saving mode, system language, time format). In addition, general system information can be found here.

- Press **5** <System>.
- Settings – System opens.

Settings menu – System

### 8.2.7.1 Switch the ambient light on/off

The illuminated ambient light ring (see chapter 4.4) can be enabled or disabled here.

- Switch the **1** <Ambient light> on or off via << or >>.
- The value is saved.

### 8.2.7.2 Enable / disable energy saving mode

Energy saving mode can be enabled or disabled here. When energy saving mode is enabled, the B.Box Evo is switched to stand-by mode after two minutes of inactivity.

- Switch **2** <Energy saving mode> on or off via << or >>.
- The value is saved.
- When <Energy saving mode> is enabled, an icon **6** appears in the status bar.

### 8.2.7.3 Select system language

The system language can be selected from 17 available languages.

- Select the **3** <System language> via << or >>.
- The value is saved.

### 8.2.7.4 Set the time format


You can choose between the 12 h or 24 h display format.

- Select the **4** <Time format> via << or >>.
- The value is saved.
- When the 12 h <Time format> is selected, the morning (AM) or afternoon (PM) indicator appears in the status bar **6** next to the current time.

### 8.2.7.5 System information

Here you will find information about the system version, software version, firmware version, kernel version, hardware version as well as the total operating time per applicator port (A1 and A2) and the region code. This information is important when communicating with customer service.


- Press **5** <Info>.
- System information opens.

System information

## 8.3 Main menu

In this main menu, the four different applications or therapies can be selected. The following application modules can be connected per main menu item:


B.Box Evo – controls

- 1** B.Body Evo, B.Bed Evo, B.Sit Evo, B.Spot Evo, B.Pad Evo
- 2** B.Body Evo, B.Bed Evo, B.Sit Evo, B.Spot Evo, B.Pad Evo
- 3** B.Body Evo, B.Bed Evo
- 4** B.Light Clear Evo, B.Light Restore Evo

- Choose the application module.
- Select port A1/A2.
- Connect the applicator or the application module (see chapter 7.3).
- Select the appropriate program (see chapter 6).


### 8.3.1 Intensity mode

After selecting the “Intensity” therapy type, you can choose between three different intensity levels and then start a therapy.

Choose the intensity level for the connected applicator

- Press **<Intensity>**.
  - If necessary, switch the connected port to A1/A2.
  - Select the intensity level. **1 - 3**
- **Intensity** - <Intensity level> opens.
- The time sequence for the therapy is shown.
- The connected application module is shown.

User interface for the -Low- intensity mode with “Plus” and “Music” controls

Before starting the treatment:

- Select **<plus Signal>** on/off.
- Select **<Music>** on/off.
- Press **<Start>** to start the therapy.
- Press **<Stop>** to end the therapy prematurely.


### 8.3.2 Program mode

After selecting the “Program” therapy type, you can choose between three different intensity levels and then start the therapy. The three different programs differ in terms of the treatment duration and intensity.

Select program P1, P2 or P3

- Press **<Programs>**.
  - If necessary, switch the connected port to A1/A2.
  - Select the program level. **1 - 3**
- **Program selection** - <Program name> opens.
- The time sequence for the therapy and the intensity is shown.
- The connected application module is shown.

Overview of the program sequence

- Expert mode enables the selection of user-created therapy programs (for time and intensity, see chapter 8.2.6).


The Expert-mode programs are not visible until they have been created (see chapter 8.2.6).

The plus Signal is activated by default before the start each treatment and must be deactivated if not required. The plus Signal cannot be deactivated during an ongoing treatment.

### 8.3.3 Sleep program

Once the "Sleep program" therapy type has been selected, the alarm time can be set. Do not set the alarm time until you are ready to go to sleep.

- Press **<Sleep>**.
- If necessary, switch the connected port to A1/A2.
- **Sleep** - Alarm time opens.
- Alarm time **1** and treatment duration **2** are displayed.


Example: setting the alarm time

- Press **3** **<|>** to adjust the hours.
- Increase or decrease the displayed hours via **4** **<+>** and **<->**.
- Press **5** **<|>** to adjust the minutes.
- Increase or decrease the displayed minutes via **4** **<+>** and **<->**.
- Select **7** **<Alarm>** on/off.

- A green circle means the sleep program can be activated. The red circle indicates that the set sleep time exceeds the minimum or maximum permitted sleep time.

The therapy time can only be set and activated for between 6 and 10 hours of sleep. From 6 hours, a green circle appears around the application time, indicating that the application can be started.

- Press **6** **<Next>**.
- **Sleep** opens.

Starting and displaying the sleep program


- Select **1** **<Alarm>** on/off.
- Select **2** **<Music>** on/off.
- Press **3** **<Start>** to start the therapy.
- Press **<Stop>** to end the therapy prematurely.

### 8.3.4 Light therapy

Once the "Light" therapy type has been selected, it is possible to choose one of three different light intensity levels and then start the treatment.


- Press <Light>.
- If necessary, switch the connected port to A1/A2.

→ **Light** – Select intensity opens.

Select light intensity

- Select the intensity level. **1** - **3**
- **Intensity** – <Intensity level> opens.
- The time sequence for the therapy and the intensity is shown.
- The connected application module is displayed.

Ongoing light treatment

- Select **1** <Music> on/off.
- Press **2** <Start> to start the therapy.
- Press **2** <Stop> to end the therapy prematurely.

## 9. Cleaning and care

### WARNING

#### Risk of infection due to transmission of disease carriers

Shared use of the application modules can lead to the transmission of diseases.

- Between each use, clean the applicators using the cleaning and disinfecting agent recommended by the manufacturer.

### 9.1 Cleaning

#### NOTE

#### Material damage due to use of unsuitable cleaning agents or cleaning procedures

The use of detergents or cleaning procedures that have not been approved by the manufacturer may result in damage to the BEMER Therapy System Evo.

- Do not use sharp or abrasive objects. These could damage the display of the B.Box Evo or the housings of the devices included with the BEMER Therapy System Evo.
- When cleaning the glass surfaces, only use microfiber cloths that are suitable for this purpose.
- To clean the system components, use only commercially available, non-abrasive, non-aggressive detergents. Follow the manufacturer's instructions.
- Do not use sharp objects or aggressive detergents to clean the application surface of the B.Light Evo applicators. Doing so could damage the surfaces and diminish or impair their effectiveness.
- Replace the light application module if the treatment surface is damaged (scratched, cloudy, etc.).

BEMER Int. AG recommends the following procedure for private users:

- If the same device is used by multiple private users, clean the applicators before each use.
- Clean and disinfect the BEMER therapy system every four weeks if it is only used by one user.

## 9.2 Disinfection

Professional users of the BEMER therapy system must clean and disinfect the applicator surfaces that are in contact with the patient's body after each application.

- Use the product CaviWipes(TM) from the company Metrex (TM) to disinfect the BEMER Therapy System Evo.
- Allow the disinfectant to work for three minutes before the next application.
- Follow the manufacturer's safety and usage instructions.

## 10. Disposal

- This device must not be disposed of as household waste. Every consumer is obliged to hand in all electrical or electronic devices, whether they contain pollutants or not, at a collection point in his or her town/city or at a retailer, to ensure that they are disposed of in an environmentally friendly manner.


Our rechargeable B.Box Evo Battery Stand for the B.Box Evo must not be disposed of as household waste. You can return the B.Box Evo Battery Stand to a municipal collection point. As a manufacturer and distributor of batteries, our obligation to take back used batteries is limited to the B.Box Evo Battery Stand for the B.Box Evo, which we either carry or have carried in our product range.

You can either return the B.Box Evo Battery Stand to us (with sufficient postage) or hand it in directly to our dispatch warehouse at the following address free of charge:


BEMER Int. AG, Austrasse 15, LI-9495 Triesen/Liechtenstein  
EEE Registration Number: M3685

## 11. Error messages and remedies

Error message	Cause	Troubleshooting
Error Code 1	An active applicator was removed during treatment.	Reconnect the applicator and restart the treatment.
Error Code 2	The device temperature is too high.	Check the ambient conditions and allow the device to cool down. If in doubt, please contact customer service.
Error Code 3-4	An application module has reported an error.	Please contact customer service.
Error Code 5-6	Power supply malfunction	Please check possible sources of error, such as the B.Box Evo Power Supply or B.Box Evo Battery Stand (see technical data). If in doubt, please contact customer service.
Error Code 7-8	The controller has reported an error.	Please contact customer service.
Error Code 9-12	An application module has reported an error.	Please contact customer service.
Error Code 101	No applicator is connected.	Please connect an applicator.
Error Code 104	An unknown applicator was detected.	The applicator module is invalid, please contact customer service.
Error Code 105	The charge level of the B.Box Evo Battery Stand is too low.	The B.Box Evo Battery Stand requires charging, please connect the B.Box Evo Power Supply to charge it.
If a warning message or error message is displayed, it can be closed via <b>&lt;OK&gt;</b> .		

Display in case of a warning (e.g., applicator not connected to the controller)

Display in case of malfunction (e.g., device overheated)

If you have any questions, please contact BEMER customer service

## 12. Meaning of the symbols on the labels (devices and packaging)

Symbol	Meaning	Location
	Active implants. Do not use on people with active implants (e.g., pacemakers)	Device and packaging
	Note	Device
	Fragile, handle with care	Device and packaging
	Temperature limit	Packaging
	Follow instructions	Device and packaging
	Store in a dry place	Packaging
	Applied part, type BF	Device
	Humidity	Packaging
	Serial number	Packaging
	Article number	Device and packaging
	Batch number	Packaging
	Air pressure limit to which the medical product can be safely exposed.	Packaging

Symbol	Meaning	Location
	Manufacturer	Device and packaging
	Date of manufacture	Packaging
	CE identification number of the notified body (e.g., MDC)	Device and packaging
	Protection class II	Device
	SGS NORTH AMERICA CERTIFICATION MARK	Device
	WEEE symbol for electrical and electronic equipment	Device and packaging
	Medical device = medical product. Indicates that it is a medical product.	Device
	Direct current	Device
	Recycling article number	Packaging
	Protection against diagonally falling dripping water up to an angle of 15° (IEC 60601-1-11 Cl. 7.4.4.)	Device
	Protection against vertically falling dripping water	Device

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BEMER Int. AG mengesyorkan prosedur berikut untuk pengguna peribadi:
BEMER Int. AG, Austraße 15, LI-9495 Triesen/Liechtenstein
Errors and omissions excepted // © BEMER Int.-AG // LI-9495 Triesen