DEUTSCH

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ENGLISH

Read the instructions for use carefully and keep them for later use. Make them accessible to other users and note the information they contain.

INCLUDED IN DELIVERY

Check that the exterior of the cardboard delivery packaging is intact and make sure that all contents are present. Before use, ensure that there is no visible damage to the device or accessories and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Service address.

Thermometer, instruction for use, 1 x 1,5V LR 41 battery

SIGNS AND SYMBOLS

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and accessories:

⚠️ WARNUNG

Indicates a potentially impending danger. If it is not avoided, there is a risk of death or serious injury.

⚠️ CAUTION

Indicates a potentially imminent hazard. If it is not avoided, slight or minor injuries may result.

⚠️ NOTICE

Indicates a potentially harmful situation. If it is not avoided, the device or something in its vicinity may be damaged.

⚠️ Important
Safety note indicating possible damage to the device/accessory

⚠️ Observe the instructions
Read the instructions before starting work and/or operating devices or machines

⚠️ The electronic device must not be disposed with household waste

Batteries must not be disposed of with household waste

⚠️ CE labelling
This product satisfies the requirements of the applicable European and national directives.

Marking to identify the packaging material.

A = Material code, B = Material number.

1-7 = Plastics, 20-22 = Paper and cardboard

⚠️ Separate the product and packaging elements and dispose of them in accordance with local regulations.

⚠️ Protected against solid foreign objects 12.5 mm in diameter and larger, and against the effects of temporary immersion in water

⚠️ Batch code

⚠️ Type BF applied part

⚠️ Temperature limit

⚠️ Importeur

Swiss authorised representative

WARNINGS AND SAFETY NOTES

- The thermometer is intended to measure the temperature of the human body.
- The thermometer is only designed for the measuring area on the human body stated in the instructions for use.
- Children should not use the thermometer without supervision.
- Before each use, check the thermometer for signs of damage or wear. Damaged or worn thermometers should not be used.
- Always wait for the acoustic signal which indicates the minimum time required for measurement. Discuss the temperature with your family physician.
- The thermometer contains sensitive electronic parts. Therefore, protect it from blows, bending, high temperature or direct exposure to sunlight.
- Using the thermometer in strong electromagnetic fields such as next to a mobile phone, can cause a malfunction.
- When the thermometer is switched on, it carries out a self-test. Therefore the accuracy of measurement need not be tested at regular intervals.
- If the temperature is very abnormal, please contact your family doctor promptly
- This device meets the requirements of the EU Directive 93/42/EEC concerning medical devices, as well as those of the Medizinproduktegesetz (German Medical Devices Act), the European Standard EN 12470-3: Clinical thermometers – part 3: Performance of compact electrical thermometers (non-predictive and predictive) with maximum device and the European Standard EN60601-1-2 (In accordance with CISPR 11, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-8), and is subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and mobile HF communication systems may interfere with this unit.

NOTES ON ELECTROMAGNETIC COMPATIBILITY

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic fields. This could result in issues such as error messages or the failure of the device/device.
- Avoid using the device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other device must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Failure to comply with the above can impair the performance of the device.

MEASURING METHODS:

Temperature measurement in the anus (rectal)

This method of measurement is most reliable and most accurate. Guide the tip of the thermometer carefully 2 to 3 cm into the anus.

Measurement in the mouth cavity (oral)

For this purpose, guide the tip of the thermometer carefully into one of the heat of the mouth cavity. The thermometer should be held in the mouth to the left or right of the tip of the tongue.

Measurement of temperature in the armpit (axillary)

This method of measurement is relatively inaccurate. Therefore, it cannot be recommended from a medical point of view.

Mode of measurement	Recommended time required for measurement	Mean divergence to the rectal method
in the anus	until the acoustic signal	–
in the oral cavity	until the acoustic signal	-0,4°C to -1,5°C
in the armpit	minimum 5 minutes	-0,7°C to -2,0°C

USE

First remove the protective foil from the display. For switching on, briefly press the on/off key. A short acoustic signal will announce that the thermometer is switched on. The thermometer will carry out a self-test for about 2 seconds. During the selftest all segments of the display will appear. Then the reference reading (37 °C) appears. Then the measurement begins. The display shows the temperature. To prolong the life of the battery, please switch off the thermometer by briefly pressing the on/off key after the measurement. Otherwise the thermometer will automatically switch off after 10 minutes.

MEMORY

The last measured value will be stored automatically. Press the ON/OFF button to display this value. The thermometer first carries out a self test. Next the last

stored measured value is displayed for approx. 2 sec and then the reference value of 37 °C and „Lo“ is shown on the display. The thermometer is ready to take a measurement and the previously stored measurement value is automatically replaced by the new measurement value.

CHANGING THE BATTERY

The clinical thermometer contains a long-lasting alkali-manganese battery (type LR 41, 1.5 V) for approximately 2 years. The battery needs to be replaced when all digital reading appear flashing on the display screen. Release the battery compartment cover by gently pressing on the catch. Remove the old battery and replace it with a battery of the same type. The + sign on the battery must point upwards. Then close the battery compartment, making sure that the catch clicks into place.

The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below apply to batteries containing harmful substances:

Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury. These batteries contain lead (Pb) and mercury (Hg).

NOTES ON HANDLING BATTERIES

⚠️ WARNING

- Risk of explosion! Risk of fire!** Failure to comply with the following points can result in personal injury or cause overheating, leakage, venting, breakage, explosion, or fire with respect to the battery.
- This device contains non-rechargeable batteries which must not be charged.
- Do not throw batteries into a fire.
- Never charge, forcibly discharge, heat, disassemble, open, crush, deform, encapsulate, or modify batteries.
- Never short-circuit batteries or battery compartment contacts.
- Protect the batteries from direct sunlight, rain, heat, and water.
- Exposure of batteries to an environment with extremely high temperatures or an extremely low air pressure may result in explosion or leakage of flammable liquids and gases.
- Dispose of defective and discharged batteries immediately and properly (see chapter on disposal).
- Do not use modified or damaged batteries.
- Always select the correct battery type.
- Always insert the batteries correctly, taking into account the polarity (+ / -).
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- If fluid from a battery comes into contact with your skin or eyes, wash the affected areas with water and seek medical assistance.
- Choking hazard!** Keep batteries out of the reach of children. Seek medical attention immediately if swallowed. Swallowing them may cause burns, severe internal injuries, and death.
- Never allow children to replace batteries without adult supervision.

⚠️ CAUTION

- Store batteries in a well-ventilated, dry, and cool place in a non-conductive container in which the batteries cannot be short-circuited to each other or by other metal objects.
- Keep batteries clean and dry.
- Keep batteries away from water.
- If the device is not going to be used for a long period of time, remove the batteries from the battery compartment.

NOTICE

- Do not use rechargeable batteries.

CLEANING/DISINFECTION

The device must be cleaned and disinfected before and after each use. The thermometer is 100% water-proof. For the purpose of cleaning or disinfection it can be dipped into water or disinfected without any difficulties. Under no circumstances should the thermometer be boiled! Use one of the following disinfectants in accordance with the manufacturer's instructions. Other cleaning agents or methods could impair the function of the thermometer or damage it.

Name of the agent: Isopropanol 70% Pharmacy, Gigasept FF, Lysosformin, Aseptisol.

STORAGE/DISPOSAL

When you are not using the thermometer, please retain it in the original protective cover.

The accuracy of this thermometer has been carefully checked and developed with regard to a long useful life.

Using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.

For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

Type	FT 09/1
Range of measurement	32°C to 42,9°C
Accuracy of measurement	±0,1°C in a water bath with a temperature of 35,5 °C to 42,0°C ±0,2°C in a water bath with a temperature of 32,0°C to 35,4°C ±0,2°C in a water bath with a temperature of 42,1°C to 42,9°C
Ambient temperature	+10°C – +40°C, with 30% to 85% relative humidity
Storage temperature	-10°C – +60°C, with 25% to 90% relative humidity

The serial number is located on the device or in the battery compartment.

WARRANTY/SERVICE

Further information on the guarantee and guarantee conditions can be found in the guarantee leaflet supplied.

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