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Model	Item No.	Description
STORP1	30038555	Not refillable, plastic shaft with Platinum Disk
STORP2	30038553	Refillable, glass shaft with Platinum Band
3M KCl	30059256	Protection solution

Specifications:

Connection: BNC
Cable length: 1.0m
Shaft length: 120mm
Shaft diameter: 12mm
Shaft material:

Plastic (STORP1)

Glass (STORP2)

Temperature range:

0-80°C (STORP1)

0-100°C (STORP2)

Reference Electrolyte:

Gel 3M KCl (STORP1)

Solution 3M KCl (STORP2)

Zero potential value:

86mV± 15mV

Grade difference: ≥ 165mV

The ORP electrode is also called Redox electrode, measures the **Oxidation-Reduction Potential**, related to the overall availability of electrons in the media, specifically the ratio of positive and negative ions in the solution.

Normal application areas are in water quality testing for swimming pools, fish farming and aquariums etc.

ORP is expressed in millivolts (mV). A range of -1000 mV to 1000mV is common with most ORP test.

The pH value influences the ORP value significantly.

Performing ORP measurement:

1. After removing the electrode from the storage solution bottles, rinse the ORP electrode with pure water before measurement; dip the ORP electrode into the sample.
2. Connect the ORP electrode to the pH meter, set the meter at mV measurement mode.
3. If the ORP electrode is refillable (STORP2), slide down the red cover to expose the filling hole when measurement, remember to close it after measurement to avoid leakage.
4. After measurement, rinse the ORP electrode with pure water and store back in the small protection bottle with soaking solution.

ORP electrode quality check:

To check if the ORP electrode quality is good:

1. Prepare 2 standard buffers solutions for ORP: Saturated quinhydrone pH4.01 and pH7.00.
2. Dip the ORP electrode into buffer 7.00, the mV value should be **86mV ± 15mV**.
3. After rinsing in pure water, dip the ORP electrode into buffer 4.01, the difference between 7.00 and 4.01 should be **≥ 165mV**.
4. If the check result is not good, should clean the ORP electrode then check again; or purchase a new one.

Storage and maintenance:

The ORP electrode should be stored in protection solution—3M KCl solution.

When reference electrolyte is lower the half full, refill 3M KCl solution through the filling hole with a syringe.

If the platinum band or disk is contaminated, this will result in slow response time or inaccurate measurement result, should clean the electrode as follows:

1. Inorganic contamination: dip the ORP electrode into 0.1M HCl for more than 10 minutes, then rinse with pure water. Soak the ORP electrode in 3M KCl solution for more than 8 hours.
2. Organic contamination such as grease coating: clean with detergent then rinse with pure water. Soak the ORP electrode in 3M KCl solution for more than 8 hours.
3. If the platinum is contaminated heavily, use toothpaste to polish it, rinse with pure water, then soak in 3M KCl solution for more than 8 hours.

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Ohaus Corporation, all rights reserved/ todos