Carbon®

FP3D

Appendix A

Protocol

HARDWARE

Printer & Cassette: Carbon M series
Print Dropdown: FP3D
Dispenser: Cox Dual-Drive Dispenser

STORAGE


Temperature: 15 - 30° C (63 – 77° F)
Packaging: 450 mL cartridge (10:1 cartridge)
Shelf Life: 12 months

Notes: Three shades are available:
• Original Pink
• Light Reddish Pink
• Dark Veined

RESIN DISPENSING & PRINTING

A/B Ratio: By mass – 10:0:1
Speed: 1
Mixing Tip: Orange
Burn off: 10 mL new cartridge
Pot life: 8 hours
Release film: Optional

Dispensing Process:
• Check the cassette.
• Dispense required resin amount into cassette.
• After printing, reuse within 8 hour pot life or discard resin.

WASHING

Solvent: IPA
Exposure: 5 min max
Smart Part Washer: No

Notes:
• Use full blade contact with the part during platform removal.
• Handle green state parts delicately.

CURING

Cure start: Within 8 hrs of printing
Salt bake: Required

Thermal

Program 1 Pattern 1
• Ramp to 120° C for 15 minutes
• Hold at 120° C for 10 hours

UV: UV cure cycle is required after baking.

Dreve PCU LED

18 MIN at 90% intensity, under vacuum.
Program Name: Dhealth 1
• Available in Configuration File 1

Denstply Sirona Digital Cure Unit (Large Capacity)

30 minute preheat at 80° C .
1 hour cure at 80° C .

Washing Process

- Washing with orbital shaker
- Remove supports.
 - Wash 2 minutes in IPA on orbital shaker at 140 RPM.
 - Brush intaglio surface and tooth pockets with tooth-brush and IPA.
 - Wash 2 minutes in IPA on orbital shaker at 140 RPM.
 - Remove from IPA and inspect. Parts should be free of uncured resin in tooth pockets and negative features.
 - Set parts aside to dry.

Curing Process

- Fill 3 inch deep baking pan with fresh salt.
- Place parts intaglio side up in salt.
- Add salt until parts are covered and shake the tray to settle the salt.
- Cover tray with aluminum foil.
- Place pan in oven.
- Start oven.
- Wait for parts to cool before handling.
- Rinse parts with warm water to remove any residual salt.
- Set parts aside to dry.
- UV cure parts according to the instructions appropriate for your tooth

Details

COMPATIBLE HARDWARE

Printing and Post-Processing

Printing:
Carbon M-series

Resin Dispensing:
Cox Dual Drive Dispenser
Orange 10:1 Mix Tips

UV Curing:
Dreve PCU LED |
Denstply Sirona Digital Cure Unit (Large Capacity)

Thermal Oven:
Yamato DKN602 convection oven

MIXING & PRINTING

Dispensing

For new cartridges, remember to dispose of the first 10 mL dispensed. The initial amount is not properly mixed and may contain air bubbles.

Printing

No special instructions.

DIGITAL PREP

Supports

Orient the parts vertically with the front teeth area facing away from the platform.
Select off platform for supporting
Use 40deg overhang detection
Manually draw fence supports along the edges If the fence terminates in the part, redraw without that area.
Anywhere you had to skip the fence add in manual bar supports, this will typically be on clasps and overhanging tooth pockets. Please take care to not under support these areas for best accuracy results.

THERMAL & UV CURING

Salt Bake

Prepare parts for baking with salt:
• Partially fill a ~3 inch deep baking pan with salt.
• Place parts intaglio side up on the salt. Ensure parts do not touch.
• Completely cover the the parts in another layer of salt by using a strainer to gently deposit salt onto the parts. Do not push down or manually pack the salt.
• Cover the tray with aluminum foil.

Washing Parts after Salt Bake

Residual salt will adhere to the parts during baking. Remove the salt from the parts by washing.
• Wait for parts to cool before handling.
• Rinse parts with warm water to remove any residual salt.If necessary, use a toothbrush.
• Set parts aside to dry completely before UV curing.

Carbon®

FP3D

Manufactured By Keystone Industries

Product Description: FP3D is designed for additive manufacturing on Carbon 3D printers (wavelength = 385 nm). FP3D is intended to be used within a computer-aided design and manufacturing (CAD/CAM) digital dentistry system that includes a 3D scanner, design software, 3D printer, post-cure unit (wavelength = 385-405 nm), and post-cure oven for the fabrication of flexible partial dentures for partially edentulous patients. For any components that are used in conjunction with FP3D, the user should review all applicable product labeling including Instructions for use, user manuals, and other associated labeling. Strict adherence to all labeling is critical in ensuring a safe and effective printed appliance.

Indications for Use: FP3D is a dual cure resin indicated for the 3D printed fabrication of flexible partial dentures in dental laboratories.

Patient Population and Medical Treatment

FP3D resins are used to fabricate 3D printed removable denture prosthetics for partially edentulous patients to replace missing teeth. They are intended for use in the general population including adults and children.

Contraindications: Contains acrylate monomers and oligomers and amines which, although rare, may cause an allergic reaction or irritation to individuals sensitive to acrylic and amine-containing products.

CAUTION: US Federal law restricts this device to sale by, or on the order of a dental professional.

Adverse Reactions: Allergic skin reactions may occur in susceptible individuals. The resin may cause respiratory irritation.

WARNINGS & PRECAUTIONS

- Follow all compatible printer settings and curing times for safe and effective results.
- For use only with Carbon printers.
- Review to product Safety Data Sheets (SDS) for both part A and part B components prior to use.
- As per the SDS, wear proper personal protective equipment when handling the FP3D and uncured printed parts. Be careful not to splash when dispensing the resin.
- Store in a cool, dry place 15°C-30°C (59°F-86°F) and away from light. Ensure that the cartridge is sealed while not in use.
- In the unlikely event of a print failure during printing, filter the liquid resin through a mesh screen with pore sizing <200 microns. It is good practice to filter the resin in a cassette periodically to prevent print failures.
- Clean the printer platform and cassette prior to using a different shade of resin or once the resin pot life has been reached. Do not mix different shades of the same product. Leftover resin should be disposed after its 8-hour pot life has been exceeded.

Compatible Equipment: To ensure the biocompatibility of the final device, FP3D is designed for fabrication on the Carbon M-series printers and with the recommended post-cure units and thermal oven. Refer to Appendix A for the full list of compatible resin dispensing equipment, Carbon printers, and post processing equipment.

DIRECTIONS FOR USE

Dispensing FP3D Resin:

- Ensure the resin is tempered to ambient temperature (20-25°C/68-77°F) prior to printing.
- Attach the 10:1 orange mix tip securely to the end of the 10:1 resin cartridge.
- When dispensing from a cartridge for the first time: Load the cartridge into the Cox Dual-Drive Dispenser and discard the first 10 mL of resin at a speed of 1 prior to the first resin dispense from the new cartridge.
- Dispense resin directly into the printing cassette. Do not exceed a dispensing speed of 1 to ensure proper mixing of FP3D Part A and Part B.
- Ensure usage of new mixing tip for each dispense required.

Printing FP3D Dentures:

- Only use Carbon M-series printers, settings, and indicated post-processing equipment to ensure the printed product is safe and accurate for use. Refer to Appendix A.
- Once your design is completed per the CAD software manufacturer's directions for use, import the design onto the Carbon printer.
- Follow Carbon's printing instructions in Appendix A.
- Clean and post-process the printed part per the instructions below.

CLEANING INSTRUCTIONS FOR PRINTED PARTS

Printed parts should be cleaned with isopropyl alcohol (IPA) (>97% purity) within approximately eight (8) hours from completion of the print using the following steps:

- Remove supports prior to post processing.
- Using agitation, wash the printed parts in a clean IPA bath for two (2) minutes.
- Using an IPA wetted soft brush, gently brush the parts, with emphasis on tooth pockets and recessed features.
- Remove and transfer the parts to a fresh IPA bath. Using agitation, wash the printed parts again for two (2) minutes.
- Remove printed parts from the bath and allow them to dry completely (air dry or blow dry using compressed air).
- Ensure all residual resin is removed from the parts, as residual resin may trap salt during the thermal baking step (see instructions below). If resin is still present, use an IPA wetted soft brush to clean the isolated area.

Note: Printed parts should not exceed five (5) minutes of IPA exposure total, brushing included. If parts experience solvent overexposure there is risk to degradation of mechanical properties and performance.

Note: The use of denatured alcohol or ethanol for cleaning is discouraged, as they may diminish or degrade the quality of the finished parts.

DIRECTIONS FOR POST-PROCESSING

Refer to Appendix A for all post processing-compatible materials and equipment. Thermal cure must be started within eight (8) hours of printing. Ensure parts are fully cleaned and dry prior to salt baking.

- Sift 2 inch layer (approximately) of powder grain salt into a baking tray or dish using a fine mesh strainer. Place printed parts onto the salt with the intaglio surface facing up. Avoid pressing parts into salt with excess force as to not cause deformation.
- Once placed, sift an additional layer of salt on top of the parts, ensuring the printed parts are completely covered. Gently shake the baking tray to ensure the parts are settled into salt.
- Cover the baking tray with aluminum foil and bake at 120°C(248°F) for 10 hours, with a 15 minute ramp to 120°C(248°F).
- Once the oven cycle is complete, allow parts to completely cool prior to handling.
- Remove parts from salt and clean parts with a soft brush and warm water to remove any excess salt from the surface. Other dental tools, such as a water pick can help reach areas that are difficult to access. We recommend paying particular attention to engraved features, clasps, and the intaglio surface.

Follow the preferred printed denture tooth IFU for final UV cure and device assembly instructions. If the tooth bonding IFU does not include a final UV cure, refer to Appendix A for final UV curing instructions in a compatible post-cure unit.

After the above steps are followed, the medical device is cured and safe with respect to residual monomers.

Note: Replace salt after 10 baking cycles to ensure part cleanliness.

Patient Cleaning Instructions: This medical device is a single-patient, customized, multi-use oral appliance that should be cleaned between usages. The patient should clean the appliance with soap and warm water, or any over-the-counter mild cleaning agents indicated for oral devices. Dry partial dentures thoroughly overnight if not wearing and do not leave submerged in the cleaner.

Repairing and Relining: If repairing or relining is required, a new flexible partial framework should be printed following the instructions provided above.

Disposal Considerations: Unused resin materials may be considered hazardous waste and should be disposed of in accordance with all local, regional, national, and international regulations.

Keystone Industries
52 West King St.
Myerstown, PA 17067 USA
800.333.3131
keystoneindustries.com

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