


# Data Handling: End Point

- Any Mass Testing methodology must match throughput to data handling/scoring
- KRAKEN has been designed as the only dedicated Mass Testing LIMS and has processed Billions of data points on over 350 genomes, on millions of amplicon assays.
- End point allows simplified data scoring reporting as a YES/NO
- Million ++ data points per day is all in house, 1 person per day



# KRAKEN LIMS

Laboratory Information Management System  
Every aspect of a Mass testing lab "Nose to tail"


SNP.Design	#	Project	Status
053-0484.2	5	074.105	Opened
053-0484.2	4	074.105	Opened
053-0484.2	3	074.105	Genotyped
053-0484.2	2	074.105	Pass

Liquid dispensing



Data reporting

All data flow

# Data calling: End Point



Data is automatically or manually scored with definable levels of user security



# Manual Auto calling data

## KRAKEN Manual calling data

Scoring can also be manual: draw round the calls to select them and clicking on the appropriate colour in the legend on the left side of the plot. Scoring is usually automatic but manual editing of calls is also possible. Draw around the wells to be changed by holding down the left mouse button while moving the mouse. The wells on the plate or cluster / ratio view can be selected this way.

It is possible to create an empty result set so that scoring can be carried out on the plate with no read data. [More info](#)

Release the mouse and the wells will be selected. To select more wells without deselecting the current wells, hold down the control key at the same time. Wells may also be selected by clicking on the plate display. NTC wells may not be selected from the plot. If you wish to select NTC wells then you may do so from the plate display. Calls that have a normalisation level below the **Min normalisation threshold** are displayed as a shorted **S** and will not be selected.

Multiple selections are made by holding down the control key while dragging the mouse.

The unselected wells appear flat while the selected wells show up in 3D. The selected wells on the plate show up with a blue background. When well locations are overlaid the 3d appearance intensifies.

Click on the legend to select the new call class. Alternatively on the keyboard press a number: 1=Red, 2=Green, 3=Blue, 4=Black, 5=Pink.

The calls will now be updated. They will be stored in Kraken when the window is closed or another read is selected. Until the window is closed you may Undo your changes. You may decide to [Submit this Assay for export](#) by clicking on the Submit icon. Reports are exported only from Submitted results.

## KRAKEN Auto calling data

Cluster calling

Display shows the scope of the Vectors in white and Excluding magnitude in the back curve. The Nearest Angle is not shown.

Click on Reader group in the diagram above, for more information

**Nearest Angle**  
This control is designed to link points that fall off in lines away from the bottom left corner. Any points that are within the area defined by a segment starting from the bottom left corner or the average control location\* to a point  $\alpha$ , the angle are linked subject to the exclusion parameters.

**Excluding vector**  
This control limits the scope of the Nearest neighbour angle to prevent picking up of out layers. The value must always be greater than the Nearest neighbour vector otherwise it will override the vector and angle.

**Excluding magnitude**

KRAKEN has been and is still used to call Millions of data points, capable of Auto calling 1.5 million or 10,000,000++ Data points per day. Its with out doubt the most significant part of any mass testing program.

# KRAKEN LIMBS

## Laboratory Information Management System

Troubleshooting  
FAQ

Kraken Overview  
Kraken Quick Guides >>  
Starting Kraken

Project types & set up >>  
Working Project management  
Import functions >>  
Design & order primers >>  
Add Samples to working plates >>  
Make & add (Assays) >>  
PCR mix

Set Incubator status  
Read plates (Import)  
Viewing & Scoring results >>  
Sequencing >>  
Reporting results >>

Global Project management  
Event log  
Searching >>  
Definitions  
Lab book  
Managing users

Advanced >>  
Installing server & clients  
Database tables

Troubleshooting  
FAQ

### Kraken

Multiple users log on to password protected Kraken. Each User has their own privileges. User names are recorded in logs, as users interact with Kraken.

Kraken has a single database and a shared fileserver area, allowing users to share data in real time.

*Kraken data is subdivided by Customer. All entities (other than Kraken users) belong to a Customer or are first created for a Customer.*

### Customer

Customers are added to Kraken. A Customer holds source information for Assays, oligonucleotide Primers, DNA Sample tracking, default QC criteria, organism, notes and contact details. A Customer may also be subdivided by Organisation.

Each Customer has one or more Projects associated with it.

### Customer

Customers are added to Kraken. A Customer holds source information for Assays, oligonucleotide Primers, DNA Sample tracking, default QC criteria, organism, notes and contact details. A Customer may also be subdivided by Organisation.

Each Customer has one or more Projects associated with it.

### Project

*{Customers} contain Projects*

Every Project is assigned to a Customer. Each Project has Master plates containing Samples provided for the Customer and Assays used on working Sample plates. It links to Assays and Primers from any Project as needed.

Projects are displayed as a navigation tree, the Project tree. From the Project tree, Assays from this or other Customers may be viewed and imported into the Project (to assist in locating the Assay Aliquot physical tubes). The Project tree contains its own working Sample plates, which may be dispensed from dilutions of the Master DNA (sisters) and daughter plates which are stamped from sisters. Assay results are linked to the Project so that experimental set up and data analysis is performed from the Project tree.

*Each Project is a self-contained information set defining everything from incoming Master sample information, to the results and the final report.*

Items stored with each Project are: Master plates containing DNA with their storage location. To each Project, Master plates may be, replicated, subsets copied and Assays added virtually, modeling the experimental workflow in the laboratory. Location tracking information is stored in real time so that you can see where a Sample plate is at any given moment. Files on the computer system, which are destinations for Reports, are made available in the Project tree.

Any item contained on a plate or in a multi-well box may be reviewed on a picture of that plate. Each spot in a picture is fully interactive, so that you can drill for information about a single well (shown in the image, right).

*Many items are stored with a Project but only a subset of the scored plate read data is required for the final report and stored as Results.*

Notes, Files, DNA samples,  
Working Sample copies, Raw  
data

### Result

Results are automatically collated from plate readers into the correct Project based on the barcode (or name) assigned to each Sample plate. This feature allows multiple Projects to be bundled in the same experimental batch, because each set of results is immediately associated with its own Project/DNA/Assay.

After data has been read it may be automatically called, and made available for display within its Project, with its identifiers for Assays and DNA Samples already assigned.


When reviewing data, spots on a cluster plot interact with a plate diagram and normalisation data. Calls may be automatic, but are also assigned manually, using these interfaces. Each spot may be drilled for information (shown in the image, right).

*Automatic reports generated from Result set*


### Report

Reports may be generated for general use. Results may be exported for use in the free SNPviewer tool. The SNPviewer can then export subsets of the data to other systems.

Reports and exported data are sent to Files available both from Kraken's Project tree and from a general computer navigation system (such as an Explorer window).



- SNPs
- 0001 GHGln90
- 0002 GHArg51Gln
- 0003 GHLeu72Met
- 0004 BDNFCC270T



# KRAKEN LIMS

## Laboratory Information Management System

Not Secure — results.lgcgenomics.com

Saddle ranch strava TR3 club Mille Miglia TR3 Sky Osprey popmax Reserchgate Kindle Kinetico Jazz Albums Kanuku Mountain tutor2u Netflix Sebring Messenger Hotmail LinkedIn BBC Smugmug facebook >>

https://testingmethods.crowdcity.com/post/3296740 Pathology Alliance webinar 12/05/2020 - YouTube **Kraken Manual** JC review.pdf - OneDrive +

### Kraken Topics

[How to use the Manual](#)  
[Website map of Manual](#)  
[New in Kraken](#)

---

[Kraken Overview](#)  
[Kraken Quick Guides >>](#)  
[Starting Kraken](#)

---

[Project types & set up >>](#)  
[Working Project management](#)  
[Import functions >>](#)  
[Design & order primers >>](#)  
[Add Samples to working plates >>](#)  
[Make & add {Assays} >>](#)  
[PCR mix](#)

---

[Set Incubator status](#)  
[Read plates \(Import\)](#)  
[Viewing & Scoring results >>](#)  
[Sequencing >>](#)  
[Reporting results >>](#)

---

[Global Project management](#)  
[Event log](#)  
[Searching >>](#)  
[Definitions](#)  
[Lab book](#)  
[Managing users](#)

---

[Advanced >>](#)  
[Installing server & clients](#)  
[Database tables](#)

The flowchart below demonstrates how a normal user would interact with the Kraken software to log in, add a Customer, assign a Project to that Customer. Within the Project, Assays are added, by adding a DNA sequence; oligonucleotides may be designed using the **Primer picker** tool, and ordered through Kraken, which streamlines production of oligos and record keeping.

Kraken reinforces workflow. When oligos arrive, they are imported into Kraken via the Oligonucleotide synthesis report. Oligos are recorded in Primer plates in Kraken (see the Primer icon in the Project tree. They are automatically associated with Assays. A Primer plate represents a physical plate or box of unmixed oligonucleotides.

Assay Aliquot plates are created to receive oligonucleotide mixes ready for PCR. Assay aliquot tube positions are recorded and start of with a status of **Empty tube**. As oligonucleotides are mixed into new Assay Aliquot tubes they are arranged in these tubes in the Assay plates (or Assay Aliquot plates) in Kraken (a plate represents a physical box). As they are dispensed into Assay mixes, the Assay tube status is changed, and as the Assays are tested (validated) then the Assay status is changed by the person assessing or scoring the results. An Assay Aliquot may be given the status **In use** or **Do not use** and others.

Similarly, {Samples} are received in Master plates and their location recorded by plate and physical storage location (e.g. walk-in freezer) in Kraken. As the {Samples} are diluted and dispensed into working Sample plates (daughter or sister or replicate plates) each dispense is monitored in Kraken. There are icons in the Project tree for the original Master plates are located and another for the working Sample plates. Each Project must have its own Master and working Sample plates, and Kraken makes it easy to transfer {Samples} between plates.

Into working Sample plates (created blank) must be dispensed, virtually in Kraken:

- first, the DNA from the Master or diluted working Sample plates (sisters) created from them.
- second: Assay mixes are dispensed from Assay aliquot tubes.
- third: plate reader files may be assigned to each plate - these are organised under the Assays icon in the Project tree.

• *Note:* This order is strictly controlled in Kraken: an Assay cannot be dispensed into a plate with no DNA samples in it and a plate read file cannot be assigned to a plate where an Assay has not been dispensed.

The Assays icon in the Project tree holds all of the results, and there you will find for each Assay, the Assay designs, and for each design, the read or reads (plate reader file data) which follow thermal or hydrocycling. From this part of the tree each set of results is rapidly accessed as a cluster or other plot. Normalisation data may be viewed.

At this point the results are scored, and the preferred, most representative plot 'Submitted' for the final report. **Submitted** is one of many statuses which may be assigned to a plate read. **\*\*\*However, more than one results set may be used to decide on the call for an individual sample. Plates are frequently re-cycled (submitted to a further 5 rounds of PCR hydrocycling) and the results compared: direction of travel of the relative intensities (data point) for a given Sample may be a critical indicator for a genotype call.**

Submitted results are copied to the Results icon in the Project tree and are seen in the final report or SNPviewer file. Note the statement above\*\*\*. In this way, plots used for validating Assays or where most of the data is not representative of the final calls, remain in Kraken but are not in a final report.

Kraken provides rapid reporting, drilling of Master plates to see % of calls (and determine whether Samples are still good), Project overviews for accounting. User privileges, and user management is available. Database maintenance tools are provided and the MySQL database tables are documented in this manual.

Kraken may be used by many users or only one, in a factory setting or in a smaller laboratory. It has XML import functions and can interface with upstream and downstream systems.

# KRAKEN LIMS

## Laboratory Information Management System

- The flowchart below demonstrates how a normal user would interact with the Kraken software to log in, add a Customer, assign a Project to that Customer. Within the Project, Assays are added, by adding a DNA sequence; oligonucleotides may be designed using the **Primer picker** tool, and ordered through Kraken, which streamlines production of oligos and record keeping.
- Kraken reinforces workflow. When oligos arrive, they are imported into Kraken via the Oligonucleotide synthesis report. Oligos are recorded in Primer plates in Kraken (see the Primer icon in the Project tree. They are automatically associated with Assays. A Primer plate represents a physical plate or box of unmixed oligonucleotides.
- Assay Aliquot plates are created to receive oligonucleotide mixes ready for PCR. Assay aliquot tube positions are recorded and start of with a status of **Empty tube**. As oligonucleotides are mixed into new Assay Aliquot tubes they are arranged in these tubes in the Assay plates (or Assay Aliquot plates) in Kraken (a plate represents a physical box). As they are dispensed into Assay mixes, the Assay tube status is changed, and as the Assays are tested (validated) then the Assay status is changed by the person assessing or scoring the results. An Assay Aliquot may be given the status **In use** or **Do not use** and others.
- Similarly, {Samples} are received in Master plates and their location recorded by plate and physical storage location (e.g. walk-in freezer) in Kraken. As the {Samples} are diluted and dispensed into working Sample plates (daughter or sister or replicate plates) each dispense is monitored in Kraken. There are icons in the Project tree for the original Master plates are located and another for the working Sample plates. Each Project must have its own Master and working Sample plates, and Kraken makes it easy to transfer {Samples} between plates.
- Into working Sample plates (created blank) must be dispensed, virtually in Kraken:
  - first, the DNA from the Master or diluted working Sample plates (sisters) created from them.
  - second: Assay mixes are dispensed from Assay aliquot tubes.
  - third: plate reader files may be assigned to each plate - these are organised under the Assays icon in the Project tree.
- *Note:* This order is strictly controlled in Kraken: an Assay cannot be dispensed into a plate with no DNA samples in it and a plate read file cannot be assigned to a plate where an Assay has not been dispensed.
- The Assays icon in the Project tree holds all of the results, and there you will find for each Assay, the Assay designs, and for each design, the read or reads (plate reader file data) which follow thermal or hydrocycling. From this part of the tree each set of results is rapidly accessed as a cluster or other plot. Normalisation data may be viewed.
- At this point the results are scored, and the preferred, most representative plot 'Submitted' for the final report. **Submitted** is one of many statuses which may be assigned to a plate read. \*\*\*However, more than one results set may be used to decide on the call for an individual sample. Plates are frequently re-cycled (submitted to a further 5 rounds of PCR hydrocycling) and the results compared: direction of travel of the relative intensities (data point) for a given Sample may be a critical indicator for a genotype call.
- Submitted results are copied to the Results icon in the Project tree and are seen in the final report or SNPviewer file. Note the statement above\*\*\*. In this way, plots used for validating Assays or where most of the data is not representative of the final calls, remain in Kraken but are not in a final report.
- Kraken provides rapid reporting, drilling of Master plates to see % of calls (and determine whether Samples are still good), Project overviews for accounting. User privileges, and user management is available. Database maintenance tools are provided and the MySQL database tables are documented in this manual.
- Kraken may be used by many users or only one, in a factory setting or in a smaller laboratory. It has XML import functions and can interface with upstream and downstream systems.
-