## INSTRUCTIONS FOR CHANGING FILTERS

### ENGLISH

#### HESPA AIR FILTERS FOR ELLIPE SAFETY RESPIRATOR

The HESPA air filter is approved as a particle filter in accordance with European Standard EN143:2000 / A1.

These filters are designed to fit the Ellipse Safety Respirator only, which is approved as a half-mask in accordance with European Standard EN140:1998.

Note: only when the above items are used together is the device approved for occupational applications.

#### WARNINGS & LIMITATIONS

Filters can only provide protection when used with a suitable, fully serviceable half mask or full face mask and in atmospheres where the type and level of hazard is known.

DO NOT USE filters where the hazard is unknown, or may be immediately dangerous to health. Filter respirators must not be used in oxygen deficient (<19.5%) or oxygen-enriched atmospheres unless the gas analysis and concentrations are not known or exceed the acceptable values for the use of filter, or explosive atmospheres.

DO NOT USE filters that have exceeded their shelf life expiry date (see packaging).

DO NOT USE in an explosive atmosphere or where high levels of heat persist.

DO NOT USE damaged filters.

DO NOT USE with beards or facial hair that prevents direct contact between the face and mask.

DO NOT USE to replace valves.

DOES NOT provide protection against hazards of oxygen and oxygen-enriched air.

ALWAYS conduct a pressure face fit check prior to each use, if the check fails the mask may be faulty or require replacement depending on age, frequency of use and replacement of filters.

Read this User Instruction Leaflet prior to using a filter.

#### PARTICULATE Filters

The main particulate filter types are shown below together with the appropriate standard:

Gas Filter Type	Colour Code	Main Applications
P3 R	White	Particulate filters

Filter Class	Main Applications
P3 Sodium Chloride * <0.05% penetration	For use against hazardous particulates including both solid and liquid particles. Dust, mist and fume of all types, and micro-organisms eg. Bacteria and viruses.

\* Note: Full details of penetration tests for both liquid and solid particles are given in EN 143:2000.

**STORAGE**  
The mask should be stored out of direct sunlight, away from sources of high temperature and in an uncontaminated environment.

Storage under conditions other than those specified may affect shelf life. Filters must be stored in a cool dry atmosphere at a temperature of between -10°C and 40°C. Shelf life is 5 years for P3 filters. The month and year of expiry is marked on the neck of the mask of packaging.

**SELECTION OF A FILTER**  
FILTERS MUST BE SELECTED BY A COMPETENT PERSON FAMILIAR WITH WORKPLACE HAZARDS AND THEIR CONCENTRATION. Particulate filters must be changed when breathing resistance increases as the filter becomes clogged.

**FITTING**  
Always check expiry date is not exceeded before unpacking replacement filters.

1. Looking inside mask, undo turnbuckle fixing (1/4 turn anti-clockwise). Turnbuckles located behind face-seal.  
2. Fold strap supports back to expose filters.  
3. Flex mask to remove both used filters.

4. Slide a new filter into position with all text facing inwards towards the mask, and the narrow end of the filter facing the "TOP" opposite "TOP" text on the mask.

5. Carefully align the filter into the mask using the channel on the filter. Make sure that the mask seal is fully located into the channel on the outside of the filter to ensure a good seal is made.


6. Fold strap support back into position and re-fit the turnbuckle (1/4 turn clockwise). Make sure that the mask rim is clear of the turnbuckle head on the inside of the mask.

7. Repeat step 4 to 6 for the other filter. Always replace both filters at the same time, and keep a record of the expiry date.

After use, dispose of the filters according to the filtered substance in accordance with current waste treatment regulations.

The manufacturer guarantees the saving of basic consumer properties during service life of the half-mask and filters.

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