

Micro 10[®]+



Economical, aldehyde-free concentrated solution for pre-sterilisation disinfection and cleaning of all dental and surgical instruments, including burs, root canal instruments, forceps, probes, mirrors, etc. Micro 10[®]+ has a broad spectrum of efficacy and can be used in a soaking tank and ultrasonic cleaner.

Key product features

- Very economical 2% dilution
- Contact time: 15 minutes
- Suitable for all instruments
- Contains corrosion inhibitors to protect instruments against rusting
- CE marked – Complies with Medical Devices Directive 93/42/EEC
- VAH/DGHM certificated



USF Healthcare S.A. | Rue François Perréard 18 | CH-1225 Chêne Bourg | Geneva | Switzerland
Tél + 41 22 839 79 00 | Fax + 41 22 839 79 10 | info@usfhc.com
www.usfhealthcare.com



Micro 10⁺


Protocol for use

Before use

Ensure appropriate gloves are worn.



Dilute Micro 10⁺ in a soaking tank at 2% (20ml per 1 litre of water).



Immerse the instruments in the solution and leave for 15 minutes.



Carefully remove the instruments from the solution and inspect. Brush if necessary to remove any remaining visible debris.



Rinse instruments under running water or by immersion in clean water.



Dry instruments with a clean disposable paper towel. Pack the instruments and prepare for sterilisation.


Precautions

Dangerous: respect the precautions for use (established according to the European regulations in force for classification and labelling of chemicals). Before each use, read the label and the information concerning the product. Storage between +5°C and +35°C. Reserved for professional use. Medical device Class IIb (complies with medical devices directive 93/42/EEC modified).

Recommended dilutions

- 2% (20ml in 1L water) for all dental instrumentation
- 5% (50ml in 1L water) for ultrasonic cleaning: 3 to 8 minutes. Supports heating up to 50°C.
- 10% (100ml in 1L water) for small instruments, e.g. canal instruments, diamond burs, etc.

Efficacy



Conforms to European standards

Bactericide: EN 1040, EN 13727, EN 14561

Yeasticide: EN 1275, EN 13624, EN 14562

Active on viruses according to EN 14476

Product composition

100 g of Micro 10⁺ contains 9 g of alkylbenzyltrimethylammonium chloride, amphoteric and non-ionic surfactants, complexant, corrosion inhibitor and additives.

Physical data

Colour: Pale yellow

Perfume: Without perfume

pH: 7.5 (diluted at 2%)

Shelf Life: Concentrate 3 years

Diluted, unused 28 days

Diluted, used 7 days depending on soil quantity

Container: Made of recyclable HDPE material

Packaging

1L	2,5L	5L
✓	✓	✓

Accessories

3L Soaking tank (INTANK)

Soaking tank for burs



ECOEngage
Our Ecological Engagement



Unident
USF Healthcare