Meisinger

# BONE MANAGEMENT® CATALOGUE 2021 / 22

BONE  
AUGMENTATION  
TECHNIQUES

– YOUR KEY  
TO SUCCESS

# INTERNATIONAL BONE MANAGEMENT® NETWORK

## BONE MANAGEMENT® ACADEMY

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When developing our surgical systems and instruments, we are not only concerned with the highest quality of "Made in Germany", but also with the correct application of these systems and instruments. Since the professional preparation of the bony implant site is essential for successful implant placement and management, MEISINGER has developed a multi-level training concept in cooperation with leading scientists and practitioners. With the Bone Management® Academy, both experienced users and beginners have the opportunity to deal intensively with the topics of bone augmentation, implantology, content related to practice management as well as providing the chance to further their education in this field.

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Our main motivation is on the continuously increasing added value for participants. The focus is therefore not on the systems developed by MEISINGER, but on professional, practical, and patient-oriented techniques that can be learned in the courses. High-quality and precisely manufactured products play an essential role in positive and successful patient outcomes. However, the correct and proficient application of these tools is the famous "key to success" and thus also the best way to a satisfied patient.

That is why we do not shy away from our competition, but deliberately involve other dental companies in order to enable our participants to form a balanced perspective, and to discuss and compare competing treatment methods. If the participants find certain techniques to be good, the products can firstly be intensively tested and practiced in the workshops. Additional work shadowing and supervision offer further opportunities to consolidate what has been learned and thus make better informed decisions.

The Young dental experts camps, which take place several times a year, also offer orientation aid for students, interns, and residents in their future specialization choice. Aspiring surgeons can also acquire solid basic knowledge in the various courses and workshops. The right suturing and cutting techniques are just as important as the right understanding of soft and hard tissue. Various augmentation techniques can also be learned and practiced within the curriculum. The Bone Management® Academy also includes complex techniques such as the Carrot and Shell techniques developed by Prof. Dr. Fouad Khoury or the Sausage technique according to Prof. Dr. Istvan Urban.

MEISINGER USA Headquarters, Colorado

MEISINGER Germany Headquarters, Neuss

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## 10 BONE SPREADING

---

10-11 CCR00 | Crest-Control  
12-13 CSP12 | Split-Control 12 mm  
14-15 CSPPL | Split-Control Plus  
16-17 BSPPR | Split-Control Professional

## 20 BONE MANAGEMENT® MASTER-LINE

---

22 BMCBA | Master-Core Basic  
23 BMCPR | Master-Core Professional  
24-25 BKM00 | Master-Mill  
26 BMPBA | Master-Pin-Control Basic  
27 BMP00 | Master-Pin-Control

## 28 BONE MANAGEMENT® KHOURY-LINE

---

30-31 BTE00 | Trephine Ejection Kit  
32 BMSBA | Micro Screw System Basic  
33 BMS00 | Micro Screw System  
34-35 BOSTO | Ost-Tray  
36 BSINO | Sinus-Tray  
37 BSOTO | Soft-Tissue-Tray  
38-39 BS003,BS004 | SAFESCRAPER® Twist

## 42 BONE TRANSFER & OSTEOSYNTHESIS

---

42-43 CTRBA | Transfer-Control Basic  
44-45 CTR00 | Transfer-Control  
46 CTRIO | Transfer-Ring-Control I  
47 BBR00 | Transfer-Ring-Control II  
48-49 CTRPL | Transfer-Control Plus  
50 BTX00 | Screw System TX  
51 BTXPR | Screw System TX Professional  
52 7120 | Threphine Basic Kit  
53 7121 | Threphine Kit  
54 7122 | Threphine Kit Long  
55 7140 | Punch Basic Kit  
55 7150 | Saw Basic Kit

## 56 SINUSLIFT

---

56 CCLBA | Crestal-Lift-Control Basic  
57 BCL00 | Crestal-Lift-Control  
58 BTL00 | External-Lift-Control  
59 CSL00 | Internal-Lift-Control  
60-61 BLIPL | Lift-Control Plus  
62 BSKSL | Surgical Kit 1  
63 BSK02 | Surgical Kit 2

## 64 IMPLANT SITE PREPARATION

64	BSK03	Surgical Kit 3
64	BSK04	Surgical Kit 4
65	ALV18	Alveoplasty Surgical Kit
66	IPK02	Implant Preparation Kit with Stop
67	BDS00	Drill-Stop-Control (BL + TL)
68	BDS02	Drill-Stop-Control (BLT)
69	BDSMP	Drill-Stop-Control (myplant two)
70-71	BGS00	Guided-Drill-Stop-Control
72	BNA00	3D-Navigation-Control
73	BNAPR	3D-Navigation-Control Professional

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**Note:** Not all products contained in this catalogue are available in all countries.

# Benex® -Control

Root Extraction System | developed with Dr. Benno Syfrig

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Benex®-Control is outstandingly suitable for the controlled and safe extraction of dental roots and dental root fragments, so the conventional extraction method using forceps and leverage can be avoided. Thanks to the innovative construction of the extractor, the root can be easily removed in a particularly controlled manner without damaging bone or soft tissue. To guarantee an optimal vertical transmission of force, there must be a flat area providing a supporting surface to the extractor. The quadrant support was developed specifically to provide this optimal support while also protecting existing dentition and providing support for the extractor in case of edentulous circumstances.

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Fig. BE001 017  
links | left

Fig. BE001 018  
rechts | right

<sup>1</sup> Teflon® is a registered trademark of DuPont, USA

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## AT A GLANCE

- Safe and controlled extraction of roots
- Possibility of immediate implant placement
- Maximum preservation of the vestibular bone wall due to vertical extraction
- Optimal protection of bone and soft tissues
- Applicable to small root fragments
- Osteotomy can be prevented
- High patient acceptance

Developed with  
**Dr. Benno Syfrig**  
 Luzern, Switzerland

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Crest-Control is a Horizontal Bone Splitting System which allows for a predetermined and minimally invasive widening of the horizontally resorbed alveolar ridge particularly in the distal lower jaw. With the well-proven Horizontal Spreaders, the alveolar ridge can be expanded up to a maximum of 5 mm, allowing an implant of every common system to be inserted in the widened segment.

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## AT A GLANCE

- Controlled widening of horizontally resorbed bone
- Effective in increasing the width of narrow alveolar ridges
- Preparation of the bone for all common implant systems
- Minimally invasive treatment concept
- Gentle widening of the alveolar ridge

# Split-Control 12 mm

## Bone Condensing and Bone Spreading System

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Split-Control 12 mm is a system for gentle widening of the alveolar ridge and simultaneous lateral bone condensing in the case of a horizontal bone deficit. With the aid of screw-like condensing and spreading instruments (Spreaders), it is possible to achieve a controlled and gentle dilatation of horizontally resorbed bone (Bone Spreading). Further, cancellous bone is condensed due to the special geometry of the spreaders (Bone Condensing), so that the primary stability of the inserted implants is increased.

Split-Control 12 mm provides an extensive range of instruments for a controlled dilatation with spreaders with a length of 12 mm.

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## AT A GLANCE

- Controlled spreading of horizontally resorbed bone
- Gentle bone condensing with the aid of non-cutting Archimedes screws
- Increased primary stability of the inserted implants
- Effective in increasing the width of narrow alveolar ridges
- Preparation of the bone for all common implant systems

# Split-Control Plus

## Combination of Split-Control and Crest-Control

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
Split-Control Plus is a combination of the proved Split-Control and Crest-Control systems, so it offers an extensive range of instruments for controlled bone spreading and condensing. It contains both specially designed screw-like condensing and spreading instruments (Spreader) and horizontal spreaders that allow for controlled and standardized spreading of horizontally resorbed bone. Due to the special geometry of the screw-like spreaders, cancellous bone is gently condensed while the Horizontal Spreaders allow for a dilatation up to a width of 5 mm, so the bone is optimally prepared for the insertion of every common implant.

Anwendungsvideo  
Application Video

Art.-No. CSPPL

© Dr. Bayer & Drs. Kistler, © Dr. Raymond Kimsey, D.D.S.

	Vorkörner Initial bur	Diamantierte Sägen Diamond coated saws		Osteotomieinstrumente Osteotomy instruments		Pilotbohrer Pilot burs		Erweiterer 12 mm Expansion burs 12 mm	
Fig.	186RF	231DC*	231DC*	HM331L**	HM254E**	A1001	D1001	A2003	D2003
Shank <sup>1</sup>	204	204	204	205	205	206	206	205	205
Size <sup>2</sup>	018	100	130	010	012	010	018	023	030
Length mm	12.0	0.3	0.3	5.5	6.0	18.0	18.0	12.0	12.0
⚠	-	-	-	-	-	-	-	1.30	2.03
* □ *	1.8	10.0	13.0	1.0	1.2	1.0	1.8	2.3	3.0

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## AT A GLANCE

- Combined System of Crest- and Split-Control
- Possibility of combining bone splitting and spreading techniques
- Horizontal spreaders for gentle widening of the alveolar ridge by up to 5 mm
- Non-cutting Archimedes screws for gentle bone condensing
- Preparation of the bone for insertion of all common implant systems

# Split-Control Professional

Bone Condensing and Bone Spreading System | used by Dr. Dr. Martin Bonsmann

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Split-Control Professional is a system for gentle widening of the alveolar ridge and simultaneous lateral bone condensing in the case of a horizontal bone deficit. With the aid of screw-like condensing and spreading instruments (Spreaders), it is possible to achieve a controlled and gentle dilatation of horizontally resorbed bone (Bone Spreading). Further, cancellous bone is condensed due to the special geometry of the spreaders (Bone Condensing), so that the primary stability of the inserted implants is increased.

To be able to respond more precisely in specific indications, the Split-Control Professional contains a broad range of spreader sizes.

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Mitra Sadrameli, DMD, MS

# BONE MANAGEMENT® MASTER-LINE

Developed by  
Prof. Istvan Urban, DMD, MD, PHD

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
# BONE MANAGEMENT® COMPETENCE CENTER

Hungary, Budapest  
**Dr. Istvan Urban DMD, MD, PHD**  
 URBAN REGENERATION INSTITUTE  
 Urbán Dental Center Kft.  
 1025 Budapest  
 Pitypang utca 7., HUNGARY  
[www.implant.hu](http://www.implant.hu)

France, Paris  
**Dr. Hélène Arnal**  
 CABINET D'IMPLANTOLOGIE  
 VILLIERS  
 100 Avenue de Villiers  
 75017 Paris  
[www.selarl-cabinetdudrarnal.chirurgiens-dentistes.fr](http://www.selarl-cabinetdudrarnal.chirurgiens-dentistes.fr)

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The Master-Core Basic is used for the simple and safe extraction of bone cylinders in the context of bone augmentation. For this purpose, the system contains trepanes with different lengths and diameters. Due to the different lengths, the trepanes have an automatic depth stop, which ensures the protection of anatomical structures and thus allows for a particularly safe and controlled bone extraction. In addition, depth markings on the trepanes' working parts enable additional depth control during the extraction. The black coating of the trepanes allows for glare-free work and the optimal visibility of the depth markings.

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## AT A GLANCE

- Simple and safe extraction of bone cylinders
- Optimal protection of anatomical structures thanks to depth-limited trepanes
- Depth markings on the trepanes for additional depth control
- Extraction of bone cylinders of different diameters and lengths
- Black coating on the trepanes for glare-free work


Fig.	DV229	DV229	BV025
Shank <sup>1</sup>	204	204	205
Size <sup>2</sup>	050	070	-
Length mm	2.0	2.0	25.5
□	5.0	7.0	-
* □ *	5.7	7.7	-

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The Master-Core Professional is used for the simple and safe extraction of bone cylinders in the context of bone augmentation. For this purpose, the system contains trepanes with different lengths and diameters. Due to the different lengths, the trepanes have an automatic depth stop, which ensures the protection of anatomical structures and thus allows for a particularly safe and controlled bone extraction. In addition, depth markings on the trepanes' working parts enable additional depth control during the extraction. The black coating of the trepanes allows for glare-free work and the optimal visibility of the depth markings. To be able to respond more precisely in specific indications, the Master-Core Professional contains a broad range of instruments.

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# Master-Mill NEW

## Surgical Bone Mill – used by Prof. Istvan Urban

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The Master-Mill is a surgical bone mill used for crushing autologous bone blocks and thus preserving small autologous bone chips. With the help of these bone chips, bone defects of all shapes can be optimally filled. The individual components of the Master-Mill are made of high quality, surgical stainless steel. All parts feature an optimal fit and a stable geometry, allowing the Master-Mill to be easy to handle with ease of assembly and intuitive use.

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## AT A GLANCE

- Crushing of autologous bone blocks
- Augmentation of bone defects of all shapes
- Individual components made of high quality, surgical stainless steel
- Stable geometry
- Easy and intuitive handling

# Master-Pin-Control Basic

Pin System for Membrane Fixation | developed with Dr. Istvan Urban

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Master-Pin-Control is designed especially for the fixation of resorbable, non-resorbable and titanium reinforced membranes. The included pins consist of titanium alloy Grade 5 and have a groove which increases the surface of the pins. Due to the groove and a unique head design, the pins can be loosened easily and safely which is especially important after a successful healing period. With their sharp tip and very stable shank, they can easily and precisely be inserted into dense cortical bone.

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## AT A GLANCE

- No bending, no breakage - Safe Fixation of Membranes with Pins made of Titanium Alloy Grade 5
- Sharp tip and high stability of the pins allow for precise placement of the pins also in cortical bone and at an angle
- Easy Removal from the tray and placement of the pins due to the specially developed Pin-Holder
- Instruments for the placement of Bleeding Points for optimal healing conditions


Fig.	186RF	203RF	203RF	MP10*
Shank <sup>1</sup>	204	205	205	-
Size <sup>2</sup>	018	006	008	-
Length mm	12.0	7.0	7.0	3.65
Minimal diameter	-	-	-	0.87
External diameter	1.8	0.6	0.8	0.95
Head-Ø mm	-	-	-	2.51

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Minimal diameter External diameter

### Nachfüll-Pins / Refill Pins :

Pins à 10 (MP10\*)

Pins à 50 (MP50\*)

Pins à 100 (MP100\*)


Fig. MP11  
illustrated 1:1 (Lenght mm: 35.0)

Fig. MP12  
illustrated 1:2 (Lenght mm: 170.0)

Fig. MP14  
illustrated 1:2 (Lenght mm: 165.0)

# Master-Pin-Control

Pin System for Membrane Fixation | developed with Dr. Istvan Urban

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Holistic system composition for the extraction, fixation and stabilization of cortical bone grafts according to Prof. Dr. Khoury - optimally complemented by selected manual instrument assortments for treatment of hard and soft tissue

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Dr. Michael Berthold | München, Germany


# Trephine Ejection Kit


Trephine System for Extraction of Cylindrical Bone Grafts | developed by Prof. Dr. Fouad Khoury

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The Trephine Ejection Kit was developed together with Prof. Dr. Fouad Khoury. It contains several unique ideas which turn this kit into something special: Two-part, internally cooled trephines allow for a safe and minimal invasive extraction of local bone cylinders. The inner cooling helps to prevent a damage of the bone tissue due to overheating, while the removable trephine working parts allow for an especially easy sampling of the cut out bone cylinders. Either with the aid of the ejection needle or the ejection instruments, which perfectly match with the diameter of the trephine working parts, the bone cylinders can be pushed out of the trephines. Furthermore, the included pre-drilling trephines allow for a safe application of the trephines. Due to the two-part trephines, the variety of different diameters and the optimal matching of the instruments, the Trephine Ejection Kit offers highest efficiency, flexibility and precision to the user for the extraction of bone cylinders.

Art.-No. BTE00


Bone core are harvested from the implant site minimal invasive with the pre- and main trephine in different diameter.

Bone core augmentation stabilized with micro screws

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

## AT A GLANCE

- Two-part trephines for a safe and minimal invasive extraction of bone cylinders
- Protection against overheating due to internally cooled trephines
- Innovative application
- Standardized diameters and precise laser markings of the trephines allow for a predictable extraction of defined bone cylinders
- Perfectly matched instruments care for highest precision and flexibility

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The Bone Management® Micro Screw Systems are used for the safe fixation and stabilisation of cortical bone grafts. The special feature of these systems are the osteosynthesis screws, which have a diameter of 1.0 mm and 1.2 mm. These Micro Screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter.

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
# Ost-Tray NEW

Selection of Manual Instruments for Oral Surgery | developed by Prof. Dr. Fouad Khoury

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The Bone Management® Ost-Tray offers a perfect combination of specially configured manual surgical instruments. In addition to a well-thought-out selection of manual instruments, such as high-quality periodontal probes, curettes, scalpel blade holders, mouth mirrors, graduated mouth wedges, a curved irrigation cannula with bulb end and various scissors, needle holders, tweezers, clamps, forceps, retractors and tongue depressors, the kit contains an individually designed, strong Khoury-probe as well as raspatories and sharp spoons according to Prof. Khoury.

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**Inklusive:** Irrigation Needle with Bulb End (Art.-No. MI165)

Handinstrumente  
Manual Instruments

Fig.	MI056	MI057	MI070	MI071	MI077	MI139	MI156/165	MI160	MI163
Name	Scharfer Löffel 'Khoury', 2,0 mm Bone Curette 'Khoury', 2,0 mm	Scharfer Löffel 'Khoury', 3,0 mm Bone Curette 'Khoury', 3,0 mm	Raspatorium 'Molt' Periosteal Elevator 'Molt'	Raspatorium 'Khoury' Periosteal Elevator 'Khoury'	Skalpellklingenhalter Scalpel Handle	Kürette 'Columbia' Curette 'Columbia'	Mundspiegel Mouth Mirror	Parodontalsonde Periodontal Probe	Sonde, kräftig Probe, strong
Length mm	175.0	175.0	180.0	183.0	152.0	166.0	165.5	160.0	158.0

Handinstrumente  
Manual Instruments


Fig.	MI083	MI101	MI102	MI104	MI105	MI106
Name	Präparierschere 'Metzenbaum', TC Dissecting Scissors 'Metzenbaum', TC	Wundhaken 'Langenbeck' 60x11 mm Tissue Retractor 'Langenbeck' 60x11 mm	Wundhaken 'Langenbeck' 42x10 mm Tissue Retractor 'Langenbeck' 42x10 mm	Arterienklemme 'Halstead-Mosquito' Haemostatic Forceps 'Halstead-Mosquito'	Tuchklemme 'Tohoku' Towel Clip 'Tohoku'	Hohlmeißelzange 'Blumenthal' Bone Rongeur Forceps 'Blumenthal'
Length mm	180.0	220.0	220.0	145.0	105.0	155.0

Handinstrumente  
Manual Instruments

Fig.	MI041	MI048	MI050	MI158	MI089	MI098
Name	Nadelhalter 'Crile Wood', TC Needle Holder 'Crile Wood', TC	Pinzette Tweezers	Pinzette, TC Tweezers, TC	Pinzette 'College' Tweezers 'College'	Chirurgische Schere 'Kelly' Surgical Scissors 'Kelly'	Zungenspatel 'Brüning's' Tongue Depressor 'Brüning's'
Length mm	150.0	175.0	175.0	175.0	160.0	200.0


# Sinus-Tray NEW

Selection of Manual Instruments for external Sinus Lift | developed by Prof. Dr. Fouad Khoury

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The new Bone Management® Sinus-Tray is ideal for performing the external sinus lift. The manual instruments designed with Prof. Khoury support the user in gently lifting the Schneiderian membrane and inserting the bone substitute material.

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**Inklusive:** Irrigation Needle with Bulb End (Art.-No. MI165)

Handinstrumente  
Manual Instruments



Fig.	MI135	MI136	MI140	MI141	MI142
Name	Kondenser 'Khoury' Condenser 'Khoury'	Knochenmateriallöffel 'Khoury' Bone Material Spoon 'Khoury'	Sinus-Lift Elevator 'Khoury'	Sinus-Lift Elevator 'Khoury'	Sinus-Lift Elevator 'Khoury'
Length mm	168.0	177.0	185.0	190.0	185.0


# Soft-Tissue-Tray NEW

Selection of Manual Instruments for Soft-Tissue Management | developed by Prof. Dr. Fouad Khoury

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The new Bone Management® Soft-Tissue-Tray was developed with Prof. Khoury and is perfectly applicable in the context of soft tissue management. A special feature is the scalpel blade holder with swivel head included in the set, with which the scalpel angle can be freely adjusted depending on the processing area. This simplifies handling, especially in areas that are difficult to access.

Art.-No. BSOT0

Handinstrumente  
Manual instruments

Fig.	MI042	MI049	MI060	MI077	MI078	MI079	MI089	MI167
Name	Nadelhalter 'Crile Wood', RH Needle Holder 'Crile Wood', RH	Pinzette 'Cooley' Tweezers 'Cooley'	Tunnelierungsinstrument Tunneling Instrument	Skalpell- klingenhalter Scalpel Handle	Skalpell- klingenhalter Scalpel Handle	Skalpell- klingenhalter mit Schwenkkopf Scalpel Handle with Swivel Head	Chirurgische Schere 'Kelly' Surgical Scissors 'Kelly'	Papillenelevator Papilla Elevator
Length mm	150.0	175.0	176.0	173.0	135.0	149.0	160.0	170.0



# SAFESCRAPER® Twist **NEW**

## Ready-To-Use Cortical Bone Collector

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With the Safescraper® Twist, autologous cortical bone chips can be obtained easily and safely. These are safely collected in the transparent chamber while the Safescraper® is in use - for optimal control of the amount of bone harvested. The bone material obtained can be used immediately thanks to the immediate mixture with blood and is suitable for filling up defects of all kinds.

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## AT A GLANCE

- Ready to use
- Excellent cutting performance for a minimally invasive and patient-friendly procedure
- Straight and angled design for every situation
- Half-round blade with 160° cutting area for maximum efficiency on all bone surfaces
- Capacity of the chambers of 3.0 cm<sup>3</sup> (straight version) or 2.5 cm<sup>3</sup> (angled version)
- Safety mechanism prevents the chambers from opening during bone harvesting

# HACIS

EUROPE

MEISINGER'S

12<sup>th</sup> Annual High Altitude Comprehensive Implant Symposium

*Meisinger*

Vail, Colorado USA

Zürs/Lech AUSTRIA

24.01 - 28.01

2023

SAVE  
THE DATE!  
Robinson Club,  
Zürs

Dr. Hélène Arnal

Dr. Valentin Loriod

and many more...

# Transfer-Control Basic

## Horizontal and Vertical Bone Replacing System

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Transfer-Control Basic permits precise and standardized transplantation of bone cylinders for horizontal and vertical bone augmentation. The instruments are available in three diameters. They are perfectly adapted to each other so that the outer diameter of the abrasive bur and the size of the wheel cutter correlates with the inner diameter of the trephine. So a press fit can be produced between the cylindrical bone and the prepared recipient site. The press fit can then be strengthened with a fixation screw, which is separately available (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound healing. So a transplantable bone site is reached after only 3-4 months.

Transfer-Control Basic is an introductory system with reduced instrument set. It is suitable for the extraction of bone cylinders with diameters of 4.0 mm, 5.0 mm and 6.0 mm.

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
## AT A GLANCE

- Easy and safe transplantation of bone cylinders
- Intelligently coordinated instruments
- Perfectly fitting bone cylinders
- Precise work due to depth marks
- Optimally suitable for vertical and horizontal augmentations
- Creation of refreshed contact surfaces for rapid vitalization and wound healing

# Transfer-Control

## Horizontal and Vertical Bone Replacing System

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Transfer-Control permits precise and standardized transplantation of bone cylinders for horizontal and vertical bone augmentation. The instruments are available in five diameters. They are perfectly adapted to each other so that the outer diameter of the ablative bur and the size of the wheel cutter correlates with the inner diameter of the trephine. So a press fit can be produced between the cylindrical bone and the prepared recipient site. The press fit can then be strengthened with a fixation screw, which is separately available (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound healing. So a transplantable bone site is reached after only 3-4 months.

Transfer-Control is suitable for the extraction of bone cylinders with diameters of 4.0 mm, 5.0 mm, 6.0 mm, 7.0 mm and 8.0 mm.

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The most important condition for a safe healing of a bone graft is a congruent and fresh contact surface of the implant area. The Transfer-Ring-Control I System allows for this precondition for the healing of vertical bone grafts in a simple and controlled manner. This system offers a variety of various instrument sizes and provides high flexibility for the extraction of bone rings. Depending on requirements, various sizes of Initial Bur Trephines, Trephines and Ablative Burs are available. With the aid of the additional tungsten carbide bur, sharp bone edges can be smoothed precisely.

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# Transfer-Ring-Control II

## One-Stage Bone Ring Augmentation System | Technology by Dr. Bernd Giesenhagen

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
Transfer-Ring-Control II is a system for the bone ring technique, developed by Dr. Giesenhagen. It allows for one-stage proceeding of augmentation and implantation in case of resorbed jaw bone, if implant insertion alone is not possible because of inadequate spatial conditions. Due to the perfectly matched instruments, there is created a congruent bone and implant site, which allows for the treatment of three-dimensional defects with the help of prefabricated, allogenic bone rings. Because of the optimal fit between jaw bone, bone ring and inserted implant and the possibility of the one-stage proceeding, high primary stability is accomplished and healing time is significantly reduced.

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Technology by  
**Dr. Bernd Giesenhagen**  
Kassel, Germany

# Transfer-Control Plus

## Combined System for Horizontal and Vertical Bone Augmentation

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
Transfer-Control Plus permits precise and standardized transplantation of bone cylinders and rings for horizontal and vertical bone augmentation. Different types of bone burs allow for high flexibility for the extraction of bone grafts. The instruments are perfectly adapted to each other often producing a press fit between the cylindrical bone or ring and the prepared recipient site. This press fit can then be strengthened with a fixation screw which is sold separately (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound-healing. So a transplantable bone site is reached after only 3-4 months.

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
## AT A GLANCE

- Combined system for easy and safe transplantation of bone cylinders and rings
- Intelligently coordinated instruments
- Optimally suitable for vertical and horizontal augmentations
- Precise work due to depth marks
- Creation of refreshed contact surfaces for rapid vitalization and wound healing

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
The Osteosynthesis System Screw System TX enables the universal use of fixation for autologous bone cylinders, bone blocks and bone shells for the purpose of preimplantological augmentation. The Torx connection guarantees reliable accommodation of the screws in the driving tool and, at the same time, high tensile forces can be transferred. For each 0.9 mm screw available, the kit contains a larger diameter rescue screw in case the transplant does not allow high tensile fixation. The small diameter and the low head size of the screws allow for a comfortable application for patient and user.

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
The Osteosynthesis System Screw System TX Professional enables the universal use of fixation for autologous bone cylinders, bone blocks and bone shells for the purpose of preimplantological augmentation. The Torx connection guarantees reliable accommodation of the screws in the driving tool and, at the same time, high tensile forces can be transferred. For each 0.9 mm screw available, the kit contains a larger diameter rescue screw in case the transplant does not allow high tensile fixation. The small diameter and the low head size of the screws allow for a comfortable application for patient and user.

The Screw System TX Professional contains additional screws in the lengths of 13.0 mm and 16.0 mm.

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
The Trephine Basic Kit provides a selected set of high-quality trephine drills in different sizes. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 10 mm and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

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### AT A GLANCE

- Selection of high quality trephines for the extraction of bone cylinders
- Various trephine sizes for every application
- Precise depth markings for optimal control
- High cutting performance for good handling

# Trephine Kit

## Trephine System for Extraction of Cylindrical Bone Grafts

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The Trephine Kit provides a broad set of high-quality trephine drills with fine size graduations. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 10 mm and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

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## AT A GLANCE

- Broad selection of high quality trephines for the extraction of bone cylinders
- Various trephine sizes for every application
- Precise depth markings for optimal control
- High cutting performance for good handling

# Trephine Kit Long

## Trephine System for Extraction of Cylindrical Bone Grafts

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The Trephine Kit provides a broad set of high-quality trephine drills with fine size graduations. These are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 14 mm\* and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

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# Punch Basic Kit

## Biopsy Punch Set for Extraction of Tissue Grafts

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The Punch Basic Kit contains Biopsy Punches in five different diameters for the performance of tissue removal. They serve for the application in immediate implant placement, for taking mucosa grafts in the scope of soft tissue management or biopsy as well as for uncovering the gingiva after implantation. The Biopsy Punches create clean cut surfaces and care for minimal tissue traumatization.

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# Saw Basic Kit

## Selection of Various Diamond Coated Saws for Oral Surgery

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The Saw Basic Kit provides a selection of diamond saws in the three most important diameters. They are used for the application in bone spreading or splitting operations. Due to the various sizes, the set offers the clinician the possibility of having a proper selection of these saws during surgery.

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© Dr. Bayer & Drs. Kistler

# Crestal-Lift-Control Basic

## Crestal Approach Sinus Lift System

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Crestal-Lift-Control Basic is a system for the performance of an easy and safe internal sinus lift. The elevation of the Schneiderian membrane occurs during the transcrestal drilling process. The specially atraumatic design of the Crestal Drill with its four cutting edges and concave head allows for a safe forming of a conical bone flap and is more ideally suited for collecting bone chips. Additionally, the stop sleeve system, which is tailored to the special instruments, prevents the membrane from being injured or punctured.

Crestal-Lift-Control Basic is an introductory system with reduced instrument set.

**CAUTION:** To prevent injury of any anatomical structure, the apical extra lengths of the CL001 of 0.58 mm must be considered.

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Crestal-Lift-Control is a system for the performance of an easy and safe internal sinus lift. The elevation of the Schneiderian membrane occurs during the transcrestal drilling process. The specially atraumatic design of the Crestal Drill with its four cutting edges and concave head allows for a safe forming of a conical bone flap and is more ideally suited for collecting bone chips. Additionally, the stop sleeve system, which is tailored to the special instruments, prevents the membrane from being injured or punctured.

Crestal-Lift-Control provides an extensive range of instruments for the performance of a safe internal sinus lift.

**CAUTION:** To prevent injury of any anatomical structure, the apical extra lengths of the CL001 of 0.58 mm must be considered.

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


External-Lift-Control is a system for the easy performance of an external sinus lift. After lateral access through the facial wall of the maxillary sinus, the Schneiderian membrane is detached using the various manual instruments and lifted cranially from the floor of the sinus. Then, a suitable augmentation material can be inserted under the membrane. Due to the special CARBOCER®-coating, the hand-held instruments provide a reflection-free working environment.

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
Internal-Lift-Control is a system for the performance of an internal sinus lift. With the aid of the elevators, the bone inside the implant site is displaced in the direction of the maxillary sinus. This procedure is suitable for the augmentation of small defects with a residual bone height of more than 5 mm.

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# Lift-Control Plus

## Combined System for Internal and External Sinus Lift

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The combined Lift-Control Plus system provides perfectly matched instruments for the performance of internal as well as external sinus floor elevation. For the internal sinus lift, the bone inside the implant site is displaced in the direction of the maxillary sinus. This procedure is suitable for the augmentation of small defects with a residual bone height of more than 5 mm. After lateral access through the facial wall of the maxillary sinus for the performance of the external sinus lift, the Schneiderian membrane is detached using the various manual instruments and lifted cranially from the floor of the sinus. After this, a suitable augmentation material can be inserted under the membrane.

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The Surgical kit 1 contains instruments with handpiece shank for the lateral approach in an external sinus lift. Carbide cutters in three diameters and diamond instruments in three grains allow for the gradual thinning of the facial bone wall of the maxillary sinus down to the Schneiderian membrane. The large diameter of the instruments protects the membrane from unintentional perforation.

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## AT A GLANCE

- Compilation of handpiece instruments for the lateral approach during the external sinus lift
- Round instruments protect the Schneiderian membrane from unintentional perforation
- Ideal for gradual, safe thinning of the fascial bone wall
- Round diamonds and carbide burs in three different diameters

# Surgical Kit 2

## Lateral Approach System for External Sinus Elevation

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The Surgical Kit 2 contains instruments with right-angled shank for the lateral approach in an external sinus lift. Carbide cutters in three diameters and diamond instruments in three grains allow for the gradual thinning of the facial bone wall of the maxillary sinus down to the Schneiderian membrane. Especially the large diameter of the instruments protects the membrane from unintentional perforation.

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## AT A GLANCE

- Compilation of contra-angle instruments for the lateral approach during the external sinus lift
- Round instruments protect the Schneiderian membrane from unintentional perforation
- Ideal for gradual, safe thinning of the fascial bone wall
- Round diamonds and carbide burs in three different diameters

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The Surgical Kit 3 contains conical Tungsten Carbide Burs with friction grip shank of various sizes and toothing which allow for a broad field of application. The particularly easy cutting and efficient cutters have an extra-long shank for an optimal view on the working area and sharp tips for the safe application in bone.

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The Surgical Kit 4 contains a selection of conical Tungsten Carbide Cutters with handpiece shank for a broad field of application in oral surgery. The cutters have a particularly easy cutting and efficient saw toothing and a specifically designed twist angle. Further, the narrow working parts allow for an optimal view on the working area.

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# Alveoplasty Surgical Kit

## Surgical Contouring and Modification Kit

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Art.-No. ALV18

The Alveoplasty Surgical Kit contains the most important instruments for surgical contouring and modification of the alveolar ridge when unwanted unevenness of the jawbone occurs after tooth extraction or bone degeneration. In this case, the kit offers the best way to smooth the jawbone. The Alveoplasty Kit contains two initial burs and two Lindemann burs, which allow for the performance of precise, initial osteotomies, even in hard cortical bone. To perform the alveoplasty itself, five different carbide cutters and a round diamond are included. The round working parts optimally protect the soft tissue against injuries. Thus, the alveolar ridge is prepared carefully and precisely for implant insertion or denture restoration.

	Vorkörner (2 Stk.) Initial burs (2 pcs.)	Chirurgiefräser Surgical cutters	Hartmetallinstrumente Tungsten Carbide Instruments				Diamantinstrument Diamond instrument		
Fig.	186RF	HM166**	165RF	HM71**	HM71**	HM72GX**	HM77GX**	HM251GX**	801*
Shank <sup>1</sup>	204	205	205	104	104	104	104	104	204
Size <sup>2</sup>	018	021	023	040	050	060	060	040	033
Length mm	12.0	11.0	7.0	-	-	12.0	11.5	8.3	-
* □ *	1.8	2.1	2.3	4.0	5.0	6.0	6.0	4.0	3.3

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# Implant Preparation Kit with Stop

## Initial Bur System with Stop for Safe Implant Site Preparation

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The Implant Preparation Kit provides the easy and safe preparation of the implant site. The stop sleeves allow for a defined depth stop during the initial and pilot drilling. For the application without stop sleeves, the instruments are equipped with laser markings.

**CAUTION:** To prevent injury of any anatomical structure, the apical extra length of the pilot drill of 0.8 mm must be considered.

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## AT A GLANCE

- Simple and safe preparation of the implant site
- Defined depth stop during initial and extension drilling
- Colored marking of the stop sleeves for easy differentiation
- Additional depth-marked instruments for optional use without stop sleeves

# Drill-Stop-Control (BL + TL)

## Drill System for Safe Implant Site Preparation

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The Drill-Stop-Control system provides an easy and safe preparation of the implant site. The increasing diameters of the twist drills are perfectly coordinated with each other. The stop sleeves provide the defined depth stop during the pilot drilling and during the expansion of the implant site. For the application without stop sleeves, the instruments are equipped with laser markings.

**CAUTION:** To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.8 mm must be considered.

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


The Drill-Stop-Control BLT provides an easy and safe preparation of the implant site for Bone Level Tapered implants of the Straumann® or Dental Ratio® implant system. The increasing diameters of the twist drills are perfectly coordinated with each other. The stop sleeves provide the defined depth stop during the pilot drilling and during the expansion of the implant site. For the application without stop sleeves, the instruments are equipped with laser markings.

**CAUTION:** To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.4 mm must be considered.

Art.-No. BDS02

# Drill-Stop-Control (myplant two)

## Drill System for Safe Implant Site Preparation

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With the specially for myplant two developed Drill-Stop-Control system, the preparation of the implant site is particularly easy and safe, since the diameters of the twist drills are perfectly matched to the myplant two implant system. The stop sleeves are used both during the pilot and the extension drilling for the defined depth stop. The lengths of the stop sleeves are ideal for subcrestal placement of the implant. For application without stop sleeves, the instruments are equipped with laser markings.

**CAUTION:** To prevent injury of any anatomical structure, the apical extra lengths of the instruments must be considered.

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# Guided-Drill-Stop-Control


Guided Drill System for Safe Implant Site Preparation | developed with Dr. Henriette Lerner

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
The Guided-Drill-Stop system is used for the easy and safe preparation of the implant site. The special feature of the system are the specially developed twist drills with guiding tip, with which an intermediate step is added prior to the actual extension drilling. With this step, the upper part of the implant site is expanded to the desired diameter, while the lower part is used to guide the drill. This avoids deviation from the desired axis. In addition, the stop sleeves provide the defined depth stop during all drillings.

**CAUTION:** To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.8 mm must be considered.

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HL-DENTCLINIC  
www.hl-dentclinic.de

Developed with  
**Dr. Henriette Lerner**

Baden-Baden, Germany


# 3D-Navigation-Control

Guided Sleeve System | developed with Dr. Hans-Joachim Nickenig

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3D-Navigation-Control is a guiding sleeve system that is suitable for all commonly used implant systems and planning programs. Due to the included closed standard sleeves and the open guiding sleeves for the molar region, the system permits the use of surgical guide templates in the anterior region and even in tight molar situations. In addition, the system provides special tool requirements of flapless surgery. The 3D-Navigation-Control system contains a separate 3D-Laboratory Kit so the dental technician is able to install the sleeves in the template precisely.

Anwendungsvideo  
Application Video

Art.-No. BNA00

© Dr. Hans-Joachim Nickenig

Inhalt | content: 3D-Navigation-Control & 3D-Navigation-Control Professional

Fig.	HN010	HN011	HN012	HN013	HN014	HN015	HN016	HN017	HN018
Shank <sup>1</sup>	205	206	206	206	206	206	206	206	206
Size <sup>2</sup>	020	020	025	030	035	040	020	025	030
Length mm	16.0	19.0	19.0	19.0	19.0	19.0	22.0	22.0	22.0
□ □	2.0	2.0	2.5	3.0	3.5	4.0	2.0	2.5	3.0

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
# 3D-Navigation-Control Professional

Guided Sleeve System | developed with Dr. Hans-Joachim Nickenig

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3D-Navigation-Control Professional is a guiding sleeve system that is suitable for all commonly used implant systems and planning programs. Due to the included closed standard sleeves and the open guiding sleeves for the molar region, the system permits the use of surgical guide templates in the anterior region and even in tight molar situations. In addition, the system provides special tool requirements of flapless surgery. The 3D-Navigation-Control Professional system contains a separate 3D-Laboratory Kit so the dental technician is able to install the sleeves in the template precisely.

© Dr. Hans-Joachim Nickenig

The Professional system offers an additional, separate 3D-Fixation Kit, which allows for the (optional) fixation of the surgical guide template.

## Inhalt | content: 3D-Navigation-Control & 3D-Navigation-Control Professional

	Standard Führungshülsen, geschlossen (je 4 Stk.) Standard guiding sleeves, closed (each 4 pcs.)					Molare Führungshülsen, offen (je 4 Stk.) Molar guiding sleeves, open (each 4 pcs.)				Schablonenbohrer Template drills	
Fig.	HN001	HN002	HN003	HN004	HN005*	HN006	HN007	HN008	HN009*	HN021	HN022
Shank <sup>1</sup>	-	-	-	-	-	-	-	-	-	205	205
Size <sup>2</sup>	-	-	-	-	-	-	-	-	-	045/055	060/065
Length mm	6.0	6.0	6.0	6.0	6.0	5.5	5.5	5.5	6.0	13.0	13.0
□	2.0	2.5	3.0	3.5	4.0	2.0	3.0	4.0	5.0	-	-
□	4.0	4.0	4.0	4.0	4.5	5.0	5.0	5.0	6.0	5.5	6.5
△	-	-	-	-	-	-	-	-	-	4.5	6.0

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## Zusätzlicher Inhalt | Additionally content: 3D-Navigation-Control Professional

	Kugelhkopfancker (je 10 Stück) Ball-head anchors (each 10 pcs.)		Schraubendreher Screw driver		Spiralbohrer Twist drill
Fig.	HN023*	HN024*	TC0SW	UG00W	HN019
Shank <sup>1</sup>	-	-	204	-	205
Size <sup>2</sup>	-	-	-	-	010
Length mm	8.0	11.0	13.0	12.0	12.0
△	1.0	1.0	-	-	-
□	1.3	1.3	-	-	1.0
Head-Ø mm	2.25	2.25	-	-	-

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Developed with  
**Dr. Hans-Joachim  
Nickenig**  
Köln, Germany

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The Periimplantitis Kit contains eight perfectly matched tungsten carbide finishing burs in egg and flame shapes for the idea peri-implantitis treatment. Both shapes come in two sizes and each in standard and ultra-fine toothing. The shapes, sizes and toothings of the working parts are perfectly suited for the intraoral treatment of titanium and for different shapes of implant neck and shoulder. All finishers have an extra-long FG shank (total length 32 mm) which allows particularly the treatment of difficult-to-access, deep areas. For the peri-implantitis treatment, the finishers are used clockwise and are guided around the exposed implant counter-clockwise so the implant surface is smoothed.

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## AT A GLANCE

- Perfectly matched tungsten carbide finishing burs for the effective, intraoral treatment of titanium
- XXL FG shanks facilitate the treatment of difficult to access implants
- Different shapes and sizes for every individual patient case
- Clean and smooth titanium surfaces

# Degranulation Kit **NEW**

## Diamond System for effective Bone Debridement

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The Degranulation Kit contains four extra coarse, round diamond instruments with a diameters of 1.0, 2.5, 3.0 and 3.5 mm. The instruments have been specially developed to enable quick and easy bone debridement. With the help of these diamonds, degranulation tissue adhering to the bone can be removed quickly and easily without damaging the surrounding bone, so that subsequently inserted grafting material comes into direct contact with healthy bone. This makes the Degranulation Kit a crucial addition to the degranulation process.

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## AT A GLANCE

- Extra coarse, round diamonds for quick and easy debridement
- Reliable removal of degranulation tissue without damaging the surrounding bone
- Four different diameters for use in all situations
- Extra long instrument necks for an optimal view of the operating area
- Optimal addition to the degranulation process

# MEISINGER

## Porcine Collagen Membrane


### Collagen Membrane from porcine Pericardium

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MEISINGER's porcine collagen membrane is a natural membrane made from porcine pericardium. The special manufacturing process ensures that the microporous, 3-layered architecture is preserved and no cross-linking chemicals are needed. This extraordinary structure leads to a very flexible membrane with particularly high tensile strength.

The porcine collagen membrane is characterized by its smooth handling, high flexibility and ideal adaptability. Once hydrated the membrane loses its original shape and becomes flexible, so that it optimally adapts to all defect surfaces. As the surface of the membrane is identical on both sides, both sides can be placed against the defect.

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### AT A GLANCE

- Made of porcine pericardium
- Special, 3-layered architecture for high tensile strength
- Smooth handling
- Especially high flexibility and ideal adaption to the defect
- Both sides of the membrane can be placed against the defect
- Resorption time of about 26 weeks

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
MEISINGER's bovine collagen membrane is made from highly purified type 1 collagen derived from bovine achilles tendon. A great advantage of the bovine membranes is their ideal handling. They combine high flexibility with optimum stability. This makes them flexible, adaptable to any shape and easy to place. At the same time they are stable enough to hold their shape and not to collapse into the defect. The unique fiber orientation ensures maximum tensile strength. Thus, the membranes can be safely fixed by sutures or pins without tearing. Due to the long and predictable resorption time of 26 to 38 weeks, there is no risk of particle loss due to early resorption.

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### AT A GLANCE

- Long resorption time of 26 to 38 weeks
- Ideal handling
- High flexibility for easy ridge preservation
- Optimum dimensional stability
- High tensile strength for safe fixing


# MEISINGER dPTFE Membrane

## Micro-textured, high-density PTFE Membrane

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MEISINGER's dPTFE membranes are made of high-density PTFE with a pore size of less than  $0.3 \mu\text{m}$ . Thus, the non-absorbable membrane is insensitive to bacteria and allows a maximum of predictability and aesthetics. When using the dPTFE membranes, no primary wound closure is necessary. The exposed membrane can be removed with tweezers after at least 21 days without surgical intervention, no additional anesthesia is required. The hexagonal depressions are designed to increase the surface area, which leads to an increased stability of the membrane and favors cell adhesion without increasing the porosity. For best results, the side with the depressions should be placed in the direction of the soft tissue. The Meisinger dPTFE membranes can be cut to size as needed.

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## AT A GLANCE

- Non-absorbable membrane made of high-density PTFE
- Maximum of predictability and aesthetics
- No primary wound closure necessary
- Non-surgical removal of the membrane
- Hexagonal depressions of the surface for high stability
- No disruption of wound healing due to early resorption

# MEISINGER Titanium-Reinforced dPTFE Membrane

## Titanium-reinforced, high-density PTFE membrane

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MEISINGER's Titanium-Reinforced membranes are made of high density PTFE. In combination with titanium reinforcement, the membranes offer exceptional dimensional stability and optimum flexibility at the same time. Therefore, the membranes are particularly suitable for the treatment of 3- to 4-wall defects in the course of the alveolar ridge augmentation. The wide variety of membrane shapes and sizes covers a variety of different defects and minimizes the need for trimming. If necessary, however, the membranes can be optimally trimmed to the required shape. The titanium-reinforced membranes are available in thicknesses of 150  $\mu\text{m}$  and 250  $\mu\text{m}$ .

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## AT A GLANCE


- Non-absorbable membrane made of high-density PTFE
- Fine Titanium Reinforcement from Titanium Grade I
- Optimum combination of flexibility and dimensional stability
- Treatment of 3- to 4-wall defects
- Wide range of different membrane shapes
- Availability of two different membrane thicknesses

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
Anwendung bei großen bukkalen Defekten.  
Designed for large buccal defects.

Art.-No. MTi2502536-1    Art.-No. MTi2502536-2  
Art.-No. MTi1502536-1    Art.-No. MTi1502536-2  
25 mm x 36 mm            25 mm x 36 mm  
(1 membrane / box)      (2 membranes / box)

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
The unique circular macroporous design of the Reinforced PTFE Meshes allows for direct contact between the bone graft and periosteum, allowing naturally occurring revascularization and infiltration of cells into the bone graft. The fine titanium frame maintains space essential for horizontal and vertical ridge augmentation. The PTFE Meshes have a thickness of 200 µm and easily conform to tissue contours.

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## AT A GLANCE

- Non-absorbable PTFE Mesh
- Circular macropores for direct contact between bone graft and periosteum
- Fine titanium frame maintains space for horizontal and vertical ridge augmentation
- Maximum of predictability and aesthetics
- No disruption of wound healing due to early resorption

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
Art.-No. RPM200PDML  
(Unterkiefer links / Mandible left)  
38 mm x 38 mm  
(1 membrane / box)

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# Cytoplast™ RTMPlug, RTMFoam, RTMTape

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## RTMFOAM

Art.-No. RTMFOAM10  
20 mm x 40 mm  
(Dicke/Thickness: 3 mm, 10 / box)

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## APPLICATIONS:

- Surgical wounds
- Periodontal surgical wounds
- Extraction sites
- Dental sores
- Oral ulcers (non-infected or viral)
- Suture sites
- Burns
- Traumatic wounds

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
The high quality of the MEISINGER PTFE suture allows for atraumatic wound closure after dental procedures. The high-quality needles are made of stainless steel and have a very high bending strength with a minimum of wear. Sharp needles are the result, even after repeated piercing. The constant cutting performance allows for a particularly tissue-saving work. The suture is made of 100% medical grade bioinert PTFE. Thanks to its smooth surface, the monofilament prevents the adhesion of bacteria to the thread and ensures a resistance-free gliding through the tissue. The PTFE is a soft, non-resorbable material with consistent tensile strength, allowing the user to work very precisely, planning with a defined healing time, and patient-friendly care.

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The Pro-Fix™ Precision Fixation System is manufactured to precise tolerances to ensure easy pick-up of screws, stable transfer to the surgical site and quick engagement in cortical bone. All kit components are labeled, organized, stored and sterilized together to simplify things for the surgeon and surgical team.

## Tecapro™ Storage Tray, Tecapro™ Storage Tray

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The autoclavable Tecapro™ Storage Tray is cleverly designed to store not only fixation instruments, but also a variety of membrane fixation, bone fixation, and tenting screws needed in dental bone grafting.

## Schrauben-Organizer, Screw Organizer Dial

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Up to 100 fixation screws – including bone fixation, membrane fixation, and tenting – may be stored in a built-in screw organizer dial designed and labeled for easy identification, easy storage and simple re-ordering.

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The stainless steel driver handle and cruciform driver blade fit together to make up a drive system that ensures easy pick-up of screws and stable transfer to the surgical site.

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Blades work universally with all Pro-Fix™ screws.

Art-No.	Inhalt/Content	
PFT	Autoklavierbares Tecapro™ Storage Tray mit Schrauben-Organizer Autoclavable Tecapro™ Storage Tray with screw organizer dial	
PFDH	Schraubendreher-Griff aus rostfreiem Stahl Stainless steel driver handle	
PFDB	76 mm Kreuzschraubendreher 76 mm cruciform driver blade	
PFDB56	56 mm Kreuzschraubendreher 56 mm cruciform driver blade	
BI1001	1,2 mm Durchmesser RA-Pilotbohrer 1.2 mm diameter latch type pilot drill	
PFDBCA	RA-Schraubendreher (24 mm lang, 10 mm exponierte distale Länge) Contra angle blade (24 mm long, 10 mm exposed distal length)	

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Pro-Fix™ Membrane Fixation Screws are designed as an attractive alternative to using tacks for membrane stabilization. Easy pick-up, solid stability of the screw during transfer to the surgical site, and easy placement make membrane fixation fast and easy.

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### Further self-drilling Membrane Fixation Screws

Art-No.	verpackt à packed à	Inhalt/Content	
PFMF-5	5	1,5 mm x 3,0 mm selbstschneidende Schrauben aus Titanlegierung zur Membranfixierung 1.5 mm x 3.0 mm self-drilling titanium alloy membrane fixation screws	

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Pro-Fix™ Tenting Screws are designed with a self-drilling tip, polished neck and broader head to maintain space under resorbable and non-resorbable membranes in horizontal and vertical bone regeneration procedures.

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### Further self-drilling Tenting Screws

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Pro-Fix™ Bone Fixation Screws are designed with finer pitched, self-tapping threads that give the screws greater clamping force while using less driver torque. The screws' threads are equipped with a cutting flute that allows for easier insertion into harder bone. The screws are placed into a 1.2 mm pre-drilled pilot hole.

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### Further self-tapping Bone Fixation Screws

Art-No.	verpackt à packed à	Inhalt Content	
PFB8	1	1,5 mm x 8,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 8.0 mm self-tapping Bone Fixation Screw	
PFB8-5	5	1,5 mm x 8,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 8.0 mm self-tapping Bone Fixation Screws	
PFB10	1	1,5 mm x 10,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 10.0 mm self-tapping Bone Fixation Screw	
PFB10-5	5	1,5 mm x 10,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 10.0 mm self-tapping Bone Fixation Screws	
PFB12	1	1,5 mm x 12,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 12.0 mm self-tapping Bone Fixation Screw	
PFB12-5	5	1,5 mm x 12,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 12.0 mm self-tapping Bone Fixation Screws	
PFB14	1	1,5 mm x 14,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 14.0 mm self-tapping Bone Fixation Screw	
PFB14-5	5	1,5 mm x 14,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 14.0 mm self-tapping Bone Fixation Screws	

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### Synthetic bone grafting material

NanoBone® contains absolutely no human or animal ingredients so no material-related contamination risks and alleviates ethical concerns for patient consent. High-tech production processes ensure consistent and excellent quality.

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### Complete remodelling

Thanks to its special structure, NanoBone® can be used for natural remodelling. Osteoclasts resorb the material as osteoblasts form new bone. Material resorption and bone formation are linked. This means predictable results and a reliable bone basis.

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### Controlled osteoinduction

NanoBone® is osteoconductive and osteoinductive. In other words, it is not simply a scaffold: it actively promotes bone formation.

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**NanoBone® | granulate** is available in either fine or coarse granules to suit the indication.

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The granulate can be mixed with the patient's blood or with a sterile saline solution. The manufacturer recommends mixing with blood. The material is extremely hydrophilic and fully absorbs the blood or saline solution quickly.

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When mixed with blood, **NanoBone® | granulate** takes on a paste-like consistency and can be easily applied with a spatula or augmentation spoon.

## NanoBone® | block

Synthetic Alternative to autologous Bone Blocks

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**NanoBone® | block** was developed in partnership with users as an alternative to the autologous bone block. The block comes in a set with two screws for attachment.

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# NanoBone® | QD

Ideal bone healing in combination with perfect handling

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
## AT A GLANCE

### The applicator

- "Ready to use" direct application from the applicator
- Easy, intuitive use
- Precise placement of the **NanoBone® QD** material
- Controlled and exact portioning
- Allows for small approaches to the application site (outer diameter of the applicator just 8 mm)
- Allows for retrograde filling of bone defects
- Difficult to access bone defects treated due to the long applicator

### The material

- Shapeable and adaptable to the defect
- Strong cohesion of the material
- Sticks to the bone (adhesion)
- Stable also in aqueous environment and profusely bleeding wounds
- Controlled osteoinduction

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# Easy-Clean Scissors NEW

## Surgical Scissors for easy and safe processing

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The MEISINGER Easy-Clean Scissors made of high-quality stainless steel can be easily and completely disassembled and put back together again. Residues below the hinge can be completely removed – for a residue-free processing\*. For dismantling, the scissors only need to be opened at a right angle. Then the upper part of the scissors is lifted briefly. In addition, the scissors have a black coating that is applied using PVD procedure (vacuum-based coating). Thanks to this black coating the products reach high hardness and durability. Furthermore, disturbing light reflections are reduced, and glare-free work is possible.

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## AT A GLANCE

- Surgical scissors made of high-quality stainless steel
- Easy and complete dismantling for residue-free processing\*
- Black coating for high hardness and durability as well as glare-free work
- Scissors in two different versions

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This ratchet was designed to transfer turning moments up to 70 Ncm. May be disassembled for simple and safe cleaning.

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
Torque ratchet for working with defined torque. In blocking position, torques for both the insertion and the loosening of screws can be transmitted, e.g. for the insertion of implants.

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
Ratchet for working in torque function and in blocked function. In torque function the application takes place using the bending rod, from which torques between 10 and 70 Ncm can be read. In blocked function, the ratchet enables a greater torque to be transmitted, for example when inserting implants or loosening connections.

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MEISINGER offers one of the most extensive programs in the area of the Surgical Round Drills, Surgical Cutters, and Trephines. The optimized special construction of the instruments guarantee a safe application and fast achievement for the desired working result.

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Surgical Cutters of Tungsten Carbide

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Surgical round drill of tungsten carbide

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
Surgical round drill of stainless Steel

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



\* im Saw Basic Kit 7150 enthalten \* contained in the Saw Basic Kit 7150

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



Innengekühlter Trepan-Schaft | Internally cooled trepan shank

### 229KH\*

Schaft / Shank	Fig.	229KH	229KH	229KH	229KH
unmontiert / unmounted		021	025	029	033
L mm		15,5	15,5	15,5	15,5
		3,1	3,5	3,9	4,3
		2,1	2,5	2,9	3,3
		1	1	1	1


Innengekühlte Trepan-Arbeitsteile | Internally cooled trepan work parts

! Technology by Dr. Bernd Giesenhagen

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
The MEISINGER gingiva trimmers are intended for use in dental mucosal surgery. They have been specially developed for the versatile and gentle treatment of gingival tissue. The working part consists of particularly high-quality and stable zirconia. The main application focus of the MEISINGER gingiva trimmers is the safe and gentle removal of unwanted granulation tissue and hyperplastic gingiva (papillectomy) - even in hard-to-reach areas. Furthermore, the trimmers are used in the enlargement of the sulcus after crown preparation, for the best possible exposition of the preparation margin and the creation of an exact impression. The gingiva trimmers are also particularly suitable for the tissue-saving exposure of covered healed implants and for use in periodontal treatment. The application takes place without cooling, so that the resulting rotational energy can be used for tissue modulation. Thanks to heat coagulation, bleeding tendency is reduced.

#### Instructions for use:

- Application without cooling
- Optimal speed: 300.000 - 450.000 rpm

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Stainless Steel Bur Blocks, non-corrosive, sterilizable

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MEISINGER drill stands are ideally suited for the effective cleaning, sterilization and organization of the required instruments.

Art.-No. BS441      B 60 x H 40 x T 45 mm  
Für 8 Winkelstück oder 8 FG-Instrumente  
For 8 RA instruments or 8 FG instruments

Art.-No. BS541      B 120 x H 40 x T 45 mm  
Für 18 Winkelstück oder 18 FG-Instrumente  
For 18 RA instruments or 18 FG instruments

Art.-No. BS440      B 60 x H 40 x T 45 mm  
Für 12 Winkelstück oder 12 FG-Instrumente  
For 12 RA instruments or 12 FG instruments

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


**FG:**

Optimal: 80.000 min<sup>-1</sup>  
 Maximal: 100.000 - 120.000 min<sup>-1</sup>

The following reference values for rotation speeds apply to surgery in general:

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Please also observe adapting the rotation speed in relation to the diameter of the instrument as well as the prevailing indication and enough cooling. As a general rule, the larger the working part of an instrument, the lower the speed should be set.

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Please also pay attention to the special preparation specifications for products made of tool steel.

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Please also note the recommended and maximum rotary speeds for the individual instruments on the product packaging and in the corresponding instructions for use.

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The symbols give merely suggestions for the possible implementation of the products. The user decides and takes full responsibility about the precise deployment according to existing indications. Please follow general application and safety instructions for MEISINGER products in the medical and dental area and also the advice for processing. Details can be found on the internet under [www.meisinger.de](http://www.meisinger.de) or you can request one by mail.

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Cleaning/Disinfection/Sterilization

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3723

Contains hazardous substances: Cobalt, CAS: 7440-48-4

This product contains cobalt in more than 0.1 percent by mass and is therefore, in accordance with current regulations, to be labeled as CMR material class 1B (carcinogenic, mutagenic and/or reprotoxic). It has been proven that there is no increased risk of cancer or adverse reproductive or genetic effects when the product is used according to the directions for its intended use.

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Meisinger stands for high quality medical devices since 1888. The quality management system of a company which manufactures medical devices must meet specific special requirements. These extremely high requirements are defined in ISO 13485 and meticulously complied with by our company. A MDSAP certificate according to ISO 13485:2016 confirms compliance with the requirements of international authorities in the USA (FDA), Canada (Health Canada), Australia (TGR), Japan (MHLW) and Brazil (ANVISA). All medical devices which you purchase from us as customer, comply with all applicable requirements of the Medical Device Directive 93/42/EEC. Our company is certified by an independent Notified Body and certification is performed according to the specifications of standards. Current certificates can be found on our homepage [www.meisinger.de](http://www.meisinger.de)

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All rights reserved.

Reproduction, also by extract and reproducing photos, are only permitted with written authorization of Hager & Meisinger GmbH, Neuss.

All offers, orders, and deliveries are subject to Meisinger's "General Sales and Delivery Terms." We reserve the right to modify our range of products and their design as well as to deviate from the illustrations and data shown.

registered trademarks

Hager & Meisinger GmbH | Hansemannstr. 10 | 41468 Neuss | Germany  
Phone: +49 2131 2012-0 | Fax: +49 2131 2012-222 | E-Mail: [info@meisinger.de](mailto:info@meisinger.de) | Internet: [www.meisinger.de](http://www.meisinger.de)

Meisinger USA. L.L.C. | 10150 E. Easter Avenue | Centennial, Colorado 80112 | USA  
Phone: +1 (303) 268-5400 | Toll free: +1 (866) 634-7464 | Fax: +1 (303) 268-5407 | E-Mail: [info@meisingerusa.com](mailto:info@meisingerusa.com) | Internet: [www.meisingerusa.com](http://www.meisingerusa.com)

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Cytoplast™ and Pro-Fix™ are registered trademarks of Osteogenics Biomedical, Inc..
NanoBone® ist ein eingetragenes Warenzeichen der ARTOSS GmbH. NanoBone® is a registered trademark of ARTOSS GmbH.
Cytoplast™ und Pro-Fix™ sind eingetragene Warenzeichen der Osteogenics Biomedical, Inc..
ARTOSS GmbH