

VITASSAY

H. pylori -Turbidimetric Assay-

Rapid test for the quantitative detection of Helicobacter pylori in human stool samples.

IUE-7115003 Ed00 April 2019



For professional *in vitro* diagnostic use only. Professional trained in Turbidimetric techniques.

INTENDED USE

Vitassay H. pylori –Turbidimetric Assay- is a rapid Turbidimetric assay for the quantitative detection of Helicobacter pylori (H. pylori) in human stool samples.

Simple, non-invasive and highly sensitive assay that helps in the diagnosis of Helicobacter pylori infection. This product has been optimized for several automated analyzer.

INTRODUCTION

Helicobacter pylori (H. pylori) is a spiral gram-negative, microaerobic human pathogen. H. pylori infection is strongly related with many gastroduodenal diseases, atrophic gastritis, mucosa associated lymphoid tissue (MALT) lymphoma and noncardia gastric cancer.

H. pylori colonize approximately 50% of world's population, but the prevalence of H. pylori is of high quality in developing countries. Risk factors for H. pylori infection varies widely by geographic area, age, race, and socioeconomic status.

PRINCIPLE

Vitassay H. pylori –Turbidimetric Assay- is a quantitative turbidimetric assay for the detection of Helicobacter pylori in human solid stool samples.

H. pylori Turbidimetric Assay is based on antigen-antibody agglutination reactions between the antigen contained in the sample and the antibodies anti-antigen coated on polystyrene latex particles. Such agglutination is measured as an increase in absorbance proportional to the quantity of antigen contained in the sample.

The use of two external controls, Control 1 and Control 2, is used to verify that the test is working properly.

PRECAUTIONS

- For professional *in vitro* use only.
- Do not use after expiration date.
- Do not use the test if its primary containers are damaged.
- Specimens should be considered as potentially hazardous and handle in the same manner as an infectious agent. Avoid contamination errors, follow proper work procedure.
- The reagents after use should be discarded in a proper biohazard container after testing.

- Reagents contain preservatives. Avoid any contact with the skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use with the **Vitassay H. pylori -Turbidimetric Assay-**. Do not use any other commercial component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not use any other commercial kit component.
- If measure range is exceeding, use the sample diluent to dilute the sample and repeat the assay again.
- Read and follow up the instructions for use provided in the kit.
- Prepare and adjust the analyzer before starting measurements.

STORAGE AND STABILITY

Store as packaged in the original primary container, the reagents should be preserved at refrigerated temperature (2-8°C/35.6-46.4°F), the sample diluent could be preserved refrigerated or at room temperature (2-30°C/35.6-86°F).

The product is stable until the expiration date printed on the label, if they have been preserved under the recommended conditions.

Do not freeze and keep away from the sunlight.

MATERIALS

MATERIAL PROVIDED	MATERIAL REQUIRED BUT NOT PROVIDED
<ul style="list-style-type: none">▪ Reagent R1▪ Reagent R2▪ Calibrator 0▪ Calibrator 1▪ Calibrator 2▪ Calibrator 3▪ Calibrator 4▪ Calibrator 5▪ Control 1▪ Control 2▪ Vials with diluent for the sample dilution.▪ Instruction for use.	<ul style="list-style-type: none">▪ Specimen collection container.▪ Disposable gloves.▪ Automated analyzer.▪ Vortex.▪ Microtubes (analyser vial).

SPECIMEN COLLECTION




Collect sufficient quantity of feces: 1-2 g or mL for liquid samples. Stool should be collected in clean and dry containers.

Samples can be stored in the refrigerator (2-8°C/36.5-46.4°F) for 1-2 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen at -20°C (-4°F). Samples must be brought to room temperature before testing. Homogenize the sample as thoroughly as possible prior to preparation.

SPECIMEN PREPARATION

1. Label the vial with the patients' identification.
2. Open the cap of the vial with diluent for sample dilution (figure 1).
3. Use the stick to collect sufficient sample quantity. For solid stool, insert the stick in 4 different areas of the stool sample (figure 2), and add it into the vial with diluent for sample dilution. For liquid stool, take approx. 20 µL of the sample using a micropipette and transfer it into the vial for the sample dilution.
4. Close the cap of the vial with the diluent and sample, shake the vial in order to assure good sample dispersion, using a vortex during 1 minute (figure 3). The sample dilution vial with dilute sample can be stored in a range of temperatures (2-8°C / 35.6-46.4°F) for 7 days in the refrigerator prior the testing.
5. Take the specimen collection vial with the sample diluted, cut the end of the cap and dispense 20 drops of the sample diluted (figure 4) into a analyser vial (microtube). This vial must be compatible with the analyzer.

Note: Do not use the sample vials directly in the analyzer.



Put the sample into the vial, close the vial and shake for a good dispersion of the samples (vortex until the sample is fully resuspended max.1 minute).

Dispense 20 drops into a analyser vial (microtube).

PROCEDURE

Allow reagents and solid stool to reach room temperature (15-30°C) prior the testing. Reagent R1 y Reagent R2 are ready to use.

Calibration curve

For calibration only use H. pylori Calibrator: Cal0, Cal1, Cal2, Cal3, Cal4 and Cal5. Contain Helicobacter pylori antigen at diferent concentrations indicated on the label of each of the vials.

	Reference	Calibrator 1	Calibrator 2	Calibrator 3	Calibrator 4	Calibrator 5
Conc.	0 ng/mL	2.5 ng/mL	5 ng/mL	10 ng/mL	20 ng/mL	40 ng/mL
Vol.	1000 µL	1000 µL	1000 µL	1000 µL	1000 µL	1000 µL

Reagents are ready to use. The frequency in the realization of the calibration curve must be established by the end user in the function of the criteria established for the clinical laboratory.

Note: See section quality control.

Analytical procedure

Measure range: 1 – 40 ng/mL.

Procedure	Steps	
R1 addition	305 µL	0 s
Sample addition	15 µL	10 s
R2 addition	25 µL	300 s
Blank measure	450 nm – 800 nm	310 s
Mainly measure	450 nm – 800 nm	610 s

*Data obtained by the analyzer Biolis i24 (Tokio Boeki)

INTERPRETATION OF RESULTS

Negative results values: Lower to the limit of detection of 0.5 ng of H. pylori antigen/mL.

Negative results determine the absence of H. pylori antigen in human stool samples.

Grey area results: Between 0.5 and 1 ng H. pylori antigen/mL.

Result between 0.5 and 1 ng/mL are between the detection and quantification limit. A second sampling is recommended. In case, the second result also shows a value in this range, the recommendation is to follow up the patient some time later.

Positive results values: Higher or equal to the limit of quantification of 1 ng H. pylori antigen/mL.

Positive results determine the presence of H. pylori antigen in human stool samples.

QUALITY CONTROL

H. pylori C1 & C2 Controls are ready to use. Allow controls to reach room temperature (15-30°C) prior to testing.

H. pylori Control 1: is liquid control at a certain concentration of Helicobacter pylori (lower than Control 2). Concentration is indicated on the label of the vial.

H. pylori Control 2: is liquid control at a certain concentration of Helicobacter pylori (higher than Control 1). Concentration is indicated on the label of the vial.

The use of controls at two different concentrations is recommended to verify test precision.

If the obtain results are out of the tolerance range, the equipment, the reagents or the technique must be reviewed. If the problems persists, stop using the reagents and contact your distributor.

LIMITATIONS

- **Vitassay H. pylori -Turbidimetric Assay** should be only used for the detection of Helicobacter pylori in human stool samples.
- The quality of **Vitassay H. pylori -Turbidimetric Assay** depends on the quality of the sample; Proper fecal specimens must be obtained.
- Positive results determine the presence of Helicobacter pylori in fecal samples; a positive result should be followed up with additional invasive techniques (endoscopy) to confirm the results.

- Positive results can be obtained even when the symptoms had disappeared.
- Negative results should not be considered as conclusive, it is possible that the concentration of antigens is lower than the detection limit. If symptoms or situation still persist, a *Helicobacter pylori* determination should be carried out, on a sample from an enrichment culture or using an invasive technique.

EXPECTED VALUES

H. pylori infection affects more than half of the adult population worldwide and is strongly related with many peptic ulcer diseases and gastric cancer. The prevalence of *H. pylori* infection varies widely by geographic area, age, race and socio-economics status. Usually, the prevalence of *H. pylori* increases with age in most countries, being able to cause more than 80% of peptic ulcer diseases and approximately 75% of gastric cancer. Peptic ulcer and gastric cancer together cause more than a million deaths per year in the world.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity:

Detection limit: 0.5 ng/mL.

Prozone:

Lower concentrations of 0.2 mg/mL of stool do not show prozone effect and no false negative results have been observed.

Within-Run Precision:

	Low (1 ng/mL)	Media (10 ng/mL)	High (40 ng/mL)
N	20	20	20
Media (µg/g)	1.08	10.23	39.76
DS (µg/g)	0.12	0.79	2.01
CV (%)	11	8	5

*Data obtained by the analyzer Biolis i24 (Tokio Boeki)

Cross reactivity:








No cross reactivity was detected against other fecal markers that are occasionally present in feces:

<i>Campylobacter Coli</i>	<i>Salmonella enteritidis</i>	<i>Shigella boydii</i>
<i>Campylobacter jejuni</i>	<i>Salmonella paratyphi</i>	<i>Shigella dysenteriae</i>
<i>C. difficile Toxin B</i>	<i>Salmonella Typhi</i>	<i>Shigella flexneri</i>
<i>Escherichia coli O157</i>	<i>Salmonella typhimurium</i>	<i>Yersinia enterocolitica O:3</i>
<i>Listeria monocytogenes</i>	<i>Salmonella Typhi</i>	<i>Yersinia enterocolitica O:9</i>

REFERENCES

1. Urrego J., Otero W., Gómez M. A review of *Helicobacter Pylori* and Colon Cancer. *Rev Col Gastroenterol/29* (3) 2014.
2. Wang Y-K, Kuo F-C, Liu C-J., Wu M-C., Shih H-Y., Wang S-S., Wu J-Y., Kuo C-H., Huang Y-K., Wu D-C. Diagnosis of *Helicobacter pylori* infection: Current options and development . *World. J. Gastroenterol.* 2015 October 28; 21(40):11221-11235.
3. Lee Y-C., Chiu H-M., Chiang T-H. Accuracy of faecal occult blood test and *Helicobacter pylori* stool antigen test for detection of upper gastrointestinal lesions. *BMJ Open* 2013; 3:e003989.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

IVD	<i>in vitro</i> diagnostic device		Keep dry
	Consult instructions for use		Temperature limitation
	Use by		Manufacturer
LOT	Batch code		Contains sufficient for <n> test
DIL	Sample diluent	REF	Catalogue number
	Keep out of the sunlight		

ADAPTED EQUIPMENT

- Biolis i24 (Tokio Boeki)
- BS200 (Mindray)
- ChemwellT (Awareness)



Vitassay
For information purposes only

