

FIDA Technology

In-solution technology for analysis
of complex molecular interactions

Now with in-solution kinetics

No restrictions on sample matrix or buffer.
Absolute measurements.
Nanoliter sample volumes.

One technology to measure:

- Affinity (K_D)
- Kinetics (k_{on} & k_{off})
- Quantity
- Sample Quality Control



FIDA IN A NUTSHELL

01 FIRST PRINCIPLE THINKING

FIDA technology is a “1st Principle” technology.

This means that FIDA does not depend on a priori assumption or on empirical calibration. It uses first principles of physics and fluid mechanics to analyse the movement of particles in a fluid. This brings simplicity and robustness straight into the users' lab.




Independently of the biology being investigated, each data point has a range of built-in QC parameters included. Thanks to that, data interpretation is straightforward, and R&D iterations can be performed instantly, which speeds up users' workflows.

02 HOW DOES IT WORK? SIMPLIFIED.

FIDA measures fluorescence of particles in the laminar flow and analyses their dispersion over time, which allows for calculation of the hydrodynamic radius of a particle of interest. The two basic principles used are Taylor Dispersion and Laminar Flow.

The sample of interest is passed through a thin capillary. Due to the difference in velocity between the walls and centre of the capillary, the sample shapes into a parabolic profile. Molecules diffuse radially, away from the flow axis. The fluorescence emitted by the molecules is acquired as a Gaussian signal by a high sensitivity detection system and is plotted against time. The size of the molecules in the sample determines their radial diffusivity, which in turn defines the extent of sample's dispersion.



$$\text{Diffusivity} = \frac{a^2}{24 \sigma^2} t_R$$

$$\text{Hydrodynamic Radius} = \frac{k_b T}{6 \pi n D}$$

FIDA can detect size changes smaller than 5%.

Scan to see
how it works!




03 IN-SOLUTION KINETICS




FIDA KINETICS TAKE ADVANTAGE OF THE TIME DEPENDENT ASPECT OF DIFFUSIVITY - THE VARIABILITY OF THE PRESSURE DRIVEN FLOW

THIS IS HOW IT LOOKS LIKE IN PRINCIPLE:


MIXING PRINCIPLES INSIDE OF THE CAPILLARY:



EQUILIBRIUM BINDING CURVE



SHIFTED KINETICS CURVE




The figure above presents equilibrium binding curves and kinetic binding curves. The top figure describes the mixing principles inside the capillary while the bottom figure describes the equilibrium binding curves and the shifted kinetics curve. The samples already prepared for the equilibrium affinity determination can be reused to measure the kinetics binding curve, minimising sample consumption.

FIDA Measurements



Size, Hydrodynamic radius (R_h)

- First Principle measurement of hydrodynamic radius in nanometers.
- Range: 0,5 nm – 500 nm Rh.
- Can detect size changes smaller than 5%.



In-solution kinetics

- k_{on} & k_{off}
- No non-specific binding issues
- No steric hindrance to high density immobilised ligands
- Seconds to minutes



Fluorescence detection

- Multiple wavelengths available.
- Labelled assays with: 480 nm, 640 nm or UV detector.



Binding Related Intensity Change (BRIC)

- In parallel with size measurement, an orthogonal BRIC signal is obtained for all measurements.



Sample Quality Control Module & Reporting Tool

- 8 Quality Control Parameters for each sample
- Custom-made reports that fit your workflow
- Export data in multiple formats



NO ENVIRONMENTAL RESTRICTIONS

Seamlessly operate in complex matrices including fermentation media, plasma or serum. With FIDA there is no need for purification, which allows You to save your sample material and time.

Robustness & Versatility



Speed & efficiency

- **Minimal sample consumption**
 - nL – uL range.
- **Walk-away automation**
 - 2x 50 vials or 2x 96 well plates.
- Up to 24 data points in <30 min.
- Up to 8 binding curves with 8 data points in triplicates, each within 4 hours of unattended operations.



Easiness of use & versatility

- Simple software interface.
- Data delivered as tables, plots, and real-time monitoring of the signal.
- Purified or un-purified samples.
- Built-in assay control.
- No limits on buffers, detergents or ionic strength.



Measurement specifications

- **Dissociation constant :**
 - pM to nM
- **Size detection**
 - Hydrodynamic radius between 0.5-500nm
- Three interchangeable detectors: 480, 640 & UV
- Detection signal-to-noise ratio of >30 (3-fold increase compared to current state-of-art detectors)



Temperature control

- **Measurement chamber:**
 - 15-45°C
 - 59-113°F
- **Autosampler**
 - 5-50°C
 - 41-122°F
- Can be controlled individually in each of the two tray holders and in the capillary compartment.

**WORK WITH
ANY BUFFER COMPOSITION
THAT FITS YOUR GOALS**

No restrictions related to detergent usage, ionic strength, temperature or sample pH.

This eliminates common research constraints, minimises assay development time and expand the scope of biological systems You can characterise.

HOW CAN YOU BENEFIT FROM IN-SOLUTION KINETICS

06

No environmental restrictions



Seamlessly operate in **complex matrices including fermentation media, plasma or serum** - no purification required.

Avoid non-specific binding




No steric hindrance to high density immobilised ligands
No non-specific binding issues
No risk of re-binding

No restrictions on detergents, ionic strengths, temperature, pH etc.



Minimise assay development time
Expand the scope of biological systems you can characterise
Increase environmental relevance

No need for regeneration



With FIDA there is no surface chemistry involved.
Eliminate the risk of denaturing immobilised protein
Rapidly determine slow off rates for high affinity interactions

Detect Strong & Weak Binders



FIDA is capable of measuring kinetics of both strong and weak interactions in-solution.

Retrospectively, it is evident how many great discoveries in science were heavily dependent on the development of new technologies - do not let technological constraints keep you away from your next discovery.

Fida Neo - No restrictions.

07

MULTIPLE SIMULTANEOUS & INSTANT READOUTS

8 Quality Control Parameters

Fida Software automatically includes Quality Control to each measurement taken. This increases the level of transparency and supports troubleshooting and assay optimisation. Our QC Module allows for building and exporting custom-made reports, which facilitates smooth workflows and data storage.




Structural integrity

- Size measured as hydrodynamic radius (Rh).
- Validate your protein stability
- Get insight into folding/unfolding and conformational changes.



Functionality/Binding

- Automated binding curves and equilibrium Kd's are obtained by loading the autosampler with your titrations.



Aggregation

- Protein/particle aggregates are clearly detectable and quantifiable whilst still leaving the core signal useful for standard measurement.



Heterogeneity (PDI)

- PDI Index allows for checking the heterogeneity of your sample.



Labelling efficiency

- Option of measuring size of up to 3 species in solution.
- Can e.g. reveal the percentage of free vs. conjugated fluorophore in your sample when you choose to use Fida 1 for labelled assays.



Stickiness

- The shape of the core signal will reveal any stickiness of your binding partners or your binding complexes.
- The core signal is useful for standard measurement despite of the stickiness.



Viscosity

- Every measurement you take provides viscosity data.
- Viscosity compensation



Sample Loss

- Transparently exposed
- Troubleshoot efficiently



PDB Correlator


- Use the absolute size as a firm reference point.
- Compatible with Protein Data Bank, Pymol or AlphaFold.



Sample Quality Control Module

CUSTOMISE & EXPORT REPORTS

Reports that meet your requirements



Easy to implement in your workflow.

Configure report

Report Structure

- Multiple Species summary table
- PDI summary table
- Generic summary table
- Datafile per page block
- Configuration appendix

Datafile per page block

- Graphics
- Generic table
- Multiple Species table
- PDI table


Report format

- PDB table
- Note
- PDF
- .txt
- .bmp graphics

Summary tables configuration

| Multiple Species | PDI | Generic |
|-------------------------------------|-------------------------------------|---|
| <input type="checkbox"/> Rh | <input type="checkbox"/> Rh | <input type="checkbox"/> Taylorgram S/N |
| <input type="checkbox"/> Area | <input type="checkbox"/> PDI | <input type="checkbox"/> Spike counter |
| <input type="checkbox"/> R2 | <input type="checkbox"/> R2 | <input type="checkbox"/> Residence time |
| <input type="checkbox"/> Statistics | <input type="checkbox"/> Statistics | <input type="checkbox"/> Viscosity |
| | | <input type="checkbox"/> Statistics |

With Fidabio Quality Control Module you can custom make and export Quality Control reports of your samples. The data can be exported as a PDF report file with graphs included, or a .txt file, which is easily processed by any data analysis software.



08

APPLICATIONS LANDSCAPE

In-solution validation of structural information.

PDB correlator

The Fida 1 readout (Rh) provides an absolute size measure (in nanometres), which is directly linked to protein structure and function.

Fida 1 software, includes a PDB correlator, which directly compares measurements obtained in solution or under native conditions (Rh) with structures obtained from X-ray, Cryo-EM, NMR, or AlphaFold.

COMPLEX INTERACTIONS

Bi-/Multi-specifics

- Assessment of **multiple Kd's** in complex binding events.
- Option of **deconvoluting multi-complex** binding events.
- Cooperativity factor.

Stoichiometry/ Oligomeric state

- **Absolute** read-out in **nanometres** facilitates reliable insights into stoichiometry of binding events.

Protacs/Molecular glues

- Ternary complex formation **in buffer and in cell lysate**.
- Assessment of **individual Kd's** in the ternary complex.
- Assessment of “**fraction bound**” of your protein
- **Cooperativity** factor.
- **Ubiquitination** assay.


STRUCTURAL BIOLOGY

- Structural biology.
- Cryo-EM and x-ray sample QC **at high concentrations**.
- PDB correlator.
- Nanobody binding.

MEMBRANE PROTEINS

Characterization & screening

- Detergent screening uses minimal sample amount in a rapid automated assay.
- Binding assays with unpurified membrane preps.
- Particularly for crude membrane preps, it is significant that Fida 1 enables you to keep your samples at 5°C in the autosampler whilst performing your actual measurement at room-temperature.



QUANTIFICATION

- **Biomarkers** in plasma/serum/cell lysate.
- **Auto-antibodies** in plasma/serum.
- **Bioprocessing** - Therapeutic proteins in fermentation media.
- **Nanoparticles** in plasma/serum/fermentation media.

PARTICLES

- Particle size.
- Polydispersity index.
- Aggregation.
- Binding characterization.

CELL & GENE THERAPY

- gRNA-nuclease complex formation.
- gRNA-nuclease **K_d** determination.
- **Shield protein** characterisation on drug delivery vesicles.

CLONE SELECTION

- Clone selection using FIDA is based on **both titre and affinity**.
- Both parameters can be obtained **simultaneously**, without purification of the cell supernatant/cell lysate.



AMYLOID FIBRILS

- In-depth fibril aggregation kinetics
- Thermodynamic characterisation of **amyloid polymorphism**
- Robust sizing and binding

LIQUID-LIQUID PHASE SEPARATION

- Phase diagrams.
- Dilute-phase concentration.
- K_d of LLPS modulating components.
- Relative droplet size distribution.
- Droplet to amyloid transition.

**Droplet count
& relative size**



09

TECHNICAL SPECIFICATIONS

and instrument characteristics

| | |
|---|---|
| Detection technology | Fluorescence - multiple wavelengths available: UV (label free), 480, 640 |
| Size accuracy | 5% |
| Kinetics | sec-hrs |
| Dissociation constant (K_d): | pM - mM |
| Size detection | Rh of 0.5 - 500nm |
| Signal-to-noise ratio | > 30 |
| Assay control | Built-in Quality Control parameters |
| Sample capacity per run | Up to 2 x 96 samples |
| Pressure range | 1 - 3500 mBar |
| Autosampler temperature control | 5°-50°C (41°-122°F) |
| Capillary chamber temperature control | 15°-45°C (59°-113°F) |
| Capillary types | Fused silica; dynamic coatings or permanently coated |
| Power | 120-240VAC, 50/60Hz |
| Operating system | Windows |



(+45) 53 72 78 70
info@fidabio.com
Generatorvej 6A
2860 Søborg
Denmark

IT'S NOT JUST AN INSTRUMENT

Become a part of FIDA users community

FIDA users receive in-person training, ad-hoc support, and advice, right to participate in EU and US user group meetings and engagement on collaborative projects, publications and conference attendance.

Book a discovery call

