

Laboratory

Water Distillation

DS 2000

DS 4000

DS 8000


DS 8008

DS 8012

DS 8025

 DragLab

Connect & Go



Laboratory Water Still DS 2000

A professional, fully automatic water distiller, the DragLab water still is constructed from high-quality stainless steel. It is designed to produce 2 liters per hour of high-purity distillate with minimal personnel supervision.

Effortless Water Purity With DragLab Distillation Solutions

DragLab Distillation Units are the most effective and reliable devices for producing pure water, requiring only a connection to raw water and a power source. This process involves transforming raw water (tap water) into steam and then back into liquid, effectively separating water from impurities with higher boiling points. The resulting steam, once condensed, produces distilled water with a purity of approximately 99.5%, removing salts, organic substances, microorganisms, pyrogens, and bacteria. DragLab water distillation units require minimal maintenance needing only regular cleaning. When carbon dioxide from the air dissolves in distilled water, its pH becomes slightly acidic, 5 to 5.6. Distilled water should be used promptly or stored in airtight conditions to maintain purity.



Application:

- Biotechnology, cell culture and media preparation
- Testing, analysis and research laboratories
- Microbiology and Bacteriology sample preparation
- Food and beverage production
- Sterilization, cleaning and Medical equipment.
- Cosmetics, skincare products production.
- Electronics and Automotive.
- Pharmaceutical manufacturing and drug formulation.

A vertical bar on the left side of the page, divided into three colored segments: black at the top, red in the middle, and yellow at the bottom.

Quality and Expertise



DragLab is certified to ISO 9001:2015, reflecting our unwavering commitment to quality and excellence. We continuously optimize our production processes, ensuring that every step is efficient and effective. Through ongoing development and innovation, we strive to meet and exceed industry standards. Our customer-centric approach prioritizes the needs and satisfaction of our clients, driving us to maintain rigorous quality assurance protocols. This dedication ensures that laboratory equipment manufactured by DragLab consistently meets the highest standards of quality and reliability.

DragLab Laboratory Devices are meticulously designed and manufactured from high-quality materials to meet the most rigorous requirements of our customers. Our commitment to excellence and innovation drives continuous development and improvement, informed by daily usage and ongoing technological advancements. 'Quality and Expertise' is not just about creating impressive product designs; it reflects our company's policy of prioritizing high preparedness and customer-centric service as our fundamental objective. This dedication extends to all DragLab products, which comply with current European standards and bear the CE mark.



Laboratory Water Still DS 2000

A professional, fully automatic water distiller, the DragLab water still is constructed from high-quality stainless steel. It is designed to produce 2 liters per hour of high-purity distillate with minimal personnel supervision.



Features:

Distillation Capacity:

Produces 2 L/h of distilled water.

Purity:

High-purity distilled water with conductivity of approximately 2.5 $\mu\text{S}/\text{cm}$ at 25 °C.

Water Quality:

Pyrogen-free, low-gas, and bacteria-free water.

Protection:

Automatic thermostatic low-water cut-off to protect the heating element in case of low water.

Operation:

Automatic water switching for continuous functionality.

Maintenance:

Easily accessed evaporator tank for effortless cleaning and maintenance.

Material:

All water contact materials are made of stainless steel (material no.: EN 1.4301).

Heating Element:

Made of stainless steel (material no.: EN 1.4301).

Housing:

Galvanized steel with electrostatic powder coating.


Installation:

Suitable for both bench and wall mounting.

Accessories:

All parts and accessories for installation included.

Dimensions:



Laboratory Water Still DS 2000

Specifications:

Capacity

Capacity	2 litter / hour
Consumption	40 litter / hour
Conductivity	Approximately 2,5 µS/cm at 25 C

Control technology

Fully Automatic	Automatic thermostatic power cut-off to protect the heating element in case of low water
Fully Automatic	Automatic power activation when the water level reaches the safety line in the evaporator tank to protect the heating element.
Water supply	Built-in solenoid valve connection for the evaporator and cooler, with automatic switching on/off based on the water level in the evaporator tank.

Standard equipment

Internal material	The tank is made of stainless steel (material no. EN 1.4301)
Internal material	The condenser (cooler) is made of stainless steel (material no. EN 1.4876).
Internal material	The heating element is made of stainless steel (material no. EN 1.4876).
External material	The housing is made from galvanized steel and coated with electrostatic powder.
Degassing system	The unit includes a top vent for effective degassing, ensuring pyrogen-free, low-gas, and bacteria-free water.
Accessing interior	The evaporator is easily accessible for cleaning and maintenance, facilitating efficient cleaning procedures.
Installation	All necessary parts and accessories for installation are included.
Distillate water withdrawal	Distillate water can be continuously withdrawn by opening the drain tap.
Main switch	The main power switch button located on the front of the unit
Led light	Two LED lights on the front: a green light indicating standby operation and a red-light indicating faults.
Mounting	Suitable for both benchtop and wall mounting.

Safety

Water	Low water level sensor with automatic power cut-off to protect the heating element.
Heating	Overheating cut-off
Led light	Red LED alarm indicator for faults & errors
Degassing	pyrogen-free, low-gas, and bacteria-free water.
Power connection	Power cable with German shock-proof type (Schuko)) plug

Laboratory Water Still DS 2000

Specifications:

Connection

Water inlet	The cooling water inlet is located on the right-hand side of the unit, with a diameter of 3/4 inch (Ø 19.05 mm).
Water pressure range	From 2 bars to 10 bars
Water outlet	The cooling water outlet is located on the right-hand side of the unit, with a diameter of 2/5 inch (inner Ø 10 mm).
Drainage tap	The drainage tap is located on the right-hand side of the unit, with a diameter of 2/5 inch (inner Ø 10 mm).
Distillate water	Distillate water can be continuously withdrawn by opening the drain tap.
Power	The main connection cable on the left-hand side of the unit

Electrical data

Voltage	230 volts, 50/60Hz
Watt	2000 watts

Shipping information

Customs tariff number	8419 4090
Packaging Dimension	53 W X 50 D X 75 H (CM)
Net weight	22 Kg
Gross weight carton	27 Kg

Order Information:

Description	Model	Article number
Water still 2L/h, 220V/50-60Hz	DS 2000	5022.000

