

CheckDigit™ | Handwritten

Designed to work with the Cerner™ Bridge Transfusion Administration product, the CheckDigit™ | Handwritten Blood Band helps ensure the right blood goes to the right patient.



1. **VERIFY** the patient's identity.

WRITE the patient's information over the prompts on the band using a ballpoint pen on a flat, hard surface.

NOTE: Do not apply the shield over handwritten information at this time.



2. **WRAP** the band tail around the patient's extremity to size.

SNAP the button closure to secure.

CUT the excess band tail material with scissors.



3. **DRAW** the pre-transfusion specimen per your protocol.

PEEL the handwritten patient information label from band.

PLACE the handwritten patient information label on the specimen tube lengthwise.



4. **PEEL** the liner from the shield.

COVER the copy of the handwritten patient information with the shield.

SEAL the edges down using your fingers.



5. **LABEL** any forms at the patient's bedside as required by your protocol.


REMOVE the tail of stickers from the band at the perforation.



6. **PEEL** the adhesive tab on the underside of extra stickers.

SECURE the extra stickers to the band.

REMOVE the extra stickers with scissor if not required by protocol.



7. **ATTACH** the tail of stickers to the specimen tube, if desired.

TRANSPORT the specimen to the blood bank.

SCAN the Blood Bank ID (BBID) into the Laboratory Information System (LIS) or enter it manually.



8. **VERIFY** the patient's identity per your facility's protocol.

SCAN the permanent BBID on the patient's blood band. If a different barcode is scanned, the system will alert that the barcode is incorrect.

COMPLETE all other required steps per your facility protocol prior to transfusion.

Visit www.typenexmedical.com or call us at 866.897.3639 to learn more.