



REF M6018

Blood Transfer Device

INDICATIONS FOR USE: The Blood Transfer Device may be used for the transfer of venous blood from a luer lock syringe to blood collection tube(s) and blood culture bottle(s).

DIRECTIONS: Use aseptic technique.

1. Peel open package at arrow and remove device.
2. Insert the syringe tip into the hub of the Blood Transfer Device. Rotate the syringe clockwise until it fits securely on the hub.
3. With the syringe tip held facing down, center the Blood Collection Tube or Blood Culture Bottle over the holder portion of the Blood Transfer Device and push it in.
4. After removing the last Tube or Blood Culture Bottle, discard the entire assembly in a biohazard Sharps container approved for this disposal. **DO NOT** disassemble and reuse Blood Transfer Device.

WARNINGS AND CAUTIONS:

- Handle all biological samples and blood collection sharps in accordance with your facility's policies and procedures.
- Discard all blood collection sharps in biohazard containers approved for their disposal.
- Obtain appropriate medical treatment in the event of any exposure to biological samples (e.g. through puncture injury) since samples may transmit infectious diseases.
- Allow natural vacuum of tube to draw specimen into tube. Filling is complete when vacuum no longer continues to draw.

Sterile EO

Made in China

For single use only



IMPORTANT
DO NOT
RE-STERILIZE



DEHP
NOT MADE
WITH DEHP



NOT MADE
WITH NATURAL
RUBBER LATEX

⚠ WARNING: THIS PRODUCT CAN EXPOSE YOU TO CHEMICALS INCLUDING DIISONONYL PHTHALATE (DINP), WHICH IS KNOWN TO THE STATE OF CALIFORNIA TO CAUSE CANCER AND BIRTH DEFECTS OR OTHER REPRODUCTIVE HARM. FOR MORE INFORMATION GO TO www.P65Warnings.ca.gov.