AN
INNOVIA
MEDICAL
COMPANY

TISSU-TRANS® FILTRON®
100cc, 250cc, 500cc, 1000cc & 2000cc

Sterile Fat Collection and Filtration Instructions for Use

3-TT-FILTRON 100, 250, 500, 1000 & 2000

DESCRIPTION

The Tissu-Trans FILTRON devices are designed to collect and process autologous fat tissue in a sterile, closed environment. The Tissu-Trans® FILTRON devices are manufactured in 5 different sizes to facilitate physician preference for treatment. Product part number is provided on pouch label and indicates size. Product is packaged in pouches and is provided sterile. DO NOT re-sterilize.

Figure 1

STERILE R Sterile Items Included:

- Harvest Canister with Filter and Waste Extraction Port
- Harvest Tubing 6 feet / 1.83m
- Transfer Tube 6 inches / 152mm with Luer Lock tip

Standard Items Also Needed:

- Universal Bracket on tall canister stand or IV pole clamp
- Sterile mayo-stand drape or equivalent to cover canister stand
- Standard suction tubing (1/4" / 6.35mm or 5/16" / 7.94mm ID)
- Waste canister
- Suction source (set at approx. 15" Hg. (380mmHg) / 51 kPa or 50% power)
- Cannula (3mm cannula recommended)
- Luer Lock syringes for transfer and re-injection

Figure 2

The FILTRON is placed in-line between the cannula handle and waste canister

INTENDED USE

The Tissu-Trans® FILTRON® is intended to be used with vacuum and/or cleared pumps, tubing and cannulas, for the collection of aspirated fat, for aesthetic body contouring. If the fat is untreated, it may be re-injected via a cleared injection apparatus.

Tissu-Trans FILTRON is used in the aspiration, harvesting, filtering and transferring of autologous tissue.

Tissu-Trans FILTRON is intended for use in the following surgical specialties when the aspiration of soft tissue is desired:

- Plastic and Reconstructive Surgery
- Neurosurgery
- Gastrointestinal and Affiliated Organ Surgery

- Urological Surgery
- General Surgery
- Orthopedic Surgery
- Gynecological Surgery
- Thoracic Surgery
- Laparoscopic Surgery

INDICATION FOR USE

The Tissu-Trans FILTRON is intended to be used with vacuum and/or cleared pumps, tubing and cannulas, for the collection of aspirated fat, for aesthetic body contouring. If the fat is untreated, it may be re-injected via a cleared injection apparatus.

CLINICAL BENEFIT

The indirect clinical benefit for the Tissu-Trans Filtron devices is:

- Enables the aspiration, filtration, and autologous transfer of fatty tissue.
- Enables aesthetic body contouring via the injection of aspirated fat.

INTENDED USER & PATIENT TARGET GROUP

The Tissu-Trans FILTRON is intended to be used by a physician who has achieved proficiency in suction lipoplasty. Other healthcare professionals such as nurses, physician assistants, and surgical assistants may assist the physician in using the Tissu-Trans FILTRON device.

The intended patient population is adults.

INTENDED LIFETIME OF THE DEVICE

The Tissu-Trans FILTRON devices have a shelf-life of three years. Use this product before the expiration date indicated on the label.

NOTES BEFORE USE

1. The Tissu-Trans FILTRON is sterile and disposable. Single use.
⊗ ⊗ Do not re-sterilize or re-use.
2. Product is ready for use. Do not open canister or lid until fat collection is complete. Opening the canister or lid prior to fat collection may affect sterility and product performance.
3. Be certain that the tissue is well hydrated - as in superwet tumescent anesthesia.
4. **Must use** a standard waste canister between the Tissu-Trans FILTRON and suction source. During harvest, waste is filtered out to waste canister.
5. **Recommended to use** a 3mm Harvest Cannula.
6. **Recommended to use** the medical grade tubing provided to assure non-toxic pathway.
7. Prior to use, check the product expiration date and verify the sterile packaging has not been compromised. DO NOT use if the package is opened or damaged or the date of expiration has passed.
8. If the Tissu-Trans FILTRON does not hold a vacuum, use a new device and call Customer Service.

WARNINGS

1. This device will not, in and of itself, produce significant weight reduction.
2. This device should be used with extreme caution in patients with chronic medical conditions.
3. The volume of blood loss and endogenous body fluid loss may adversely affect intra and/or postoperative hemodynamic stability and patient safety. The capability of providing

adequate, timely replacement is essential for patient safety.

4. This device should not be exposed to suction pressure higher than 25 "Hg / 85 kPa.
5. Overfilling the device can result in harvest tubing blockage which may lead to extended operation times.
6. Do not attempt to re-use or re-sterilize the device, re-use may compromise the sterility of the device, and adversely affect the device performance. Re-use may lead to cross contamination and patient infection.
7. DO NOT USE if the device or its packaging has been damaged.

PRECAUTIONS

1. This device is used to collect aspirated fat for autologous fat transfer. It should be used with appropriate suction devices, tubing and cannulas which remove localized deposits of excess fat through small incisions.
2. Use of this device is limited to those physicians who have attained proficiency in suction lipoplasty. Users of this device may also include healthcare professionals who assist the physician.
3. Results of this procedure will vary depending on patient age, surgical site, and experience of the physician.
4. Results of this procedure may or may not be permanent.
5. The amount of fat removed should be limited to that necessary to achieve a desired medical effect.
6. All disposable components must be discarded.
7. Use the appropriate size canister for harvesting the required amount of fat, overfilling may result in sub-optimal filtration, and a poorer quality fat.

Do not re-use or re-sterilize. Single patient use.

DISPOSAL

Dispose in standard medical biohazard waste receptacle. Do not re-sterilize or re-use.


POSSIBLE ADVERSE EFFECTS

Autologous Fat Grafting procedures may result in the development of fat necrosis, cyst formation, infection, inflammation, donor or recipient site cellulitis, palpable abnormality and contour irregularities.

SERIOUS INCIDENT REPORTING

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the FDA/competent authority of the Member State in which the user and/or patient is established.

Figure 3

INSTRUCTIONS



1. Drape canister stand with sterile Mayo stand cover (or equivalent).
2. Remove sterile components and place on sterile surface. Remove and discard protective packaging.
3. Secure Tissu-Trans FILTRON to canister stand in an upright position by inserting mounting bracket on canister into universal slots on stand - being careful to maintain sterility.
4. Connect standard suction tubing to the VAC port on the canister lid. Connect the other end to waste canister and suction source.
5. Connect L-Connector end of included harvest tubing to PATIENT port on the canister lid. Connect other end of tubing to cannula.
6. Turn on suction. Check all connections. **Adjust vacuum pressure to 15"Hg (380mmHg) / 51 kPa or 50% power.** System may be tested with sterile saline to ensure proper flow.
7. Harvest desired volume of fat. Virtually pure* fat will remain in filter and waste will be suctioned out to waste canister. Fat remaining in filter will be ready for re-injection (if applicable).

Filtered fat may still have some liquid, which is necessary to ensure proper flow during transfer. After fat is transferred to syringes, it may be decanted further and excess liquid expelled.


*Demonstrated in a study to produce virtually pure fat containing trace amounts (1%) of residual oil and fluid

Figure 4

REMOVAL / TRANSFER**OPTION 1: DEVICE IS INVERTED.**

1. Turn off suction. Disconnect tubes from Tissu-Trans FILTRON lid, maintaining sterility.
2. On lid, close VAC port with attached cap. All ports except PATIENT should now be closed.
3. Attach the Transfer Tube to PATIENT transfer port. **See Figure 5.** Close the white pinch clamp on the Transfer Tube. **All ports are now closed.**
4. Carefully invert the Tissu-Trans FILTRON upside down.
5. Reattach inverted canister to the canister stand using the same bracket.
6. Attach Luer Lock syringe to Transfer Tube. Open pinch clamp and draw out fat to be re-injected. **See Figure 6.** Repeat with additional syringes.
7. When desired volume has been transferred, discard entire system.
8. Attach re-injection cannula. Fat is ready to be re-injected to recipient site (if applicable).

Figure 5**Figure 6****REMOVAL / TRANSFER****OPTION 2 (500CC, 1000CC, & 2000CC ONLY):****DEVICE REMAINS UPRIGHT.**

1. Turn off suction. Disconnect tubes from Tissu-Trans FILTRON lid, maintaining sterility.
2. On lid, close both VAC and PATIENT ports with attached caps. All ports should now be closed.
3. Attach transfer or re-injection cannula to luer lock syringe.
4. Remove cap from large CAP1 Port.
5. Insert cannula into CAP1 Port and draw out fat to be re-injected. Repeat with additional syringes.
6. When desired volume has been transferred, discard entire system.
7. Attach re-injection cannula. Fat is ready to be re-injected to recipient site.

Figure 7**PERFORMANCE CHARACTERISTICS**

- The Tissu-Trans FILTRON devices filter lipoaspirate to provide virtually pure fat* for autologous fat grafting
- Approximate volume of fat which can be filtered within the device
 - 100cm³ (3-TT-FILTRON 100)
 - 250cm³ (3-TT-FILTRON 250)
 - 500cm³ (3-TT-FILTRON 500)
 - 1000cm³ (3-TT-FILTRON 1000)
 - 2000cm³ (3-TT-FILTRON 2000)
- Maximum suction pressure: 25"Hg / 85 kPa

*Demonstrated in a study to produce virtually pure fat containing trace amounts (1%) of residual oil and fluid

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