

# INFINITY<sup>®</sup>

MASSAGE CHAIRS



**User Manual**

EVOLUTION MAX 4D



Welcome to

# INFINITY®

THE ULTIMATE MESSAGE EXPERIENCE

Thank you and congratulations on purchasing an Infinity Massage Chair - you've made a great choice.

With over a decade of experience designing, building and servicing massage chairs, Infinity takes pride in every chair sold, including yours. We are here for you. The US-based Customer Care team can be reached at (603) 910-5000 or visit us online at [www.infinitymassagechairs.com](http://www.infinitymassagechairs.com) and find valuable self service resources under the SUPPORT tab, including electronic copies of all documents and video tutorials.

Inside this manual, you will find everything you need to setup and safely enjoy your chair. Please review the entire manual, including all safety information, before assembling or using your chair. If you are assembling the chair yourself, please see the Assembly Guide next and reference the online assembly video.

For convenience, we've also created a Quick Start Guide. We suggest you keep it in the remote pocket as a tool to get to know all the features of your new chair.

**Get ready to sit back, relax and enjoy the Ultimate Massage Experience.**

## Warranty Registration

Register your product at [www.infinitymassagechairs.com/warranty](http://www.infinitymassagechairs.com/warranty) to ensure it is in our system. This will simplify the warranty process should you ever need it.

Alternatively, **scan the QR code** to go directly to the URL above.



SCAN CODE

## Assembly Video

Before assembling your chair, watch an assembly video at:  
[www.infinitymassagechairs.com/assembly-tutorials](http://www.infinitymassagechairs.com/assembly-tutorials)

Alternatively, **scan the QR code** to go directly to the URL above.



SCAN CODE



We know you can't wait to sit in your new Infinity Massage chair! Take a picture of you enjoying the ultimate massage experience and make sure to tag us on social media with the hashtag [#infinitymassagechairs](https://www.instagram.com/explore/tags/infinitymassagechairs)



# Table of Contents

## Your Chair

Chair Features	6
External Structure	7

## Before Getting Started

Installation Site	8
Clearance Space for the Chair	8
Floor Protection	8
Move Methods	9
Connect to Power	9

## Chair Setup

Getting Comfortable	10
Main Power Switch	10
Adjust Pillow	10
Adjustable Shoulder Airbags	10
Reduce the Intensity of Your Massage - Manually	11
Truefit™ Body Scanning	11
Getting Out Of The Chair	11
USB Charging Station	12

## Chair Entry & Exit

Chair Entry and Exit Hazards	13
------------------------------	----

## Taking Control

Remote Control Layout	14
Remote Control Operation	15
Auto Programs	16
Navigating the Remote	16
4D Intensity Control	16
Program Descriptions	17
Airbag Massage	18
Airbag Region	18
Airbag Intensity	18
Back Heat	18
Foot Reflexology Rollers & Total Calf™ Kneading and Oscillation	19
Joystick Control and Wireless Charging/USB Port	19

## Set Your Environment

Air Ionizer	20
Timer	20
App Download & Installation	20
Activating Voice Control	21
Intelligent Voice	
Command+Control	21
Turn Off Listening	21
Bluetooth® Audio	22
Sounds	22
Volume	23
Remote Control Beep	23
Demo Voice	23
Chromotherapy Programs	24
Light Themes	24
Light Colors	24
Light Brightness	25

## Settings

Language	26
Voice Control	26
Remote Control Screen Brightness	26
Button Lights	27
Remote Control Sound	27
Bluetooth® Audio	27
Sleep Mode	27

## Build Your Own

Manual Programs	28
Back Massage Technique Mode	28
Massage Part (Back Region)	28
Back Node Width	28
Massage Speed	29
Foot Reflexology Rollers & Total Calf™ Kneading and Oscillation	29
Weightless Zero Gravity Recline	29
Remote Control In-Use Screen	29

## Staying Safe

30

## Fixing Issues


Trouble Shooting	33
Specifications	33
Frequently Asked Questions	34




Warning! Read important safety information found on page 30 before using or operating chair.

# Your Chair

## Chair Features

<p>Zero Wall Fit™ Space-Saving Technology</p>  <p>See page 8</p>	<p>Adjustable Shoulder Fit</p>  <p>See page 10</p>	<p>Truefit™ Body Scanning</p>  <p>See page 11</p>	<p>4D Back Massage Mechanism</p>  <p>See page 16</p>	<p>Back Heat</p>  <p>See page 18</p>
<p>Complete Calf™ Kneading + Oscillation</p>  <p>See page 19 &amp; 29</p>	<p>Foot Reflexology Rollers</p>  <p>See page 19 &amp; 29</p>	<p>Joystick Control and Wireless Charging/USB Port</p>  <p>See page 19</p>	<p>Air Ionizer</p>  <p>See page 20</p>	<p>Apple &amp; Android App Functionality</p>  <p>See page 20</p>
<p>Intelligent Voice Command + Control</p>  <p>See page 21</p>	<p>Nature Sounds</p>  <p>See page 22</p>	<p>Bluetooth® Speakers</p>  <p>See page 22 &amp; 27</p>	<p>Chromotherapy Programs</p>  <p>See page 24</p>	<p>Manual Massage Customization</p>  <p>See page 28</p>
<p>Weightless Zero Gravity Recline</p>  <p>See page 29</p>				

# External Structure




# Before Getting Started

## Installation Site

### Clearance Space for the Chair

- 4" clearance for the recline of the backrest; 24" clearance for the extension of footrest and sliding forward of the chair.
- Keep it 3 feet away from the TV, radio or other audio & video appliance to avoid signal interference.




- Keep the distance between the power socket and power switch within 60" to avoid damage of power cord while massage chair sliding forward.
- Straighten the power cord while in use to avoid tangle or damage.



- To avoid serious injury or death, before sliding chair forward or backward check to make sure area is clear of children, pets, or obstacles.

## Floor Protection

We recommend to place a piece of thin carpet or soft mat underneath it to avoid damaging the floor.




- Do not use the massage chair in high moisture environment such as near a swimming pool or bathroom to prevent electric shock accident.
- Please keep the massage chair away from heated surfaces to prevent fire hazard or damage to the product.



## Move Methods

Apply pressure to the back of the massage chair from behind to pivot on the Castors. Have a second person up front to guide the Footrest until in the correct position. Maintain clearances as mentioned above. Gently lower into place. Do not pull on the side arms.



- Do not move the chair when there is user inside.
- Chair is heavy. Use caution and get assistance when lifting to avoid injury. Do NOT lift chair by footrest or armrests to avoid damage to product.

## Connect to Power

Plug into a 110V grounded outlet.



### Grounding Instructions

**WARNING**


To avoid risk of electrical shock, turn power switch off before plugging unit into electrical outlet. Use only with properly grounded outlet (3 prong). Do not plug into ungrounded outlet.

**WARNING** – Improper connection of the equipment-grounding conductor can result in a risk of electric shock. Check with a qualified electrician or serviceman if you are in doubt as to whether the product is properly grounded. Do not modify the plug provided with the product – if it will not fit the outlet, have a proper outlet installed by a qualified electrician.

This product is for use on a nominal 120-V circuit and has a grounding plug that looks like the plug illustrated. Make sure that the product is connected to an outlet having the same configuration as the plug. No adapter should be used with this product.

### Earthing

This massage chair is a Class 1 electrical appliance. Make sure to use a grounded outlet to the grounded power supply circuit and a properly installed earth grounding rod with ground connection, to avoid electric shock.




# Chair Setup

## Getting Comfortable

### Main Power Switch


To power up the chair, flip the main power switch located on the side of the base of the chair to **ON**. Ensure the electrical cord is connected to a 110V grounded outlet and plugged into the power cord port located on the back of the base of the chair.



- To avoid possible damage to product or injury to unsupervised children, always turn off power switch and unplug massage chair from electrical outlet immediately after use.
- When not in use for a long period of time, we recommend that the chair be covered, the power cord be coiled and that the unit be placed in an environment free of dust and moisture.


### Adjust Pillow

Take a seat in the chair and take note of where the neck pillow rests on your neck. Adjust the pillow so that it rests comfortably in the arch of your neck. Do this by removing it from the velcro, aligning it, and then reattaching to the velcro once it is in the desired position.



### Adjustable Shoulder Airbags

The shoulder airbags on your massage chair can be adjusted to better your shoulder height. Should you need to adjust the shoulder airbags up or down, simply grab the top of the airbag and pull down. To bring the airbag back up, reverse this process by grabbing the bottom edge and pulling up.




**Smart tip:** To increase intensity in the neck and shoulder region, flip the pillow or the pillow and pillow pad up so you lay directly on the seat back pad.

## Reduce the Intensity of your Massage - Manually


To make the massage less intense, you can **add the back pad** to the seat back of your chair, or you can also **fold down the neck rest** to provide more cushioning between you and the massage mechanism in the seat back. To make the foot rollers less intense, insert the foot pads that came with your chair into the foot bed of the footrest.

## Truefit™ Body Scanning

Select any auto program and press the **OK** button. The chair will automatically begin scanning your body, mapping your back, shoulders, and neck. During this process it is important to remain seated back into the chair without lifting your head. The scan can be viewed on the remote control screen. If fine tuning is needed just after your body scan, wait for the **shoulder adjust prompt** to appear and then use the **Back Adjustment** buttons to move the massage rollers up or down for the perfect shoulder kneading position.



Body Scanning



Shoulder Adjust Prompt

## Getting out of the Chair

To easily exit the chair at any time, press the **ON/OFF** power button on the remote control to restore the chair to the upright position.

**WARNING!** Check to make sure children and pets are clear of chair **before** pressing on/off button to return chair to upright position. See next page for detailed hazard warnings.




# Chair Setup

## USB Charging Station

Your massage chair includes a USB-A port to charge your phone or tablet while using the chair. This feature is especially helpful when using the app to control the chair.


The port is located on the left arm panel on the **Joystick Control**.




# Chair Entry & Exit

## Chair Entry and Exit Hazards


**WARNING!** The areas under the footrest and behind the back rest are potentially dangerous for a child or pet. The chair back is heavy and reclines with force. A small child or pet located behind the chair could be crushed as chair is reclined resulting in serious injury or death. When chair is reclined, small children or pets who have crawled under the leg rest may not be visible to the user and could be crushed when leg rest is lowered or when chair is returned to upright position. **To avoid serious injury or death, check these areas carefully before reclining, lowering leg rest, or returning chair to upright position.**



Make sure to check for children or pets **between the chair back and the base of the chair** before entering zero gravity recline.




Make sure to check for children or pets **underneath the footwell** before exiting zero gravity recline and returning to upright position.



Make sure to check for children or pets **between the footrest and the chair body** before lowering the legrest.

# Taking Control

## Remote Control Layout



## Remote Control Operation



**On/Off** – Powers up the remote or returns the chair to the exit position and powers down remote.



**Pause** – Press pause to immediately stop all massage chair activity, press again to resume.



**Home** – Shortcut to the home page of all menu screens.



**Feet Adjust** – Press and hold up arrow to bring foot rest up, press and hold down arrow to bring foot rest down.



**4D** – Press to increase or decrease 4D intensity of back mechanism. In-use screen displays intensity dial.



**Auto** – Shortcut to the Auto Programs menu screen.



**Back Adjust** – Back mechanism position adjustment; press and hold up to raise, press and hold down to lower. Use during initial body scan fine tuning or in manual point or Partial zone mode.



**Foot Rollers** – Press this button to toggle the foot reflexology rollers on or off.



**OK** – After navigating to the command desired on the display, push OK to select. During the massage, press OK button to switch to zero gravity.



**Directional Button** – Press up, down, left or right to navigate the display screen. During the massage, press up or down buttons to adjust massage speed. Press left or right button to adjust airbag intensity.



**Seat Back Position** – Press and hold down arrow to recline the seat back, press and hold the up arrow to upright the seat back.



**Zero Gravity** – Press for Zero Gravity 1, press again for Zero Gravity 2, press again to restore recline to upright position.



**Leg Length** – Press and hold right arrow to extend the leg length adjustment, press and hold the left arrow to retract the leg length adjustment.

# Taking Control

## Auto Programs

When massage chair is operating, press the **AUTO** button on the remote control to enter Auto Program Menu. Choose from **Extension, Release, Full Body, Thai, Chinese, Post Workout, Vigor, Night, Calm, Neck & Shoulder, Waist & Glutes, or Back & Spine**. Press **OK** to enter the program you desire. Selected program will show on the home screen.



## Navigating the remote

*To show the Auto Programs list:* Use the directional buttons to highlight **AUTO** from Home display, press **OK** button to select. Press up or down button to navigate amongst the 12 auto programs.

*To select an auto program:* Use the up and down directional buttons to highlight the program you wish to run. Press **OK** button to select and begin program.



## 4D Intensity Control

While the in-use screen is showing, use the **4D** button on your remote control to increase or decrease intensity. Intensity dial shows on in-use display area.






Program	What for?	Description
<b>Sports Refresh</b>	Soothe muscles after sports activities and exercise.	Performs a strong massage with deep, penetrating massage techniques to alleviate your aches and pains, improve circulation, and recuperate the body after intense physical activity.
<b>Extension</b>	This program is the best choice to relax tired muscles and increase physical vigor.	Inspired by Thai massage. The strong kneading and stretching functions effectively decompress the entire body.
<b>Rest &amp; Sleep</b>	Try this program as you prepare your body to rest for the night.	A progression of mild massage techniques and chair movements brings harmony to your body. This program helps you unwind from a long day, and prepares you for a restful night's sleep.
<b>Working Relief</b>	This is a great program for people who work at a desk or travel a lot.	Performs a strong massage to ease muscles and recuperate the body.
<b>Neck &amp; Shoulder</b>	Soothing pain and tension in the neck and shoulders.	A powerful combination of kneading, grabbing, and pressing massage techniques target the neck, shoulder, and upper back to help ease unwanted pain and stiffness.
<b>Waist &amp; Spine</b>	Helps to relax, recover, and revitalize the main back muscles along your spine.	Designed to be an incredible injury prevention program, this massage focuses on the lumbar vertebra with targeted massage to loosen rigid muscles, relieve tension, improve circulation, and revitalize your energy levels.
<b>Deep Shiatsu</b>	This program is recommended to relieve lower back pain.	Deep shiatsu on the back for 10 minutes.
<b>Healthy Breath</b>	Designed for morning or midday relaxing.	Rhythmic five-minute massage to aid with deep breathing.
<b>Massage Extend</b>	Reduce built up tension carried in the body from everyday stress.	Designed for morning or midday relaxing, this massage helps to boost the vitality of muscles.
<b>Relaxation</b>	Best for winding down and re-energizing yourself after a long day.	Assists in relieving lower and upper back pain through a series of massage techniques.
<b>All Air</b>	Many find this program beneficial to relax and ease stress when winding down.	Full-body airbag-only massage without any kneading and tapping.
<b>Demo</b>	A quick way to feel all the massage techniques and chair movements in a brief 8 minute session.	A quick show of massage chair's features and functions. That can be accompanied by an explanatory voice over. The voice can be turned off and on in SOUNDS.

# Taking Control

## Airbag Massage

Select **AIRBAG** from Home display, press up or down button to navigate among: **PART** and **FORCE**. Select by pressing **OK**.

**PART** allows you to select which airbags to use. Select Full for all airbags. Otherwise pick and choose amongst the other Parts to turn them on or off. You may select more than one Part. To turn off all Airbags, select any lit parts and hit **OK** to turn off. When no Parts are lit, all airbags are off.



## Airbag Region

Navigate to **PART**, press right button, then use the up or down button to select your desired airbag region: Full, Arm & Shoulder, Leg & Foot, Buttock. Press **OK** button to select the airbag region.



## Airbag Intensity



Navigate to **FORCE**, press right button, then use the up or down button to select your desired intensity: 1 = light, 5 = firm, or Off. Press **OK** button to select the airbag intensity.

While the in-use screen is displayed during a massage program, press the **left** or **right directional button** to adjust the airbag intensity (this feature works in both Auto and Manual programs).



## Back Heat

Select **HEAT** from Home display, press right button, then press up or down button to navigate to desired heat setting: 1 = low, 2, 3 = high or Off. Press **OK** button to select the heat setting.



**Smart tip:** While the in-use screen is displayed during a massage program, press the **left** or **right directional button** to adjust the airbag intensity (this feature works in both Auto and Manual programs).

## Foot Reflexology Rollers & Total Calf™ Kneading and Oscillation

Select **MANUAL** from Home display, then press up or down button to navigate to **ROLLER**. Press right button, then use the up and down buttons to select your desired foot roller speed: 1 = slow, 2, 3 = fast, or Off (deactivate foot rollers). Press **OK** button to select the foot roller setting.

Alternatively, press the **Foot Roller button** on your remote control. If the foot rollers are on, pressing this button will turn the foot rollers off. Press again to turn on. Your current foot roller speed is displayed in the in-use screen in the bottom right corner.



## Joystick Control and Wireless Charging/USB Port


As well as having a corded remote control, your chair also comes with a built in joystick controller which is located on the top of the left arm panel. The joystick controller allows you to turn the chair on and off, select from six of the auto program's, adjust the recline angle of the seat back, and wirelessly charge your mobile device.

**Wireless Charging Pad:** Place any wireless charging-enabled phone directly on to pad to charge.

**Button:** Press and hold to turn chair ON. Press again to cycle through six of the auto program's. The program name is announced through the speakers and starts immediately. Press and hold to turn chair OFF.

**Dial:** Use dial to recline and upright the Evolution Max 4D massage chair by turning left or right.

**USB Charging port:** Attach a standard USB-A cord to charge devices.



# Set your Environment


## Air Ionizer

Select **AIR IONIZER** from Home display, press right button, then press up or down button to select your desired air ionizer setting: On, or Off. Press **OK** button to select the setting.



## Timer

Select **TIME** from Home display, press right button, then press up or down button to select your desired timer: **10 minutes, 20 minutes, or 30 minutes**. Press **OK** to select your desired massage timer length. 20 minutes is the default setting.



## App Download & Installation

### Android System Software (Android System V2.2 or later)

- Open the Google play store. In the search bar, search for '**Infinity Message Chair Evolution Max**'. Click on the Evolution Max app tile.
- Click "Install" to download and install the Evolution Max **APK**, then press Settings button on your computer or Android device > press Wireless and Web Setting > Bluetooth® > Scan Bluetooth® device > Evolution Max, then pair.
- Start the Evolution Max APK, press Settings menu > Bluetooth® > and pair with the Evolution Max to link the APK


### Apple iOS: (iPhone or iPad devices)

- Open the App Store. In the search bar, search for '**Evolution Max by Infinity**' to find the app. Click "Get" at the top right corner to download and install the software.
- In your device settings, open Bluetooth® and connect it to the massage chair.
- *Please note: Only Apple equipment compatible with Bluetooth® 4.0 can connect to the chair's Bluetooth®.*

## Activating Voice Control

Navigate to **VOICE**, press right button, then press up or down button to select your desired voice control setting: On, or Off. Press **OK** button to select the setting.

To activate voice control function without going into settings, touch the button on the left speaker or say “Switch to voice mode”.



## Intelligent Voice Command + Control

Your chair comes with built in intelligent voice control. This feature allows you to control your chair and massage experience completely hands-free by using one of the many voice commands.

To activate the intelligent voice control function, simply say '**Switch to Voice Mode**'. You may now give your chair any of the commands to the right **within 5 seconds**. Alternatively, you can touch the voice command pad on the left speaker to activate, and then speak your command within 5 seconds. Make sure to speak clearly and slowly. After 5 seconds the Intelligent Voice Command function will turn off automatically and must be awakened again.

### Turn Off Listening

You can turn off intelligent voice control listening by using the remote control. Select **SETTING** from Home display, press up or down button to navigate to **VOICE**, press right button, then press up or down button to select **Off**. Press **OK** button to select the setting. Intelligent voice control will now no longer activate until turned back on.

### Commands

- Power On
- Switch Off
- Activate Massage
- Sports Refresh
- Extension
- Rest & Sleep
- Working Relief
- Neck & Shoulder
- Waist & Spine
- Airbag Massage
- Deep Shiatsu
- Healthy Breath
- Massage Extend
- Zero Gravity
- Start Heat
- Heat Off

# Set your Environment

## Bluetooth® Audio

Your massage chair includes Bluetooth® compatibility so you can play and enjoy your favorite music. To turn on your chair's Bluetooth®, select **SETTING** from Home display, then press up or down button to navigate to **BLUETOOTH**. Press right button, then press up or down button to select your desired Bluetooth® setting: On or Off. Press **OK** to select setting.



Open Bluetooth® from audio device and pair with **Evolution Max** to play music.

## Sounds

Press **SOUNDS** from Home display, press right button then press up or down button to navigate among: **SOUND**, **VOLUME**, **BEEP**, and **DEMO VOICE**. Press **OK** button to enter the sounds setting you desire to control.




Note: Bluetooth must be turned on in settings for Sounds to work.

Sounds	What is it?
<b>White Noise</b>	A constant, even loop of pure white noise to mask out distractions
<b>Ocean Surf</b>	Soothing ocean waves rolling onto the shore
<b>Crickets</b>	Night time in the suburbs. Churping Crickets, water trickling from a pond, and faint traffic noises
<b>Rain</b>	Quiet, light, refreshing rain falling to the ground
<b>Nature</b>	Refreshing sounds of small birds peacefully singing along a tranquil stream.
<b>Air Plane</b>	A Cessna propeller airplane motor idling
<b>Car Ride</b>	Calming inner-city traffic featuring cars, trucks, and motorcycles
<b>Fireplace</b>	Crackling of fire in a fireplace
<b>Thunderstorm</b>	Heavy rainstorm with thunder
<b>Windy Night</b>	Heard from inside a cabin, a howling wind blows constantly outside

Note: Bluetooth must be **ON** for sounds to work.


## Volume

Navigate to **VOLUME**, press right button, then press up or down button to select your desired volume setting: 1 = lowest, 2, 3, 4, 5 = highest. Press **OK** button to select the volume setting.




## Remote Control Beep

Navigate to **BEEP**, press right button, then press up or down button to select your desired remote control beep setting: 1 = lowest, 2, 3, 4, 5 = highest. Press **OK** button to select the volume setting.



## Demo Voice


Navigate to **DEMO VOICE**, press right button, then press up or down button to select your desired setting: On, or Off. Press **OK** button to select the setting. Voice demo only works in the Demo program.



# Set your Environment

## Chromotherapy Programs

Your massage chair has chromotherapy lights located on the speakers near your head. Select **LIGHTS** from Home display, then press up or down button to navigate among: **THEMES**, **COLORS**, and **BRIGHTNESS**. Press **OK** button to enter the lights setting you desire to control.



### Light Themes

Navigate to **THEMES**, press right button, then use the up or down button to select your desired light theme: Solid, Breathe, Wake Up, Wind Down, Multi-Wave, or Multi-Rotate. Press **OK** button to select the setting.

**Solid:** Single color light.


**Breathe:** Single color that dims and brightens on a loop.

**Wake Up:** Blue enriched lights - gradually increasing brightness.

**Wind Down:** Colors of sunset - gradually decreasing brightness.


**Multi-wave:** All the colors of the rainbow morphing into each other.

**Multi-rotate:** Colors rotate from one to another.



### Light Colors


Navigate to **COLORS**, press right button, then use the up or down button to select your desired light color: Ice Blue, Red, Yellow, Lime, Green, Cyan, Blue, or Purple. Press **OK** button to select the setting. You may only choose colors in Solid and Breathe themes.





## Light Brightness


Navigate to **BRIGHTNESS**, press right button, then use the up or down button to select your desired light brightness: 1 = dim, 2, 3, 4, 5 = bright. Press **OK** button to select the setting.



# Settings

## Setting

Select **SETTING** from Home display, then use the up or down button to toggle between **Language, Bluetooth, Power Save** or **Voice**.



## Language


Navigate to **LANGUAGE**, press right button, then use the up or down button to select your desired language: English, Chinese, Japanese, or Korean. Press **OK** button to select language.



## Voice Control


Navigate to **VOICE**, press right button, then press up or down button to select your desired voice control setting: On, or Off. Press **OK** button to select the setting.

To activate voice control function without going into settings, touch the button on the left speaker or say "Switch to voice mode".



## Remote Control Screen Brightness


Navigate to **BRIGHTNESS**, press right button, then press up or down button to select your desired remote control screen brightness setting: High, Medium, or Low. Press **OK** button to select the setting.



**Smart tip:** In order to play music from your Bluetooth enabled audio device, **Bluetooth** must first be turned on in settings.


## Button Lights

Navigate to **BUTTONS**, press right button, then press up or down button to select your desired remote key back-light setting: On, or Off. Press **OK** button to select the setting.



## Remote Control Sound

Navigate to **BEEP**, press right button, then use the up or down button to select your desired remote control beep setting: On or Off. Press **OK** button to select setting.



## Bluetooth® Audio

To turn on your chair's Bluetooth®, Navigate to **BLUETOOTH**, press right button, then use the up or down button to select your desired Bluetooth® setting: On or Off. Press **OK** button to select setting. Open Bluetooth from audio device and pair with **Evolution Max** to play music.


In order to use the chair app, play in-built sounds from your chair, or to play music from your Bluetooth® enabled audio device, Bluetooth® must first be turned on in settings.



## Sleep Mode

Navigate to **SLEEP MODE**, press right button, then use the up or down button to select your desired sleep mode setting: On or Off. Press **OK** button to setting.

When sleep mode is **on** the chair will stop in the reclined position after Auto Programs, when sleep mode is **off** the chair will return to a fully upright position after each Auto Program.




# Build your own

## Manual Programs


In Manual mode you can operate components of your massage chair individually or combine them to create unique massage scenarios.

Press **MANUAL** from Home display, press up or down button to navigate among: **MODE**, **PART**, **WIDTH**, **SPEED**, and **ROLLER**.



## Back Massage Technique Mode


Navigate to **MODE**, press right button, then use the up or down button to select your desired mode: Kneading, Tapping, Knocking, Shiatsu, Sync, or 4D Shiatsu. Press **OK** button to select mode.



## Massage Part (Back Region)


Navigate to **PART**, press right button, then use the up or down button to select your desired region: Whole, Partial, and Point massage. Press **OK** button to select region.

**Whole** travels the entire length of the track, **Partial** travels a short distance and reverses direction for focused relief, **Point** stops on the track for targeted relief. Partial and Point location can be fine tuned using the arrow up or arrow down buttons.



## Back Node Width

Navigate to **WIDTH**, press right button, then use the up or down button to select your desired width: Wide, Medium, and Narrow roller width. Press **OK** button to select mode. (Massage Width can't be adjusted when in Kneading, or Sync Manual modes).




**Smart tip:** **Message Width** can't be adjusted when in Kneading, or Sync mode techniques. **Message Speed** can't be adjusted when in Shiatsu.

## Massage Speed

Navigate to **SPEED**, press right button, then use the up or down button to select your desired speed: 1 = slow, 6 = fast. Press **OK** button to select the speed. (Massage Speed can't be adjusted when in Shiatsu).

While the in-use screen is displayed during a massage program, press the **up** or **down directional button** to adjust the massage speed. Your current massage speed is displayed in the bottom left corner.




## Foot Reflexology Rollers & Total Calf™ Kneading and Oscillation

Navigate to **ROLLER**, press right button, then use the up or down button to select your desired foot roller and calf oscillation speed: 1 = slow, 2, 3 = fast, or OFF (deactivate foot rollers + oscillation). Press **OK** button to select the setting.

Alternatively, press the **Foot Roller button** on your remote control. If the foot rollers are on, pressing this button will turn the foot rollers off. Press again to turn on. Your current foot roller speed is displayed in the in-use screen in the bottom right corner.



## Weightless Zero Gravity Recline


Zero Gravity mode will recline the back and elevate the feet to reduce the gravitational force on your spine. To activate, simply press the  button on your remote control.



Press once for Zero Gravity 1, press again for Zero Gravity 2, press again to restore recline to upright position.



## Remote Control In-Use Screen

Whilst running an Auto Program or Manual mode, the remote control screen will return to the in-use screen (as shown below) **after 10 seconds of being idle**.

 - Do not spill liquids on remote control and avoid placing heavy objects on remote control or power cord to avoid damage or malfunction.



# Staying Safe

## SAFETY INFORMATION

To ensure safe use of the massage chair and avoid injury or unsafe conditions, please read and follow all instructions before using product.

### WARNING

#### USE NOT RECOMMENDED FOR CERTAIN PERSONS

Use of massage chair may not be safe or suitable for persons with certain medical conditions and could cause serious injury or death in some circumstances. If you have any of the following conditions, please consult a physician before using the massage chair:

- Use of implanted electronic medical devices (such as a cardiac pacemaker)
- Pregnancy
- Cancer, malignant tumors, or other growths
- Osteoporosis or spinal fracture
- Other fractures or musculoskeletal injuries
- Loss of or reduced sensitivity to heat
- Implanted medical fixation devices such as rods, pins, etc.
- Recent surgery or sutures
- Skin diseases or skin injury
- High fever
- Certain pulmonary or breathing conditions
- Impaired circulation

### WARNING

To avoid possible death or serious injury to user or others, please read and follow these instructions:

### WARNING

- **The areas under the footrest and behind the back rest are potentially dangerous for a child or pet. The chair back is heavy and reclines with force. A small child or pet located behind the chair could be crushed as chair is reclined resulting in serious injury or death. When chair is reclined, small children or pets who have crawled under the leg rest may not be visible to the user and could be crushed when leg rest is lowered or when chair is returned to upright position. To avoid serious injury or death, check these areas carefully before reclining, lowering leg rest, or returning chair to upright position.**
- **Chair should not be used by children under age 14.** Small children could be trapped or injured by airbags, heating mechanism, or moving parts. When chair is not in use, unplug unit and secure chair so that children cannot accidentally power on and operate chair.
- Children over age 14, persons with reduced physical, sensory, or mental capabilities, and elderly or frail persons should exercise caution and use massage chair only when supervised or assisted by a person over 16 years of age.
- Chair has a heated surface and could cause burns, particularly in persons with reduced sensitivity to heat. Such persons should use the chair only after consulting with a physician and should use extreme caution during any physician-approved use.

## CAUTION

To avoid the risk of minor injury or discomfort, or damage to the massage chair or other property, please read and follow these instructions for safe use of massage chair:

### CAUTION

- Do not stand or stack heavy objects on footrest.
- Do not operate heating features of chair while under a blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
- To avoid pinching or crushing, please do not place body parts, especially fingers and toes, between parts of chair and keep clear of moving parts.
- To avoid overheating, injury or discomfort from overuse, do not fall asleep while using massage chair and do not use for extended periods of time (longer than 20 minutes).
- Avoid use while under the influence of alcohol, medications or other substances that may cause drowsiness, disorientation, or impairment of judgment or motor skills.
- Do not use the massage chair right after a meal, as it may cause discomfort.
- If you feel any abnormal discomfort during the operation of the massage chair, turn off the power immediately and discontinue use.
- Do not use this appliance when ambient temperature is over 40 C (104 F).
- Do not use this appliance in humid or dusty environment, or when any part of the body is in contact with plumbing or any similar ground.
- Never operate the appliance with the air openings blocked. Keep the air openings free of lint, hair, dust and the like.
- Do not use outdoors.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- Place the product on a flat, level floor to reduce risk of tipping or overturning.
- Do not use or store near heat or open flame.
- Intended only for users who weigh below 300 lbs.
- Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Massage chair intended only for use by one person. Do not use by two or more persons at the same time to avoid damage to the product.
- Do not allow children to stand, jump, or play on or around massage chair, even when it is not powered on or in use.
- Remove personal accessories from arms and empty pockets while using the massage features of the chair.
- Check the upholstery under the cushions and other spots for damage or loose objects that may cause injury or damage chair while operating. If any damage is present, do not use the massage chair and refer to a qualified service professional for servicing.
- To avoid damage to chair, walls, or other objects, make sure there are no obstacles behind the chair and there is enough clearance before adjusting the backrest or leg rest angle.
- Please be seated properly while using the product.
- To avoid damage to chair or floors, do not drag or push the chair or roll its casters on wood floor or carpet. See instructions on p.8 for safely moving chair.
- If massage chair is not working properly, has been dropped or damaged, or has been dropped into water or exposed to flooding, return the appliance to a service center for examination and repair before attempting to use.
- Do not drop or insert any object into any openings in chair.

# Staying Safe

## WARNING

Read and follow these instructions concerning safe use of the power cord and power unit to avoid possible electric shock, fire hazard, use of chair by and injury to unattended children, or damage to chair or other property.

### WARNING

- To avoid risk of electric shock, connect this unit to a properly grounded outlet only. See Grounding Instructions on pg. 9.
- Always unplug massage chair from the electrical outlet immediately after use and before cleaning. An appliance should never be left unattended when plugged in. To disconnect, turn all controls to the off position, then remove plug from outlet.
- Do not use massage chair and unplug from outlet during thunder or lightning storms, and before putting on or taking off parts.
- To ensure safe operation, power must be supplied only at standard 110 voltage, to the specifications of the appliance. Confirm that the power voltage stated on the rating label corresponds to your local electricity supply.
- The massage chair should be used only with the power supply unit provided. Do not use with any other power source.
- Do not use a wet or damp hand to unplug the cord to avoid electric shock.
- Do not pull on the power cord to unplug. Grip by plastic end nearest outlet to unplug power cord.
- Do not squeeze, over bend, pull hard, or twist the power cord.
- Keep the power cord away from heated surfaces.
- Do not carry or drag the chair by power cord or use cord as a handle.
- If there is any smoke or unusual smell while using massage chair, stop using immediately and unplug the power cord.
- Never operate the chair if it has a damaged power cord or plug. If the power cord is damaged, it must be replaced by the manufacturer or its service agent or a similarly qualified person to avoid a hazard.

**Do not disassemble, dismantle or attempt to repair massage chair. Such actions, as well as misuse of product, may cause damage to or failure of massage chair and may void the manufacturer's warranty. For repairs, please contact an authorized dealer or repair agent at 603-910-5000. For additional details concerning terms, conditions and limitations of warranty, see <https://infinitymassagechairs.com/warranty-info>.**



# Fixing Issues

## Trouble Shooting

No.	Symptoms	Possible Causes	Possible Solution
1	Noise is heard while in use	These noises are normal and arise from friction between the massage rollers and the faux leather cover	No need to take any measures
2	Massage Chair does not start	<ul style="list-style-type: none"> <li>- The power cord is not secured into the socket</li> <li>- Not switched ON</li> <li>- The fuse is blown</li> </ul>	<ul style="list-style-type: none"> <li>- Insert the power plug properly into the power socket</li> <li>- Turn on the switch found on chair side. Then push power button on remote</li> <li>- Replace only with a fuse of same type and rating</li> </ul>
3	No music comes from the speaker	The volume is turned to the lowest setting	Adjust the speaker volume from the mobile device. Make sure Bluetooth setting on chair and phone are both ON
4	Backrest or footrest can not raise or recline	Excessive load has been pressed on the backrest or legrest	Reduce the load and try again
5	Airbag malfunction	The air hose is blocked	Straight the air hose and clearout the blockage

If problem remains, please contact Infinity for maintenance.

**Tip:** Do not attempt to repair this product yourself, our company shall not be liable for the injury or damage caused by any servicing without authorizations. Thank you for your cooperation

## Specifications

Model	Evolution Max 4D
Dimension of Main Body (in)	66 x 35 x 48 (L x W X H)
Necessary Clearance to Wall (in)	4"
Net Weight (lbs)	323
Shipping Weight (lbs)	225
Max. User Weight (lbs)	300
Voltage	110-120V
Rated Current	1.55A
Power Consumption	160W

# Fixing Issues

## Frequently Asked Questions

### **Q: I find the chair is too intense to be comfortable. What can I do?**

**A:** Your body needs to become acclimated to therapeutic massage. To be effective, the massage techniques must compress the tissue in your body. We recommend the following steps to help get started. Then over time as your body acclimates you can increase the intensity by reducing the layers between you and the chair.

1. Use body the Pillow and Pillow pad to reduce upper body intensity.
2. Go into Airbag on the main menu, scroll down to FORCE, select air intensity setting 1  
- If a specific airbag is too intense, scroll down to PART, turn off Full and turn on only the airbags you want.
3. Go into Manual, select Rollers, and turn the foot reflexology rollers off.
4. Adjust the recline on the remote. More upright will be less intense. Weightless Zero Gravity Recline will be more intense.
5. If after adjustments and attempt to acclimate, you still experience discomfort, discontinue use and consult a physician.

### **Q: I'm not a foot massage person, what should I do?**

**A:** Turn off foot reflexology rollers, navigate to Manual on the home display, arrow down and select Rollers, then switch to off. You may also want to turn off the Leg/Foot airbags, or reduce the air intensity to level 1.

### **Q: My chair is not performing as expected. What should I do?**

**A:** Your Infinity massage chair is an amazing combination of programming and robotics. On occasion the chair just may need to be rebooted. To do this power down the chair and unplug from the wall for 5 minutes then plug back in and power back up to see if it self-healed. If you are still having an issue, give Infinity a call at 603-910-5000.

### **Q: I can never seem to get in my Infinity Massage chair, someone is always in it.**

**A:** We know. It happens. Tell your friends to get their own Infinity Chair.

For more information and access to further questions, and resources, visit our website at [www.infinitymassagechairs.com](http://www.infinitymassagechairs.com)



# INFINITY<sup>®</sup>

72 Stard Road, Seabrook, NH 03874  
603-910-5000  
[www.infinitymassagechairs.com](http://www.infinitymassagechairs.com)

Version 1.1