# FL-6000

## LED Illuminator

### Operation Manual

# Model FL-6000 LED Fiber Illuminator

## User's Manual

### 1. Scope

Applicable to FL-6000

### 2. Symbol Description ; Symbolbeschreibung; Description des symboles; Descrizione del simbolo

CAUTION: Refer to Operator's Manual

### 3. Warnings and Precautions ; Warnungen und Sicherheitshinweise; Avertissements et précautions; Avvertenze e precauzioni

#### 3.1. READ ALL INSTRUCTIONS BEFORE USE

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#### 3.2. WHEN CHANGING FIBER, DEVICE SHOULD BE TURNED OFF

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During operation, some internal components can become hot. To avoid the possibility of injury, always turn the device off and allow 30 seconds for cooling prior to changing the fiber or adapter or otherwise exposing the internal components.

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### 3.3. PROVIDE ADEQUATE VENTILATION AT ALL TIMES

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Allow clearance for the ventilation slots as described in section 5.3 to ensure the device does not overheat and trigger the internal shut-off mechanism.

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### 3.4. DIRECT VIEWING OF EXITING LIGHT CAN BE HARMFUL

DER DIREKTE BLICK IN DEN AUSTRETENDEN LICHTSTRAHL KANN  
SCHÄDLICH SEIN  
IL PEUT ÊTRE DANGEREUX DE REGARDER DIRECEMENT LA LUMIÈRE DE  
SORTIE  
GUARDARE DIRETTAMENTE VERSO LA LUCE IN USCITA PUÒ ESSERE  
DANNOSO

Avoid direct viewing of the light. All LEDs and LED-based products are sources of light which can in some circumstances cause damage to the eye. The potential for any eye injury will depend upon many factor including but not limited to:

- duration of exposure and extent of pupil dilation prior to exposure
- distance from light source to eye
- drive current supplied to LEDs determined by the intensity light setting
- beam pattern (degree of focus)

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Use of other power converters can damage product, will void warranty,  
and may result in unsafe operation

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### 3.6. PLUG AC CORD ONLY INTO SUITABLY GROUNDED AC RECEPTACLES

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### 3.7. ALWAYS USE THE DEVICE IN A MANNER SPECIFIED BY THE MANUFACTURER

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If the illuminator is used in a manner not specified by the manufacture,  
the protection provided by the equipment may be impaired.

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### 3.8. DO NOT USE FOR MEDICAL APPLICATIONS

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This product is not intended for medical use.

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## 4. General Information

### 4.1. Intended use of equipment

4.2. The FL-6000 series Illuminator is a LED-based fiber optic illuminator designed for general microscopy and for industrial applications such as assembly illumination, machine vision, and boroscope inspection

## Technical specifications

Technical data	Value
Light output at maximum intensity setting	780 lumens typical*
Color temperature	5,000 to 6,000°K
Input voltage	90 to 264VAC, 1.5 A, 50 to 60Hz
Overvoltage Category	II
Dimensions	3.9" W x 5.1" H x 7.0" L (9.9cm W x 13.0cm T x 17.8cm L)
Weight	6.8 lbs/3.1kg

Environmental use conditions	Value
Ambient temperature	0°C to 40°C (32°F to 104°F)
Relative humidity	80% max at 31°C (88°F) decreasing linearly to 50% at 40°C (104°F)
Altitude	Up to 2000 m/6600 ft
Pollution degree	2

\*Theoretical value at fiber input, contact Dolan-Jenner for details.

### 4.3. Instructions for use

See section 6 "Equipment Operation"

#### 4.4 Manufacturer

Name and address of manufacturer or supplier from whom technical assistance may be obtained.

Dolan-Jenner Industries  
159 Swanson Road  
Boxborough, MA 01719  
1-800- 833-4237  
www.dolan-jenner.com

#### 4.5. Equipment components

The product includes the following components:

- FL-6000 Fiber-Optic Illuminator
- External power converter (AC power supply) with attached cord
- AC power cord

#### 5. Equipment installation


##### 5.1. Positioning the illuminator

- Choose a suitable location that provides the following:
- Adequate clearance to connect the power converter per section 5.2.
- Within reach of a suitable grounded AC outlet.
- Adequate clearance to operate the controls per sections 6.2 and 6.3.
- Clearance for ventilation per section 5.3.


##### 5.2. Connections to the supply

- 5.2.1. Insert the male plug from the external power converter into the female power jack on the back of the unit.
- 5.2.2. Plug the AC cord into the receptacle on the external power converter and into the AC outlet.
- 5.2.3. Plug AC cord only into suitably grounded AC receptacles
- 5.2.4. Use only the power converter provided with the product - use of other power converters can damage product, will void warranty, and may result in unsafe operation."

Inserting power converter plug into rear of Fiber Illuminator


AC plug connection

### 5.3. Ventilation requirements

Do not block the ventilation slots provided for the internal cooling fan. Slots are located on both sides of the unit and on the back.

**Do not block vent slots.**


## 6. Equipment operation

### 6.1. Attaching a fiber optic adapter

If fiber adapter is not already attached to Fiber Illuminator unit, orient the adapter such that the adapter code designation is viewable on the top of the adapter as shown below. Insert the adapter into the Fiber Illuminator aligning the three pins with the corresponding three slots, push and twist clockwise until the thumbscrew is oriented horizontally as shown.

### 6.2. Turning the illuminator on and off

To turn the illuminator on, push the power switch on the lower left of the front panel to the ON position. A green indicator light will turn on when the illuminator is powered on.

To turn the illuminator off, push the switch to the OFF position.

### 6.3. Changing intensity

To change the intensity of the light output, rotate the Intensity dial on the bottom right of the front panel. Rotating clockwise increases light output and rotating counterclockwise decreases intensity.

### 6.4. Instructions for cleaning

6.4.1. The exterior of the illuminator can be gently cleaned with a soft cloth dampened with water and mild detergent.

Do not use any solvents or cleaner containing chlorine as these will damage the illuminator and may void your warranty.

6.4.2. The ventilation slots can be cleaned by blowing aerosol- based dusting air into the slots.

Do not blow compressed “shop-air” into the illuminator as such air often contains oil and water mist which may degrade the internal optics.

## 7. Product Ratings

For all CE technical questions, contact Setra Systems, USA. EU customers may contact our EU representative Hengstler GmbH, Uhlandstr 49, 78554 Aldingen, Germany (Tel: +49-7424-890; Fax: +49-7424-89500).