Hotplates & Stirrers

Instructions for use

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	Model No.
Hotplates	SB300
	SB500
	CB300
	CB500
Digital Hotplates	SD300
	SD500
Magnetic Stirrer	SB301
Stirrer Hotplates	SB302
	CB302
Infrared Hotplate	CR300
Infrared Hotplate Stirrer	CR302

Before Use

Thank you for purchasing this Bibby Scientific product. To get the best performance from the equipment, and for your own safety, please read these instructions carefully before use.

If the equipment is not used in the manner described in this manual the protection provided may be impaired.

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuation not exceeding 10%
- ❖ Overvoltage category II IEC60364-4-443
- ❖ Pollution degree 2 IEC664
- ❖ Use with a minimum distance all around of 200mm from walls or other items

Electrical Installation

Before connection please read and understand this instruction manual and ensure that the line supply corresponds to that shown on the rating plate. All models require a supply rated at: 230V, 50Hz, ~, single phase.

Power consumption is:

Model	Power	Model	Power
SB300	600W	CB300	1200W
SB301	50W	CB302	1250W
SB302	650W	CB500	2250W
SB500	1500W	CR300	900W
SD300	600W	CR302	950W
SD500	1500W		

THIS INSTRUMENT MUST BE EARTHED

All models are supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a U.K. 3 pin plug and the other has a 2 pin "Shuko" plug for connection to the mains. Choose the lead appropriate for your electrical installation and discard the other.

Should neither lead be suitable, take the lead with the U.K. plug and replace the plug with a suitable alternative. This involves cutting off the moulded plug, preparing the cable and connecting to the rewirable plug in accordance with its instructions.

IT IS IMPORTANT THAT THIS OPERATION SHOULD ONLY BE UNDERTAKEN BY A QUALIFIED ELECTRICIAN

NOTE: Refer to the equipment's rating plate to ensure that the plug and fusing are suitable for the voltage and wattage stated.

The wires in the mains cable are coloured as follows:

BROWN	- LIVE
BLUE	- NEUTRAL
GREEN/YELLOW	- EARTH

The instruments are fitted with an IEC socket at the rear of the instrument for connection of the mains lead. The appropriate mains lead should be connected BEFORE connection to the mains supply.

Should the mains lead need replacement a cable of 1mm² of harmonised code H05W-F connected to an IEC320 plug should be used.

N.B. The UK mains lead is protected by a 10A fuse mounted in the plug top.

IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN

When connected to the mains electricity supply a green light on the front panel illuminates, except for the digital SD300 and SD500 models where the display will illuminate showing the word "OFF".

Safety Advice Before Use

- ❖ If the mains supply is interrupted the units will restart with its current settings on the restoration of the electricity supply.
- ❖ Do not use an instrument with a ceramic top if the top plate is scratched, chipped or shows signs of chemical etching.
- ❖ Position the instrument on a firm level surface away from any heat sensitive or flammable material.
- ❖ Do not use the units with flammable liquids or in a hazardous atmosphere.
- ❖ There is a danger of liquid spills if vessels are stirred too vigorously. Always build stirrer speed slowly and never stir faster than necessary.
- ❖ Never move or carry the instrument until it has been switched off and allowed to cool.
- ❖ Never move or carry the unit with containers on the top plate or while the instrument is connected to the mains electricity supply.

Hot Surfaces

WARNING: The top surface of the instrument may be HOT. Do not leave heaters switched on when not in use.

WARNING: When the surface becomes too hot to touch the red "HOT" warning light on the front panel will begin to flash. This will continue to flash while the plate temperature is above 50°C even if the unit is switched off as long as it is connected to the electricity supply.

WARNING: Do not disconnect or switch off the electricity supply until the "HOT" warning light has ceased flashing.

WARNING: The plate may still be hot when large masses are left on the plate, even though the "HOT" warning light has gone out.

WARNING: The top surface of the instrument may be HOT, especially in free air when a surface temperature of up to 450°C can be achieved on some of the ceramic top models and 300°C on the metal top models.

Heating

Models SB300, SB302, SB500, CB300, CB302, CB500:

The knob labelled "Heat" on the front panel controls the heat output and hence the plate temperature.

The knob is graduated from 1 - 9. Increasing the value increases the heat output.

Turn the hotplate control knob clockwise to the desired setting. The amber load indicator light will illuminate and it will begin to cycle on and off when the desired temperature is reached and the hotplate begins to control.

Models SD300 and SD500:

When connected to the mains electricity supply the digital display will illuminate and show the word OFF.

Switch the unit on by pressing the left hand one of the 2 buttons located to the left of the digital display. The display will now show the actual plate temperature in degrees Celsius.

Press and hold the "press to set" button, the right hand one of the 2 buttons located to the left of the digital display. The display will now show the set temperature.

The set temperature may be adjusted by turning the knob located to the right of the digital display whilst continuing to hold in the "press to set" button. When the desired setting is displayed release the "press to set" button. The display will revert to show the actual plate temperature and the instrument will begin to heat.

The set temperature can be observed at any time by pressing the "press to set" button.

Turning the adjusting knob without pressing the "press to set" switch will have no effect on the hotplate setting.

Models CR300 and CR302:

The knob labelled "Heat" on the front panel controls the heat output and hence the plate temperature. When the heat is switched on the amber pilot light illuminates.

The knob is graduated from 1 - 9. Increasing the value increases the heat output.

In operation only the area defined by the pattern will become fully heated. The outer edges remain cooler for safety considerations.

Vessels larger than the described area should not be used. Failure to observe this could give rise to overheating of the outer parts of the instrument.

NOTE: To prevent overheating of the infrared heater a temperature limiter is fitted which switches off the heater should the temperature exceed 580°C.

Stirring

Models SB301, SB302, CB302 and CR302:

The knob labelled "Stir" on the front panel controls the stirrer speed.

It is labelled with an arbitrary numeric scale, 1 - 9 and turning the knob to a higher setting increases stirrer speed.

When the stirrer is switched on the amber pilot light illuminates.

Cleaning and Care

Before attempting cleaning disconnect from the mains electricity supply and ensure that the top is cool.

Periodically clean the body of the instrument using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents.

Ceramic top models: The ceramic top is very easy to keep clean. A damp cloth will normally remove most types of contamination. For more difficult stains a domestic cream cleanser is recommended. Cleaning is made easier if spillages are attended to promptly. In any case, spillages of alkali, phosphoric acid and hydrofluoric acid **MUST** be removed immediately as these chemicals can attack and damage the glass ceramic. Ensure that the appropriate safety precautions are observed.

During cleaning and general operation take care not to scratch the surface of the top plate as this could result in subsequent thermal breakage.

WARNING:

A ceramic top which is scratched, chipped, chemically etched or damaged must not be used.

Metal top models: The metal top plates should be cleaned using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents. Cleaning is made easier if spillages are attended to promptly. In any case, spillages of acids and alkalis **MUST** be removed immediately as these chemicals can attack and damage the aluminium alloy. Ensure that the appropriate safety precautions are observed.

Servicing and Repair

This product range does not require any routine servicing.

Note: There are no internal user replaceable parts

In the event of product failure any repair or replacement of parts should only be undertaken by suitably qualified personnel. For advice please contact the Service Department of Bibby Scientific Ltd. quoting model and serial number.

Only spare parts supplied or specified by Bibby Scientific Ltd. or its agents should be used. Fitting of non-approved parts may affect the performance of the safety features designed into the instrument.

Note: The magnetic stirrer drive utilises strong magnets.

If in any doubt, please contact the Service Department of Bibby Scientific Ltd. or the point of sale.

Bibby Scientific Ltd.
Stone, Staffordshire ST15 0SA
United Kingdom
Tel: +44 (0)1785 810475
Fax: +44 (0)1785 810405
e-mail service@bibby-scientific.com
www.bibby-scientific.com

Warranty

Bibby Scientific Ltd warrants these instruments to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of **three (3)** years. In the event of a justified claim Bibby Scientific will replace any defective component or replace the unit free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by installation, adaptation, modification, fitting of non approved parts or repair by unauthorised personnel.

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These products meet the relevant EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in their vicinity will meet these standards

and we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Technical Service Department of Bibby Scientific Ltd.

Declaration of Conformity

Large Hotplate, Models SB300, SB301, SB302, SB500, CB300, CB302, CB500, SD300, SD500

These products comply with the requirements of the EU Directives listed below:

2004/108/EC	EMC Directive.
2006/95/EC	Low voltage Directive (LVD)
2011/65/EC	RoHS Directive

CE mark affixed '02.

Signed:  (Mr S. Marriott)

Date: 17/11/15

Authority: Technical Director
Bibby Scientific Ltd

 Bibby Scientific Ltd - Stone - Staffs - ST15 0SA - UK
Tel: +44 (0) 1785 812121 - Fax +44 (0) 1785 813748
DofC_88_Rev01

INSPECTION REPORT

MODEL _____

SERIAL No. _____


ELECTRICAL SAFETY

1. Earth continuity
2. Insulation
3. Flash test

FUNCTIONAL

1. Indicators
2. Heater/stirrer control
3. Visual acceptance

QUALITY CONTROL INSPECTOR

Bibby Scientific - France

Bâtiment le Deltaparc Silic pars Nord 2,
7 rue du Canal,
BP 55437 Villepinte,
95944 ROISSY Charles de Gaulle,
France
Tel: +33 (0) 1 48 63 78 00
Fax: +33 (0) 1 48 63 78 01
Email: contact@bibby-scientific.com
Web: www.bibby-scientific.fr

Bibby Scientific US

3 Terri Lane,
Suite 10,
Burlington,
NJ 08016, USA
Tel: +1 609 589 2560
Fax: +1 609 589 2571
Email: labproducts@techneusa.com
Web: www.bibby-scientific.com

Bibby Scientific UK (Group HQ)

Beacon Road,
Stone,
Staffordshire
United Kingdom ST15 0SA
Tel: +44 (0) 1785 812121
Fax: +44 (0) 1785 810405
Email: sales@bibby-scientific.com
Web: www.bibby-scientific.com

Bibby Scientific - Middle East

PO Box 27842
Engomi 2433,
Nicosia, Cyprus
Tel: +357 22 660 423
Fax: +357 22 660 424
Email: sales@bibby-scientific.me
Web: www.bibby-scientific.com

Bibby Scientific - Singapore

Prudential Tower,
Level 26
30 Cecil Street
Singapore 049712
Tel: + 65 6631 2976
Fax: +44 (0) 1785 810405
Email: info@bibby-scientific.com
Web: www.bibby-scientific.com