
# Rotary Evaporators

RE300 RE301 RE302

RE300P RE301P RE302P


Assembly and operating instructions  
Instructions de montage et d'utilisation  
Istruzioni per il montaggio e per l'uso  
Instrucciones de ensamblaje y funcionamiento  
Montage- und Bedienungsanleitung

# Rotary Evaporators

## RE300 RE301 RE302

## RE300P RE301P RE302P

## Assembly and operating instructions

### Before use

***If the equipment is not used in the manner described in this manual the protection provided by the equipment may be impaired.***

Bibby Scientific rotary evaporators are designed to operate under the following conditions:

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude up to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply voltage fluctuations not greater than  $\pm 10\%$  of nominal
- ❖ Over voltage category 2 IEC 60364-4-443
- ❖ Pollution degree 2 IEC664

### Location & Services

In order to get the optimum performance from your Bibby Scientific rotary evaporator careful consideration should be given to the intended location. Please study the following notes before commencing assembly.

Choose a convenient location, which will allow easy access to both rotating and receiving flask.

In order to allow the glassware to be raised and lowered, a minimum headroom of 850mm is required.

The location should have access to the following services:

- Cold water supply capable of delivering 60l/hr is required for RE300, RE300P, RE301, RE301P.
- Electricity supply
- Drain
- Vacuum line (optional)

## Electrical Supply

### THIS EQUIPMENT MUST BE EARTHED

Before connection please read and understand this instruction manual and ensure that the line supply corresponds to that shown on the rating plate. Bibby Scientific rotary evaporators require a supply rated at 220-240V, 50/60Hz, single phase, ~. The power consumption of the rotary evaporator is 50W. F3.15AL (20x5mm) fuses fitted in both live and neutral lines.

#### IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN

The unit is supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a U.K. 3 pin plug and the other has a 2-pin "Shuko" plug for connection to the mains. Choose the lead appropriate for your electrical installation and discard the other.

Should neither lead be suitable, take the lead with the U.K. plug and replace the plug with a suitable alternative. This involves cutting off the moulded plug, preparing the cable and connecting to the rewirable plug in accordance with its instructions.

#### IT IS IMPORTANT THAT THIS OPERATION SHOULD ONLY BE UNDERTAKEN BY A QUALIFIED ELECTRICIAN

The wires in the mains cable are coloured as follows:

BROWN	- LIVE
BLUE	- NEUTRAL
GREEN/YELLOW	- EARTH

The unit is fitted with an IEC socket at the rear of the instrument for connection of the mains lead. The appropriate mains lead should be connected BEFORE connection to the mains supply.

Should the mains lead need replacement a cable of 1mm<sup>2</sup> of harmonised code H05W-F connected to an IEC 320 plug should be used.

N.B. The UK mains lead is protected by a 10A fuse mounted in the plug top.

## Safety Advice before Use

- If the jacking mechanism is to be raised without the glassware fitted great care must be taken. Without the weight of the glass the jack will rise rapidly and can cause damage to the mechanism. Unlock the jack carefully while applying downward pressure to control the movement.**
- As the apparatus is likely to be used under vacuum, before operation, examine all glassware carefully for scratches or chemical etching. Use of damaged glassware under vacuum could result in an implosion. In any case it is recommended that the apparatus be operated behind a safety screen.**
- Do not use the equipment in a hazardous atmosphere or for mixing of hazardous materials.**
- Special care should be taken when using the equipment with flammable solvents. The unit is not spark or explosion proof.**
- In case of mains interruption, a fault or mechanical failure, the unit will continue to operate on removal of the fault.**
- Beware of safety issues associated with rotating glassware and when jacking glassware up and down.**


Figure 2


Figure 3

## Assembly

**If the jacking mechanism is to be raised without the glassware fitted great care must be taken. Without the weight of the glass the jack will rise rapidly and can cause damage to the mechanism. Unlock the jack carefully while applying downward pressure to control the movement.**

- Unpack the rotary evaporator and identify the following components:

Item No	Component	Catalogue Number
1	Condenser	RE100/CO or RE200/VC or RE200/CF and RE200/CFD
2	Conical joint clip	KCM29
3	Feed/vacuum release adapter	RE100/VR
4	Florentine flask 1L	FD1L/4RE
5	Receiving flask 1L	RE100/RF/1L
6	Rotary evaporator motor stand	RE300/MS
7	Spherical joint clip	JC35
8	Vapour tube	RE100/VT or RE100/VT/CF
9	Grey retaining cap	RE100/CS
10	Circular metal spring	RE100/RS
11	Vacuum seal	} RE100/VS
12	Sub-seal	
13	Extractor	
-	UK lead	All components are included in the rotary evaporator motor stand box
-	European lead	
-	Spare fuse	

2. Place the stand on a level surface. You will require easy access to a water supply, drain, electric power and vacuum line to use your rotary evaporator.
3. Pass the plain end of the vapour tube (8) through the right side of the motor assembly and push into position until it locates firmly behind the spring clip. The conical joint should be to the right of the motor when viewed from the front.
4. Place the main vacuum seal (11) and sub-seal (12) onto the plain end of the vapour tube (left) and slide them down until they locate in the bearing housing. It is important that these elements are fitted in the correct order and that the spring loaded side of the main seal faces the bearing housing (see fig.2).
5. Place the retaining clip (13) on the jointed side of the vapour tube (right) and slide until it clips in the plastic screw bush (fig 2).
6. Place the grey retaining cap (9) over the butt joint of the condenser (1) with the screwthread facing outwards and retain in position with the circular spring (10) (see figure 3).
7. Place the condenser over the vapour tube so that it butts up to the vacuum seal. Ensure the spherical joint faces downwards. Secure by screwing the retaining cap onto the bearing housing. If using a RE300 diagonal coil condenser, check that the vapour tube is central and will not foul the condenser cooling coil.
8. The condenser may now be set to the required angle:
 

<b>RE300 and RE300P</b>	Diagonal position
<b>RE301, RE301P and RE302, RE302P</b>	Vertical position
9. Support the condenser while loosening the glass angle control.
10. Position condenser and tighten the control.
11. Hold the receiving flask (5) in position and secure with the spherical joint clip (7).
12. Place the rotating flask (4) onto the end of the vapour tube and secure using the conical joint clip (2).
13. Remove the grey plastic screwcap from the narrow end of the condenser and assemble it to the feed/vacuum release adapter (3) (see fig. 4).
14. Pass the PTFE tube through the condenser and vapour tube so that the end is inside the rotating flask. Note that for the RE301 and RE302, the PTFE tube may have to be reduced in length to fit into the Florentine flask (see fig. 4).
15. Secure by tightening the screw cap.
16. **RE302 and RE302P only:** Assemble the drain stopcock (RE200/CFD) to the condenser using the grey plastic cap located on the lower side arm.
17. To connect the condenser to the water supply, drain and vacuum line, flexible hose of 9mm I.D. is required. Ensure you identify every connection correctly before starting. Note tubing suitable for use under vacuum must be used for the vacuum connection.
 

**RE300, RE300P** First remove the plastic connectors from the condenser and fit into the hose. Secure with a suitable clip, e.g. Jubilee. The connectors and hose may now be fitted to the inlet and outlet water connections on the condenser (bottom two connections linked to the coil). Proceed in the same way for the vacuum connection (top connection).

**RE301, RE301P** First remove the plastic connectors from the condenser and fit into the hose. Secure with a suitable clip, e.g. Jubilee. The connectors and hose may now be fitted to the inlet and outlet water connections on the condenser (bottom two connections linked to the coil). Proceed in the same way for the vacuum connection (top connection).

**RE302, RE302P** First remove the plastic connector at the back of the condenser and fit into the hose. Secure with a suitable clip, e.g. Jubilee. The connector and hose may now be fitted to the vacuum line.


Note: A plastic cover is provided to reduce evaporation of cooling agent.

Your rotary evaporator is now ready to be used.

## Operation

### Batch Operation

- a) Remove the rotating flask (4) from the vapour tube (8) and fill with the liquid to be evaporated. Ensure that the flask is never more than half full. Secure the flask onto the conical joint using the clip (2)
- b) Examine the feed/release adapter (3) and ensure that both the stopcock and vacuum release screwcap are closed (see fig. 4).

**Figure 4**

- c) Ensure that water is flowing through the condenser at approx. 60l/hr (RE300 and RE301) or that the cold trap is filled with ice or dry ice (RE302).
- d) Ensure the speed control is set to 0.
- e) Switch on the electricity supply.
- f) Raise the glassware assembly to its highest position by turning the handle anticlockwise. Note that the mechanism is spring loaded. Once at the right height, turn the handle clockwise to secure the position.

Note: Further height adjustment can be made by unscrewing the additional height adjustment knob at the back of the stand (see fig 1).

- g) Place a suitable water bath directly under the rotating flask.

**N.B. A water bath designed for use with the RE300 is available from Bibby Scientific, code RE300B. Please contact the sales office for details.**

- h) Lower the glassware until the rotating flask is partially immersed in the water.

**N.B. Ensure that the glassware does not touch the bath and that no water overflows as the flask is immersed.**

- i) Switch the on/off switch to the on position. The green switch will illuminate.
- j) Set the speed control to the desired speed ensuring that rotation is not so fast that water is splashed out of the bath.

- k) If evaporation at reduced pressure is required a vacuum line should be attached to the vacuum side arm on the condenser. Turn on the vacuum and set to required value.
- N.B. A vacuum pump designed for use with the RE300 is available from Bibby Scientific, code RE2022C. Please contact the sales office for details.**
- l) Set the temperature of the water bath as required.
- m) When evaporation is complete release the vacuum by slowly unscrewing the vacuum release screwcap. Set the rotation speed to 0 and raise the glassware clear of the bath.

### **Continuous Feed Operation**

When the apparatus is used under reduced pressure it is possible to replenish the liquid in the rotating flask without interrupting the evaporation or removing the flask.

The following procedure should be adopted.

- a) Attach a length of flexible tubing to the continuous feed inlet of the feed/release adapter (see fig.4).
- b) Carry out instructions b-k for batch operation.
- c) Immerse the free end of the flexible tubing into the container holding the liquid to be evaporated. Ensure there is no strain on the feed/release adapter.
- d) With the apparatus under vacuum slowly open the stopcock on the feed/release adapter. The vacuum will draw liquid into the rotating flask.
- e) Close the stopcock when sufficient liquid has entered the flask. Further quantities of liquid may be added during evaporation by reopening the stopcock.

**N.B. Ensure there is sufficient capacity in the flask to cope with the addition. Never fill the flask more than half full.**

**N.B. Ensure there is sufficient capacity in the receiving flask to cope with the condensation.**

**N.B. This evaporator should not be used with rotating flasks bigger than 2 litres or receiving flasks bigger than 1 litre capacity.**

- f) If the receiving flask should require emptying during operation the following procedure should be adopted:
  - i Release the vacuum and stop rotation
  - ii Raise the glassware clear of the water bath
  - iii Wait until boiling stops and liquid stops dripping into the receiving flask
  - iv The receiving flask may now be removed.
- g) When evaporation is complete release the vacuum by slowly unscrewing the vacuum release screwcap. Set the rotation speed to 0 and raise the glassware clear of the bath.

## Care and Maintenance

**IMPORTANT: Before commencing any maintenance operation or replacement of components parts the unit should be isolated from the electricity supply, water supply and vacuum line.**

**These operations should only be undertaken by suitably qualified personnel.**

**Note: F3.15AL (20x 5mm) fuses fitted in both live and neutral lines.**

### **Glassware**

All glassware parts should be cleaned regularly and examined for scratches, cracks and chemical etching. Replace any damaged parts.

### **Plastic coated glassware**

The plastic coated glassware should not be exposed to temperature above 80°C and is not dishwasher proof.

### **Cleaning**

The main stand and motor should be cleaned using a mild detergent solution.

### **Vacuum Seal**

It is recommended that the vacuum seal be regularly washed with distilled or deionised water and allowed to dry naturally.

### **Removal of the Vacuum Seal**

- a) Remove the rotating and receiving flasks
- b) Remove feed/release adapter
- c) Remove condenser
- d) Withdraw vapour tube from the bearing housing by completely unscrewing the plastic screw bush on the right side of the bearing house.
- e) The vacuum seal may now be removed.
- f) Reassemble following the assembly instructions in this manual.

## Servicing

It is recommended that any servicing or repair is only undertaken by suitably qualified personnel.

Only spare parts supplied by Bibby Scientific or its agent should be used. Fitting of non-approved parts may affect the performance of the safety features of the instrument.

If in doubt, contact the Technical Service Department of Bibby Scientific Ltd.

Stone, Staffordshire ST15 0SA

United Kingdom

Tel: +44 (0)1785 812121

Fax: +44 (0)1785 813748

e-mail [info@bibby-scientific.com](mailto:info@bibby-scientific.com)

[www.bibby-scientific.com](http://www.bibby-scientific.com)

## Spares

The following spare parts are available from your laboratory supplier. (For a comprehensive parts list or circuit and wiring diagrams please contact the Technical Service Department of Bibby Scientific Ltd quoting model and serial number).

Cat. No.	Description
RE100/CO	Diagonal coil condenser (feed tube not included)
RE200/VC	Vertical coil condenser (feed tube not included)
RE200/CF	Cold finger condenser (feed tube and drain not included)
RE100/COP	Diagonal coil condenser plastic coated (feed tube not included)
RE200/VCP	Vertical coil condenser plastic coated (feed tube not included)
RE200/CFP	Cold finger condenser plastic coated (feed tube and drain not included)
RE200/CFD	Drain for cold finger condenser
RE100/VR	Feed tube
RE100/VT	Vapour tube for diagonal condenser
RE100/VT/CF	Vapour tube for vertical coil and cold finger condensers
RE100/VS	Vacuum seal
FD1L/4RE	Florentine flask 1000ml
FD1L/4REP	Florentine flask 1000ml plastic coated
RE100/RF/1L	Receiving flask 1000ml
RE100/RF/1LP	Receiving flask 1000ml plastic coated
KCM29	Conical joint clip
JC35	Spherical joint clip

## Warranty

Bibby Scientific Ltd warrants this mechanical and electrical equipment to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of three (3) years. In the event of a justified claim Bibby Scientific will replace any defective component free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage by installation adaptation, modification, fitting of non-approved parts or repair by unauthorised personnel.

This warranty does not apply to the glassware or vacuum seal.


These products meet the applicable EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in its vicinity will meet these standards and so we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Technical Department of Bibby Scientific.

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These products meet the relevant EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in their vicinity will meet these

standards and we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Technical Service Department of Bibby Scientific Ltd.

## Declaration of Conformity

### Rotary Evaporator Model RE300, RE300P, RE301, RE301P, RE302, RE302P

These products comply with the requirements of the EU Directives listed below:

**2004/108/EC**      **EMC Directive.**  
**2006/95/EC**      **Low voltage Directive (LVD)**

Compliance with the requirements of these Directives is claimed by meeting the following standards:

**EN 61326-1:2006 (Electrical Equipment for Measurement, Control and Laboratory use).**  
**EN 61010-1: 2001**  
**(Safety Requirements Electrical Equipment for Measurement, Control and Laboratory use)**  
**EN 61010-2-010: 2003 (Particular Requirements for Laboratory Equipment for Heating of Materials).**  
**EN61010-2-051: 2003 (Particular Requirements for Laboratory Equipment for Mixing and Stirring)**

Compliance Certificates and Full Reports.

Ref: RETS0424/1 and RETS0424/2

The above certificates and reports, from an independent test house, are available upon request.

**CE mark affixed '02.**

Signed:  (Mr C. Warren)

Date: JUNE 2009

Authority: Technical Manager  
Bibby Scientific Ltd

 Bibby Scientific

Bibby Scientific Ltd - Stone - Staffs - ST15 0SA - UK  
Tel: +44 (0) 1785 812121 - Fax +44 (0) 1785 813748  
DofC\_88\_Rev01

**Bibby Scientific Ltd**

Beacon Road Stone  
Staffordshire ST15 0SA  
United Kingdom  
Tel: +44 (0)1785 812121  
Fax: +44 (0)1785 813748  
e-mail: [info@bibby-scientific.com](mailto:info@bibby-scientific.com)  
[www.bibby-scientific.com](http://www.bibby-scientific.com)

**Bibby Scientific France SAS**

ZI du Rocher Vert - BP 79  
77793 Nemours Cedex  
France  
Tel: +33 1 64 45 13 13  
Fax: +33 1 64 45 13 00  
e-mail: [bsf@bibby-scientific.fr](mailto:bsf@bibby-scientific.fr)  
[www.bibby-scientific.com](http://www.bibby-scientific.com)

**Bibby Scientific US Ltd**

3 Terri Lane Suite 10  
Burlington NJ 08016  
USA  
Tel: 800-225-9243  
Fax: 609-589-2571  
[www.bibby-scientific.com](http://www.bibby-scientific.com)

**Bibby Scientific Italia Srl**

Via Alcide de Gasperi 56  
20077 Riozzo di Cerro al Lambro  
Milano Italia  
Tel: +39 (0)2 98230679  
Fax: +39 (0)2 98230211  
e-mail: [marketing@bibby-scientific.it](mailto:marketing@bibby-scientific.it)  
[www.bibby-scientific.it](http://www.bibby-scientific.it)

**Bibby Scientific (Asia) Ltd**

Room 607 Yen Sheng Centre  
64 Hoi Yuen Road Kwun Tong  
Kowloon Hong Kong  
Tel: +852 3583 1581  
Fax: +852 3583 1580  
e-mail [bibby@bibby-scientificasia.com](mailto:bibby@bibby-scientificasia.com)  
[www.bibby-scientificasia.com](http://www.bibby-scientificasia.com)