# Thermo Scientific Merlin LR Recirculating Chiller

Thermo Scientific Manual P/N U01230  
Rev. 11/07/2017

**Multilingual  
Quick Starts**

**Multilingual  
Essential Safety  
Instructions**

**Installation**

**Operation**

**Maintenance**

**Troubleshooting**

Visit our Web site at:

**<http://www.thermofisher.com/tc>**  
Product Service Information, Applications  
Notes, SDS Forms, e-mail.

Voice Info: (800) 258-0830

# **Thermo Scientific Merlin LR Recirculating Chiller**

Thermo Scientific Manual P/N U01230  
Rev. 11/07/2017

**Multilingual  
Quick Starts**

**Multilingual  
Essential Safety  
Instructions**

**Installation**

**Operation**

**Maintenance**

**Troubleshooting**

Visit our Web site at:

**<http://www.thermofisher.com/tc>**  
Product Service Information, Applications  
Notes, SDS Forms, e-mail.

Voice Info: (800) 258-0830

## **Thermo Fisher Scientific**

25 Nimble Hill Road  
Newington, NH 03801  
Tel : (800) 258-0830 or  
(603) 436-9444  
Fax : (603) 436-8411  
[www.thermoscientific.com/tc](http://www.thermoscientific.com/tc)

## **Sales, Service, and Customer Support**

25 Nimble Hill Road  
Newington, NH 03801  
Tel: (800) 258-0830  
Sales: 8:00 am to 5:00 pm  
Service and Support: 8:00 am to 6:00 pm  
Monday through Friday (Eastern Time)  
Fax: (603) 436-8411  
[service.tc.us@thermofisher.com](mailto:service.tc.us@thermofisher.com)

Dieselstrasse 4  
D-76227 Karlsruhe, Germany  
Tel : +49 (0) 721 4094 444  
Fax : +49 (0) 721 4094 300  
[info.tc.de@thermofisher.com](mailto:info.tc.de@thermofisher.com)

Building 6, No. 27  
Xin Jinqiao Rd., Shanghai 201206  
Tel : +86(21) 68654588  
Fax : +86(21) 64457830  
[info.china@thermofisher.com](mailto:info.china@thermofisher.com)

### **Statement of Copyright**

Copyright © 2017 Thermo Fisher Scientific. All rights reserved.

This manual is copyrighted by Thermo Fisher Scientific.

Users are forbidden to reproduce, republish, redistribute, or resell any materials from this manual in either machine-readable form or any other form.

# Contents

<b>Preface</b>	<b>iii</b>
Compliance	iii
WEEE	iii
After-sale support	iv
Feedback	iv
Warranty	iv
Unpacking	v
Safety Warnings	vi
Quick Reference Operating Procedures	ix
<b>Chapter 1 Description</b>	<b>1</b>
Specifications	2
<b>Chapter 2 Installation and Operation</b>	<b>7</b>
Site	7
Electrical Requirements	8
Plumbing Requirements	10
Approved Fluids	11
Additional Fluid Information	11
Filling Requirements	12
Reservoir Tank Isolation Valves	12
Controller	13
Start Up/Shut Down	15
Pre-Start	15
Circuit Protector	15
Starting	15
Setpoint	15
Temp Alarms	16
Remote Serial Comm	16
External Sensor	17
Stopping	17
Setup/Tuning Loop	18
<b>Chapter 3 Options</b>	<b>21</b>
Communication Connector	21
Remote Start/Stop	22
External Temperature Sensor	23
Reservoir Level Warning/Fault	24
Low Pump Flow Warning/Fault	24
Additional Accessories	24

<b>Chapter 4</b>	<b>Maintenance.....</b>	<b>25</b>
	Error Codes.....	25
	Displaying Software Version Number.....	26
	Controller PID Values .....	27
	Cleaning .....	27
	Power Cable Mating Ends.....	28
	Leaks.....	28
	Phase Rotation Interlock .....	28
	PD Pressure Relief Valve.....	29
	Pump Strainer .....	30
	BOM Decoder .....	31
	Temperature Sensor Calibration.....	32
	Decommissioning/Disposal.....	36
	Shipment Storage.....	36
<b>Chapter 5</b>	<b>Troubleshooting.....</b>	<b>37</b>
	Checklist.....	37
	Service Assistance.....	39
<b>Appendix A</b>	<b>Water Quality Standards and Recommendations.....</b>	<b>41</b>
<b>Appendix B</b>	<b>NC Serial Communications Protocol .....</b>	<b>47</b>
<b>Declaration of Conformity</b>		
<b>Warranty</b>		

# Preface

## Compliance

The Declaration of Conformity is located in the back of this manual.

## WEEE

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2012/19/EU. It is marked with 'wheelie bin' symbol:

Thermo Fisher Scientific has contracted with one or more recycling/disposal companies in each EU Member State, dispose of or recycle this product through them. Further information on Thermo Fisher Scientific's compliance with these Directives is available at:

[www.thermoscientific.com/WEEERoHS](http://www.thermoscientific.com/WEEERoHS)

## After-sale Support

Thermo Fisher Scientific is committed to service both during and after the sale. If you have questions concerning operation, contact our Sales Department. If your chiller fails to operate properly, or if you have questions concerning spare parts or Service Contracts, contact our Service Department.

Before calling, please obtain the following information:

- chiller BOM number (see page 32)
- chiller serial number
- chiller software version ( see page 26)
- voltage of power source

The chiller's BOM and serial number are located on the name plate label on the rear of the chiller.

## Feedback

We appreciate any feedback you can give us on this manual. Please e-mail us at [tcmanuals@thermofisher.com](mailto:tcmanuals@thermofisher.com). Be sure to include the manual part number and the revision date listed on the front cover.

## Warranty

Thermo Scientific Merlin chillers have a warranty against defective parts and workmanship for 24 months (**12 months for positive displacement pumps**) from date of shipment. See back page for more details.

## Unpacking

Retain all cartons and packing material until the chiller is operated and found to be in good condition. If the chiller shows external or internal damage contact the transportation company and file a damage claim. Under ICC regulations, this is your responsibility.

### **Out of Box Failure**

An Out of Box Failure is defined as any product that fails to operate in conformance with seller's published specifications at initial power up. The chiller must be installed in accordance with manufacturer's recommended operating conditions within 30 days of shipment from the seller.

Any Temperature Control product meeting the definition of an Out of Box Failure must be packed and shipped back in the original packaging to Thermo Fisher Scientific for replacement with a new chiller; Seller to pay the cost of shipping. Customer must receive a Return Material Authorization (RMA) from Thermo Fisher prior to shipping the chiller.

## Safety Warnings

Make sure you read and understand all instructions and safety precautions listed in this manual before installing or operating your chiller. If you have any questions concerning the operation of your chiller or the information in this manual, please contact us. See inside cover for contact information.

**DANGER** indicates an imminently hazardous situation which, if not avoided, *will* result in death or serious injury.

**WARNING** indicates a potentially hazardous situation which, if not avoided, *could* result in death or serious injury.

**CAUTION** indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It is also be used to alert against unsafe practices.

The lightning flash with arrow symbol, within an equilateral triangle, is intended to alert the user to the presence of non-insulated "dangerous voltage" within the chiller's enclosure. The voltage magnitude is significant enough to constitute a risk of electrical shock.

This label indicates read the manual.

**Note** The chiller's equipment design incorporates a complete sheet metal enclosure for personnel protection from mechanical and electrical hazards. ▲

**Observe all warning labels.** ▲

**Never remove warning labels.** ▲

**Never place the chiller in a location where excessive heat, moisture, or corrosive materials are present.** ▲

**The chiller construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided.** ▲

**Never connect the process fluid inlet or outlet fittings to your building water supply or any water pressure source.** ▲

**On chillers with MD pumps never completely restrict flow to your application. Dead-heading the pump damages its coupling and requires pump replacement.** ▲

**Never use corrosive fluids. Highly distilled and deionized water may be aggressive and cause material corrosion. Please contact Thermo Fisher before subjecting this chiller to prolonged exposure to highly distilled or deionized water. ▲**

**Do not use automotive antifreeze. Commercial antifreeze contains silicates that can damage the pump seals. Use of automotive antifreeze voids the manufacturer's warranty. ▲**

In addition to the specific warnings listed on the previous page the following general warnings also apply to your chiller:

**Performance of installation, operation, or maintenance procedures other than those described in this manual may result in a hazardous situation and may void the manufacturer's warranty. ▲**

**Transport the chiller with care. Sudden jolts or drops can damage the its components. ▲**

**If the chiller is to be transported and/or stored it needs to be drained. If the chiller is to be transported and/or stored in cold temperatures it needs to be drained and then flushed with a 50/50 glycol/water mixture. ▲**

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**Observe all warning labels. ▲**

**Never remove warning labels. ▲**

**Never operate damaged or leaking equipment. ▲**

**Never operate the chiller without cooling fluid in the reservoir. ▲**

**Always turn off the chiller and disconnect the power cord from the power source before performing any service or maintenance procedures, or before moving the chiller. ▲**

**Never operate equipment with damaged power cords. ▲**

**Refer service and repairs to a qualified technician. ▲**

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## **Fluorinated Greenhouse Gases**

The following information is included to comply with REGULATION (EU) No 517/2014 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 April 2014 on fluorinated greenhouse gases:

**This product contains fluorinated greenhouse gases in a hermetically sealed system.**

**If a leak in the sealed system is detected, the operator shall repair without undue delay.**

Refer to the F-Gas Declaration of Conformity for additional information.

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**Declaration of conformity with Article 14 of Regulation (EU) No 517/2014 of the European Parliament and of the Council**

We, Thermo Fisher Scientific (Asheville) LLC, F-Gas Portal Registration Number 23643, declare under our sole responsibility that when placing on the market pre-charged equipment, which we import to or manufacture in the Union, the hydrofluorocarbons contained in that equipment are accounted for within the quota system referred to in Chapter IV of Regulation (EU) No 517/2014 as:

We hold authorisation(s) issued in accordance with Article 18(2) of Regulation (EU) No 517/2014 and registered in the registry referred to in Article 17 of that Regulation, at the time of release for free circulation to use the quota of a producer or importer of hydrofluorocarbons subject to Article 15 of Regulation (EU) No 517/2014 that cover(s) the quantity of hydrofluorocarbons contained in the equipment.

The hydrofluorocarbons contained in the equipment have been placed on the market in the Union, subsequently exported and charged into the equipment outside the Union, and the undertaking that placed the hydrofluorocarbons on the market made a declaration stating that the quantity of hydrofluorocarbons has been or will be reported as placed on the market in the Union and that it has not been and will not be reported as direct supply for export in the meaning of Article 15(2)(c) of Regulation (EU) No 517/2014 pursuant to Article 19 of Regulation (EU) No 517/2014 and Section 5C of the Annex to Commission Implementing Regulation (EU) No 1191/2014.

January 1, 2017

\_\_\_\_\_  
Date

Mark Pearson, Director Global Regulatory Affairs

\_\_\_\_\_  
[name and position of legal representative]



\_\_\_\_\_  
[signature of legal representative]

Manufacturer:

Thermo Fisher Scientific (Asheville) LLC  
275 Aiken Road  
Asheville, NC 28804  
U.S.A.

EU Only Representative:

Thermo Electron LED GmbH  
Robert-Bosch-Strasse 1  
D-63505 Langenselbold  
Germany  
**VAT ID Number: DE 812 403 137**

## Quick Reference Operating Procedures

### Installation

The chiller has an air-cooled refrigeration system. Air is drawn in the front of the chiller and discharged through rear and sides. Position the chiller so the intake and discharge are not impeded. Inadequate ventilation will reduce cooling capacity and, in extreme cases, result in compressor failure.

Avoid excessively dusty areas and institute a periodic cleaning schedule. For proper operation, the chiller needs to pull substantial amounts of air through a condenser. A build up of dust or debris on the fins of the condenser will lead to a loss of cooling capacity.

The chiller will retain its full rated capacity in ambient temperatures up to approximately +77°F (+25°C).

Make sure the voltage of the power source meets the specified voltage,  $\pm 10\%$ .

The plumbing connections are located on the rear of the chiller and are labelled  and . These connections are 1/2 inch FPT (3/4 inch FPT for chillers with CP-55 pumps). Remove the plastic protective plugs from both plumbing connections. Connect the  fitting to the inlet of your application. Connect the  fitting to the outlet of your application.

To fill the reservoir, remove the reservoir access panel and locate the reservoir cap. Remove the cap and fill the reservoir with clean cooling fluid.

For approved fluids refer to page 11.

### Operation

Before starting the chiller, double-check all electrical and plumbing connections. Make sure the circulating system is properly filled with cooling fluid. Place the circuit breaker located on the rear of the chiller to the up position.

To start the chiller, press . To turn the chiller off press  again.

The  LED indicates the status of the refrigeration system. It illuminates to indicate the refrigeration system is removing heat from the cooling fluid. As the operating temperature approaches the set point, the LED flashes.

### Temperature Adjustment

To display the temperature set point, press  on the controller. The  indicator illuminates and the display flashes the current set point value.

To adjust the temperature set point, press the arrow keys until the desired temperature set point is indicated. Press  again to confirm the change. The display rapidly flashes the new value for a short time and then returns to the recirculating fluid temperature.

### Periodic Maintenance

Periodically inspect the reservoir fluid. If cleaning is necessary, flush the reservoir with a cleaning fluid compatible with the circulating system and the cooling fluid.

Replace the cooling fluid periodically. Frequency depends on the operating environment and amount of usage.

Before changing the cooling fluid ensure it is at a safe temperature.

Periodic vacuuming of the condenser fins is necessary. The frequency of cleaning depends on the operating environment. We recommend a monthly visual inspection of the condenser after initial installation. After several months, the cleaning frequency will be established.

Chillers with PD or T1 pumps have a strainer. If debris is in the system, the strainer prevents the material from being drawn into the pump and damaging the pump vanes. A clogged strainer also causes increased pump discharge pressure.

After initial installation, the strainer may become clogged. **Clean the strainer after the first week of installation.** After this first cleaning, we recommend a monthly visual inspection. After several months, the cleaning frequency will be established. Before cleaning the strainer, disconnect the power cord from the power source and drain the chiller.

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
# Chapter 1 Description

The Thermo Scientific Merlin LR Recirculating Chiller is designed to provide a continuous supply of cooling fluid at a constant temperature and flow rate. The chiller consists of an air-cooled refrigeration system, plate heat exchanger, recirculating pump, polyethylene reservoir, and a microprocessor controller.

Throughout the manual, you will be asked to consult the chiller's name plate label for specific information. The label is located on the rear of the chiller. See below for an example.



## Specifications

	M 25	M 33	M 75	M 150
<b>Setpoint Temperature Range</b>	-15°C to +35°C	-15°C to +35°C	-15°C to +35°C	-15°C to +35°C
<b>Ambient Temperature Range</b>	+10°C to +35°C	+10°C to +35°C	+10°C to +35°C	+10°C to +35°C
<b>Temperature Stability</b>	±0.1°C	±0.1°C	±0.1°C	±0.15°C
<b>Reservoir Volume</b>				
Gallons	0.5	0.5	0.5	0.5
Liters	1.8	1.8	1.8	1.8
<b>Refrigerant</b>	R134a	R134a	R134a	R404a
<b>Weight<sup>1</sup></b>				
Pounds	120	120	176	254
Kilograms	54.5	54.5	79.8	115.2

Thermo Fisher Scientific reserves the right to change specifications without notice.

1. Chillers with PD1 and PD2 pumps. M33 chillers with T1 pumps weigh 141 pounds (64 kilograms). M75 chillers with T1 pumps weigh 200 pounds (91 kilograms). M150 chillers with TU1 weigh 310 pounds (141 kilograms) or with CP55 pumps weigh 275 pounds (124.7 kilograms).




**Cooling Capacity\***


Thermo Fisher Scientific reserves the right to change specifications without notice.

\*Specifications obtained at sea level using 50/50 ethylene glycol/water mixture as the recirculating fluid, at a +20°C process setpoint, +20°C ambient condition, at nominal operating voltage. Other fluids, process temperatures, ambient temperatures, altitude or operating voltages will affect performance. Cooling capacity based on chillers with PD1 pumps (PD2 pumps for M150) with no backpressure. Other pumps will affect cooling capacity performance.

**Pumping Capacities**

Dimensions (Inches)	M 25	M 33	M 75	M 150
Dimension A <sup>1</sup>	23 1/2	23 1/2	26 1/8	30 3/8
Dimension B <sup>1</sup>	12 5/8	12 5/8	16 1/4	21 1/4
Dimension C <sup>1</sup>	20 7/8	20 7/8	24	29 1/4
Dimension D <sup>2</sup>	2	2	2	3
Dimension E <sup>3</sup>	18 1/8	18 1/8	20 3/4	25
Dimension F <sup>4</sup>	3	3	3	3
Dimension G <sup>5</sup>	16 1/4	16 1/4	18 7/8	23
Dimension H <sup>6</sup>	1 3/8	1 3/8	1 3/8	1 3/8
Dimension I <sup>7</sup>	2 7/8	2 7/8	3 5/8	4 3/4
Dimension J <sup>8</sup>	5	5	7	7 1/2
Dimension K <sup>9</sup>	7 1/8	7 1/8	7 1/8	7 1/8
Dimension L <sup>10</sup>	1 3/4	1 3/4	1 3/4	1 3/4
Crate Dimensions (H x W x D) <sup>11</sup>	36 x 18½ x 26	36 x 18½ x 26	39 x 22¼ x 29¼	42 x 27 x 35¼

1. Dimension A is the height of the chiller. Chillers with CP-55 pumps are 6" higher. The chiller width and depth (dimensions B and C) are the case dimensions, add approximately 1/2 inch to include the plumbing connections.
2. Dimension D is the distance from the floor to the bottom of the chiller case (height of the casters).
3. Dimension E is the distance from the floor to the center of the outlet connection.
4. Dimension F is the distance between the center of the outlet and inlet connections.
5. Dimension G is the distance from the floor to the center of the pressure relief valve (chillers with PD pumps).
6. Dimension H is the distance from the chiller's side to the center of the outlet/inlet connections, and the pressure relief valve. 2" from the chiller's side to the center of the outlet/inlet connections for M-150 chillers with CP-55 pumps.
7. Dimension I is the distance from the chiller's side to the center of the drain connection. 5 1/4" for M-150 chillers with CP-55 pumps.
8. Dimension J is the distance from the floor to the center of the drain connection. 7 1/2" for M-150 chillers with CP-55 pumps.
9. Dimension K is the distance from the top of the chiller to the center of the power connection.
10. Dimension L is the distance from the chiller's side to the center of the power connection. NA for M-150 chillers rated 400V/50Hz.
11. Crate dimensions for chillers with CP-55 pumps are 6 inches higher. All packages tested to ISTA 2B.

## Chapter 2 Installation and Operation

### Site

Ambient Temperature Range*	10°C to 35°C (50°F to 94°F)
Relative Humidity Range	10% to 80% (non-condensing)
Operating Altitude	Sea Level to 6000 feet (1830 meters)
Overvoltage Category	II
Pollution Degree	2


\*The chiller retains its full rated capacity in ambient temperatures up to approximately 25°C (77°F). Reduce the cooling capacity 1% for every 0.5°C (1°F) above 25°C (77°F), up to a maximum ambient temperature of 35°C (94°F).

**Never place the chiller in a location where excessive heat, moisture, or corrosive materials are present. ▲**

The chiller has an air-cooled refrigeration system. Air is drawn through the front of the chiller and discharged through rear and side panels. The chiller must be positioned so the intake and discharge are not impeded. A minimum clearance of 3 feet (1 meter) on all vented sides is necessary for adequate ventilation. Inadequate ventilation will cause a reduction in cooling capacity and, in extreme cases, compressor failure.

Direction of Airflow

Excessively dusty areas should be avoided and a periodic cleaning schedule should be instituted (see Chapter 4, Maintenance). Optional air filters are available, contact our Service Department. See Preface, After-Sale Support.

## Electrical Requirements

The chiller construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided. ▲

The user is responsible to ensure that the power cord provided meets local electrical codes. If not, contact qualified installation personnel.

The chiller is intended for use on a dedicated outlet. The Merlin has an internal circuit protection that is equivalent (approx.) to the branch circuit rating. This is to protect the Merlin, and is not intended as a substitute for branch circuit protection.

60 Hertz chillers available power options:

Chiller	Pump	Voltage	Current Value	Line Cord Plug
M 25	MD 30	115	8.1	5-15P
	PD 1	115	9.8	5-15P
M 33	MD 30	115	11.5	5-15P
	PD 1	115	13.2	5-20P
	PD 2	115	15.5	5-20P
	PD 1	208/230	6.5	6-15P
	PD 2	208/230	7.8	6-15P
M 75	PD 1	208/230	8.8	6-15P
	PD 2	208/230	10.1	6-15P
	T1	230	10.3	IEC receptacle # 001395
M 150	PD 2	208/230	15.4	L6-30P
	T1	230	15.6	L6-30P
	CP 55	208/230	18.4	L6-30P

**NOTE** All 60 Hertz chillers are UL Listed and Certified to Canadian Standards. ▲

Amperage based on extreme operating conditions. Operating in a 20°C ambient and at nominal voltage results in a lower amperage draw.

Refer to the chiller's name plate label for specific electrical requirements.

**NOTE** Three-phase chillers are not supplied with a line cord. Instructions for wiring three-phase chillers are located behind the power box cover. To access the power box remove the chiller's top and left-side panel. Remove the screws securing the cover. ▲

50 Hertz chillers:

Chiller	Pump	Voltage	Current Value	Line Cord Plug
M 25	PD 1	230	5.1	Country Specific
M 33	PD 1	230	6.3	Country Specific
	PD 2	230	7.6	Country Specific
	T1	230	7.0	IEC receptacle # 001395
M 75	PD 1	230	8.8	Country Specific
	PD 2	230	10.1	Country Specific
	T1	230	10.1	IEC receptacle # 001395
M 150	PD 2	230	15.6	Country Specific
	PD 2	400/3Ø/5W	7.7	Country Specific
	CP 55	400/3Ø/5W	5.5	Country Specific

Amperage based on extreme operating conditions. Operating in a 20°C ambient and at nominal voltage results in a lower amperage draw.

Refer to the name plate label on the rear of the chiller for specific electrical requirements.

**NOTE** Three-phase chillers are not supplied with a line cord. Instructions for wiring three-phase chillers are located behind the power box cover. To access the power box remove the chiller's top and left-side panel. Remove the screws securing the cover. ▲


## Plumbing Requirements

The plumbing connections are located on the rear of the chiller and are labeled  (supply) and  (return). The connections are 1/2 inch Female Pipe Thread. For chillers with CP-55 pumps, the connections are 3/4 inch Female Pipe Thread. Chillers with 1/2 inch fittings are supplied with 3/8 inch and 1/2 inch barbed adapters, chillers with 3/4 inch fittings are supplied with 1/2 inch and 3/4 inch barbed adapters.

Remove the plastic protective plugs from both plumbing connections. Install the barbed adapters to these connections.

Connect  to the hose feeding the inlet of your application. Connect  to the hose from the outlet of your application. Clamp all connections.

Plumbing Connections

**Never connect the fittings to your building water supply or any water pressure source. ▲**

**NOTE** On chillers equipped with PD pumps, ensure your plumbing is rated to withstand 110 psi at the highest operating temperature. ▲

**On chillers with MD pumps never completely restrict flow to your application. Dead-heading the pump damages its coupling and requires pump replacement. ▲**

It is important to keep the distance between the chiller and the instrument being cooled as short as possible. Tubing should be straight and without bends. If diameter reductions must be made, do them at the inlet and outlet of your application, not at the chiller.

All chillers have a ball-valve reservoir drain located on the rear. Add the desired length of tubing to the drain valve. To help prevent any spillage from your application, consider installing external isolation valves on the supply and return lines.

## Approved Fluids

Never use corrosive fluids with this chiller. Do not use automotive antifreeze. Commercial antifreeze contains silicates that can damage the pump seals. Use of automotive antifreeze voids the manufacturer's warranty. ▲

Approved fluids and their normal operating temperature ranges are:

- Filtered/Single Distilled water, +10°C to +35°C
- 50/50 Uninhibited Ethylene Glycol/Water, -15°C to +35°C
- 50/50 Inhibited Ethylene Glycol/Water, -15°C to +35°C
- 50/50 Uninhibited Propylene Glycol/Water, -15°C to +35°C
- 50/50 Inhibited Propylene Glycol/Water, -15°C to +35°C
- Deionized water (1 - 3 MΩ-cm, compensated), +10°C to +35°C (For applications requiring resistivity greater than 1 MΩ-cm please call and speak to an applications engineer for additional information.)

Ethylene glycol (EG) is hygroscopic, it absorbs water from its environment. This can affect the freezing point and boiling point of the fluid over time. This may result in system failure. ▲

When using EG/water or PG/water, top-off with EG/water or PG/water. Do not top-off with plain water. Topping-off with plain water can severely affect the freezing point and boiling point of the fluid. This may result in system failure. ▲

Do not use a Deionization (DI) filter with Inhibited EG or Inhibited PG. A DI filter removes inhibitors from the solution rendering the fluid ineffective against corrosion protection. Also, inhibitors increase fluid conductivity. ▲

When using the Merlin LR chiller to circulate through aluminum, use a compatible corrosion inhibitor to prevent galvanic corrosion.


Fluid viscosity should be 50 cSt or less at the lowest temperature used.

Refer to Appendix A for additional information.

## Filling Requirements

The polyethylene reservoir is translucent for easy fluid level monitoring. Remove the reservoir access panel, see page 10. Locate and remove the reservoir cap. Fill the reservoir with clean cooling fluid.

Since the reservoir capacity may be small compared to your application, have extra cooling fluid on hand to keep the system topped off when external circulation is started.

## Reservoir Tank Isolation Valves

Chillers are equipped with two reservoir tank isolation valves. The valves are used to isolate the tank from the fluid flow. The chillers are factory shipped with the valves open.


Ensure the valves are open when circulating to a closed system. When operating at low temperatures and circulating to a closed system, should condensation appear on the reservoir, closing the upper valve *after initial start up* helps eliminate it. Leave the valve open for initial start up so process fluid flow through the reservoir can rapidly bleed air from the system. To close the upper valve, open the reservoir access panel and turn the valve 90°.

**System temperature changes with the isolation valves closed, and while circulating to a closed system, can create a vacuum and/or high pressure conditions. ▲**

## Controller

The controller controls temperature using a Proportional-Integral-Derivative (PID) algorithm. It is designed with self-diagnostic features and easy to use operator interface.

Use this key to toggle the **chiller on or off**.

Use this key to **scroll** through the controller's LEDs. It is also used to accept and save changes.

If the alarm sounds, use this key to toggle the **alarm off and on**. **NOTE:** If the cause of the alarm is cleared but then reoccurs, the alarm sounds again.

Use this key to **increase displayed numerical values**.

Use this key to **decrease displayed numerical values**.

Indicates **refrigeration system status**. It illuminates to indicate the refrigeration system is removing heat from the cooling fluid. As the operating temperature approaches the temperature setpoint, the LED flashes. The indicator is off when heat is not being removed.

**NOTE** The following indications are not visible until  is depressed. ▲

 Indicates the controller is displaying the **set point**. Adjust the display value using the arrow keys.

**NOTE** You cannot adjust the set point closer than 2°C to either of the temperature limits discussed below. ▲

 Indicates the controller is displaying the **low temperature alarm** setting. The indicator flashes and the alarm sounds if this limit is exceeded. Adjust the display value using the arrow keys.


 Indicates the controller is displaying the **high temperature alarm** setting. The indicator flashes and the alarm sounds if this limit is exceeded. Adjust the display value using the arrow keys.

 **(Optional)** Indicates the chiller is in the **remote/serial communication** mode of operation. The indicator flashes if the chiller is in the remote start/stop mode of operation. Toggle between the local and remote modes using the arrow keys.

**NOTE** For all four indicators, if neither arrow key is pressed within 10 seconds the display returns to the recirculating temperature display. Save any change using . After pressing an arrow key, if  is not pressed within 60 seconds the display returns to the recirculating temperature display and ignores any change. ▲

## Start Up/Shut Down

### Pre-Start

Before starting the chiller, double check all electrical and plumbing connections. Ensure the power cord is secured to the rear of the chiller. Have extra recirculating fluid on hand. If the chiller will not start refer to Chapter 5 Troubleshooting.

### Circuit Protector

Place the circuit protector located on the rear of the chiller to the up position, the controller flashes and the alarm momentarily sounds.

**The chiller has automatic restart. If the chiller was shutdown as a result of a power failure and power is restored, it will restart. ▲**

### Starting

Press . The controller does a self-test (the controller quickly sequences through its LEDs and momentarily sound the alarm) and then display the recirculating fluid temperature. The refrigeration system and the recirculation pump then starts. The RECIRCULATING PRESSURE gauge displays the pump operating pressure. If the pressure needs adjusting, refer to page 29.

**NOTE** If on start up the chiller's recirculating fluid is outside either temperature limit, the chiller operates but the appropriate indicator flashes until the fluid is within the limit. ▲

### Setpoint

To display/change the set point press  until  illuminates. The display flashes the current set point value. Use the arrow keys to change the value. The range is -15°C to 35°C.

**NOTE** If the arrow keys are not pressed within 10 seconds the display returns to the current reservoir temperature. The controller will not allow you to enter a set point closer than 2°C of either temperature alarm setting discussed on the next page. Trying to use a set point within 2°C causes the appropriate indicator to flash and sounds the audible alarm twice. ▲

Once the desired set point is displayed, press  again to confirm the change. The display rapidly flashes the new value for a short period and then returns to the recirculating fluid temperature.

**NOTE** If the value is not confirmed within 60 seconds the display returns to the recirculating fluid temperature and ignores any changes. ▲

## Temp Alarms

To display/change the temperature alarm setting press  until either  or  illuminates. The display flashes the current limit value. Use the arrow keys to change the value. The low-end range is -20°C to 30°C. The high-end range is -10°C to 40°C. You cannot set a low limit higher than the high limit or a high limit lower than the low limit.

**NOTE** If the arrow keys are not pressed within 10 seconds the display returns to the current reservoir temperature. You cannot set either alarm closer than 2°C of the set point. ▲

Once the desired setting is displayed, press  again to confirm the change. The display rapidly flashes the new value for a short period and then returns to the recirculating fluid temperature.

**NOTE** If the value is not confirmed within 60 seconds the display returns to the recirculating fluid temperature and ignores any changes. ▲

## Remote/Serial Comm

**(Optional)** Once your remote/serial comm cable is connected to the chiller, enable or disable remote operation or serial communication by pressing  until  illuminates. The display flashes the communication status. Use the up arrow key to enable, the down arrow to disable. See pages 21 and 22 for additional information.

**NOTE** If the arrow keys are not pressed within 10 seconds the display returns to the current reservoir temperature. ▲

Once the desired status is displayed, press  again to confirm the change. The display rapidly flashes the new status for a short period and then returns to the recirculating fluid temperature.

**NOTE** If the status is not confirmed within 60 seconds the display returns to the recirculating fluid temperature and ignores any change. ▲

**NOTE** Enabling serial communications disables the remote start/stop feature. In either mode, the set point and alarm values can not be changed using the controller but they can still be viewed. ▲

## External Sensor

**(Optional)** Once the external sensor cable is connected to the chiller, enable or disable the sensor by pressing  a fifth time. The display alternates with *rPrb* and either *on* or *off*. If desired, use the arrow keys to change the status and then press  to return to the temperature display. See page 23 for additional information.

## Stopping

To turn the chiller off, press .

When the display goes blank it is safe to place the circuit protector located on the rear of the chiller to the down position.

**Using any other means to shut the chiller down can reduce the life of the compressor. ▲**

**The circuit protector located on the rear of the chiller is not intended to be used as a disconnecting means. ▲**

**Always turn the chiller off and disconnect it from its supply voltage before moving it. ▲**

## Setup/Tuning Loop

The Setup/Tuning Loop, see illustration next page, is used to configure the controller temperature display and various operating parameters. To enter the loop you must be displaying the reservoir fluid temperature and then, while pressing and holding  $\blacktriangledown$ , press  $\odot$ . The display indicates *tunE*. Use the arrow keys to enter/bypass the loop. Once in the loop press  $\odot$  to sequence through it. Use the arrow keys to change any display.

The loop is used to determine how the chiller reacts when a fault occurs (*FLTS*) — either shut down (*on*) or continue to run (*off*). The chiller is shipped configured to run.

It is also used to configure the temperature displays to indicate to a tenth of a degree (*tEnt*).

The loop is used to select which temperature sensor the high and low temperature monitors uses (*APrb*). *rtd1* is the internal sensor, *rtd2* is the optional external sensor. *rtd2* is available only if the external sensor is properly connected to the rear of the chiller.

If the optional external/remote probe is enabled, then the loop is also used to set the remote probe temperature deviation (*rPrd*) limit. This limit is applied to the setpoint to calculate the minimum fluid temperature supplied to your application. For example, if the setpoint is 10°C and the deviation is set to 3°C, the coldest fluid temperature supplied to your application is limited to 7°C.

The deviation range is 1°C to 53°C. The factory default setting is 50°C.

The loop can also be used to set the controller's PID parameters (*COOL HEAT*). The factory preset values are on page 27.


**Thermo Fisher recommends that only a qualified technician change PID parameters. Incorrect PID values hamper chiller performance. ▲**

The heat PID is only applicable to chillers configured with a heater.

The loop is also used to adjust/verify serial comm parameters (*SEr*). This section of the loop is applicable only to chillers configured with serial communications capability.

Changes made in the Tuning/Setup Loop take affect when ▲ is pressed at the *Star* prompt. Pressing  $\blacktriangledown$  at the prompt aborts *all* changes.

## Chapter 3 Options

### Communication Connector

A female 9-pin D-connector is located on the rear of the chiller. The connector is used for RS-232 serial communication (see Appendix B).

When the chiller is configured for serial communication the  LED on the controller is illuminated. If the chiller is shut off while still configured for serial communication, the temperature display goes blank, but the  LED remains illuminated.

When the chiller is configured for serial communication it can be stopped using the local controller. To restart, send another start command.

If the chiller is running and serial communications is enabled (  LED on), and you need to operate the chiller using the keypad – the procedure is to press and hold  for approximately two seconds until the low temperature alarm icon illuminates. Press ,  to scroll to the  LED, press the down arrow and then press  again to turn communications off. The  LED will extinguish and the chiller can now be operated from the keypad.

**NOTE** If you need to start the chiller using the local controller, simultaneously depress and hold both arrow keys for approximately 10 seconds. The display shows the temperature, and the alarm sounds. Use  and scroll to the  LED and turn it OFF. The chiller can now be started and operated from the keypad. ▲

The pin out information is:


Pin	RS-232 COMM
1	No connection
2	TX
3	RX
4	No connection
5	GND = Signal ground
6-9	No connection

TX = Transmitted data from controller

RX = Received data to controller.

Hardware Mating Connector AMP Part# 745492-2 or equivalent

## Remote Start/Stop

A male 9-pin D-connector is located on the rear of the chiller. The connector is used for either 12 - 24VAC or 12 - 24VDC remote start/stop; and as a status relay dry contact.

To start the chiller from a remote location apply the correct voltage to pins 6 & 8. When the chiller is started by a remote start/stop voltage  on the controller flashes.

Once started remotely, stop the chiller by removing the applied voltage from pins 6 & 8 or by pressing  on the local controller. To restart remotely after stopping locally, remove and reapply the remote start voltage to J101.

**NOTE** Enabling serial communications disables the remote start/stop feature. ▲

Should the chiller experience a high/low temperature warning/fault, a reservoir level warning/fault, or a low pump flow warning/fault, the status relay de-energizes.

The pin out information is:

Pin #	Function	Pin #	Function
1	No Connection	6	Remote On +
2	Status, NO	7	RTD2B
3	Status, NC	8	Remote On -
4	RTD2A	9	No Connection
5	Status, Common		

Hardware Mating Connector AMP Part# 745491-2 or equivalent

## Remote Start Input

12 to 24 Volts DC  $\pm 10\%$  maximum input current, 9ma.

12 to 24 Volts AC  $\pm 10\%$  maximum input current, 5ma.

## Absolute Maximum Input Voltages

40 volts rms. continuous.

56 volts peak.

Rating for both normally open and normally closed status relay contacts are: Maximum rated load - 1.0 Amp, 24 Volts AC/DC, resistive. Minimum permissible load - 1 ma @ 5 VDC.

## External Temperature Sensor

Connect the external sensor to the male 9-pin D-connector located on the rear of the chiller, see previous page.

The external temperature sensor is enabled using the controller. With the controller displaying the temperature press  five times. The display alternates with *rPrb* and either *on* or *off*. If desired, use the arrow keys to change the status and then press  to return to the temperature display.

If the external sensor is selected the controller displays the temperature at the external sensor. Also, the controller indicates the external sensor is selected by flashing the left-most temperature decimal point.

Flashes to indicate the external sensor is enabled.

**NOTE** To display the internal temperature press and hold  for one second. The display indicate *rt di* for ½ second and then displays the internal temperature for five seconds. The display then indicate *rt d2* for ½ second and returns to the external temperature. The  key needs to be released and then pressed and held again to repeat the cycle. To display the remote probe temperature deviation use the Setup/Tuning Loop, see page 18. 

The controller continues to use the existing setpoint. When using the external sensor the controller will not allow internal chiller temperatures outside -17°C to 37°C.

The selected sensor is retained by the controller if it is turned off, or if power is removed.

If there is an external sensor error the chiller automatically returns to the internal sensor. An error code flashes and if the fault mode is on the chiller shuts down. See Error Codes on page 25. Press  to restart the chiller.

The external probe can be used to control the high and low temperature monitors using the Setup/Tuning Loop, see page 18.

**NOTE** If you also need to monitor status while using the external sensor, an optional Y-connector cable is available from Thermo Fisher. 

## Reservoir Level Warning/Fault

Should the reservoir fluid level drop below normal the mutable alarm sounds, the status relay de-energizes (if the chiller is equipped with the communications package), an error message is displayed and, depending on the chiller's configuration, it continues to run or shut down. The chiller is configured to run or shut down using the controller's Setup/Tuning Loop on page 18.

Refill the reservoir and, if necessary, restart the chiller.

**NOTE** The chiller will not start without any fluid in the reservoir regardless of how it is configured. ▲

## Low Pump Flow Warning/Fault

Should the pump flow drop below 0.5 gpm (1.0 gpm for chillers with a CP-55 pump) the mutable alarm sounds, the status relay de-energizes (if the chiller is equipped with the communications package), an error message is displayed, and, depending on the chiller's configuration, it continues to run or shuts down.

Once the cause of the low flow is identified and corrected restart the chiller, if necessary.

The chiller is configured to run or shut down using the controller's Setup/Tuning Loop, see page 18.

## Additional Accessories

### External Filtration Package

A partial flow liquid filter available in 5, 10, 25, and 40 micron filter sizes.

### External DEI Package

Partial flow DEI package will maintain between 1 and 3 M $\Omega$ cm water resistivity. A cartridge indicator lets you know when the filter needs changing.

### Condenser Air Filter

An easy to install air filter protects the condenser in dusty environments.

### External Strainer Accessory Kit

An easy to install kit to make routine strainer preventative maintenance easier.

## Chapter 4 Maintenance

This chapter covers the chiller's basic maintenance.

### Error Codes

- Er 00* ROM checksum. Clear with  key.
- Er 01* Test failure. Locks up the program.
- Er 02* Display board failure. Clears when display board is fixed.
- Er 03* Critical checksum failure. Clear with  key.
- Er 14* Synchronous comm - check connections.
- Er 15* Asynchronous comm
- Er 1b* Bad calibration data - redo calibration

### Unusual Hardware Conditions

These errors flashes on the display and cannot be cleared. These are internal controller problems. *Er 04* through *Er 13* - Interrupt errors during runtime. *Er 0F* BOM invalid.

### Functional/ Machine errors

These errors clear themselves once the problem disappears.

### Display

- L 0t* Low temp set point warning/fault.
- R 0d* Low level warning, fluid below minimum fill level, see page 12.
- H 1t* High temp set point warning/fault.
- Er 19* Lo t Variable low temp - Customer set.
- Er 21* Hi t Variable high temp - Customer set.
- Er 22* Fixed High Temp - 93°C.

*Er 23* Shorted external temperature sensor (rdt2).

*Er 24* Open external temperature sensor (rdt2).

*Er 25* Shorted internal temperature sensor (rdt1).

*Er 26* Open internal temperature sensor (rdt1).

*Er 51* ConF BOM - invalid BOM error.

*LofL* Low process flow warning/fault.

*LLF* Fluid low level fault - chiller shuts down.

Chiller reaction to warning/fault errors depend on how the it is configured, see Setup/Tuning Loop on page 18. The chiller is shipped configured to continue running. If any other code appears contact our customer service.

## Displaying Software Version Number

The controller can display the installed software version number. For example, for a chiller with software version **026950.1A**:

1. Chiller is running normally and displaying recirculating fluid temperature.
2. Press and hold  for at least 10 seconds. The display shows the first two digits, for example: **02**.
3. Press .  
The display shows the remaining digits to the left of the decimal, for example: **6950**.
4. Press .  
The display shows the decimal point and the digit to the right of the decimal point, for example, **.1**.
5. Press .  
The display shows the revision letter, as its equivalent number, for example, **1 = A**.
6. Press .  
Disregard this display.
7. Press .  
The display returns to the recirculating fluid temperature.

## Controller PID Values

The factory preset PID values are:

	P	I	D
COOL	12	0.5	0.0
HEAT	5	0.5	0.0

Thermo Fisher recommends that only a qualified technician change PID parameters. Incorrect PID values will hamper chiller performance. ▲

## Cleaning

### Reservoir

Periodically inspect the fluid inside the reservoir. If cleaning is necessary, flush the reservoir with a cleaning fluid compatible with the circulating system and the cooling fluid.

The cooling fluid should be replaced periodically. Replacement frequency depends on the operating environment and amount of usage.

Before changing the cooling fluid ensure it is at a safe handling temperature. ▲

**NOTE** It will be necessary to use a wet-vac or siphon to completely drain the reservoir. Use care to prevent damage the reservoir float switch. ▲

See page 12, Filling Requirements, for cooling fluid replacement.

### Condenser

For proper operation, the chiller needs to pull substantial amounts of air through a condenser. A build up of dust or debris on the fins of the condenser leads to a loss of cooling capacity. Optional air filters are available, contact our Service Department. See Preface, After-Sale Support.

The lower front of the chiller has a one-piece grille assembly. Using your hands gently pry the assembly off. Use care not to scratch the paint.

Periodic vacuuming of the condenser fins is necessary. The cleaning frequency depends on the operating environment. After initial installation we recommend a monthly visual inspection of the condenser. After several months, the cleaning frequency will be established.

Use care cleaning the condenser fins, they can be easily bent. ▲

## Power Cable Mating Ends

Every five years inspect the power cable mating ends for signs of corrosion and damage. Clean or replace as required.

## Leaks

A very small leak between the pump and motor is normal. If this is a new pump, allow it to run for 24 hours and reinspect. A continuous leak requires pump replacement.

## Phase Rotation Interlock

Three phase chillers have a rotation interlock which prevents it from starting if the phase rotation is wrong. If the chiller will not start, see Chapter 5, Troubleshooting. If the options in the checklist are not applicable, the problem may be phase rotation.

The interlock is located in the power box. Remove the chiller's top and left-side panel to access it. The interlock has a green light which illuminates when the phase is correct.

If the phase is incorrect, disconnect the chiller from its power source. Reverse any two line conductors on the line side of the main circuit breaker.

**Never remove the green ground wire. ▲**

Reconnect the chiller to its power source. If it will not start, contact our Service Department.

## PD Pressure Relief Valve

Chillers with PD pumps have a Pressure Relief Valve. Refer to the name plate label on the chiller and the BOM decoder on page 31 to identify the pump typer. The valve is used to adjust the chiller's fluid pressure.

**NOTE** The valve is factory preset for the most common applications and normally requires no further adjustment. The valve is factory preset not to exceed  $80 \pm 5$  psi ( $550 \pm 1$  kPa). ▲

**NOTE** PD pumps are capable of 110 psi maximum output pressure. Ensure your plumbing is rated to withstand 110 psi at the highest operating temperatures. ▲

System back pressure is installation and application dependent. Due to internal Merlin back pressure, the minimum pressure setting for a deadheaded pump is 28 psi for a PD-2 pump, and 10 psi for a PD-1 (these settings prohibit external flow from the Merlin).

Before adjusting the valve turn the chiller off.

Locate the circular relief valve opening on the rear of the chiller. Turn the threaded stem fully clockwise.

If the chiller is not plumbed to an application, install a loop of hose equipped with a shut-off valve between the supply and return fittings.

Turn the chiller on. Ensure that there is back pressure on the system. Use the pressure gauge to read the relief valve setting.

Back out the threaded stem on the relief valve counterclockwise. Continue until the gauge indicates 80 psi (550 kPa) or the desired setting.

**NOTE** The valve may drip if the threaded stem is backed out too far. ▲

**NOTE** If the chiller starts to vibrate the valve setting may be the cause. Changing the pressure setting  $\pm 5$  psi will eliminate the vibration. ▲

## Pump Strainer

Chillers with PD or T1 pumps, and chillers with CP-55 pumps that are equipped with reservoir level/low pump flow switches, have a strainer. Refer to the name plate label on the rear of the chiller and the BOM decoder on page 31 to identify the type of pump and configuration.

If debris is in the system, the strainer prevents the material from being drawn into the flow switch and the pump, damaging the pump vanes. A clogged strainer also causes increased discharge pump pressures.


After initial installation, the strainer may become clogged with debris and scale. Therefore, **the strainer must be cleaned after the first week of installation.** After this first cleaning, we recommend a monthly visual inspection. After several months, the cleaning frequency will be established.

**Before cleaning ensure the cooling fluid is at a safe handling temperature. Disconnect the power cord from the power source. ▲**

We also recommend draining the reservoir.

The strainer is in the pump return line, the exact location depends on the chiller's configuration. For most configurations it is located behind the rear panel. For chillers configured with low level and low temperature switches, it is behind the right-side panel.

**NOTE** A small amount of fluid will come out of the internal strainer housing when the fitting is removed. Have a rag handy to put under the housing before removing the fitting. ▲

Remove the panel from the chiller. Use one wrench to support the filter housing, and a second wrench to unscrew the fitting,  $1\frac{3}{8}$  inch for CP pumps, or  $\frac{3}{4}$  inch for PD and T1 pumps. Remove the screen.


**NOTE** Some chillers with  $\frac{3}{4}$  inch fittings also have a brass end cap, see illustration. ▲

Clean the screen by rinsing it with water. When the screen is clean, replace it and the cap. Tighten the fitting and replace the panel. Refer to Chapter 2 for instructions on replacing the cooling fluid.

**NOTE** An external strainer accessory kit is available to make routine strainer preventative maintenance easier. ▲

## BOM Decoder

The Bill of Material (BOM) number helps identify the configuration of your chiller. The number is on the label located on the rear of the chiller.

## Temperature Sensor Calibration

### Internal Sensor


If the temperature fluid temperature display disagrees with your reference thermometer, the internal temperature sensor (rtd1) may need calibration.

Do not pick calibration points that are outside the safe operating limits of the fluid in your application. For example with water, 35°C and 5°C would be typical calibration points.

Run the chiller to a high-end calibration point. Place a calibrated reference thermometer in the reservoir, ensure the fluid temperature is stable.


Enter Calibration by pressing and holding  and then pressing  three times. The display reads CAL. Press  and the display reads rtd1.

To calibrate the high-end temperature press  and the display reads r1H. Press  again and the display flashes between r1H and the current probe temperature.

Use  and  to adjust the r1H display. Once the display temperature matches the reference thermometer press  twice. Press  until Stor is displayed, press  to accept the new value.

Repeat for the low-end temperature r1L.

**NOTE** Both the high and low temperatures must be entered for a valid calibration. If power is lost before the procedure is complete, critical data needed for the calculation of calibration parameters is lost. After pressing  at the Stor prompt wait several seconds before proceeding to ensure that a bad calibration message does not appear. Premature use of the keypad after pressing  may cancel the bad calibration error message. 



### Optional External Sensor

When using any external sensor (rtd2) for the first time, Thermo Fisher recommends a two-point calibration. This calibration only affects the temperature read by the external sensor.

Do not pick calibration points that are outside the safe operating limits of the fluid in your application. For example with water, 35°C and 5°C would be typical calibration points.

Enter the Calibration Loop by pressing and holding  and then pressing the  three times. The display reads CAL. Press  and the display reads rtd1, then press  and the display reads rtd2.

To calibrate the high-end temperature press  and the display reads r2H. Press  again and the display flashes between r2H and the current probe temperature.

Place the external sensor and a calibrated reference thermometer in a source containing fluid that is just below the high-end temperature range of the chiller, ensure the fluid temperature is stable. Use ▲ and ▼ to adjust the temperature display to match the reference thermometer.

Once the display temperature matches the reference thermometer press  twice. Press  until Stor is displayed, pressing ▲ accepts the new value.

Repeat for the low-end temperature, r2L. Use a fluid that is just above the chiller's low-end temperature range.

**NOTE** Both the high and low temperatures must be entered for a valid calibration. If power is lost before the procedure is complete, critical data needed for the calculation of calibration parameters is lost. After pressing ▲ at the Stor prompt wait several seconds before proceeding to ensure that a bad calibration message does not appear. Premature use of the keypad after pressing ▲ may cancel the bad calibration error message. ▲

## Decommissioning/ Disposal

Decommissioning prepares equipment for safe and secure transportation.

Laboratory Grade Ethylene glycol (EG) is poisonous and flammable. Before disposing refer to the manufacturer's most current SDS for handling precautions. ▲

Decommissioning must be performed only by qualified dealer using certified equipment. All prevailing regulations must be followed. ▲

Consider decommissioning the chiller when:

- It fails to maintain desired specifications
- It no longer meets safety standards
- It is beyond repair for its age and worth

Refrigerant and compressor oil must be recovered from equipment before disposal.

Keep in mind any impact your application may have had on the chiller.

Direct questions about chiller decommissioning or disposal to our Sales, Service and Customer Support.

Handle and dispose in accordance with the manufacturers specification and/or the SDS for the material used. ▲

## Shipment Storage

Follow the manufacturer's SDS instructions if decontamination is required. ▲

Transporting and/or storing the chiller in near or below freezing temperatures requires draining, see **Draining** in this Section. Store the chiller in the temperature range of -25°C to 60°C (with packaging), and <80% relative humidity. ▲

If the chiller is stored for more than 90 days it must be flushed with clean fluid before operating. ▲

# Chapter 5 Troubleshooting

This chapter covers the basic troubleshooting procedures.

## Checklist

### 1. Chiller will not start or shuts down

- a. Check the line cord; ensure it is plugged in.
- b. Check the circuit protector on the rear of the chiller.
- c. Check the controller for error codes, see page 25.

**NOTE** On chillers with a Low Flow Switch, and that are configured to shut down with a low flow fault (**LoFL**), several starting attempts may be necessary. ▲

- d. On chillers with a Reservoir Level Switch, check the fluid level.
- e. On chillers with optional serial communications or remote start/stop, ensure the desired mode is enabled.
- f. Ensure power source voltage is within the chiller rated voltage  $\pm 10\%$ .
- g. 400 Volt, 50 Hertz chillers are equipped with a phase monitor. Ensure proper phase to the chiller. See Phase Rotation Interlock on page 28.
- h. If  is lit, the chiller was shut off using serial communications. Start the chiller using serial communications. To start the chiller using the local controller press and hold both arrow keys for 10 seconds, the controller will power. Turn serial communications off using the procedure on page 16, and then press  to turn off the controller. Press  again to start the chiller.

## 2. Chiller will not circulate fluid

- a. Check the reservoir level. Fill, if necessary.
- b. Check the position of the reservoir tank valves. See Reservoir Tank Isolation Valves/Open Tank Circulation, Chapter 2.
- c. Check for restrictions in the cooling lines to the application.
- d. Check the pump strainer (chillers with PD or T1 pumps).
- e. The pump motor may be too hot. Allow it to cool.
- f. Check the pressure gauge, adjust the relief valve as necessary (chillers with PD pumps).
- g. On chillers with PD pumps, check the pump and coupling. Ensure the power is off and remove the top and right-side panels. Locate the brass pump and grey motor. A silver clamp ring holds the two together. Remove the clamp and separate the pump and motor. First look to see if the coupling is broken or not, and then check the pump shaft by rotating. If it is seized or sticks then it has failed.
- h. On chillers with MD pumps never completely restrict flow to your application. Dead-heading the pump will damage its coupling and require pump replacement.

## 3. Inadequate temperature control

- a. Verify the setpoint.
- b. If the temperature continues to rise, make sure your application's heat load does not exceed the rated specification, see Chapter 2, Specifications.
- c. Make sure the air intake and discharge are not impeded and the ambient temperature does not exceed +35°C.
- d. Make sure the condenser is free of dust and debris, see Chapter 4, Maintenance and Cleaning.
- e. Verify the controller's PID values, see Setup/Tuning Loop.

**4. Chiller runs loudly**

- a. Check the position of the reservoir tank valves. See Reservoir Tank Isolation Valves/Open Tank Circulation, Chapter 2.

**5. Condensation on reservoir**

- a. Check the position of the upper reservoir tank valve. See Reservoir Tank Isolation Valves/Open Tank Circulation, Chapter 2.

**6. Chiller vibrates**

- a. The external pressure relief valve may need adjustment. Change the pressure setting  $\pm 5$  psi, see Chapter 4.

**7. No serial communications (optional)**

- a. All chillers are tested for serial communications before they leave the factory. Check the controller, ensure serial communication is enabled.
- b. Check all communications commands, they *must be exact*. See Appendix B.
- c. Check all wiring for proper connections or possible shorts.
- d. Software to verify serial communication is available from Thermo Fisher.

**Service Assistance**

If, after following these troubleshooting steps, your chiller fails to operate properly, contact our Service Department for assistance (see Preface, After-sale Support).

Before calling, *please* obtain the following information:

- chiller BOM number
- chiller serial number
- chiller software version (page 26)
- voltage of power source

# Appendix A Water Quality and Standards

## Contaminants

	Permissible (PPM)	Desirable (PPM)
Microbiologicals (algae, bacteria, fungi)	0	0
Inorganic Chemicals		
Calcium	<40	<0.6
Chloride	<250	<25
Copper	<1.3	<1.0
Iron	<0.3	<0.1
Lead	<0.015	0
Magnesium	<12	<0.1
Manganese	<0.05	<0.03
Nitrates\Nitrites	<10 as N	0
Potassium	<20	<0.3
Silicate	<25	<1.0
Sodium	<20	<0.3
Sulfate	<250	<50
Hardness	<17	<0.05
Total Dissolved Solids	<50	<10
Other Parameters		
pH	6.5-8.5	7-8
Resistivity	0.01*	0.05-0.1*


\* MΩcm (Compensated to 25°C)

Unfavorably high total ionized solids (TIS) can accelerate the rate of galvanic corrosion. These contaminants can function as electrolytes which increase the potential for galvanic cell corrosion and lead to localized corrosion such as pitting which can be observed at the studs and on the outside surface of cooling coils. Eventually, the pitting becomes so extensive that the coil will leak refrigerant into the water reservoir.

For example, raw water in the U.S. averages 171 ppm (of NaCl). The recommended level for use in a water system is between 0.5 to 5.0 ppm (of NaCl).

Recommendation: Initially fill the tank with distilled/deionized water. Do not use untreated tap water as the total ionized solids level may be too high.

The desired level for long time usage is 1 to 3 MΩcm (compensated to 25°C). This reduces the electrolytic potential of the water and prevent or reduce the galvanic corrosion observed.

## Compatibility with Approved Fluids

Never use corrosive fluids with this unit. Do not use automotive antifreeze. Commercial antifreeze contains silicates that can damage the pump seals. Use of automotive antifreeze will void the manufacturer's warranty. ▲

### Filtered/Single Distilled water (+5°C to +35°C)

This fluid is recommended primarily because it has all microorganisms that cause biological fouling removed through vaporizing and condensing the water. However, distilled water does not remain pure for very long when exposed to the atmosphere. Air-borne spores can contaminate the water and activate algae growth. An effective maintenance plan would include switching out the fluid with newly distilled water every six months. The particulates that have been filtered out in the process are also preventative in keeping the system "clean" of contaminants.

### 50/50 Uninhibited Ethylene Glycol/Water (-15°C to + 35°C)

Ethylene glycol is used to depress the freezing point of water as a coolant. We recommend not using the uninhibited (no corrosion additives) ethylene glycol. It is more corrosive to copper than plain water so it is not recommended unless required for the application.

### 50/50 Inhibited Ethylene Glycol/ Water (-15°C to +35°C)

Inhibited glycol can be used to increase the operating temperature range of the fluid but should not be used as a "pre-mixed anticorrosive" solution. Industry standards use a pH standard of 8 to determine when the fluid has become corrosive. Dowtherm<sup>®</sup> is an ethylene based product that contains dipotassium phosphates in a 4% concentration. The recommended use of Dowtherm<sup>®</sup> is mixing with distilled or deionized water or water that contains less than 25 ppm chloride and sulfate and less than 100 ppm total hardness of CaCO<sub>3</sub>.

The general term, inhibited glycol/water, is too close to meaning inhibited water. Inhibited water can have many types of additives including chromate that will quickly foul the cooling system. Some inhibitor additives can release the bonding agent in the carbon graphite in the PD2 pumps so they are incompatible, such as Sodium Hydroxide.

**50/50 Uninhibited Propylene Water (-15°C to +35°C)**

Although the use of propylene glycol is similar to ethylene glycol, propylene glycol is considered “safe” to use in the food industry. Propylene is less dense than ethylene and has a tendency to weep through mechanical seals.

**50/50 Inhibited Propylene/ Water (-15°C to +35°C)**

Same issues as with uninhibited propylene and inhibited ethylene glycol.

**Deionized water 1-3 MΩcm, compensated (+5°C to +35°C)**

Deionized water has had the conductive ions that cause galvanic corrosion between dissimilar metals removed.

## Wetted Materials

### **PD PUMPS (PD1 and PD2):**

Carbon Graphite (Graphitar Type B)

Brass

303 Stainless Steel

Ceramic (85% Aluminum Oxide and Balance is Considered Glass)

Fluorocarbon (Viton equivalent)

Ultem (relief Valve Material)

### **T1 PUMP:**

Stainless Steel AISI 304

Bronze ASTM B62

Bronze ASTM B16

Buna N

Buna/Ceramic

Buna/Carbon

### **CP55 PUMP:**

Stainless Steel

Viton

### **Evaporator**

AISI 316 Stainless Steel

Copper Brazing

### **Tank Body**

Polyethylene (LLPDE, Mobile NRP=135)

### **Hose Assemblies/Plumbing:**

Buna-N (Parker Push-Lok O-ring)

Nitrile based synthetic rubber (Parker Push-Lok hoses)

Brass

Copper

### **Ball Valve:**

Viton

Chrome plated brass body and ball

PTFE

### **Sealants:**

95/5 solder

Bcup 5 brazing

Anaerobic Methacrylate Ester

### **Low Level/Low Flow Switches (Optional):**

Polysulphone Low Level

Brass Low Flow

Continued on next page

**Tank Fittings:**

Brass

**Strainer**

Bronze

304 Stainless Steel Mesh

**Tank Cap:**

Nylon

# Appendix B RS-232 NC Serial Communications Protocol


**NOTE** This appendix assumes you have a basic understanding of communications protocols. ▲

All data is sent and received in binary form, do not use ASCII. In the following pages the binary data is represented in hexadecimal (hex) format.

The NC Serial Communications Protocol is based on a master-slave model. The master is a host computer, while the slave is the chiller's controller. Only the master can initiate a communications transaction (half-duplex). The slave ends the transaction by responding to the master's query. The protocol uses an RS-232 serial interface with the default parameters: 9600 baud, 8 data bits, 1 stop bit, no parity and no handshaking.

**NOTE** Before the chiller will communicate, RS-232 must be turned on using the controller. Ensure the  indicator is illuminated. ▲

The chiller can be controlled through your computer's serial port by using the chiller's standard 9-pin RS-232 connection. Data read of the serial port connects to the data transmit (pin 2) of the chiller. Data transmit of the serial port connects to data read (pin 3) of the chiller.

Communication cables are available from Thermo Fisher. Contact our sales department for additional information.

All commands must be entered in the exact format shown in the tables on the following pages. The tables show all commands available, their format

and responses. Controller responses are either the requested data or an error message. The controller response *must* be received before the host sends the next command.

The host sends a command embedded in a single communications packet, then waits for the controller's response. If the command is not understood or the checksums do not agree, the controller responds with an error command. Otherwise, the controller responds with the requested data. If the controller fails to respond within 1 second, resend the command.

**NOTE** All byte values are shown in hex, hex represents the binary values that must be sent to the chiller. **Do not use ASCII.** ▲

The framing of the communications packet in both directions is:

Checksum region								
Lead char CA	Addr-MSB 0	Addr-LSB 1	Command	n d-bytes	d-byte 1	...	d-byte n	Checksum
<i>Lead char</i>			CA (hex).					
<i>Addr-msb</i>			Device address is 1					
<i>Addr-lsb</i>			Most significant byte of device address is 0.					
<i>Command</i>			Least significant byte of device address is 1.					
<i>n d-bytes</i>			Command byte (see Table 1).					
<i>d-byte 1</i>			Number of data bytes to follow (00 to 03 hex).					
...			1 <sup>st</sup> data byte (the qualifier byte is considered a data byte).					
<i>d-byte n</i>			...					
<i>Checksum</i>			n <sup>th</sup> data byte.					
			Bitwise inversion of the 1 byte sum of bytes beginning with the most significant address byte and ending with the byte preceding the checksum. (To perform a bitwise inversion, "exclusive OR" the one byte sum with FF hex.)					


The master requests information by sending one of the Read Functions as shown in Table 1. Since no data is sent to the chiller during a read request, the master uses 00 for the number of data bytes following the command byte.

The chiller will respond to a Read Function by echoing the lead character, address, and command byte, followed by the requested data and checksum. When the chiller sends data, a qualifier byte is sent first, followed by a two byte signed integer (16 bit, MSB sent first). The qualifier byte indicates the precision and units of measure for the requested data as detailed in Table 2.

As an example, the master requests to read internal temperature by sending:


If the temperature is -12°C, the chiller would reply:


The master sets parameters in the chiller by sending one of the Set Functions as shown in Table 1. The master does not send a qualifier byte in the data field. The master should be preprogrammed to send the correct precision and units (it could also read the parameter of interest first to decode the correct precision and units needed).

For example, if the master wants to set the setpoint to 30°C, it would send :

The chiller responds:

**Table 1** Commands (All bytes are in hex)

FUNCTION	MASTER SENDS	CHILLER RESPONDS
Read Acknowledge	CA 00 01 <b>00</b> 00 FE	CA 00 01 <b>00</b> 02(v1)(v2)(cs)
Read Status (see Table 3)	CA 00 01 <b>09</b> 00 F5	CA 00 01 <b>09</b> 02(d1)(d2)(cs)
Error		CA 00 01 <b>0F</b> 02(en)(ed)(cs)
Read Internal Temperature (RTD1)	CA 00 01 <b>20</b> 00 DE	CA 00 01 <b>20</b> 03(qb)(d1)(d2)(cs)
Read External Temperature(RTD2)	CA 00 01 <b>21</b> 00 DD	CA 00 01 <b>21</b> 03(qb)(d1)(d2)(cs)
Read Setpoint (control point)	CA 00 01 <b>70</b> 00 8E	CA 00 01 <b>70</b> 03(qb)(d1)(d2)(cs)
Read Low Temperature Limit	CA 00 01 <b>40</b> 00 BE	CA 00 01 <b>40</b> 03(qb)(d1)(d2)(cs)
Read High Temperature Limit	CA 00 01 <b>60</b> 00 9E	CA 00 01 <b>60</b> 03(qb)(d1)(d2)(cs)
Read Cool Proportional Band (P)	CA 00 01 <b>74</b> 00 8A	CA 00 01 <b>74</b> 03(qb)(d1)(d2)(cs)
Read Cool Integral (I)	CA 00 01 <b>75</b> 00 89	CA 00 01 <b>75</b> 03(qb)(d1)(d2)(cs)
Read Cool Derivative (D)	CA 00 01 <b>76</b> 00 88	CA 00 01 <b>76</b> 03(qb)(d1)(d2)(cs)
Read Heat Proportional Band (P)	CA 00 01 <b>71</b> 00 8D	CA 00 01 <b>71</b> 03(qb)(d1)(d2)(cs)
Read Heat Integral (I)	CA 00 01 <b>72</b> 00 8C	CA 00 01 <b>72</b> 03(qb)(d1)(d2)(cs)
Read Heat Derivative (D)	CA 00 01 <b>73</b> 00 8B	CA 00 01 <b>73</b> 03(qb)(d1)(d2)(cs)

command bytes shown in **bold**  
 qb = qualifier byte, see Table 2  
 d1,d2 = 16 bit signed integer of the value being sent or received  
 cs = the checksum of the string (see text)

v1. v2 = protocol version  
 ed Bad command byte gets echoed)  
 en (error number) 01 = Bad command,  
 02 = Bad checksum  
 xx = no valid data, include in checksum

**Table 1** Commands (continued) (All bytes are in hex)

FUNCTION	MASTER SENDS	CHILLER RESPONDS
Set Setpoint (control point)*	CA 00 01 <b>F0</b> 02(d1)(d2)(cs)	CA 00 01 <b>F0</b> 03(qb)(d1)(d2)(cs)
Set Low Temp	CA 00 01 <b>C0</b> 02(d1)(d2)(cs)	CA 00 01 <b>C0</b> 03(qb)(d1)(d2)(cs)
Set High Temp	CA 00 01 <b>E0</b> 02(d1)(d2)(cs)	CA 00 01 <b>E0</b> 03(qb)(d1)(d2)(cs)
Set Cool Proportional Band (P = 0.1 - 99.9)	CA 00 01 <b>F4</b> 02(d1)(d2)(cs)	CA 00 01 <b>F4</b> 03(qb)(d1)(d2)(cs)
Set Cool Integral (I = 0 - 9.99)	CA 00 01 <b>F5</b> 02(d1)(d2)(cs)	CA 00 01 <b>F5</b> 03(qb)(d1)(d2)(cs)
Set Cool Derivative (D = 0 - 5.0)	CA 00 01 <b>F6</b> 02(d1)(d2)(cs)	CA 00 01 <b>F6</b> 03(qb)(d1)(d2)(cs)
Set Heat Proportional Band (P = 0.1 - 99.9)	CA 00 01 <b>F1</b> 02(d1)(d2)(cs)	CA 00 01 <b>F1</b> 03(qb)(d1)(d2)(cs)
Set Heat Integral (I = 0 - 9.99)	CA 00 01 <b>F2</b> 02(d1)(d2)(cs)	CA 00 01 <b>F2</b> 03(qb)(d1)(d2)(cs)
Set Heat Derivative (D = 0 - 5.0)	CA 00 01 <b>F3</b> 02(d1)(d2)(cs)	CA 00 01 <b>F3</b> 03(qb)(d1)(d2)(cs)
Set On/Off Array		

The Set On/Off Array command is used to set and/or request chiller settings. The array consists of 5 data bytes, d1 - d5. Each byte represents the state of one setting, see \*\* below. The master can send the array with a value of 0 (turn off/disable), 1 (turn on/enable), or 2 (do not change). 0 or 1 causes the chiller to change the setting, 2 only requests the state of a setting. The slave returns values of 0 (off/disabled) or 1 (on/enabled).

Sample - turn chiller on (02950.1J-K) CA 00 01 **81** 05 01 02 02 02 02 6F

\* = limited to the range of the chiller

command bytes shown in **bold**

qb = qualifier byte, see Table 2

d1,d2 = 16 bit signed integer of the value being sent or received

cs = the checksum of the string (see text)

\*\*d1 - chiller off = 0, on =1

d2 - external sensor disabled = 0, enabled = 1

d3 - fault mode disabled = 0, enabled = 1

d4 - tenths display disabled = 0, enabled = 1

d5 - Alarms on internal = 0, or external sensor = 1, (1 is invalid when d2 = 0, external sensor disabled)

**Table 2****QUALIFIER BYTE**

00	0 precision, no units of measure
01	0 precision, °C
10	1 precision, no units of measure
11	1 precision, °C
20	2 precision, no units of measure

Example: The temperature of 45.6°C would be represented by the qualifier 11 hex, followed by the 2 bytes 01 C8 hex (456 decimal).

**Table 3** READ STATUS

<b>BIT</b>	<b>d1</b>	<b>d2<sup>t</sup></b>
b.7 = 1	External Temp Sensor Enabled	Reserved
b.6 = 1	Level 1 Warning/Auto Refill	Reserved
b.5 = 1	Low Flow Warning	Internal Temp Sensor Fault
b.4 = 1	Low Level 2 Warning	External Temp Sensor Fault
b.3 = 1	High or Low Temp Warning	High Temp/Fixed High Temp Fault
b.2 = 1	High or Low Temp Bypass	Low Temp/Fixed Low Temp Fault
b.1 = 1	Chiller Faulted*	Low Flow Fault
b.0 = 1	Chiller Running*	Low Level 2 Fault

<sup>t</sup> d2 is the detail of fault d1-b.1

\* mutually exclusive

## DECLARATION OF CONFORMITY

Manufacturer: Thermo Fisher Scientific  
Address: 25 Nimble Hill Road  
Newington, NH 03801

We declare that the equipment named below has been designed to comply with the relevant sections of the below referenced specifications and is in accordance with the requirements of the indicated directives.

Product: Merlin Recirculating Chiller

Models: Units have a 12 digit part number consisting of UUU T VV PP XXXX defined as follows:

UUU = Unit type can be:

262 = M25; 263 = M33; 264 = M75; 265 = M100; 266 = M150

T = Temperature Range can be 1-3 inclusive, where:

1 = Std temp; 2 = Low Temp; 3 = High Temp

V = Unit voltage rating:

12 = 115 / 60 100 / 50; 16 = 200 / 230 / 60 200 / 50; 25 = 230 / 50; 26 = 400 / 50 5-Wire

P = Pump type

02 = CP-55; 03 = PD-1; 04 = PD-2; 15 = MD-30; T1 = 99

X = Any digit 0-9 inclusive, used as sequential numbering only.

Directives and Standards:

2014/30/EC ± Electromagnetic Compatibility Directive (EMC)

- EN 61326-1: 2013 Electrical equipment for measurement, control, and laboratory use - EMC requirements. General requirements

2014/35/EC - Low Voltage Directive (LVD):

- EN 61010-1: 2010 Safety requirements for electrical equipment for measurement, control, and laboratory use -- Part 1: General requirements.
- EN 61010-2-010: 2003 Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials


2011/65/EU – Restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS)

- EN 50581: 2012 - Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

Authorised representative in the EC:

Name: Thermo Fisher Scientific  
Address: Dieselstrasse 4  
76227 Karlsruhe Germany

On behalf of the manufacturer:

  
John Parker  
Electrical Engineering Manager  
Temperature Control  
Thermo Fisher Scientific  
Newington, NH, USA

Date:

15 August 2016

## WARRANTY

Thermo Fisher Scientific warrants for 24 months (**excluding Positive Displacement pumps which are warranted for 12 months**) from date of shipment the Thermo Scientific Merlin chiller according to the following terms.

Any part of the chiller manufactured or supplied by Thermo Fisher Scientific and found in the reasonable judgment of Thermo Fisher to be defective in material or workmanship will be repaired at an authorized Thermo Fisher Repair Depot without charge for parts or labor. The chiller, including any defective part must be returned to an authorized Thermo Fisher Repair Depot within the warranty period. The expense of returning the chiller to the authorized Thermo Fisher Repair Depot for warranty service will be paid for by the buyer. Our responsibility in respect to warranty claims is limited to performing the required repairs or replacements, and no claim of breach of warranty shall be cause for cancellation or rescission of the contract of sales of any chiller. With respect to chillers that qualify for field service repairs, Thermo Fisher Scientific's responsibility is limited to the component parts necessary for the repair and the labor that is required on site to perform the repair. Any travel labor or mileage charges are the financial responsibility of the buyer.

The buyer shall be responsible for any evaluation or warranty service call (including labor charges) if no defects are found with the Thermo Scientific product.

This warranty does not cover any chiller that has been subject to misuse, neglect, or accident. This warranty does not apply to any damage to the chiller that is the result of improper installation or maintenance, or to any chiller that has been operated or maintained in any way contrary to the operating or maintenance instructions specified in this Instruction and Operation Manual. This warranty does not cover any chiller that has been altered or modified so as to change its intended use.

In addition, this warranty does not extend to repairs made by the use of parts, accessories, or fluids which are either incompatible with the chiller or adversely affect its operation, performance, or durability.

Thermo Fisher Scientific reserves the right to change or improve the design of any chiller without assuming any obligation to modify any chiller previously manufactured.

THE FOREGOING EXPRESS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO WARRANTIES OR MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

OUR OBLIGATION UNDER THIS WARRANTY IS STRICTLY AND EXCLUSIVELY LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE COMPONENT PARTS AND Thermo Fisher Scientific DOES NOT ASSUME OR AUTHORIZE ANYONE TO ASSUME FOR IT ANY OTHER OBLIGATION.

Thermo Fisher Scientific ASSUMES NO RESPONSIBILITY FOR INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES INCLUDING, BUT NOT LIMITED TO LOSS OR DAMAGE TO PROPERTY, LOSS OF PROFITS OR REVENUE, LOSS OF THE CHILLER, LOSS OF TIME, OR INCONVENIENCE.

This warranty applies to chillers sold by Thermo Fisher Scientific. (Refer to the warranty for chillers sold by the affiliated marketing company of Thermo Fisher Scientific for any additional terms.) This warranty and all matters arising pursuant to it shall be governed by the law of the State of New Hampshire, United States. All legal actions brought in relation hereto shall be filed in the appropriate state or federal courts in New Hampshire, unless waived by Thermo Fisher Scientific.

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