# **Thermo Scientific Merlin M75 Recirculating Chiller**

Thermo Scientific Manual P/N U01116  
Rev. 09/21/2017

**Multilingual  
Quick Starts**

**Multilingual  
Essential Safety  
Instructions**

**Installation**

**Operation**

**Maintenance**

**Troubleshooting**

Visit our Web site at:

**<http://www.thermofischer.com/tc>**  
Product Service Information, Applications  
Notes, SDS Forms, e-mail.

Voice Info: (800) 258-0830

# **Thermo Scientific Merlin M75 Recirculating Chiller**

Thermo Scientific Manual P/N U01116  
Rev. 09/21/2017

**Multilingual  
Quick Starts**

**Multilingual  
Essential Safety  
Instructions**

**Installation**

**Operation**

**Maintenance**

**Troubleshooting**

Visit our Web site at:

**<http://www.thermofischer.com/tc>**  
Product Service Information, Applications  
Notes, SDS Forms, e-mail.

Voice Info: (800) 258-0830

## **Thermo Fisher Scientific**

25 Nimble Hill Road  
Newington, NH 03801  
Tel : (800) 258-0830 or  
(603) 436-9444  
Fax : (603) 436-8411  
[www.thermoscientific.com/tc](http://www.thermoscientific.com/tc)

## **Sales, Service, and Customer Support**

25 Nimble Hill Road  
Newington, NH 03801  
Tel: (800) 258-0830  
Sales: 8:00 am to 5:00 pm  
Service and Support: 8:00 am to 6:00 pm  
Monday through Friday (Eastern Time)  
Fax: (603) 436-8411  
[service.tc.us@thermofisher.com](mailto:service.tc.us@thermofisher.com)

Dieselstrasse 4  
D-76227 Karlsruhe, Germany  
Tel : +49 (0) 721 4094 444  
Fax : +49 (0) 721 4094 300  
[info.tc.de@thermofisher.com](mailto:info.tc.de@thermofisher.com)

Building 6, No. 27  
Xin Jinqiao Rd., Shanghai 201206  
Tel : +86(21) 68654588  
Fax : +86(21) 64457830  
[info.china@thermofisher.com](mailto:info.china@thermofisher.com)

### **Statement of Copyright**

Copyright © 2017 Thermo Fisher Scientific. All rights reserved.

This manual is copyrighted by Thermo Fisher Scientific.

Users are forbidden to reproduce, republish, redistribute, or resell any materials from this manual in either machine-readable form or any other form.

# Contents

<b>Preface</b>	<b>iii</b>
Compliance .....	iii
WEEE .....	iii
After-sale support .....	iv
Unpacking .....	iv
Warranty .....	iv
Safety Warnings .....	v
Quick Reference Operating Procedures .....	viii
<b>Chapter 1 Description</b> .....	<b>1-1</b>
Specifications .....	1-1
Wetted Parts .....	1-2
<b>Chapter 2 Installation and Operation</b> .....	<b>2-1</b>
Site .....	2-1
Electrical Requirements .....	2-2
Plumbing Requirements .....	2-2
Filling Requirements .....	2-4
Reservoir Tank Isolation Valves .....	2-4
Controller .....	2-5
Start Up/Shut Down .....	2-7
Setup/Tuning Loop .....	2-10
<b>Chapter 3 Options</b> .....	<b>3-1</b>
Communication Connector .....	3-1
Remote Start/Stop .....	3-2
Reservoir Level Warning/Fault .....	3-4
Low Pump Flow Warning/Fault .....	3-4
<b>Chapter 4 Maintenance</b> .....	<b>4-1</b>
Error Codes .....	4-1
Displaying Software Version Number .....	4-2
Controller PID Values .....	4-3
Cleaning .....	4-3
Power Cable Mating Ends .....	4-4
Leaks .....	4-4
PD Pressure Relief Valve .....	4-4
Pump Strainer .....	4-5
Temperature Sensor Calibration .....	4-6
Decommissioning/Disposal .....	4-8
Shipment/Storage .....	4-8

<b>Chapter 5</b>	<b>Troubleshooting.....</b>	<b>5-1</b>
	Checklist.....	5-1
	Service Assistance.....	5-3
<b>Appendix</b>	<b>NC Serial Communications Protocol .....</b>	<b>A-1</b>
<b>Declaration of Conformity</b>		
<b>WARRANTY</b>		

# Preface

**Compliance** The Declaration of Conformity is located in the back of this manual.

**WEEE** This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2012/19/EU. It is marked with 'wheelie bin' symbol:

Thermo Fisher Scientific has contracted with one or more recycling/disposal companies in each EU Member State, dispose of or recycle this product through them. Further information on Thermo Fisher Scientific's compliance with these Directives is available at:

[www.thermoscientific.com/WEEERoHS](http://www.thermoscientific.com/WEEERoHS)


## After-sale Support

Thermo Fisher Scientific is committed to service both during and after the sale. If you have questions concerning operation, contact our Sales Department. If your chiller fails to operate properly, or if you have questions concerning spare parts or Service Contracts, contact our Service Department.

Before calling, please obtain the following information:

- BOM number
- serial number
- software version
- voltage of power source

The chiller's BOM and serial number are located on the name plate label on the rear of the chiller.

## Unpacking

Retain all cartons and packing material until the chiller is operated and found to be in good condition. If the chiller shows external or internal damage contact the transportation company and file a damage claim. Under ICC regulations, this is your responsibility.

## Warranty

Thermo Scientific Merlin chillers have a warranty against defective parts and workmanship for 12 months from date of shipment.

## Safety Warnings

Make sure you read and understand all instructions and safety precautions listed in this manual before installing or operating your chiller. If you have any questions concerning the operation of your chiller or the information in this manual, please contact us. See inside cover for contact information.

**DANGER** indicates an imminently hazardous situation which, if not avoided, *will* result in death or serious injury.

**WARNING** indicates a potentially hazardous situation which, if not avoided, *could* result in death or serious injury.

**CAUTION** indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It is also be used to alert against unsafe practices.

The lightning flash with arrow symbol, within an equilateral triangle, is intended to alert the user to the presence of non-insulated "dangerous voltage" within the chiller's enclosure. The voltage magnitude is significant enough to constitute a risk of electrical shock.

This label indicates read the manual.

**Note** The chiller's equipment design incorporates a complete sheet metal enclosure for personnel protection from mechanical and electrical hazards. ▲

**Observe all warning labels.** ▲

**Never remove warning labels.** ▲

**Never place the chiller in a location where excessive heat, moisture, or corrosive materials are present.** ▲

**The chiller construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided.** ▲

**Never connect the process fluid inlet or outlet fittings to your building water supply or any water pressure source.** ▲

**Never use flammable or corrosive fluids with this chiller. Highly distilled and deionized water may be aggressive and cause material corrosion.** ▲

**Please contact Thermo Fisher before subjecting this chiller to prolonged exposure to highly distilled or deionized water. ▲**

**Do not use automotive antifreeze. Commercial antifreeze contains silicates that can damage the pump seals. Use of automotive antifreeze will void the manufacturer's warranty. ▲**

In addition to the specific warnings listed on the previous page the following general warnings also apply to your chiller:

**Performance of installation, operation, or maintenance procedures other than those described in this manual may result in a hazardous situation and may void the manufacturer's warranty. ▲**

**Transport the chiller with care. Sudden jolts or drops can damage the chiller's components. ▲**

**Drain the chiller before it is transported and/or stored. ▲**

**The circuit protector located on the rear of the chiller is not intended to be used as a disconnecting means. ▲**

**Observe all warning labels. ▲**

**Never remove warning labels. ▲**

**Never operate damaged or leaking equipment. ▲**

**Never operate the chiller without cooling fluid in the reservoir. ▲**

**Always turn off the chiller and disconnect the power cord from the power source before performing any service or maintenance procedures, or before moving the chiller. ▲**

**Never operate equipment with damaged power cords. ▲**

**Refer service and repairs to a qualified technician. ▲**

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## **Fluorinated Greenhouse Gases**

The following information is included to comply with REGULATION (EU) No 517/2014 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 April 2014 on fluorinated greenhouse gases:

**This product contains fluorinated greenhouse gases in a hermetically sealed system.**

**If a leak in the sealed system is detected, the operator shall repair without undue delay.**

Refer to the F-Gas Declaration of Conformity for additional information.

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**Declaration of conformity with Article 14 of Regulation (EU) No 517/2014 of the European Parliament and of the Council**

We, Thermo Fisher Scientific (Asheville) LLC, F-Gas Portal Registration Number 23643, declare under our sole responsibility that when placing on the market pre-charged equipment, which we import to or manufacture in the Union, the hydrofluorocarbons contained in that equipment are accounted for within the quota system referred to in Chapter IV of Regulation (EU) No 517/2014 as:

We hold authorisation(s) issued in accordance with Article 18(2) of Regulation (EU) No 517/2014 and registered in the registry referred to in Article 17 of that Regulation, at the time of release for free circulation to use the quota of a producer or importer of hydrofluorocarbons subject to Article 15 of Regulation (EU) No 517/2014 that cover(s) the quantity of hydrofluorocarbons contained in the equipment.

The hydrofluorocarbons contained in the equipment have been placed on the market in the Union, subsequently exported and charged into the equipment outside the Union, and the undertaking that placed the hydrofluorocarbons on the market made a declaration stating that the quantity of hydrofluorocarbons has been or will be reported as placed on the market in the Union and that it has not been and will not be reported as direct supply for export in the meaning of Article 15(2)(c) of Regulation (EU) No 517/2014 pursuant to Article 19 of Regulation (EU) No 517/2014 and Section 5C of the Annex to Commission Implementing Regulation (EU) No 1191/2014.

January 1, 2017

\_\_\_\_\_  
Date

Mark Pearson, Director Global Regulatory Affairs

\_\_\_\_\_  
[name and position of legal representative]

\_\_\_\_\_  
[signature of legal representative]



Manufacturer:

Thermo Fisher Scientific (Asheville) LLC  
275 Aiken Road  
Asheville, NC 28804  
U.S.A.

EU Only Representative:

Thermo Electron LED GmbH  
Robert-Bosch-Strasse 1  
D-63505 Langenselbold  
Germany  
**VAT ID Number: DE 812 403 137**

## Quick Reference Operating Procedures

### Installation

The chiller has an air-cooled refrigeration system. Air is drawn in the front of the chiller and discharged through rear and sides. Position the chiller so the intake and discharge are not impeded. Inadequate ventilation will cause a reduction in cooling capacity and, in extreme cases, compressor failure.

Excessively dusty areas should be avoided and a periodic cleaning schedule should be instituted. For proper operation, the chiller needs to pull substantial amounts of air through a condenser. A build up of dust or debris on the fins of the condenser will lead to a loss of cooling capacity.

The chiller will retain its full rated capacity in ambient temperatures up to approximately +77°F (+25°C).

Make sure the voltage of the power source meets the specified voltage,  $\pm 10\%$ .

The plumbing connections are located on the rear of the chiller and are labelled  and . These connections are 1/2 inch FPT. Remove the plastic protective plugs from both plumbing connections. Connect the  fitting to the inlet of your application. Connect the  fitting to the outlet of your application.

To fill the reservoir, remove the reservoir access panel and locate the reservoir cap. Remove the cap and fill the reservoir with clean cooling fluid.

The chiller is designed to use HFE-7100, HFE-7200, FC-3283 and FC-77.

### Operation

Before starting the chiller, double-check all electrical and plumbing connections. Make sure the circulating system is properly filled with cooling fluid. Place the circuit protector located on the rear of the chiller to the up position.

To start the chiller, press . To turn the chiller off press  again.

The  LED indicates the status of the refrigeration system. It illuminates to indicate the refrigeration system is removing heat from the cooling fluid. As the operating temperature approaches the setpoint, the LED will flash.

### Temperature Adjustment

To display the temperature setpoint, press  on the controller. The  indicator will illuminate and the display will flash the current setpoint value. To adjust the temperature setpoint, press the arrow keys until the desired temperature setpoint is indicated. Press  again to confirm the change. The display will rapidly flash the new value for a short time and then return to the recirculating fluid temperature.

### Periodic Maintenance

Periodically inspect the reservoir fluid. If cleaning is necessary, flush the reservoir with a cleaning fluid compatible with the circulating system and the cooling fluid.

The cooling fluid should be replaced periodically. Frequency depends on the operating environment and amount of usage.

Before changing the cooling fluid ensure it is at a safe temperature.

Periodic vacuuming of the condenser fins is necessary. The frequency of cleaning depends on the operating environment. We recommend a monthly visual inspection of the condenser after initial installation. After several months, the cleaning frequency will be established.

Chillers have a pump strainer. If debris is in the system, the strainer will prevent the material from being drawn into the pump and damaging the pump vanes. A clogged strainer will also cause increased pump discharge pressure.



After initial installation, the strainer may become clogged. **Clean the strainer after the first week of installation.** After this first cleaning, we recommend a monthly visual inspection. After several months, the cleaning frequency will be established. Before cleaning, disconnect the power cord from the power source and drain the chiller.

# Chapter 1 Description


## Description

The Thermo Scientific Merlin Recirculating Chiller is designed to provide a continuous supply of cooling fluid at a constant temperature and flow. The chiller is hard-plumbed and has an air-cooled refrigeration system, plate heat exchanger, recirculating pump, polyethylene reservoir, and a microprocessor temperature controller.

## Specifications

<b>Temperature Range<sup>1</sup></b>	-15°C to +35°C
<b>Temperature Stability<sup>1,2</sup></b>	±0.15°C
<b>Reservoir Volume</b> <i>Gallons/Liters</i>	0.5/1.9
<b>Refrigerant</b>	R134a
<b>Cooling Capacity (Watts)<sup>1</sup></b>	
<b>Pump Capacity</b>	
<b>Shipping Weight</b> <i>(Pounds/Kilograms)</i>	~200/~91

1. Specifications subject to change.
2. At the chiller. Display resolution 1.0°C. 0.1°C display resolution is available using the controller's Setup/Tuning Loop, see Chapter 2.

Chillers are shipped in a crate measuring 39 x 22 1/4 x 29 1/4.

### Wetted Parts

Material	Locations
Copper	Plumbing lines, plate exchanger, fittings
Brass	Pump housing, fittings, relief valve
Carbon/Graphite	Pump vanes, pump seal
Ceramic	Pump seal
Stainless steel (303,304,316)	Plate exchanger, pump parts
Soft Solder material (tin)	All solder joints
Phosphor Bronze	Pressure gage bourdon tube
Nitrile rubber	Pump o-rings, seals
Polyethylene (HDPE)	Tank
Bronze	Relief valve
Nickel (plating)	Tank isolation valve

## Chapter 2 Installation and Operation

<b>Site</b>	Ambient Temperature Range*	10°C to 35°C (50°F to 94°F)
	Relative Humidity Range	10% to 80% (non-condensing)
	Operating Altitude	Sea Level to 6000 feet (1830 meters)
	Overvoltage Category	II
	Pollution Degree	2

---


\*The chiller retains its full rated capacity in ambient temperatures up to approximately 25°C (77°F). Reduce the cooling capacity 1% for every 0.5°C (1°F) above 25°C (77°F), up to a maximum ambient temperature of 35°C (94°F).

**Never place the chiller in a location where excessive heat, moisture, or corrosive materials are present. ▲**

The chiller has an air-cooled refrigeration system. Air is drawn through the front and discharged through rear and side panels. The chiller must be positioned so the intake and discharge are not impeded. A minimum clearance of 3 feet (1 meter) on all vented sides is necessary for adequate ventilation. Inadequate ventilation will cause a reduction in cooling capacity and, in extreme cases, compressor failure.

Direction of Airflow

Excessively dusty areas should be avoided and a periodic cleaning schedule should be instituted (see Chapter 4, Maintenance). Optional air filters are available, contact our Service Department. See Preface, After-Sale Support.

## Electrical Requirements

The chiller construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided. ▲

Amperage based on extreme operating conditions. Operating in a 20°C ambient and at nominal voltage will result in a lower amperage draw.

Refer to the name plate label on the rear of the chiller for specific electrical requirements


Chiller	Pump	Voltage	Current Value	Line Cord Plug
M 75	PD 2	208/230	10.1	6-15P

## Plumbing Requirements

The plumbing connections are located on the rear of the chiller and are labeled  (supply) and  (return). The connections are 1/2 inch Female Pipe Thread. Chillers are supplied with 3/8 inch and 1/2 inch barbed adapters.

Remove the plastic protective plugs from both plumbing connections. If requires, install the barbed adapters to these connections.

Do not install the adapters or your fittings dry. Use two wraps of Teflon® tape on the male thread wound in a clockwise direction while facing the thread, or a suitable Teflon® sealant/lubricant. The 1/2" male connection can be tightened 2 full turns  $\pm 1/4$  turn to seal. Overtightening will damage the female thread.

Connect  to the hose feeding the inlet of your application. Connect  to the hose from the outlet of your application. Clamp all connections.

**Never connect the fittings to your building water supply or any water pressure source. ▲**

**NOTE** Ensure your plumbing is rated to withstand 110 psi at the highest operating temperature. ▲

It is important to keep the distance between the chiller and the instrument being cooled as short as possible. Tubing should be straight and without bends. If diameter reductions must be made, they should be made at the inlet and outlet of your application, not at the chiller.

All chillers have a ball-valve reservoir drain located on the rear. Add the desired length of tubing to the drain valve. To help prevent any spillage from your application, consider installing external isolation valves on the supply and return lines.

**Never use flammable or corrosive fluids with this chiller. Do not use automotive antifreeze. Commercial antifreeze contains silicates that can damage the pump seals. Use of automotive antifreeze will void the manufacturer's warranty. ▲**

The chiller is designed to use:

HFE-7100

HFE-7200


FC-3283

FC-77.

## Filling Requirements

The polyethylene reservoir is translucent for easy fluid level monitoring. Remove the reservoir access panel. Locate and remove the reservoir cap. Fill the reservoir with clean cooling fluid.

Since the reservoir capacity may be small compared to your application, have extra cooling fluid on hand to keep the system topped off when external circulation is started.

## Reservoir Tank Isolation Valves

Chillers are equipped with two reservoir tank isolation valves. The valves are used to isolate the tank from the fluid flow. The chillers are factory shipped with the valves open.


Ensure the valves are open when circulating to a closed system. When operating at low temperatures and circulating to a closed system, should condensation appear on the reservoir, closing the upper valve *after initial start up* will help eliminate it. Leave the valve open for initial start up so process fluid flow through the reservoir can rapidly bleed air from the system. To close the upper valve, open the reservoir access panel and turn the valve 90°.

**System temperature changes with the isolation valves closed, and while circulating to a closed system, can create a vacuum and/or high pressure conditions. ▲**

## Controller

The controller controls temperature using a Proportional-Integral-Derivative (PID) algorithm. It is designed with self-diagnostic features and easy to use operator interface.

Use this key to toggle the **chiller on or off**.

Use this key to **scroll** through the controller's LEDs. It is also used to accept and save changes.

If the alarm sounds, use this key to toggle the **alarm off and on**. **NOTE:** If the cause of the alarm is cleared but then reoccurs, the alarm will sound again.

Use this key to **increase displayed numerical values**.

Use this key to **decrease displayed numerical values**.

Indicates **refrigeration system status**. It illuminates to indicate the refrigeration system is removing heat from the cooling fluid. As the operating temperature approaches the temperature setpoint, the LED will flash. The indicator is off when heat is not being removed.

**NOTE** The following indications are not visible until  is depressed. ▲

 Indicates the controller is displaying the **setpoint**. Adjust the display value using the arrow keys.

**NOTE** You cannot adjust the setpoint closer than 2°C to either of the temperature limits discussed below. ▲


 Indicates the controller is displaying the **low temperature alarm** setting. The indicator flashes and the alarm sounds if this limit is exceeded. Adjust the display value using the arrow keys.

 Indicates the controller is displaying the **high temperature alarm** setting. The indicator flashes and the alarm sounds if this limit is exceeded. Adjust the display value using the arrow keys.

 **(Optional)** Indicates the unit is in the **remote/serial communication** mode of operation. The indicator flashes if the unit is in the remote start/stop mode of operation. Toggle between the local and remote modes using the arrow keys.

**NOTE** For all four indicators, if neither arrow key is pressed within 10 seconds the display will return to the recirculating temperature display. Save any change using . After pressing an arrow key, if  is not pressed within 60 seconds the display will return to the recirculating temperature display and ignore any change. ▲

## Start Up/Shut Down

### Pre-Start

Before starting the chiller, double check all electrical and plumbing connections. Ensure the power cord is secured to the rear of the chiller. Have extra recirculating fluid on hand. If the chiller will not start refer to Chapter 5 Troubleshooting.

### Circuit Protector

Place the circuit protector located on the rear of the chiller to the up position, the controller flashes and the alarm momentarily sounds.

**The chiller has automatic restart. If the chiller was shutdown as a result of a power failure and power is restored, it will restart. ▲**

### Starting

Press . The controller will do a self-test (the controller will quickly sequence through its LEDs and momentarily sound the alarm) and then display the recirculating fluid temperature. The refrigeration system and the recirculation pump will then start. The RECIRCULATING PRESSURE gauge will display the pump operating pressure. If the pressure needs adjusting, refer to Chapter 4.

**NOTE** If on start up the chiller's recirculating fluid is outside either temperature limit, the chiller will operate but the appropriate indicator will flash until the fluid is within the limit. ▲

### Setpoint

To display/change the setpoint press  until  illuminates. The display will flash the current setpoint value. Use the arrow keys to change the value. The range is -15°C to 35°C.

**NOTE** If the arrow keys are not pressed within 10 seconds the display returns to the current reservoir temperature. The controller will not allow you to enter a setpoint closer than 2°C of either temperature alarm setting discussed on the next page. Trying to use a setpoint within 2°C causes the appropriate indicator to flash and sounds the audible alarm twice. ▲

Once the desired setpoint is displayed, press  again to confirm the change. The display rapidly flashes the new value for a short period and then return to the recirculating fluid temperature.

**NOTE** If the value is not confirmed within 60 seconds the display will return to the recirculating fluid temperature and ignore any changes. ▲

## Temp Alarms

To display/change the temperature alarm setting press  until either  or  illuminates. The display flashes the current limit value. Use the arrow keys to change the value. The low-end range is -20°C to 30°C. The high-end range is -10°C to 40°C. You cannot set a low limit higher than the high limit or a high limit lower than the low limit.

**NOTE** If the arrow keys are not pressed within 10 seconds the display returns to the current reservoir temperature. You cannot set either alarm closer than 2°C of the setpoint. ▲

Once the desired setting is displayed, press  again to confirm the change. The display will rapidly flash the new value for a short period and then return to the recirculating fluid temperature.

**NOTE** If the value is not confirmed within 60 seconds the display returns to the recirculating fluid temperature and ignore any changes. ▲

## Remote/Serial Comm

**(Optional)** Once your remote/serial comm cable is connected to the chiller, enable or disable remote operation or serial communication by pressing  until  illuminates. The display flashes the communication status. Use the up arrow key to enable, the down arrow to disable. See Chapter 3 for additional information.

**NOTE** If the arrow keys are not pressed within 10 seconds the display returns to the current reservoir temperature. ▲

Once the desired status is displayed, press  again to confirm the change. The display rapidly flashes the new status for a short period and then returns to the recirculating fluid temperature.

**NOTE** If the status is not confirmed within 60 seconds the display returns to the recirculating fluid temperature and ignore any change. ▲

**NOTE** Enabling serial communications disables the remote start/stop feature. In either mode, the setpoint and alarm values can not be changed using the controller but they can still be viewed. ▲

## External Sensor

**(Optional)** Once the external sensor cable is connected to the chiller, enable or disable the sensor by pressing  a fifth time. The display alternates with *rPrb* and either *on* or *off*. If desired, use the arrow keys to change the status and then press  to return to the temperature display. See Chapter 3 for additional information.

## Stopping

To turn the chiller off, press .

When the display goes blank it is safe to place the circuit protector located on the rear of the chiller to the down position.

Using any other means to shut the chiller down can reduce the life of the compressor. ▲

The circuit protector located on the rear of the chiller is not intended to be used as a disconnecting means. ▲

Always turn the chiller off and disconnect it from its supply voltage before moving it. ▲

## Setup/Tuning Loop

The Setup/Tuning Loop, see illustration next page, is used to configure the controller temperature display and various operating parameters. To enter the loop you must be displaying the reservoir fluid temperature and then, while pressing and holding  $\blacktriangledown$ , press  $\odot$ . The display indicates *tunE*. Use the arrow keys to enter/bypass the loop. Once in the loop press  $\odot$  to sequence through it. Use the arrow keys to change any display.

The loop is used to determine how the chiller reacts when a fault occurs *FLtS*— either shut down *on* or continue to run *off*. The chiller is shipped configured to run.

It is also used to configure the temperature displays to indicate to a tenth of a degree *tEnt*.

The loop is used to select which temperature sensor the high and low temperature monitors uses *APrb*. *rtd1* is the internal sensor, *rtd2* is the optional external sensor. *rtd2* is available only if the external sensor is properly connected to the rear of the chiller.

If the optional remote probe is enabled, the loop is used to set the remote probe temperature deviation *rPrd* limit. This limit is applied to the setpoint to calculate the minimum fluid temperature supplied to your application. For example, if the setpoint is 10°C and the deviation is set to 3°C, the coldest fluid temperature supplied to your application is limited to 7°C.

The deviation can be set from 1°C to 53°C. The factory default setting is 50°C.

The loop can also be used to set the controller's PID parameters. The factory preset values are in Chapter 4.


**Thermo Fisher recommends that only a qualified technician change PID parameters. Incorrect PID values will hamper chiller performance. ▲**

The heat PID is only applicable to chillers configured with a heater.

The loop is also used to adjust/verify serial communication parameters *SEr*. This section of the loop is applicable only to chillers configured for serial communications capability.

Changes made in the Tuning/Setup Loop take affect when  $\blacktriangle$  is pressed at the *Stor* prompt. Pressing  $\blacktriangledown$  at the prompt aborts *all* changes.


## Chapter 3 Options

### Communication Connector

A female 9-pin D-connector is located on the rear of the chiller. The connector is used for RS-232 serial communication (see Appendix B).

When the chiller is configured for serial communication the  LED on the controller is illuminated. If the chiller is shut off while still configured for serial communication, the temperature display goes blank, but the  LED remains illuminated.

When the chiller is configured for serial communication it can be stopped using the local controller. To restart, send another start command.

If the chiller *is* running and serial communications is enabled  (LED on), and you need to operate the chiller using the keypad – the procedure is to press and hold  for approximately two seconds until the low temperature alarm icon illuminates. Press  or , to scroll to the  LED, press the down arrow and then press  again to turn communications off. The  LED will extinguish and the chiller can now be operated from the keypad.

**NOTE** If serial communications is enabled, and the chiller *is not* running, ensure the circuit protector is on and then simultaneously depress and hold both the up and down arrow keys for approximately 10 seconds. The display will show the temperature and the alarm, if enabled, will sound. Then use the above procedure. ▲

The pin out information is:


Pin	RS-232 COMM
1	No connection
2	TX
3	RX
4	No connection
5	GND = Signal ground
6 - 9	No connection

TX = Transmitted data from controller

RX = Received data to controller.

Hardware Mating Connector AMP Part# 745492-2 or equivalent

## Remote Start/Stop

A male 9-pin D-connector is located on the rear of the chiller. The connector is used for either 12 - 24VAC or 12 - 24VDC remote start/stop; and as a status relay dry contact.

To start the chiller from a remote location apply the correct voltage to pins 6 & 8. When the chiller is started by a remote start/stop voltage  on the controller will flash.

Once started remotely, stop the chiller by removing the applied voltage from pins 6 & 8 or by pressing  on the local controller. To restart remotely after stopping locally, remove and reapply the remote start voltage to J101.

**NOTE** Enabling serial communications disables the remote start/stop feature. ▲

Should the chiller experience a high/low temperature warning/fault, a reservoir level warning/fault, or a low pump flow warning/fault, the status relay will de-energize.

The pin out information is:

Pin #	Function	Pin #	Function
1	No Connection	6	Remote On +
2	Status, NO	7	RTD2B
3	Status, NC	8	Remote On -
4	RTD2A	9	No Connection
5	Status, Common		

Hardware Mating Connector AMP Part# 745491-2 or equivalent

### Remote Start Input

12 to 24 Volts DC  $\pm 10\%$  maximum input current, 9ma.

12 to 24 Volts AC  $\pm 10\%$  maximum input current, 5ma.

### Absolute Maximum Input Voltages

40 volts rms. continuous.

56 volts peak.

Rating for both normally open and normally closed status relay contacts are: Maximum rated load - 1.0 Amp, 24 Volts AC/DC, resistive. Minimum permissible load - 1 ma @ 5 VDC.

## External Temperature Sensor

Connect the external sensor to the male 9-pin D-connector located on the rear of the chiller, see previous page.

The external temperature sensor is enabled using the controller. With the controller displaying the temperature press  five times. The display alternates with *rPrb* and either *on* or *off*. If desired, use the arrow keys to change the status and then press  to return to the temperature display.

If the external sensor is selected the controller displays the temperature at the external sensor. Also, the controller indicates the external sensor is selected by flashing the left-most temperature decimal point.

Flashes to indicate the external sensor is enabled.

**NOTE** To display the internal temperature press and hold  for one second. The display indicates *rtd1* for ½ second and then displays the internal temperature for five seconds. The display then indicates *rtd2* for ½ second and returns to the external temperature. The  key needs to be released and then pressed and held again to repeat the cycle. To display the remote probe temperature deviation use the Setup/Tuning Loop, see Chapter 2. 

The controller continues to use the existing setpoint. When using the external sensor the controller will not allow internal chiller temperatures outside -17°C to 37°C.

The selected sensor is retained by the controller if it is turned off, or if power is removed.

If there is an external sensor error the chiller automatically returns to the internal sensor. An error code flashes and if the fault mode is on the chiller shuts down. See Error Codes in Chapter 4. Press  to restart the chiller.

The external probe can be used to control the high and low temperature monitors using the Setup/Tuning Loop, see Chapter 2.

**NOTE** If you also need to monitor status while using the external sensor, an optional Y-connector cable is available from Thermo Fisher. 

## Reservoir Level Warning/Fault

Should the reservoir fluid level drop below normal the mutable alarm will sound, the status relay will de-energize (if the chiller is equipped with the communications package), an error message will be displayed and, depending on the chiller's configuration, the chiller will continue to run or shut down. The chiller is configured to run or shut down using the controller's Setup/Tuning Loop in Chapter 2.

Refill the reservoir and, if necessary, restart the chiller.

**NOTE** The chiller will not start without any fluid in the reservoir regardless of how it is configured. ▲

## Low Pump Flow Warning/Fault

Should the pump flow drop below 0.5 gpm the mutable alarm will sound, the status relay will de-energize (if the chiller is equipped with the communications package), an error message will be displayed, and, depending on the chiller's configuration, the chiller will continue to run or shut down.

Once the cause of the low flow is identified and corrected restart the chiller, if necessary.

The chiller is configured to run or shut down using the controller's Setup/Tuning Loop, see Chapter 3.

## Chapter 4 Maintenance

This chapter covers the chiller's basic maintenance.

### Error Codes

*Er00* ROM checksum. Clear with  key.

*Er01* Test failure. Locks up the program.

*Er02* Display board failure. Clears when display board is fixed.

*Er03* Critical checksum failure. Clear with  key.

*Er14* Synchronous comm - check connections.

*Er15* Asynchronous comm

*Er16* Bad calibration data - redo calibration

### Unusual Hardware Conditions

These errors will flash on the display and cannot be cleared. These are internal controller problems.

*Er04* through *Er13* - Interrupt errors during runtime.

*Conf* BOM invalid

### Functional/ Machine errors

These errors will clear themselves once the problem disappears.

Display	Indication
---------	------------

<i>Lot</i>	Low temp setpoint warning/fault
------------	---------------------------------

<i>Add</i>	Low level warning, fluid below minimum fill level, see Chapter 2
------------	------------------------------------------------------------------

<i>H it</i>	High temp setpoint warning/fault
-------------	----------------------------------

<i>Er25</i>	Shorted internal temperature sensor (rtd1)
-------------	--------------------------------------------

<i>Er26</i>	Open internal temperature sensor (rtd1)
-------------	-----------------------------------------

**LoFL** Low process flow warning/fault

**LLF** Fluid low level fault - chiller shuts down

Chiller reaction to warning/fault errors depend on how the chiller is configured, see Setup/Tuning Loop in Chapter 2. The chiller is shipped configured to continue running. If any other code appears contact our customer service.

## Displaying Software Version Number

The controller can display the installed software version number. For example, for a chiller with software version **026950.1A**:

1. Chiller is running normally and displaying recirculating fluid temperature.

2. Press and hold  for at least 10 seconds. The display will show the first two digits, for example: **02**.

3. Press .  
The display will show the remaining digits to the left of the decimal, for example: **6950**.

4. Press .  
The display will show the decimal point and the digit to the right of the decimal point, for example, **.1**.

5. Press .  
The display will show the revision letter, as its equivalent number, for example, **1 = A**.

6. Press .  
Disregard this display.

7. Press .

The display returns to the recirculating fluid temperature.

# Controller PID Values

The factory preset PID values are:

P	I	D
12	0.5	0.0

Thermo Fisher recommends that only a qualified technician change PID parameters. Incorrect PID values will hamper chiller performance. ▲

## Cleaning

### Reservoir

Periodically inspect the fluid inside the reservoir. If cleaning is necessary, flush the reservoir with a cleaning fluid compatible with the circulating system and the cooling fluid.

The cooling fluid should be replaced periodically. Replacement frequency depends on the operating environment and amount of usage.

Before changing the cooling fluid ensure it is at a safe handling temperature. ▲

**NOTE** It will be necessary to use a wet-vac or siphon to completely drain the reservoir. Use care so as not to damage the reservoir float switch. ▲

See Chapter 2, Filling Requirements, for instructions on cooling fluid replacement.

### Condenser

For proper operation, the chiller needs to pull substantial amounts of air through a condenser. A build up of dust or debris on the fins of the condenser will lead to a loss of cooling capacity. Optional air filters are available, contact our Service Department. See Preface, After-Sale Support.

The lower front of the chiller has a one-piece grille assembly. Using your hands gently pry the assembly off. Use care not to scratch the paint.

Periodic vacuuming of the condenser fins is necessary. The cleaning frequency depends on the operating environment. After initial installation we recommend a monthly visual inspection of the condenser. After several months, the cleaning frequency will be established.

Use care cleaning the condenser fins, they can be easily bent. ▲

## Power Cable Mating Ends

Every five years inspect the power cable mating ends for signs of corrosion and damage. Clean or replace as required.

## Leaks

A very small leak between the pump and motor is normal. If this is a new pump, allow it to run for 24 hours and reinspect. A continuous leak will require pump replacement.

## PD Pressure Relief Valve

**NOTE** The valve is factory preset for the most common applications and normally requires no further adjustment. The valve is factory preset not to exceed  $80 \pm 5$  psi ( $550 \pm 1$  kPa). ▲

**NOTE** PD pumps are capable of 110 psi maximum output pressure. Ensure your plumbing is rated to withstand 110 psi at the highest operating temperatures. ▲

System back pressure is installation and application dependent. Due to internal Merlin back pressure, the minimum pressure setting for a deadheaded pump is 28 psi for a PD-2 pump (these settings prohibit external flow from the Merlin).

Before adjusting the valve turn the chiller off. Locate the circular relief valve opening on the rear of the chiller. Turn the threaded stem fully clockwise.

If the chiller is not plumbed to an application, install a loop of hose equipped with a shut-off valve between the supply and return fittings.

Turn the chiller on. Ensure that there is back pressure on the system. Use the pressure gauge to read the relief valve setting.

Back out the threaded stem on the relief valve counterclockwise. Continue until the gauge indicates 80 psi (550 kPa) or the desired setting.

**NOTE** The relief valve may drip if the threaded stem is backed out too far. ▲

**NOTE** Should the chiller start to vibrate the valve setting may be the cause. Changing the pressure setting  $\pm 5$  psi will eliminate the vibration. ▲

## Pump Strainer

If debris is in the system, the strainer will prevent the material from being drawn into the flow switch and the pump, damaging the pump vanes. A clogged strainer will also cause increased discharge pump pressures.

After initial installation, the strainer may become clogged with debris and scale. Therefore, **the strainer must be cleaned after the first week of installation.** After this first cleaning, we recommend a monthly visual inspection. After several months, the cleaning frequency will be established.

**Before cleaning the strainer ensure the cooling fluid is at a safe handling temperature. Disconnect the power cord from the power source. ▲**


We also recommend draining the reservoir.

The strainer is behind the right-side panel.

**NOTE** A small amount of fluid will come out of the internal strainer housing when the fitting is removed. Have a rag handy to put under the housing before removing the fitting. ▲

Remove the panel from the chiller. Use one wrench to support the filter housing, and a second wrench to unscrew the  $\frac{3}{4}$  inch fitting. Remove the screen.

**NOTE** Some chillers with  $\frac{3}{4}$  inch fittings also have a brass end cap, see illustration below. ▲

Clean the screen by rinsing it with water. When the screen is clean, replace it and the cap. Tighten the fitting and replace the panel. Refer to Section III for instructions on replacing the cooling fluid.

## Temperature Sensor Calibration

### Internal Sensor

If the temperature fluid temperature display disagrees with your reference thermometer, the internal temperature sensor (rtd1) may need calibration.

Do not pick calibration points that are outside the safe operating limits of the fluid in your application. For example with water, 35°C and 5°C would be typical calibration points.


Run the chiller to a high-end calibration point. Place a calibrated reference thermometer in the reservoir, ensure the fluid temperature is stable.

Enter Calibration by pressing and holding  and then pressing  three times. The display will read CAL. Press  and the display will read rtd1.

To calibrate the high-end temperature press  and the display will read r1H. Press  again and the display will flash between r1H and the current probe temperature.

Use  and  to adjust the r1H display. Once the display temperature matches the reference thermometer press  twice. Press  until Stor is displayed, press  to accept the new value.

Repeat for the low-end temperature r1L.

**NOTE** Both the high and low temperatures must be entered for a valid calibration. If power is lost before the procedure is complete, critical data needed for the calculation of calibration parameters will be lost. After pressing ▲ at the Stor prompt wait several seconds before proceeding to ensure that a bad calibration message does not appear. Premature use of the keypad after pressing ▲ may cancel the bad calibration error message. ▲

## Decommissioning/ Disposal

Decommissioning prepares equipment for safe and secure transportation.

Laboratory Grade Ethylene glycol (EG) is poisonous and flammable. Before disposing refer to the manufacturer's most current SDS for handling precautions. ▲

Decommissioning must be performed only by qualified dealer using certified equipment. All prevailing regulations must be followed. ▲

Consider decommissioning the chiller when:

- It fails to maintain desired specifications
- It no longer meets safety standards
- It is beyond repair for its age and worth

Refrigerant and compressor oil must be recovered from equipment before disposal.

**Note** Keep in mind any impact your application may have had on the chiller.

Direct questions about chiller decommissioning or disposal to our Sales, Service and Customer Support.

Handling and disposal should be done in accordance with the manufacturers specification and/or the SDS for the material used. ▲

## Shipment Storage

Follow the manufacturer's SDS instructions if decontamination is required. ▲

Transporting and/or storing the chiller in near or below freezing temperatures requires draining, see **Draining** in this Section. Store the chiller in the temperature range of -25°C to 60°C (with packaging), and <80% relative humidity. ▲

If the chiller is stored for more than 90 days it must be flushed with clean fluid before operating. ▲

# Chapter 5 Troubleshooting

This chapter covers the basic troubleshooting procedures.

## Checklist

### 1. Chiller will not start or shuts down

- a. Check the line cord; ensure it is plugged in.
- b. Check the position of the circuit protector on the rear of the chiller.
- c. Check the controller for error codes, see Chapter 4.

**NOTE** On chillers with a Low Flow Switch, and that are configured to shut down with a low flow fault (LFFL), several starting attempts may be necessary. ▲

- d. On chillers with a Reservoir Level Switch, check the reservoir level.
- e. On chillers with optional serial communications or remote start/stop, ensure the desired mode is enabled.
- f. Make sure power source voltage is within the chiller rated voltage  $\pm 10\%$ .
- g. If  is lit, the chiller was shut off using serial communications. Start the chiller using serial communications. To start the chiller using the local controller press and hold both arrow keys for 10 seconds, the controller will power. Turn serial communications off using the procedure in Chapter 2, and then press  to turn off the controller. Press  again to start the chiller.

## 2. Chiller will not circulate fluid

- a. Check the reservoir level. Fill, if necessary.
- b. Check the position of the reservoir tank valves. See Reservoir Tank Isolation Valves/Open Tank Circulation, Chapter 2.
- c. Check the instrument being cooled for restrictions in the cooling line.
- d. Check the pump strainer.
- e. The pump motor may be too hot. Allow it to cool.
- f. Check the pressure gauge, adjust the relief valve as necessary.
- g. Check the pump and coupling. Make sure the power is off and remove the top and right-side panels. Locate the brass pump and grey motor. A silver clamp ring holds the two together. Remove the clamp and separate the pump and motor. First look to see if the coupling is broken or not, and then check the pump shaft by rotating. If it is seized or sticks then it has failed.

## 3. Inadequate temperature control

- a. Verify the setpoint.
- b. If the temperature continues to rise, make sure your application's heat load does not exceed the rated specification, see Chapter 2, Specifications.
- c. Make sure the air intake and discharge are not impeded and the ambient temperature does not exceed +35°C.
- d. Make sure the condenser is free of dust and debris, see Chapter 4, Maintenance and Cleaning.
- e. Verify the controller's PID values, see Setup/Tuning Loop.

#### 4. Chiller runs loudly

- a. Check the position of the reservoir tank valves. See Reservoir Tank Isolation Valves/Open Tank Circulation, Chapter 2.

#### 5. Condensation on reservoir

- a. Check the position of the upper reservoir tank valve. See Reservoir Tank Isolation Valves/Open Tank Circulation, Chapter 2.

#### 6. Chiller vibrates

- a. The external pressure relief valve may need adjustment. Change the pressure setting  $\pm 5$ psi, see Chapter 4.

#### 7. No serial communications (optional)

- a. All chillers are tested for serial communications before they leave the factory. Check the controller, ensure serial communication is enabled.
- b. Check all communications commands, they *must be exact*. See Appendix.
- c. Check all wiring for proper connections or possible shorts.
- d. Software to verify serial communication is available from Thermo Fisher.

## Service Assistance

If, after following these troubleshooting steps, your chiller fails to operate properly, contact our Service Department for assistance (see Preface, After-sale Support). Before calling, *please* obtain the following information:

Before calling, *please* obtain the following information:

- chiller BOM number
- chiller serial number
- chiller software version
- voltage of power source

# Appendix NC Serial Communications Protocol

**NOTE** This appendix assumes you have a basic understanding of communications protocols. ▲

All data is sent and received in binary form, do not use ASCII. In the following pages the binary data is represented in hexadecimal (hex) format.

The NC Serial Communications Protocol is based on a master-slave model. The master is a host computer, while the slave is the chiller's controller. Only the master can initiate a communications transaction (half-duplex). The slave ends the transaction by responding to the master's query. The protocol uses an RS-232 serial interface with the default parameters: 9600 baud, 8 data bits, 1 stop bit, no parity and no handshaking.

**NOTE** Before the chiller will communicate, RS-232 must be turned on using the controller. Ensure the  indicator is illuminated. ▲

The chiller can be controlled through your computer's serial port by using the chiller's standard 9-pin RS-232 connection. Data read of the serial port connects to the data transmit (pin 2) of the chiller. Data transmit of the serial port connects to data read (pin 3) of the chiller.


Communication cables are available from Thermo Fisher. Contact our sales department for additional information.

All commands must be entered in the exact format shown in the tables on the following pages. The tables show all commands available, their format and responses. Controller responses are either the requested data or an error message. The controller response *must* be received before the host sends the next command.

The host sends a command embedded in a single communications packet, then waits for the controller's response. If the command is not understood or the checksums do not agree, the controller responds with an error command. Otherwise, the controller responds with the requested data. If the controller fails to respond within 1 second, the host should resend the command.

**NOTE** All byte values are shown in hex, hex represents the binary values that must be sent to the chiller.  
**Do not use ASCII. ▲**

The framing of the communications packet in both directions is:


<i>Lead char</i>	CA (hex).
<i>Addr-msb</i>	Device address is 1
<i>Addr-lsb</i>	Most significant byte of device address is 0.
<i>Command</i>	Least significant byte of device address is 1.
<i>n d-bytes</i>	Command byte (see Table 1).
<i>d-byte 1</i>	Number of data bytes to follow (00 to 03 hex).
<i>...</i>	1 <sup>st</sup> data byte (the qualifier byte is considered a data byte).
<i>d-byte n</i>	...
<i>Checksum</i>	n <sup>th</sup> data byte.
	Bitwise inversion of the 1 byte sum of bytes beginning with the most significant address byte and ending with the byte preceding the checksum. (To perform a bitwise inversion, "exclusive OR" the one byte sum with FF hex.)


The master requests information by sending one of the Read Functions as shown in Table 1. Since no data is sent to the chiller during a read request, the master uses 00 for the number of data bytes following the command byte.

The chiller will respond to a Read Function by echoing the lead character, address, and command byte, followed by the requested data and checksum. When the chiller sends data, a qualifier byte is sent first, followed by a two byte signed integer (16 bit, MSB sent first). The qualifier byte indicates the precision and units of measure for the requested data as detailed in Table 2.

As an example, the master requests to read internal temperature by sending:

If the temperature is -12°C, the chiller would reply:


The master sets parameters in the chiller by sending one of the Set Functions as shown in Table 1. The master does not send a qualifier byte in the data field. The master should be preprogrammed to send the correct precision and units (it could also read the parameter of interest first to decode the correct precision and units needed).

For example, if the master wants to set the setpoint to 30°C, it would send :

The chiller responds:

**Table 1** Commands (All bytes are in hex)

FUNCTION	MASTER SENDS	CHILLER RESPONDS
Read Acknowledge	CA 00 01 <b>00</b> 00 FE	CA 00 01 <b>00</b> 02(v1)(v2)(cs)
Read Status (see Table 3)	CA 00 01 <b>09</b> 00 F5	CA 00 01 <b>09</b> 02(d1)(d2)(cs)
Error		CA 00 01 <b>0F</b> 02(en)(ed)(cs)
Read Internal Temperature (RTD1)	CA 00 01 <b>20</b> 00 DE	CA 00 01 <b>20</b> 03(qb)(d1)(d2)(cs)
Read Setpoint (control point)	CA 00 01 <b>70</b> 00 8E	CA 00 01 <b>70</b> 03(qb)(d1)(d2)(cs)
Read Low Temperature Limit	CA 00 01 <b>40</b> 00 BE	CA 00 01 <b>40</b> 03(qb)(d1)(d2)(cs)
Read High Temperature Limit	CA 00 01 <b>60</b> 00 9E	CA 00 01 <b>60</b> 03(qb)(d1)(d2)(cs)
Read Cool Proportional Band (P)	CA 00 01 <b>74</b> 00 8A	CA 00 01 <b>74</b> 03(qb)(d1)(d2)(cs)
Read Cool Integral (I)	CA 00 01 <b>75</b> 00 89	CA 00 01 <b>75</b> 03(qb)(d1)(d2)(cs)
Read Cool Derivative (D)	CA 00 01 <b>76</b> 00 88	CA 00 01 <b>76</b> 03(qb)(d1)(d2)(cs)
Set Setpoint (control point)*	CA 00 01 <b>F0</b> 02(d1)(d2)(cs)	CA 00 01 <b>F0</b> 03(qb)(d1)(d2)(cs)
Set Low Temp	CA 00 01 <b>C0</b> 02(d1)(d2)(cs)	CA 00 01 <b>C0</b> 03(qb)(d1)(d2)(cs)
Set High Temp	CA 00 01 <b>E0</b> 02(d1)(d2)(cs)	CA 00 01 <b>E0</b> 03(qb)(d1)(d2)(cs)
Set Cool Proportional Band (P = 0.1 - 99.9)	CA 00 01 <b>F4</b> 02(d1)(d2)(cs)	CA 00 01 <b>F4</b> 03(qb)(d1)(d2)(cs)
Set Cool Integral (I = 0 - 9.99)	CA 00 01 <b>F5</b> 02(d1)(d2)(cs)	CA 00 01 <b>F5</b> 03(qb)(d1)(d2)(cs)
Set Cool Derivative (D = 0 - 5.0)	CA 00 01 <b>F6</b> 02(d1)(d2)(cs)	CA 00 01 <b>F6</b> 03(qb)(d1)(d2)(cs)

command bytes shown in **bold**  
qb = qualifier byte, see Table 2  
d1,d2 = 16 bit signed integer of the value being sent or received  
cs = the checksum of the string (see text)

v1. v2 = protocol version  
ed Bad command byte gets echoed)  
en (error number) 01 = Bad command,  
02 = Bad checksum  
xx = no valid data, include in checksum

**Table 1** Commands (continued) (All bytes are in hex)

FUNCTION	MASTER SENDS	CHILLER RESPONDS
Set On/Off Array		
The Set On/Off Array command is used to set and/or request chiller settings. Depending on the software version, the array consists of up to 5 data bytes, d1 - d5. Each byte represents the state of one setting, see ** below. The master can send the array with a value of 0 (turn off/disable), 1 (turn on/enable), or 2 (do not change). 0 or 1 causes the chiller to change the setting, 2 only requests the state of a setting. The slave returns values of 0 (off/disabled) or 1 (on/enabled).		
(SW Version 02950.1A-E)	CA 00 01 <b>81</b> 01(d1)(cs)	CA 00 01 <b>81</b> 01(d1)(cs)
(SW Version 02950.1F-H)	CA 00 01 <b>81</b> 04(d1)... (d4)(cs)	CA 00 01 <b>81</b> 04(d1)... (d4)(cs)
(SW Version 02950.1J-K)	CA 00 01 <b>81</b> 05(d1)... (d5)(cs)	CA 00 01 <b>81</b> 05(d1)... (d5)(cs)
(SW Version 02950.1L- )	Any of the above formats	
Sample - turn chiller on (02950.1J-K)	CA 00 01 <b>81</b> 05 01 02 02 02 02 6F	

\* = limited to the range of the chiller command bytes shown in **bold**

qb = qualifier byte, see Table 2

d1,d2 = 16 bit signed integer of the value being sent or received

cs = the checksum of the string (see text)

**Table 2****QUALIFIER BYTE**

00	0 precision, no units of measure
01	0 precision, °C
10	1 precision, no units of measure
11	1 precision, °C
20	2 precision, no units of measure

\*\*d1 - chiller off = 0, on =1

d2 - external sensor disabled = 0, enabled = 1

d3 - fault mode disabled = 0, enabled = 1

d4 - tenths display disabled = 0, enabled = 1

d5 - Alarms on internal = 0, or external sensor = 1, (1 is invalid when d2 = 0, external sensor disabled)

Example: The temperature of 45.6°C would be represented by the qualifier 11 hex, followed by the 2 bytes 01 C8 hex (456 decimal).

**Table 3** READ STATUS

BIT	d1	d2 <sup>t</sup>
b.7 = 1	External Temp Sensor Enabled	Reserved
b.6 = 1	Level 1 Warning/Auto Refill	Reserved
b.5 = 1	Low Flow Warning	Internal Temp Sensor Fault
b.4 = 1	Low Level 2 Warning	External Temp Sensor Fault
b.3 = 1	High or Low Temp Warning	High Temp/Fixed High Temp Fault
b.2 = 1	High or Low Temp Bypass	Low Temp/Fixed Low Temp Fault
b.1 = 1	Chiller Faulted*	Low Flow Fault
b.0 = 1	ChillerRunning*	Low Level 2 Fault

<sup>t</sup> d2 is the detail of fault d1-b.1

\* mutually exclusive

## DECLARATION OF CONFORMITY

Manufacturer: Thermo Fisher Scientific  
Address: 25 Nimble Hill Road  
Newington, NH 03801

We declare that the equipment named below has been designed to comply with the relevant sections of the below referenced specifications and is in accordance with the requirements of the indicated directives.

Product: Merlin Recirculating Chiller

Models: Units have a 12 digit part number consisting of UUU T VV PP XXXX defined as follows:

UUU = Unit type can be:

262 = M25; 263 = M33; 264 = M75; 265 = M100; 266 = M150

T = Temperature Range can be 1-3 inclusive, where:

1 = Std temp; 2 = Low Temp; 3 = High Temp

V = Unit voltage rating:

12 = 115 / 60 100 / 50; 16 = 200 / 230 / 60 200 / 50; 25 = 230 / 50; 26 = 400 / 50 5-Wire

P = Pump type

02 = CP-55; 03 = PD-1; 04 = PD-2; 15 = MD-30; T1 = 99

X = Any digit 0-9 inclusive, used as sequential numbering only.

Directives and Standards:

2014/30/EC ± Electromagnetic Compatibility Directive (EMC)

- EN 61326-1: 2013 Electrical equipment for measurement, control, and laboratory use - EMC requirements. General requirements

2014/35/EC - Low Voltage Directive (LVD):

- EN 61010-1: 2010 Safety requirements for electrical equipment for measurement, control, and laboratory use -- Part 1: General requirements.
- EN 61010-2-010: 2003 Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials


2011/65/EU – Restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS)

- EN 50581: 2012 - Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

Authorised representative in the EC:

Name: Thermo Fisher Scientific  
Address: Dieselstrasse 4  
76227 Karlsruhe Germany

On behalf of the manufacturer:

  
John Parker  
Electrical Engineering Manager  
Temperature Control  
Thermo Fisher Scientific  
Newington, NH, USA

Date:

15 August 2016

## WARRANTY

Thermo Fisher Scientific warrants for 24 months (**excluding Positive Displacement pumps which are warranted for 12 months**) from date of shipment the Thermo Scientific Merlin chiller according to the following terms.

Any part of the chiller manufactured or supplied by Thermo Fisher Scientific and found in the reasonable judgment of Thermo Fisher to be defective in material or workmanship will be repaired at an authorized Thermo Fisher Repair Depot without charge for parts or labor. The chiller, including any defective part must be returned to an authorized Thermo Fisher Repair Depot within the warranty period. The expense of returning the chiller to the authorized Thermo Fisher Repair Depot for warranty service will be paid for by the buyer. Our responsibility in respect to warranty claims is limited to performing the required repairs or replacements, and no claim of breach of warranty shall be cause for cancellation or rescission of the contract of sales of any chiller. With respect to chillers that qualify for field service repairs, Thermo Fisher Scientific's responsibility is limited to the component parts necessary for the repair and the labor that is required on site to perform the repair. Any travel labor or mileage charges are the financial responsibility of the buyer.

The buyer shall be responsible for any evaluation or warranty service call (including labor charges) if no defects are found with the Thermo Scientific product.

This warranty does not cover any chiller that has been subject to misuse, neglect, or accident. This warranty does not apply to any damage to the chiller that is the result of improper installation or maintenance, or to any chiller that has been operated or maintained in any way contrary to the operating or maintenance instructions specified in this Instruction and Operation Manual. This warranty does not cover any chiller that has been altered or modified so as to change its intended use.

In addition, this warranty does not extend to repairs made by the use of parts, accessories, or fluids which are either incompatible with the chiller or adversely affect its operation, performance, or durability.

Thermo Fisher Scientific reserves the right to change or improve the design of any chiller without assuming any obligation to modify any chiller previously manufactured.

THE FOREGOING EXPRESS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO WARRANTIES OR MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

OUR OBLIGATION UNDER THIS WARRANTY IS STRICTLY AND EXCLUSIVELY LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE COMPONENT PARTS AND Thermo Fisher Scientific DOES NOT ASSUME OR AUTHORIZE ANYONE TO ASSUME FOR IT ANY OTHER OBLIGATION.

Thermo Fisher Scientific ASSUMES NO RESPONSIBILITY FOR INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES INCLUDING, BUT NOT LIMITED TO LOSS OR DAMAGE TO PROPERTY, LOSS OF PROFITS OR REVENUE, LOSS OF THE CHILLER, LOSS OF TIME, OR INCONVENIENCE.

This warranty applies to chillers sold by Thermo Fisher Scientific. (Refer to the warranty for chillers sold by the affiliated marketing company of Thermo Fisher Scientific for any additional terms.) This warranty and all matters arising pursuant to it shall be governed by the law of the State of New Hampshire, United States. All legal actions brought in relation hereto shall be filed in the appropriate state or federal courts in New Hampshire, unless waived by Thermo Fisher Scientific.

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