# **BioCane™** Type CK50900

## **Cryobiological Storage Vessel**

Operation Manual and Parts List LT509X17 (7000807) Rev. 9

Visit us online to register your warranty  
[www.thermoscientific.com/labwarranty](http://www.thermoscientific.com/labwarranty)

**Thermo**  
SCIENTIFIC

Catalog Number	Model Number	Name
CK509X2	807	BioCane 20
CK509X3	808	BioCane 34
CK509X4	809	BioCane 47
CK509X6	820	BioCane 73

**Manual Number LT509X17 (7000807)**

9	41864	11/07/17	Removed CE mark, number	bpg
8	41238	12/20/16	Combined languages to one file	slb
7	40578	2/23/16	Clarified Do Not Roll warning	ccs
6	--	10/27/15	Updated symbols and labels, added cautions	ccs
5	30559	11/12/13	Added warnings and cautions to pg 1-2, environmental conditions to pg 2-1	ccs
4	27797	9/21/11	Added manufacturer icon and euro rep	ccs
3	27547	7/28/11	Updated specs - BioCane 73 vials from 2880 to 2016	ccs
2	27457/SI-10618	7/21/11	Added CE number to pg 2-1	ccs
1	--	2/24/11	Added part numbers for replacement canisters	ccs

**Important** Read this instruction manual. Failure to read, understand and follow instructions in this manual may result in damage to the unit, injury to operating personnel, and poor equipment performance. ▲

**Caution** All internal adjustments and maintenance must be performed by qualified service personnel. ▲

**Caution** It is recommended to have this vessel tested by the manufacturer or qualified cryovessel service technician every 7-10 years, regardless of any problems (or lack thereof) you may have had in the past. This will help insure your samples against sudden loss of liquid nitrogen due to vacuum failure. ▲

**Warning** Nitrogen gas can cause suffocation without warning. Store and use liquid nitrogen only in a well-ventilated place. As the liquid evaporates, the resulting gas displaces the normal air in the area. (The cloudy vapor that appears when liquid nitrogen is exposed to the air is condensed moisture, not the gas itself. The issuing gas is invisible.) In closed areas, excessive amounts of nitrogen gas reduces the concentration of oxygen and can result in asphyxiation. Because nitrogen gas is colorless, odorless, and tasteless, it cannot be detected by the human senses. Breathing an atmosphere that contains less than 19.5% oxygen can cause dizziness and quickly result in unconsciousness and death. Therefore, the use of oxygen monitoring equipment is strongly recommended. ▲

**Warning** Contact of liquid nitrogen with the skin or eyes may cause serious (freezing) injury. ▲

**Warning** Do not roll. Always keep liquid nitrogen vessel in an upright position. Exercise caution to prevent spilling and splashing. ▲

Material in this manual is for information purposes only. The contents and the product it describes are subject to change without notice. Thermo Fisher Scientific makes no representations or warranties with respect to this manual. In no event shall Thermo be held liable for any damages, direct or incidental, arising out of or related to the use of this manual.

**Intended Use:**

The Thermo Fisher Scientific cryogenic storage vessels are intended to store biological samples at cryogenic temperatures. These products are comprised of a vessel and a lid, with or without a level monitor.

**Unintended Use:**

These products are not intended for use within the patient environment.

When this manual is translated into other languages, the US English version is binding.

© 2010 Thermo Fisher Scientific. All rights reserved. Printed in U.S.A.

Important operating and/or maintenance instructions. Read the accompanying text carefully.

Extreme temperature hazards. Only qualified persons should perform procedures associated with this symbol.

Potential biological hazards. Proper protective equipment and procedures must be used.

WEEE Compliance: Thermo Fisher Scientific has contracted with companies for recycling/disposal in each EU Member State. For further information, send an email to [wEEE.recycle@thermofisher.com](mailto:wEEE.recycle@thermofisher.com).

- ✓ Always use the proper protective equipment (clothing, gloves, goggles, etc.)
- ✓ Always dissipate extreme cold or heat and wear protective clothing.
- ✓ Always follow good hygiene practices.
- ✓ Each individual is responsible for his or her own safety.

## Do You Need Information or Assistance on Thermo Scientific Products?

If you do, please contact us 8:00 a.m. to 6:00 p.m. (Eastern Time) at:

1-740-373-4763

1-800-438-4851

1-877-213-8051

<http://www.thermofisher.com>

[service.led.marietta@thermofisher.com](mailto:service.led.marietta@thermofisher.com)

[www.unitylabservices.com](http://www.unitylabservices.com)

Direct

Toll Free, U.S. and Canada

FAX

Internet Worldwide Web Home Page

Tech Support Email Address

Certified Service Web Page

Our **Sales Support** staff can provide information on pricing and give you quotations. We can take your order and provide delivery information on major equipment items or make arrangements to have your local sales representative contact you. Our products are listed on the Internet and we can be contacted through our Internet home page.

Our **Service Support** staff can supply technical information about proper setup, operation or troubleshooting of your equipment. We can fill your needs for spare or replacement parts or provide you with on-site service. We can also provide you with a quotation on our Extended Warranty for your Thermo Scientific products.

Whatever Thermo Scientific products you need or use, we will be happy to discuss your applications. If you are experiencing technical problems, working together, we will help you locate the problem and, chances are, correct it yourself...over the telephone without a service call.

When more extensive service is necessary, we will assist you with direct factory trained technicians or a qualified service organization for on-the-spot repair. If your service need is covered by the warranty, we will arrange for the unit to be repaired at our expense and to your satisfaction.

Regardless of your needs, our professional telephone technicians are available to assist you Monday through Friday from 8:00 a.m. to 6:00 p.m. Eastern Time. Please contact us by telephone or fax. If you wish to write, our mailing address is:

Thermo Fisher Scientific (Asheville) LLC  
401 Millcreek Road, Box 649  
Marietta, OH 45750

International customers, please contact your local Thermo Scientific distributor.

# Table of Contents

<b>Section 1</b>	<b>Safety Information</b> .....	<b>1-1</b>
	Please Read .....	1-1
<b>Section 2</b>	<b>General Specifications</b> .....	<b>2-1</b>
	Unpacking .....	2-1
<b>Section 3</b>	<b>Operation</b> .....	<b>3-1</b>
	Measuring Liquid Nitrogen Quantity .....	3-2
<b>Section 4</b>	<b>Replacement Parts</b> .....	<b>4-1</b>
<b>Section 5</b>	<b>Warranty Information</b> .....	<b>5-1</b>

## Section 1 Safety Information

This manual contains important operating and safety information. The user must carefully read and understand the contents of this manual prior to the use of this equipment.

Your Thermo Scientific cryogenic vessel has been designed with function, reliability, and safety in mind. It is the user's responsibility to install it in conformance with local electrical codes. For safe operation, pay attention to the warnings and cautions throughout the manual.

### **Please Read**


All cryopreservation vessels, regardless of who manufactures them, use a vacuum to provide for the super insulative properties needed to keep nitrogen in a liquid form. They are very similar, in fact, to the thermos vessels you may have used to store coffee, soup or milk.

Remember how fragile they were? Since no vacuum is perfect or will last forever, we suggest that you monitor the consumption of liquid nitrogen used by your vessel on a regular basis, i.e. every 1-2 days. By monitoring, you may be able to anticipate subtle changes in consumption and possibly react to problems before they arise. You may opt to use a level monitor for this purpose.

Also, if your samples are super-critical, one-of-a-kind, or irreplaceable, consider allocating your samples into 2 separate vessels.

Remember the old adage about placing all your eggs into one basket...

Finally, because no vacuum will last forever, consider having your vessel tested by a qualified cryovessel manufacturer every 7-10 years regardless of any problems (or lack thereof) you may have had in the past. This will help insure your samples against sudden loss of nitrogen due to vacuum failure.

## Warnings

- This unit is intended to be used by trained personnel in a well-ventilated environment.
- Liquid nitrogen is extremely cold; it boils at  $-196^{\circ}\text{C}$ .
- Follow all necessary precautions and warnings relevant to applicable hazards such as frostbite, contact burn, asphyxiation, etc. in cryogenic areas.
- To avoid injury due to frostbite or ruptured vials, use extreme care whenever handling liquid nitrogen, liquid nitrogen storage or transfer vessels or any objects which have come in contact with liquid nitrogen.
- Leave no areas of skin exposed.
- Wear adequate PPE where necessary: face shield, cryogenic gloves, cryogenic apron.
- Only use the power supply provided with the unit.
- Do not tightly seal liquid nitrogen containers or prevent nitrogen gas from escaping.
- Always handle vials with tongs.
- Use extreme care to prevent spilling and splashing liquid nitrogen during transfer and removal of storage contents and holders.
- Immediately remove any clothing or safety attire on which liquid nitrogen has been spilled or splashed .
- Get immediate medical attention for any frostbite injuries due to liquid nitrogen.
- Never overfill liquid nitrogen vessels.

**When inserting or removing cans, be careful not to come in contact with the neck tube area of the vessel. Remove or insert cans slowly, in a vertical manner. Scratches on the neck tube can cause premature vacuum failure.**

**Do not tamper with, or remove vacuum port (covered by black plastic cap on side of vessel); this will void warranty.**

Do not spill liquid nitrogen on vacuum port - this can cause vacuum failure.

Always keep liquid nitrogen in an upright position.

## Section 2 General Specifications

	BioCane 20	BioCane 34	BioCane 47	BioCane 73
Height (Inches)	25.70"	26.60"	26.50"	27.00"
Diameter (Inches)	14.50"	18.20"	20.00"	22.00"
LN2 Capacity	20.50L	34.80L	47.40L	73.0L
Static Holding Time (Days)	227	193	121	121
Static LN2 Evaporation Rate*	0.1L/day	0.18L/day	0.40L/day	0.60L/day
Liquid Phase Capacity 1.2 ml	180 vials	720 vials	1512 vials	2016 vials

*\* Actual working time may vary due to current atmospheric conditions, container history, manufacturing tolerances, and any individual patterns of use.*

### Environmental Conditions

Operating temperature: 17°C - 27°C

20% to 80% relative humidity, non-condensing.

Installation Category II (overvoltage) in accordance with IEC 664.

Pollution Degree 2 in accordance with IEC 664.

Altitude limit: 2,000 meters.

Storage: -25°C to 65°C; 10% to 85% relative humidity.

Locator and Locator Plus units are intended for indoor use only.

## Unpacking

Before using your new vessel, carefully inspect it. Check for signs of damage which may have occurred in shipment. It is advisable to fill (see Filling Instructions) all new units with liquid nitrogen and watch liquid nitrogen loss rate for a few days. If there are any problems, call Technical Services as soon as possible.

The on-line warranty registration address is listed on the front of this manual. The unit should be registered within 30 days to make warranty effective. This information must include the serial number which is located on the handle of the vessel.

**Note** The most prevalent cause of failure of liquid nitrogen storage vessels is mechanical. The vessel neck tube supports the full weight of the inner shell and all liquid nitrogen it contains. A side or corner blow to the vessel causes the inner shell to swing in a pendulum motion causing the neck tube to be damaged. Any storage vessel which has been exposed to an accident, has been dropped or lowered to hit on one corner will tend to fail more rapidly than one which has not.

**Caution** Exercise caution when moving your BioCane™ vessel. BioCane cryobiological storage systems are not transportation vessels. Transport carts are designed for mobility within the lab or lab to lab only. Moving full vessels long distances, over cracks in floor, thresholds, on inclined ramps or in elevators can cause premature vacuum failure.

**Note** If samples must be transported under cryogenic conditions, consider the Thermo Scientific Arctic Express line of Dry Shippers.

## Section 3 Operation

To avoid damage to your cryogenic storage vessel which may result in premature vacuum loss, it is important that the following procedure be used during the addition of liquid nitrogen to a warm vessel:

1. Add only a small amount of liquid nitrogen (2-10 liters) to new or warm vessels.
2. Allow this small amount of liquid nitrogen to sit in the covered vessel for a minimum of 2 hours. This will limit stress caused by the sudden temperature change associated with adding liquid nitrogen to a warm vessel.
3. Add an additional 2-10 liters LN2 to vessel.
4. Allow vessel to sit for 48 hours and monitor liquid nitrogen consumption.
5. Fill vessel as desired. Remember to allow for displacement of liquid nitrogen when cannisters and canes are inserted.
6. Insert and remove cannisters slowly. Allow liquid nitrogen to run out of cannisters.

**Caution** Never overfill liquid nitrogen vessels.

**When inserting or removing cans, be careful not to come in contact with the neck tube area of the vessel. Remove or insert cans slowly, in a vertical manner. Scratches on the neck tube area can cause premature vacuum failure.**

Do not spill liquid nitrogen on vacuum port - this can cause vacuum failure. ▲

## **Measuring Liquid Nitrogen Quantity**

1. Use a wooden yardstick to measure liquid nitrogen level. Never use a hollow tube or plastic dipstick to measure liquid nitrogen level.
2. Level will be indicated by frost line which develops when dipstick is removed and waved in a back and forth motion away from the user.

## Section 4 Replacement Parts

Part Number	Description
AY509X13	Nalge Cryovials 1.2ml, 500/case
AY509X12	Nalge Cryovials 2.0ml, 500/case
AY509X20	CryoCane (aluminum) - each cane holds 5 vials for storage in BioCane cannisters (48/case)
AY509X21	CryoSleeve enclosures - Clear plastic sleeve fits over CryoCane for added security during handling or storage (100/case)
CS509X17	Nalgene System 100-ampule box accommodates most internally threaded 1.2 and 2.0 ml vials
CS509X18	Nalgene System 100-ampule box accommodates most internally threaded 1.2 and 2.0 ml vials
AY509X32	Nalgene System 100 Cryogenic Vials - 1 ml (500/case)
AY509X33	Nalgene System 100 Cryogenic Vials - 1.5 ml (500/case)
AY509X34	Nalgene System 100 Cryogenic Benchtop Vial Holder (1/package)
AY509X35	Nalgene System 100 Cryogenic Benchtop Vial Holder (1/package)
<b>CryoCane coders - Aluminum identification tabs that slip over the end of a CryoCane to identify the cane written on for further identification (100/case)</b>	
AY509X22	White
AY509X23	Yellow
AY509X24	Blue
AY509X25	Green
AY509X26	Red
<b>Cryovial closure color coders - Plastic disks that fit into the top of Nalge Cryovials to identify individual (100/package)</b>	
AY509X27	White
AY509X28	Yellow
AY509X29	Blue
AY509X30	Green
AY509X31	Red

**Section 4**  
Replacement Parts

<b>Part Number</b>	<b>Description</b>
HR509X25	BioCane 20 (CK509X2) [807]) replacement canister
HR509X26	BioCane 34 (CK509X3 [808]) replacement canister
HR509X27	BioCane 47 (CK509X4 [809]) replacement canister
HR509X28	BioCane 73 (CK509X6 [820]) replacement canister

## Ordering Procedures

Refer to the Specification Plate for the complete model number, serial number, and series number when requesting service, replacement parts or in any correspondence concerning this unit.

All parts listed herein may be ordered from the Thermo Scientific dealer from whom you purchased this unit or can be obtained promptly from the factory. When service or replacement parts are needed, check first with your dealer. If the dealer cannot process your request, then contact our Technical Services Department.

Prior to returning any materials, contact our Technical Services Department for a “Return Materials Authorization” number (RMA). Material returned without an RMA number will be refused.

## **THERMO FISHER SCIENTIFIC STANDARD PRODUCT WARRANTY (LN<sub>2</sub> Vacuum)**

The Warranty Period starts two weeks from the date your equipment is shipped from our facility. This allows for shipping time so the warranty will go into effect at approximately the same time your equipment is delivered. The warranty protection extends to any subsequent owner during the first year warranty period.

During the first year, component parts proven to be non-conforming in materials or workmanship will be repaired or replaced at Thermo's expense, labor included. *LN<sub>2</sub> Vacuum Integrity is covered for five years.* Installation and calibration are not covered by this warranty agreement. The Technical Services Department must be contacted for warranty determination and direction prior to performance of any repairs. Expendable items, glass, filters and gaskets are excluded from this warranty.

Replacement or repair of components parts or equipment under this warranty shall not extend the warranty to either the equipment or to the component part beyond the original warranty period. The Technical Services Department must give prior approval for return of any components or equipment. At Thermo's option, all non-conforming parts must be returned to Thermo postage paid and replacement parts are shipped FOB destination.

**THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, WHETHER WRITTEN, ORAL OR IMPLIED. NO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE SHALL APPLY.** Thermo shall not be liable for any indirect or consequential damages including, without limitation, damages relating to lost profits or loss of products.

Your local Thermo Sales Office is ready to help with comprehensive site preparation information before your equipment arrives. Printed instruction manuals carefully detail equipment installation, operation and preventive maintenance.

If equipment service is required, please call your Technical Services Department at 1-800-438-4851 (USA and Canada) or 1-740-373-4763. We're ready to answer your questions on equipment warranty, operation, maintenance, service and special application. Outside the USA, contact your local distributor for warranty information.

Rev. 0 5/10

**thermoscientific.com**

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

---

**Thermo Fisher Scientific (Asheville) LLC**  
401 Millcreek Road  
Marietta, Ohio 45750  
United States

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