***Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册***

ECO 6 Thermoreactor

F10100120, F10110120

**General Information / Informazioni Generali / Informations Générales / Información General /
Allgemeine Hinweise / 一般信息**

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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IEC/EN 61010-1

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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COD Analysis

It is the amount of oxygen from chemical oxidants required for the complete oxidation of substances dissolved or suspended in water. By COD determination it is possible to evaluate quantitatively substances which can be oxidized by strong oxidants like potassium dichromate in highly acidic solutions by concentrated sulfuric acid. The method allows the determination of substances which are not oxidized by the corresponding biological method (BOD, Biological Oxygen Demand), e.g. cellulose. Piridine, benzene and toluene are not oxidized. Chemical oxidation is more complete than biological oxidation and therefore COD values for a given sample are generally higher than BOD values. According to the oxidant used and to its concentration, to the concentration of acid, to the catalyst, to the temperature and time, different substances are more or less completely oxidized and analytical conditions must be strictly standardized. The determination of COD is performed by "return titration" of excess dichromate, added in known amount, after the completion of the oxidation of the solution containing chemically oxidable substances. The amount of oxidable substances which are present in a given sample is directly proportional to the amount of dichromate consumed. The addition of silver sulfate as a catalyst increases the oxidation speed of alcohols and acids but not that of aromatic hydrocarbons. The presence of high concentrations of iodides, bromides or chlorides can result in results higher than real. This inconvenient is overcome by adding mercury sulfate which forms insoluble complexes with halogens.

Reagents

Dichromate standard solution 0.250 N - 15-20 g of potassium dichromate analytical grade ($K_2Cr_2O_7$) are oven dried for two hours at 150 °C and then left to cool in a desiccator. 12.259 g are weighed and dissolved in 400-500 ml distilled water in a volumetric flask of 1000 ml volume. The flask is then filled to the mark with distilled water.

Sulphuric acid reagent - 5.4 g of silver sulphate analytical grade (Ag_2SO_4) are added to 1 kg of concentrated sulphuric acid 96% analytical grade (H_2SO_4 d = 1.835) corresponding to 545 ml. Complete solution requires two days. The reagent is preserved indefinitely in a well stoppered and dark glass bottle.

Ammonium ferrous sulfate standard solution 0.125 N - 49.02 g of ammonium ferrous sulfate exahydrate analytical grade ($Fe(NH_4)_2(SO_4)_2 \cdot 6H_2O$), are weighed and dissolved in half a liter about of distilled water. 20 ml of concentrated sulphuric acid analytical grade (H_2SO_4 96% d = 1.835) are slowly added to the ammonium and iron solution while shaking.

Do not add the solution to the concentrated acid in order to avoid dangerous splashes. Fill to 1000 ml in a graduated flask with distilled water. Keep in a dark bottle. The titre decreases with time and is to be controlled every time an analysis is performed by comparison with a standard solution of potassium dichromate.

Solution of Ferroin indicator - 1.48 g of O-phenanthroline analytical grade (1.60 g of monohydrate) and 0.695 g of ferrous sulfate heptahydrate analytical grade ($FeSO_4 \cdot 7H_2O$) are weighed and dissolved together in 50 ml about of distilled water. The final volume of 100 ml is obtained by adding distilled water.

Mercury (II) sulphate - Small crystals of mercury (II) sulphate analytical grade ($HgSO_4$).

Potassium hydrogen phthalate standard solution - Using potassium hydrogen phthalate RPE ($H_5C_8O_4K$) dried to constant weight at 120 °C, a solution containing 425 mg in 1000 ml is prepared. The COD of this solution is 500 mg O_2/l . If stored in a refrigerator at + 4 °C its shelf life is 3-4 months.

Procedure

Sample digestion

The reagents and the sample to be tested are introduced by calibrated pipettes in perfectly clean digestion vials in the order and in the amounts shown in the following table.

Reactor	Mercuric sulphate mg	Dichromate 0,25 N ml	Sulphuric acid ml	Sample ml
ECO 6	400	10	30	20

The amounts of mercury sulphate, the only solid reagent, can be measured by a small spoon purposely calibrated, taking into account that 100 mg of mercury sulphate are able to complex 10 mg of chloride ions. If the water sample contains chloride ions in higher concentration, the amount of mercury sulphate added is to be increased on the basis of mercury sulphate and chloride ions weight ratio of 10 to 1. The appearance of a precipitate after the addition of reagents does not influence the results. Place the air condensers with antisplash bells upon the vials. Blank prepared by heating tend to overheat and then to start boiling in a violent manner. To avoid this inconvenient it is advisable to add to the vials with blank solutions some bumping stones (small spheres or capillaries or chips of glass; small pieces of pumice). The vials are introduced into the VELP reactor. After the end of the heating period of 120 min. at 150 °C (the boiling point of the used 50% sulphuric acid solution), which is shown by a buzzer, the vials are extracted and left to cool.

Sample titration

The content of each vial is poured in a wide mouth Erlenmeyer flask together with the distilled water used for rinsing 3-4 times. 5-6 drops of ferroin solution are then added. After cooling, the flask content is titrated with ammonium iron sulfate standard solution until the color changes from blue-green to orange.

Control of the titer of ammonium ferrous sulfate standard solution 0.125 N

The titer of ammonium ferrous sulfate solution is to be controlled during each day because it decreases with time.

The variation of titer must be considered when computing COD.

This is obtained by introducing in an Erlenmeyer flask of suitable capacity.

Dichromate 0,25 N ml	Concentrated sulfuric acid ml
10	30

bringing to 100 ml volume using distilled water. After cooling the titration is performed with ammonium ferrous sulfate solution till turning to orange. The used volume in ml is used to compute the correction factor. If the titer of solution is exact, 20 ml are necessary to titre the used amount of dichromate while a larger volume is necessary with lower normality. Accordingly the correction factor is given by $F = 20/n$ where "n" is the number of ml actually used.

Evaluation of blank

The same procedure used at 2.3 is used to measure the amount of ammonium ferrous sulfate standard solution necessary to turn the colour in a volume of distilled water.

Calculation of COD

The C.O.D. value of the examined sample is obtained by:

$$\text{C.O.D.mg/l} = (b - a) \cdot \frac{8000 \cdot N}{X} \cdot F_c$$


where:

b = ml of ammonium ferrous sulfate solution used to titer blank.

a = ml of ammonium ferrous sulfate solution used to titer the sample.

N = normality of ammonium ferrous sulfate. **X** = ml of sample under examination. **F_c** = correction factor.

The reactions interested in the analysis are the following:

Bibliography

- Water Research Institute - National Research Council of Italy (1993) Method n° 5110. Chemical Oxygen Demand (C.O.D.) pages 285-289 (in Italian) Graphical Institute and Mint of the State. Rome, Italy.
- APHA-AWWA-WPCF (1985) Standard Methods for the examination of water and wastewater. Method n° 508 - Oxygen demand (chemical). pages 532-538. 16th Ed. Washington DC USA.

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Thermo reaction can be carried out quickly and safely by placing the tubes into the thermoreactor heating block. Set the working temperature from room to 200 °C and the working time up to 199 minutes, or continuous. ECO 6 offers the possibility to thermoreact 6 tubes (Ø 42 mm).

2. Assembly and installation

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box includes:

- ECO 6
- Power supply cord
- Instruction manual

2.1 Electrical connection

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and easy to reach.

2.2 Start-up

Use the on-off switch on the front of the unit to turn the unit on (position “-”). When switched on the unit all the leds blink for few seconds. When the leds are off the unit is ready to be programmed at different work temperatures and times.

3. Operating controls

Shift (A): it moves the messenger under all the parameters to be set during programming.

Increase (B): it shifts the value of the parameter corresponding to the messenger position.

Enter (C): it stores the values shown by the display.

Start (D): it sets going the work cycle.

Main switch (E): it turns on and off the instrument.


Even 20 ml of sample can be tested. The following glassware is to be used:

Equipment

Two notch calibrated pipettes
Wide mouth Erlenmeyer flasks
Volumetric flasks
Burette
Graduated cylinders

ECO 6 - sample volume 20 ml

5-10-20-30-100 ml
500 ml
100-1000 ml
10 ml (subdivision 1/20)
50-250 ml

3.1 Programming

Turn on the main switch on control panel (E).

The unit performs in a few second a "Self-Test", after which: VELP SCIENTIFICA 20 °C *READY*

Pressing the Shift push-button the display shows two values contained in the memory (temperature and time):

Set Temp. Time
100 °C 30'

Where:

Temp. (100 °C) is the temperature set. Range from room temperature to 200 °C.

Time (30') is the time of staying at the temperature set. Range from 0 to 199 minutes.

A continuous operation can be set by a 200, 300... value of time. A symbol === is displayed in this case.

The parameters set are stored by pressing Enter (C). When the unit is turned off the set values are retained in the memory.

3.2 Work cycle

It is set going by pressing the Start (D) and the display shows: Working 20 °C 60'

Where:

Working Shows the operation of the work cycle.

20 °C Is the present temperature of heating block.

60' Is the time of staying at the set temperature.

When the set value of temperature is reached the value of time displayed starts to decrease with a resolution of one second.

At the end of work cycle a beeper warns and the display shows: *END*

When a new work cycle should start with a temperature of the heating block higher than the set value, the display shows: VELP SCIENTIFICA 160 °C NOT READY

Where the 160 °C refers to the temperature of the heating block and NOT READY refers to the instrument operation.

If the Start push-button is pressed the cycle does not start and the display shows: NOT ACTIVE Temp. too high

Should a work cycle be interrupted the push buttons (C) and (B) must be pressed contemporaneously.

3.3 Heating block temperature verification

The temperature is controlled electronically which ensures stability without temperature oscillations around the set point. The Pt 100 probe provides high precision and does not require calibration because the internal software performs a self-calibration every time the instrument is turned on. If, in spite of everything, a control is necessary, it is possible to measure the temperature of the heating block with an external 3mm probe, introducing the 3mm probe adapter into the hole in the upper part of the block

3.4 Power interruptions

The unit will register a black out or a lack of power supply as a manual stop. The work cycle will have to be restarted.

3.5 Safety features

The unit is protected against overheating by a thermostat that turns off the power supply when the temperature exceeds the safety temperature. Temperature setting is continuously controlled by electronic equipment that stops thermoregulation if anomalous measured are sensed. In case, all the leds start blinking and an acoustical warning signal sounds.

3.6 End of work cycle

At the end of the count-down the work cycle will stop automatically. To stop in every moment the work cycle turn off the general switch. When the work session is completed, should the unit remain connected to the power supply, switch it off using the on-off switch on the frontal panel.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Technical data

Power	700 W
Dimensions (WxHxD)	198x132x319 mm
Weight	5,6 Kg
Holes and diameter	6 pos. (Ø 42 mm)
Construction material	Epoxy painted aluminum cast
Thermoregulation	Pid microprocessor
Selectable temperatures	room - 200 °C
Selectable times	0 - 199 min. - infinite
Set temperature reached	Visual signal
Count-down	Visual signal
End of cycle	Acoustic and visual signal
Heating block temperature stability and homogeneity	± 0.5 °C
Temperature precision	± 1.0 °C
Enviromental ambient temperature	5 - 40 °C (41 - 104°F)
Storage temperature range	-10... + 60 °C (14 - 140°F)
Max humidity	85%
Pollution degree	2

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61326-1
2015/863/EU (RoHS III)

2012/19/EU (WEEE)

and satisfies the essential requirements of the following directives:

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9. Declaration of conformity ^{UK}CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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VELP Scientifica Srl