Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

OV5 Homogenizer

R20900010

**General Information / Informazioni Generali / Informations Générales / Información General /
Allgemeine Hinweise / 一般信息**

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de funcionamiento.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch
在使用本装置之前，请仔细阅读以下使用说明书。

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.
No tirar el aparato en los desechos urbanos, como exige la Directiva 2002/96/CE.
Dieses Gerät unterliegt der Richtlinie 2002/96/EG und darf nicht mit dem normalen Hausmüll entsorgt werden.
根据 EEC 指令 2002/96/CE，请不要将本设备作为城市垃圾处理。

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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IEC/EN 61010-1

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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
1. Introduction

To homogenize or disperse the sample can be carried out quickly and safely by using the dispersing tools inside the beaker containing the sample. The OV5 homogenizer is composed by a smooth motor (1) with integrated electronic speed control and permits to work in a very homogeneous way up to 30,000 rpm.

Smooth start, thanks to general switch (2), prevents unnecessary spills and an additional feature of automatic overload protection increases the life span of the motor. The OV5 homogenizer has only one shaft model (3) suitable for all the applications and configurable with different models of rotors and stators (4) in a very safe and comfortable way.

- | | | |
|-----------------|-------------------|----------------------------|
| 1. Motor | 2. General switch | 3. Dispersing tool (shaft) |
| 4. Rotor/Stator | 5. Strap clamp | 6. Homogenizer support rod |
| 7. Bosshead | 8. H-stand | 9. Speed controller |

2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- OV5 Homogenizer
- Instruction manual

2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

2.2 Start-up

Screw the support rod into the flange of the drive motor and tighten it. For correct use the OV5 must be secured to a stand with a bosshead. For safety reason the container with the sample should be fixed with a strap clamp.

Insert the dispersing tool into the drive flange by pressing upwards and turning clockwise until it locks into place. To remove the dispersing tool from the flange, hold the dispersing tool and turn anticlockwise and remove the dispersing tool downwards. The distance between the bottom of the vessel and the shaft should not be less than 10 mm. Check that the knob of the speed control is in position "A". Use the general switch on the front of the unit to turn the unit on (position "I").

3. Operating controls

With the unit switched on by the general switch, use the knob on the top of the unit to select the desired speed. Thanks to the different positions, identified by numbers, it is possible to select the following speed:

1 = 10,000 rpm 2 = 14,000 rpm 3 = 18,000 rpm 4 = 22,000 rpm 5 = 26,000 rpm max = 30,000 rpm

For reaching the best effects, the rotor speed of OV5 should be adapted to the dispersion problem.

NOTE: The maximum rotor speed depends on the kind of sample and its viscosity. The speed is reduced by its viscosity in order to protect the motor from overload. If the viscosity is too high, the motor will stop automatically to prevent damage to the equipment. The electronic control system prevents damage to the motor.

NOTE: The instrument cannot be used without put the dispersion tool inside a sample.

NOTE: The OV5 is designed for continues operation, however normally the ultimate fineness will be reached within few minutes. Longer application times bring no fineness improvement, but introduce unnecessary heat into the medium and could even cause the damage of the dispersing tool.

3.1 Overload protection

When the motor overloads and the power supply cuts off, switch off the general switch, disconnect the dispersing tool and re-start the instrument as fast as possible, without the load in order to accelerate cooling by the built-in fan.

3.2 Disassembling the dispersing tool

1. Hold the dispersing tool with a cloth (to avoid cutting of your hand due to sharp edges at the stator). At the same time turn the shaft tube clockwise. Screw the shaft tube of the stator and pull it over the axle and put it aside.
2. Now hold the rotor with a cloth and turn the axle anticlockwise. Slowly screw the axle from the rotor and put it aside.

3. Remove the PTFE bearing from the axle. After this disassembly you should have the following parts: shaft tube, axle, rotor, stator, and PTFE sheat.

To assemble the dispersing tool please follow above 3 steps in reverse order. (3, 2 and 1).

NOTE: During the assembling of the dispersing tool, make sure that the PTFE bearing is between the flange and the nut to avoid scratching and rubbing marks on the flange. Constantly control the PTFE bearing.

The optimal immersion depth of the dispersing shaft is approximately 2/3 below the liquid surface and 1/3 above the bottom of the beaker. When inserting the dispersing shaft slightly slanted (approximately at an angle 15 degrees), it improves the efficiency of the disperser. The dispersing shaft should not be immersed more than 30mm below the flange.

4. End of work session

When the work session is completed, should the unit remain connected to the power supply, switch it off using the general switch on the front.

4.1 Cleaning the dispersing tool

Immediately after finishing use the OV5 homogenizer, the dispersing element must be cleaned so that substance residues do not stick to the threads. In fact small bacterial cultures can form in the small grooves, and this could contaminate the samples that will be treated subsequently. For this reason, and to clean the dispersing tool effectively, it should be run in a solvent that dissolves the residue, but does not harm the Teflon bearing and the steel (i.e. n-hexane). Then the dispersing tool must be disassembled. Check the status of PTFE bearing and replace if necessary. Proper care and cleaning of the equipment will ensure a longer and better use of the equipment.

4.2 Sterilization

There are three different procedures of sterilization:

- Chemical process: germicidal solutions (formalin, phenol, alcohol, etc.) can disinfect in most cases. However, residues of the disinfectant must subsequently be removed with sterilized water.
- Humid heat: this means sterilizing with steam at a pressure of 2 bar and a temperature of 120 °C.
- Hot air: hot air sterilization is normally carried out in a range of temperature between 160 °C – 190 °C.

5. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

5.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

6. Technical data

Power	500 W
Dimensions (WxHxD)	70x255x70 mm
Weight	1.3 Kg
Speed range	10,000 ÷ 30,000 rpm
Operating mode	Continuous
Construction material	Technopolymer
Environmental ambient temperature	5 - 40 °C
Max humidity	85%

7. Accessories / Spare parts

A00000045	H-stand with strap clamp, bosshead clamp	A00001301	Double clamp
A00001300	Support rod and base	A00000044	Strap clamp
A00000279	Stator rotor release kit	10002666	PTFE bearing

The OV5 homogenizer is normally fixed to either H-stand code A00000045 or to the support rod and base code A00001300. It is always important to safely fix the unit even when not using the original stand for it supplied by VELP Scientifica, as centre of gravity is relatively high. The H-stand code A00000045 consists of a double rod, two bossheads, a safety ring and a strap clamp to secure the vessel containing the sample. Make sure that the motor drive is fixed safely and cannot move downwards at any time as this may cause damage to the instrument, to the user or to the sample.

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MULTI

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-2-051

EN 61326-1

and satisfies the essential requirements of the following directives:

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9. Declaration of conformity ^{UK}CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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