

Quantitative real-time PCR Standard *Brachyspira hyodysenteriae*


Product Name: Quantitative real-time PCR Standard *Brachyspira hyodysenteriae*

Order No.: PC400

Unit: 25 µl


Description: The Quantitative real-time PCR Standard *Brachyspira hyodysenteriae* is a quantitative real-time PCR standard for use with BactoReal® Kit *Brachyspira hyodysenteriae* (order no. DVEB01111, DVEB01113, DVEB01151, or DVEB01153). The standard contains a fragment of the *Brachyspira hyodysenteriae* DNA at a concentration of 10⁷ copies/µl. The standard curve has been determined using ABI PRISM® 7500 (Fast).

Use: Add 1 µl of a 1:10 dilution series to 20 µl PCR reaction. Please refer to the corresponding BactoReal® Kit *Brachyspira hyodysenteriae* manual for further instructions.



ABI Prism® 7500:

1:10 serial dilution of *Brachyspira hyodysenteriae* standard (FAM channel, 530 nm)



ABI Prism® 7500:

Standard curve generated with *Brachyspira hyodysenteriae* standard (FAM channel, 530 nm)

For veterinary use only

For further information on our products please visit our homepage (www.ingenetix.com)