







# Oocyte Collection Needle

1. Company profile		Update : November 2016 Issue date : August 2015
1.1	<b>Name :</b> Laboratoire CCD	
1.2	<b>Address :</b> 48, rue des Petites Ecuries 75010 Paris FRANCE	<b>Tel :</b> + 33 (0)1 44 95 14 95 <b>E-mail :</b> contact@ccdlab.com <b>Site internet :</b> www.ccd-lab.com
1.3	<b>Medical Device Vigilance contact :</b> Emilie FEZELOT, Q/A manager	<b>Tel :</b> +33 (0)3 44 26 89 56 <b>E-mail :</b> efezelot@ccdlab.com

2. Product information																							
2.1	<b>Common designation :</b> Oocyte collection needle																						
2.2	<b>Trade name :</b> Oocyte collection needle																						
2.3	<b>Medical device classification :</b> Class IIa <b>Relevant EU Directive :</b> 93/42/CEE <b>Certified Body identification number :</b> CE 0459 <b>Manufacturer :</b> Laboratoire CCD																						
2.4	<b>Product features :</b> Simple flow, triple bevel needle, echogenic on the last 2 cm, then on 5 mm as of the 3 <sup>rd</sup> distal cm. The needle is fitted with a Luer cone to connect a syringe. An embossed arrow on the base indicates the bevel orientation.  2 models: <table border="1" data-bbox="239 1317 1449 1590"> <thead> <tr> <th>Designation</th><th>Oocyte Collection Needle 30 cm</th><th>Oocyte Collection Needle 35 cm</th></tr> </thead> <tbody> <tr> <td>Product code</td><td>1403001</td><td>1311101</td></tr> <tr> <td>Usable length</td><td>30 cm</td><td>35 cm</td></tr> <tr> <td>Outer diameter</td><td>1.5 mm – G 17</td><td>1.6 mm – G 16</td></tr> <tr> <td>Inner diameter</td><td>1.2 mm</td><td>1.3 mm</td></tr> <tr> <td>Echogenic zone</td><td>Ground surface on 2 cm and 5 mm</td><td>Ground surface on 2 cm and 5 mm</td></tr> <tr> <td>Connection</td><td>Luer-lock</td><td>Luer-lock</td></tr> </tbody> </table>		Designation	Oocyte Collection Needle 30 cm	Oocyte Collection Needle 35 cm	Product code	1403001	1311101	Usable length	30 cm	35 cm	Outer diameter	1.5 mm – G 17	1.6 mm – G 16	Inner diameter	1.2 mm	1.3 mm	Echogenic zone	Ground surface on 2 cm and 5 mm	Ground surface on 2 cm and 5 mm	Connection	Luer-lock	Luer-lock
Designation	Oocyte Collection Needle 30 cm	Oocyte Collection Needle 35 cm																					
Product code	1403001	1311101																					
Usable length	30 cm	35 cm																					
Outer diameter	1.5 mm – G 17	1.6 mm – G 16																					
Inner diameter	1.2 mm	1.3 mm																					
Echogenic zone	Ground surface on 2 cm and 5 mm	Ground surface on 2 cm and 5 mm																					
Connection	Luer-lock	Luer-lock																					
2.5	<b>Catalogue reference :</b> <b>PRODUCT CODE :</b> 1403001      Oocyte Collection Needle - 30 cm G17 1311101      Oocyte Collection Needle - 35 cm G16 <b>Packaging :</b> Unit packaging <b>Ordering unit :</b> Unit <b>Minimum order :</b> box of 10 units  <b>Labelling :</b>																						

	<div> <div> <b>10</b> REF.1403001 </div> <div> <b>Alguille Ponction Ovocytaire 30cm G17</b>  <b>Oocyte Collection Needle 30cm / 1.50mm</b> </div> <div>  0459   </div> <div>   11111111  11111111 </div> <div> USA caution : Only for sale by or on the order of a physician </div> </div> <div> Laboratoire CCD  48 Rue des Petites Ecuries  75010 Paris - France  Tel : 33 1 44 96 14 95 - Fax : 33 1 44 96 14 99 </div>
<b>2.6</b>	<b>Components and accessories :</b> Steel <b>Active substances :</b> Latex free Phtalates free Contains no component of biological or animal origin
<b>2.7</b>	<b>Indications :</b> <ul style="list-style-type: none"> <li>▪ Oocyte collection</li> <li>▪ Ovarian cyst aspiration</li> </ul>

### 3. Sterilization method

	<b>Sterile medical device :</b> Yes <b>Sterilization method :</b> Irradiation
--	--

### 4. Storage and shelf life

	<b>Current conservation and storage conditions :</b> Store in a dry place, away from light <b>Specific precautions :</b> Do not use if unit packing is damaged <b>Product shelf life :</b> 5 years
--	---

### 5. Operational safety

<b>5.1</b>	<b>Technical safety :</b> Single use medical device. Do not re-sterilize. Do not re-use. Read carefully the instructions before use.
<b>5.2</b>	<b>Biological safety :</b> N/A

## 6. Instructions for use

<b>6.1</b>	<b>Operating instructions :</b> <ul style="list-style-type: none"> <li>Check the integrity of the device (bevel sharpness, connection part).</li> <li>Rinse the needle and syringe set with a biological medium suited to fine needle aspiration.</li> <li>Depending on the approach chosen, perform the fine needle aspiration using ultrasound guidance or coeloscopy.</li> </ul>
<b>6.2</b>	<b>Indications :</b> <ul style="list-style-type: none"> <li>Oocyte collection</li> <li>Ovarian cyst aspiration</li> </ul>
<b>6.3</b>	<b>Precautions for use :</b> <ul style="list-style-type: none"> <li>Check that the packaging is not damaged</li> <li>Before use, check the integrity of the medical device</li> <li>Prior to use, rinse the device with appropriate medium according to state of the art practices</li> </ul>
<b>6.4</b>	<b>Contra-indications :</b> <ul style="list-style-type: none"> <li>None</li> </ul>

## 7. Additional information on the product

	N/A
--	-----

## 8. Appendices

	N/A
--	-----

## 9. Pictures

	 <p>Classic needle</p>
--	--