

# **Blood DNA Isolation Mini Kit**

Norgen's Blood DNA Isolation Mini Kit is designed for the rapid preparation of genomic DNA from up to 200 mL of whole blood. Purification is based on spin column chromatography as the separation matrix. Norgen's column binds DNA under optimized salt concentrations and releases the bound DNA under low salt and slightly alkali conditions. The purified genomic DNA is fully digestible with all restriction enzymes tested, and is completely compatible with downstream applications including real-time PCR and southern blot analysis.

Norgen's Blood DNA Isolation Mini Kit allows for the isolation of genomic DNA from the blood of various



species, including humans. The genomic DNA is preferentially purified from other cellular proteinaceous components. Typical yields of genomic DNA will vary depending on the cell density of the blood sample. Preparation time for a single sample is less than 30 minutes, and each kit contains sufficient materials for 50 preparations.

Kit Specifications				
Maximum Blood Input	200 µL	Average Yield (from 200 µL)	4 - 12 μg*	
Column Binding Capacity	> 50 µg	Time to Complete 10 Purifications	30 minutes	

\* Yeild will vary depending on the type of blood processed

Blood Genomic DNA Isolation Kit Benefits				
Fast and easy processing	Rapid spin-column format allows for the processing of multiple samples in 30 minutes.			
High quality DNA	Isolated DNA is of the highest quality and free from RNA contamination.			
Recovered genomic DNA is suitable for downstream applications	Purified genomic DNA is fully compatible with restriction enzyme digestions, PCR, and Southern Blot analysis.			



## **Blood DNA Isolation Mini Kit**




Figure 1. High Yields of Genomic DNA Isolated from 20µL to 200µL of Whole Blood. Genomic DNA was isolated from 20, 50, 100 and 200 µL of whole blood using Norgen's Blood DNA Isolation Mini Kit and a leading competitor's kit. Following isolation, 15 µL from each 200 µL elution was loaded on 1% TAE agarose gel. Norgen's Blood Genomic DNA Isolation Mini Kit demonstrated a better DNA yield than the leading competitor's kit. Lane M: Norgen's UltraRanger 1kb DNA Ladder.

### **Blood DNA Isolation Mini Kit Contents**

- 1. Lysis Buffer B
- 2. Solution WN
- 3. Wash Solution A
- 4. Elution Buffer B
- 5. Proteinase K
- 6. Spin Columns
- 7. Collection Tubes
- 8. Elution Tubes
- 9. Product Insert

### **Storage Conditions**

All solutions should be kept tightly sealed and stored at room temperature. These reagents should remain stable for at least 1 year in their unopened containers. The kit contains a ready-to-use Proteinase K solution, which is dissolved in a specially prepared storage buffer. The Proteinase K is stable for up to 1 year after delivery when stored at room temperature. To prolong the lifetime of Proteinase K, storage at 2–8°C is recommended.




Figure 2. Purified DNA Can be Amplified in a Realtime PCR (TagMan) Reaction. Genomic DNA was isolated from 20, 50, 100 and 200 µL of whole human blood using Norgen's Blood DNA Isolation Mini Kit (Blue) and a leading competitor's kit (Red). Nine µL of the DNA from each 200 µL elution was used in a real-time PCR reaction (total reaction volume of 20 µL) with GAPDH TaqMan probe and primers. The real-time PCR was successful in amplifying the GAPDH gene, indicating that the DNA is of a high quality and can be used in sensitive downstream applications. Furthermore, Norgen-isolated DNA was amplified with a lower Ct value from all DNA isolated from the different blood input volumes, indicating the higher yield and purity of DNA isolated using Norgen's kit.

Customer-Supplied Reagents and Equipment

- Benchtop microcentrifuge
- 1.5 mL microcentrifuge tubes
- 55°C water bath or heating block
- 96 100% ethanol
- Vortex

#### **Shipping Conditions**

The Blood Genomic DNA Isolation Mini Kit is shipped at room temperature.

(	Cat # Description		Quantity
	46300	Blood DNA Isolation Mini Kit	50 preps

3430 Schmon Parkway, Thorold, ON, Canada L2V 4Y6 www.norgenbiotek.com

Phone: 905-227-8848 Fax: 905-227-1061 North American Toll Free: 1-866-667-4362

© Norgen Biotek Corp. 2014