PRO-FLEX® PIVOT

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RÉFÉRENCES

1. Hunter, David J., Deborah Schofield, and Emily Callander. "The individual and socioeconomic impact of osteoarthritis." *Nature Reviews Rheumatology* 10.7 (2014): 437-441.
2. Nho, Shane J., Steven M. Kymes, John J. Callaghan, and David T. Felson 2013, The Burden of Hip Osteoarthritis in the United States: Epidemiologic and Economic Considerations. *The Journal of the American Academy of Orthopaedic Surgeons* 21 Suppl 1: S1–6.
3. HCUP: Healthcare Cost and Utilization Project, June 2015. Agency for Healthcare Research and Quality, Rockville, MD, U.S. Department of Health and Human Services. <http://hcupnet.ahrq.gov/Hcupnet.jsp> Agency for Healthcare Research and Quality
4. Ziegler-Graham, Kathryn, et al. "Estimating the prevalence of limb loss in the United States: 2005 to 2050." *Archives of physical medicine and rehabilitation* 89.3 (2008): 422-429.
5. Miller, William C., et al. "The influence of falling, fear of falling, and balance confidence on prosthetic mobility and social activity among individuals with a lower extremity amputation." *Archives of physical medicine and rehabilitation* 82.9 (2001): 1238-1244.
6. Burger, Helena, C. R. T. Marincek, and Eli Isakov. "Mobility of persons after traumatic lower limb amputation." *Disability & Rehabilitation* 19.7 (1997): 272-277.
7. Geertzen JH, Bosmans JC, Van der Schans CP. Claimed walking distance of lower limb amputees. *Disabil Rehabil* 2005;27:101-4.
8. Nolan L, Wit A, Dudzinski K, Lees A, Lake M, Wychowanski M. Adjustments in gait symmetry with walking speed in trans-femoral and trans-tibial amputees. *Gait Posture*.2003;17(2):142–51 prosthetic knee. *Arch Phys Med Rehabil* 2007;88:207-17.
9. Burke MJ, Roman V, Wright V. Bone and joint changes in lower limb amputees. *Ann Rheum Dis*. 1978;37(3): 252–54.
10. Struyf, Pieter A., et al. "The prevalence of osteoarthritis of the intact hip and knee among traumatic leg amputees." *Archives of physical medicine and rehabilitation* 90.3 (2009): 440-446.
11. Gailey R, Allen K, Castles J, Kucharik J, Roeder M. Review of secondary physical conditions associated with lower-limb amputation and long-term prosthesis use. *J Rehabil Res Dev* 2008;45(1):15–29.
12. Hurley GR, Mckenney R, Robinson M, Zdravac M, Pierrynowski, MR. The role of the contralateral limb in below knee amputee gait. *Prosthet Orthot Int*. 1990;14(1):33–42.
13. Struyf, Pieter A., et al. "The prevalence of osteoarthritis of the intact hip and knee among traumatic leg amputees." *Archives of physical medicine and rehabilitation* 90.3 (2009): 440-446.
14. Kulkarni J, Adams J, Thomas E, Silman A. Association between amputation, arthritis and osteopenia in British male war veterans with major lower limb amputations. *Clin. Rehabil.*, 12 (4) (1998), pp. 348–353
15. Centers for Disease Control and Prevention (CDC). Prevalence and impact of chronic joint symptoms—seven states, 1996. *MMWR Morb. Mortal. Wkly Rep*. 47, 345–351 (1998).
16. Dunlop, D. D., Manheim, L. M., Song, J. & Chang, R. W. Arthritis prevalence and activity limitations in older adults. *Arthritis Rheum*. 44, 212–221 (2001).
17. Hunter, D. J. Lower extremity osteoarthritis management needs a paradigm shift. *Br. J. Sports Med*. 45, 283–288 (2011).
18. Hootman, J. M. & Helmick, C. G. Projections of US prevalence of arthritis and associated activity limitations. *Arthritis Rheum*. 54, 226–229 (2006).
19. Perruccio, A. V., Power, J. D. & Badley, E. M. Revisiting arthritis prevalence projections—it's more than just the aging of the population. *J. Rheumatol*. 33, 1856–1862 (2006).
20. Nolan L, Wit A, Dudzinski K, Lees A, Lake M, Wychowanski M. Adjustments in gait symmetry with walking speed in trans-femoral and trans-tibial amputees. *Gait Posture*.2003;17(2):142–51 prosthetic knee. *Arch Phys Med Rehabil* 2007;88:207-17.
21. Hungerford D, Cockin J. Fate of the retained lower limb joints in World War II amputees. *J. Bone Jt. Surg.*, 57 (1975), p. 111
22. Melzer I, Yekutieli M, Sukenik S. Comparative study of osteoarthritis of the contralateral knee joint of male amputees who do and do not play volleyball. *J. Rheumatol.*, 28 (1) (2001), pp. 169–172
23. Norvell DC, Czerniecki JM, Reiber GE, Maynard C, Pecoraro JA, Weiss NS. The prevalence of knee pain and symptomatic knee osteoarthritis among veteran traumatic amputees and nonamputees. *Arch Phys Med Rehabil* 2005;86(3):487–93.
24. Lemaire ED, Fisher FR. Osteoarthritis elderly amputee gait. *Arch Phys Med Rehabil* 1994;75(10):1094–9.
25. Breakley J. Gait of unilateral trans-tibial amputees. *Orthot Prosthet*. 1976;30:17–24.
26. Murray MP, Mollinger LA, Sepic SB, Gardner GM, Linder MT. Gait patterns in above-knee amputee patients: Hydraulic swing control vs constant-friction knee components. *Arch Phys Med Rehabil*. 1983;64(8):339–45.
27. Engsborg JR, Lee AG, Tedford KG, Harder JA. Normative ground reaction force data for able-bodied and below knee amputee children during walking. *J Pediatr Orthop*. 1993;13(2):169–73.
28. Zernicke RF, Hoy MG, Whiting WC. Ground reaction forces and center of pressure patterns in the gait of children with amputation: Preliminary report. *Arch Phys Med Rehabil*. 1985;66(11):736–41.
29. Menard MR, McBride ME, Sanderson DJ, Murray D. Comparative biomechanical analysis of energy-storing prosthetic feet. *Arch Phys Med Rehabil*. 1992;73(5):451–58.
30. Schneider K, Hart T, Zernicke RF, Setoguchi Y, Oppenheim W. Dynamics of below-knee amputee child gait: SACH foot versus Flex foot. *J Biomech*. 1993;26(10): 1191–1204.
31. Suzuki K. Force plate study on the artificial limb gait. *J Jpn Orthop Assoc*. 1972;46:503–16.
32. Engsborg JR, Lee AG, Patterson JL, Harder JA. External loading comparisons between able-bodied and below knee amputee children during walking. *Arch Phys Med Rehabil*. 1991;72(9):657–61
33. E. Levy, A. Ferme, D. Perocheau, and I. Bono, "Socioeconomic costs of osteoarthritis in France," *Revue du Rhumatisme*, vol. 60, no. 6, pp. 63S–67S, 1993.
34. Chen, A., et al. "The global economic cost of osteoarthritis: how the UK compares." *Arthritis* 2012 (2012).
- 34a. Centers for Disease Control and Prevention (CDC). Racial disparities in total knee replacement among Medicare enrollees— United States, 2000-2006. *MMWR Morb Mortal Wkly Rep* 2009;58(6):1338.
35. Murray, C. J. et al. Disability-adjusted life years (DALYs) for 291 diseases and injuries in 21 regions, 1990–2010: a systematic analysis for the Global Burden of Disease Study 2010. *Lancet* 380, 2197–2223 (2013).
36. Hunter, David J., Deborah Schofield, and Emily Callander. "The individual and socioeconomic impact of osteoarthritis." *Nature Reviews Rheumatology* 10.7 (2014): 437-441.
37. Snyder, R.D., et al., The effect of five prosthetic feet on the gait and loading of the sound limb in dysvascular below-knee amputees. *J Rehabil Res Dev*, 1995. 32(4): p. 309-15.
38. Lehmann JF, Price R, Boswell-Bessette S, Dralle A, Questad K. Comprehensive analysis of dynamic elastic response feet: Seattle Ankle/Lite Foot versus SACH foot. *Archives of Physical Medicine & Rehabilitation* 1993;74(8):853-61.
39. Heitzmann DWW. et al; Evaluation of a novel prosthetic foot while walking on level ground, ascending and descending a ramp; *Gait & Posture* 42 (2015): S94-S95. Abstract, Oral Presentation at the ESMAC 24th annual Meeting Heidelberg, Germany, September 10-12, 2015; E-mail: daniel.heitzmann@med.uni-heidelberg.de
40. Morgenroth, David C., et al. "The effect of prosthetic foot push-off on mechanical loading associated with knee osteoarthritis in lower extremity amputees." *Gait & posture* 34.4 (2011): 502-507.
41. Heitzmann, D. W. W., et al. "Evaluation of a novel prosthetic foot while walking on level ground, ascending and descending a ramp." *Gait & Posture* 42 (2015): S94-S95.
42. Powers, Christopher M., et al. "Influence of prosthetic foot design on sound limb loading in adults with unilateral below-knee amputations." *Archives of physical medicine and rehabilitation* 75.7 (1994): 825-829.
43. Segal, Ava D., et al. "The effects of a controlled energy storage and return prototype prosthetic foot on transtibial amputee ambulation." *Human movement science* 31.4 (2012): 918-931.
44. Kuo, Arthur D. "The six determinants of gait and the inverted pendulum analogy: A dynamic walking perspective." *Human movement science* 26.4 (2007): 617-656.
45. Kuo, Arthur D., J. Maxwell Donelan, and Andy Ruina. "Energetic consequences of walking like an inverted pendulum: step-to-step transitions." *Exercise and sport sciences reviews* 33.2 (2005): 88-97.
46. Osteoarthritis kneebracing – A health economic evaluation – USA, 2012. On file at Össur