
PROPRIO FOOT®

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PROPRIO FOOT®

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1. Miller, William C., Mark Speechley, and Barry Deathe. "The prevalence and risk factors of falling and fear of falling among lower extremity amputees." *Archives of physical medicine and rehabilitation* 82.8 (2001): 1031-1037.
2. Pauley T, Devlin M, Heslin K. Falls sustained during inpatient rehabilitation after lower limb amputation: prevalence and predictors. *Am J Phys Med Rehabil.* 2006; 85:521-532; quiz 533-535.
3. Miller, William C., et al. "The influence of falling, fear of falling, and balance confidence on prosthetic mobility and social activity among individuals with a lower extremity amputation." *Archives of physical medicine and rehabilitation* 82.9 (2001): 1238-1244.
4. Dite, Wayne, Helen J. Connor, and Heather C. Curtis. "Clinical identification of multiple fall risk early after unilateral trans-tibial amputation." *Archives of physical medicine and rehabilitation* 88.1 (2007): 109-114.
5. Asano, Miho, et al. "Predictors of quality of life among individuals who have a lower limb amputation." *Prosthetics and orthotics international* 32.2 (2008): 231-243.
6. A. A. Bohl, P. A. Fishman, M. A. Ciol, B. Williams, J. LoGerfo, and E. A. Phelan, "A Longitudinal Analysis of Total 3-Year Healthcare Costs for Older Adults Who Experience a Fall Requiring Medical Care: Longitudinal costs of older adult fallers," *Journal of the American Geriatrics Society*, vol. 58, no. 5, pp. 853-860, May 2010.
7. Kaufman, K. "Risk factors and costs associated with accidental falls among adults with above-knee amputations: a population-based study," *American Orthotic and Prosthetic Association* 2016. (Mayo Clinic). <http://www.aopanet.org/resources/research/>
8. E. R. Burns, J. A. Stevens, and R. Lee, "The direct costs of fatal and non-fatal falls among older adults — United States," *Journal of Safety Research*, vol. 58, pp. 99-103, Sep. 2016.
9. B. Mundell, H. Maradit Kremers, S. Visscher, K. Hoppe, and K. Kaufman, "Direct medical costs of accidental falls for adults with transfemoral amputations," *Prosthet Orthot Int*, p. 0309364617704804, Jun. 2017.
10. Rosenblatt, Noah J., et al. "Active dorsiflexing prostheses may reduce trip-related fall risk in people with transtibial amputation." *J Rehabil Res Dev* 51.8 (2014): 1229-1242.
11. Ludviksdottir A, Gruben K, Gunnsteinsson K, Ingvarsson Th, Nicholls M. Effects on user mobility and safety when changing from a carbon fiber prosthetic foot to a bionic prosthetic foot. Presented at Orthopadie&Reha-Technik Congress, Leipzig, May 2012.
12. Esquenazi, Alberto, and Robert DiGiacomo. "Rehabilitation after amputation." *Journal of the American Podiatric Medical Association* 91.1 (2001): 13-22.
13. Paysant, Jean, et al. "Influence of terrain on metabolic and temporal gait characteristics of unilateral trans-tibial amputees." *Journal of rehabilitation research and development* 43.2 (2006): 153.
14. Agrawal, Vibhor, et al. "Symmetry in external work (SEW): A novel method of quantifying gait differences between prosthetic feet." *Prosthetics and orthotics international* 33.2 (2009): 148-156.
15. Delussu, Anna Sofia, et al. Assessment of the effects of carbon fiber and bionic foot during overground and treadmill walking in trans-tibial amputees. *Gait & posture*, 2013, 38. Jg., Nr. 4, S. 876-882.
16. Alimusaj M, Fradet L, Braatz F, Gerner HJ, Wolf SI. Kinematics and kinetics with an adaptive ankle foot system during stair ambulation of trans-tibial amputees. *Gait & Posture*. 2009; 30:3:356-363.
17. Fradet L, Alimusaj M, Braatz F, Wolf SI. Biomechanical analysis of ramp ambulation of trans-tibial amputees with an adaptive ankle foot system. *Gait & Posture*. 2010; 32(2): 191 - 198.
18. Wolf, S.I, Alimusaj M, Fradet L, Siegel J, Braatz F. Pressure characteristics at the stump/socket interface in trans-tibial amputees using an adaptive prosthetic foot. *Clinical Biomechanics*. 2009; 24(10), 860-5.

P-962031ES