## Genium X3 3B5-3 / 3B5-3=ST

<b>DE</b> Gebrauchsanweisung (Benutzer) .....	3
<b>EN</b> Instructions for use (user) .....	41
<b>FR</b> Instructions d'utilisation (utilisateur) .....	77
<b>IT</b> Istruzioni per l'uso (utilizzatore) .....	117
<b>ES</b> Instrucciones de uso (usuario) .....	153
<b>PT</b> Instruções de utilização (utilizador) .....	191
<b>NL</b> Gebruiksaanwijzing (gebruiker) .....	229
<b>HR</b> Upute za uporabu (korisnik) .....	267
<b>TR</b> Kullanım kılavuzu (kullanıcı) .....	303
<b>JA</b> 取扱説明書(装着者用) .....	339

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<b>1</b>	<b>Foreword .....</b>	<b>43</b>
<b>2</b>	<b>Product description .....</b>	<b>43</b>
2.1	Design .....	43
2.2	Function .....	43
<b>3</b>	<b>Application .....</b>	<b>44</b>
3.1	Indications for use .....	44
3.2	Conditions of use .....	44
3.3	Indications .....	44
3.4	Contraindications .....	45
3.4.1	Absolute Contraindications .....	45
3.5	Qualification .....	45
<b>4</b>	<b>Safety .....</b>	<b>45</b>
4.1	Explanation of warning symbols .....	45
4.2	Structure of the safety instructions .....	45
4.3	General safety instructions .....	45
4.4	Information on the Power Supply/Battery Charging .....	48
4.5	Battery charger information .....	48
4.6	Information on Proximity to Certain Areas .....	49
4.7	Information on Use .....	50
4.8	Notes on the safety modes .....	52
4.9	Instructions for use with an osseointegrated implant system .....	53
<b>5</b>	<b>Scope of Delivery and Accessories .....</b>	<b>53</b>
5.1	Scope of delivery .....	53
5.2	Accessories .....	53
<b>6</b>	<b>Charging the battery .....</b>	<b>54</b>
6.1	Connecting the power supply and battery charger .....	54
6.2	Charging the prosthesis battery .....	55
6.3	Display of the current charge level .....	55
6.3.1	Display of battery charge level without additional devices .....	55
6.3.2	Display of the current charge level during the charging process .....	56
<b>7</b>	<b>Use .....</b>	<b>56</b>
7.1	Standing .....	56
7.1.1	Stance function .....	56
7.2	Walking .....	57
7.3	Running short distances ("walk-to-run" function) .....	57
7.4	Sitting down .....	58
7.5	Sitting .....	58
7.5.1	Sitting function .....	58
7.6	Standing up .....	58
7.7	Walking up stairs step-over-step .....	58
7.8	Overcoming obstacles .....	59
7.9	Walking down stairs .....	60

7.10	Walking down a ramp .....	60
7.11	Bluetooth.....	60
7.11.1	Establishing the Bluetooth connection.....	60
7.12	Mute mode (silent mode) .....	60
7.13	Deep sleep mode .....	60
<b>8</b>	<b>MyModes .....</b>	<b>61</b>
8.1	Running function as configured MyMode.....	61
8.2	Switching MyModes using motion patterns.....	61
8.3	Switching from a MyMode back to basic mode.....	62
<b>9</b>	<b>Additional operating states (modes) .....</b>	<b>63</b>
9.1	Empty battery mode .....	63
9.2	Mode for charging the prosthesis .....	63
9.3	Safety mode.....	63
9.4	Overheating mode .....	64
<b>10</b>	<b>Storage and bleeding .....</b>	<b>64</b>
<b>11</b>	<b>Cleaning .....</b>	<b>64</b>
<b>12</b>	<b>Maintenance .....</b>	<b>64</b>
<b>13</b>	<b>Legal information .....</b>	<b>65</b>
13.1	Liability .....	65
13.2	Trademarks .....	65
13.3	CE conformity .....	65
13.4	Local Legal Information .....	65
<b>14</b>	<b>Technical data .....</b>	<b>66</b>
<b>15</b>	<b>Appendices .....</b>	<b>68</b>
15.1	Symbols Used.....	68
15.2	Operating states/error signals.....	69
15.2.1	Signals for operating states .....	69
15.2.2	Warnings/error signals.....	70
15.2.3	Status signals .....	72
15.3	Directives and manufacturer's declaration .....	73
15.3.1	Electromagnetic environment .....	73

# 1 Foreword

## INFORMATION

Date of last update: 2025-03-20

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The product "3B5-3, 3B5-3=ST Genium X3" is referred to below as the product/prosthesis/knee joint/component.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Design

The product consists of the following components:

1. Proximal pyramid
2. Optional flexion stop
3. Battery
4. Hydraulic unit
5. LED (blue) as indicator for the Bluetooth connection
6. Receiver of the inductive charging unit

### 2.2 Function

This product features microprocessor control of the stance and swing phase.

The microprocessor uses the measurements of an integrated sensor system as a basis to control a hydraulic unit that influences the damping behaviour of the product.

These sensor data are updated and evaluated 100 times per second. As a result, the behaviour of the product is adapted to the current motion situation (gait phase) dynamically and in real time.

The product can be customised to your individual requirements using an adjustment software/adjustment app.

The product features MyModes for special motion types (e.g. cycling ...). These are pre-configured by the O&P professional and can be activated using special movement patterns (see page 61) and the "Cockpit" user app (see the section "Accessories" (see page 53)).

In case of a product malfunction, safety mode makes restricted operation possible. Resistance parameters that are predefined by the product are configured for this purpose (see page 63).

Empty battery mode permits safe walking when the battery is drained. Resistance parameters that are predefined by the product are configured for this purpose (see page 63).

### **The microprocessor-controlled hydraulic unit offers the following advantages**

- Approximation of the physiological gait pattern
- Stability while standing and walking
- Adaptation of product characteristics to various surfaces, inclines, gait situations and walking speeds

### **Essential performance of the product**

- Stability in the stance phase
- Adjustable swing phase extension resistance

## **3 Application**

### **3.1 Indications for use**

The product is to be used **solely** for lower limb exoprosthesis fittings.

### **3.2 Conditions of use**

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.).

Permissible ambient conditions are described in the technical data (see page 66).

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

The MOBIS classification describes the mobility grade and body weight, and makes it easy to identify compatible components.

The product is recommended for mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands). Approved for a max. body weight of **150 kg (330 lbs)**.

### **3.3 Indications**

- For users with knee disarticulation, transfemoral amputation or hip disarticulation.
- For unilateral or bilateral amputation
- Dysmelia patients with residual limb characteristics corresponding to knee disarticulation, transfemoral amputation or hip disarticulation
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

## 3.4 Contraindications

### 3.4.1 Absolute Contraindications

- Body weight over 150 kg

### 3.5 Qualification

The product may be fitted only by qualified personnel authorised by Ottobock after completing the corresponding training.

If the product is to be connected to an osseointegrated implant system, the qualified personnel must also be authorised for the connection to the osseointegrated implant system.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

 <b>WARNING</b>
<b>The heading describes the source and/or the type of hazard</b>
The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:
> E.g.: Consequence 1 in the event of failure to observe the hazard
> E.g.: Consequence 2 in the event of failure to observe the hazard
▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

### 4.3 General safety instructions

 <b>WARNING</b>
<b>Non-observance of safety notices</b>
Personal injury/damage to the product due to using the product in certain situations.
▶ Observe the safety notices and the stated precautions in this accompanying document.

 <b>WARNING</b>
<b>Using the prosthesis while operating a vehicle</b>
Accident due to unexpected behaviour of the prosthesis because of changed damping behaviour.
▶ All users are required to observe their country's national and state driving laws when operating vehicles with a prosthesis. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre.
▶ Observe national legal regulations for retrofitting your vehicle in accordance with the type of fitting.
▶ The leg on which the prosthesis is worn may not be used to control the vehicle or its peripheral components (e.g. clutch pedal, brake pedal, gas pedal, etc.).

**⚠ WARNING**

**Use of damaged power supply unit, adapter plug or battery charger**

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

**⚠ CAUTION**

**Failure to observe warning/error signals**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see page 70) and corresponding change in damping settings must be observed.

**⚠ CAUTION**

**Failure to observe activated mute mode (silent mode)**

Falling due to unexpected product behaviour because of changed damping behaviour.

The following feedback signals are deactivated when mute mode is activated:

- > Long vibration signal if the hydraulic unit overheats.
- > Beep and vibration signal to confirm that the movement pattern has been recognised (switching to a MyMode/basic mode with movement pattern).
- > Beep and vibration signal to indicate successful switching to a MyMode/basic mode.
- > Beep and vibration signal upon successfully switching to deep sleep mode.
- ▶ Before activating mute mode, note that these feedback signals will be deactivated. For more information about mute mode, see the section "Mute mode" (see page 60).
- ▶ After switching to a MyMode/basic mode, verify the changed damping settings.
- ▶ Ensure that you stand securely during all switching processes.
- ▶ To deactivate mute mode, connect and then disconnect the battery charger.

**⚠ CAUTION**

**Independent manipulation of the product and the components**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by authorised, qualified Ottobock personnel (no replacement by the user).
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

**⚠ CAUTION**

**Mechanical stress on the product**

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

**⚠ CAUTION**

**Use of the product when battery charge level is too low**

Falling due to unexpected behaviour of the prosthesis because of changed damping behaviour.

- ▶ Check the current charge level before use and charge the prosthesis if required.
- ▶ Note that the operating time of the product may be reduced at low ambient temperatures or due to ageing of the battery.

**⚠ CAUTION**

**Risk of pinching in the joint flexion area**

Injuries due to pinching of body parts.

- ▶ Ensure that fingers/body parts or soft tissue of the residual limb are not in this area when bending the joint.

**⚠ CAUTION**

**Penetration of dirt and humidity into the product**

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ Ensure that no solid particles or foreign objects can penetrate into the product.
- ▶ The knee joint and the AXON tube adapter are waterproof, corrosion-resistant and protected against the penetration of jets of water. The knee joint and AXON tube adapter can be operated in fresh and salt water. Do not use the knee joint under extreme conditions like diving or jumping into water. The knee joint and AXON tube adapter are designed for underwater use (for the maximum duration and water depth, see the section "Technical data" (see page 66).
- ▶ After contact with water, hold the prosthesis with the sole of the foot facing up until the water has drained from the knee joint and AXON tube adapter.
- ▶ After using the knee joint in a salt water environment, remove the Protective Cover and rinse knee joint, AXON tube adapter and Protective Cover with fresh water. Dry the knee joint and components with a lint-free cloth and allow the components to fully air dry.
- ▶ Should the knee joint or AXON tube adapter come into contact with **solutions other than fresh or salt water, promptly** remove the Protective Cover and clean the **knee joint**. In order to do so, rinse knee joint, AXON tube adapter and Protective Cover with fresh water and let them dry.
- ▶ In case of a malfunction after drying, the knee joint and AXON tube adapter must be inspected by an authorised Ottobock Service Center.
- ▶ The knee joint and the AXON tube adapter are not protected against the penetration of steam.

**⚠ CAUTION**

**Use of the product without protector or with damaged protector**

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ If the protector has been removed, ensure prior to the next use of the product that the protector has been properly installed.
- ▶ Using the product with a damaged protector or without the protector is not permitted.
- ▶ Using the product with a foam cover is not possible, since the protector would have to be removed for this purpose.

### CAUTION

#### **Signs of wear and tear on the product components**

Falling due to damage or malfunction of the product.

- ▶ Regular service inspections (maintenance) are mandatory in the interest of your own safety and in order to maintain operating reliability and protect the warranty.

### CAUTION

#### **Use of unapproved accessories**

> Falling due to product malfunction as a result of reduced interference resistance.

- > Interference of other electronic devices due to increased emissions.
- ▶ Use the product only in combination with the accessories, signal converters and cables listed in the sections "Scope of delivery" (see page 53) and "Accessories" (see page 53).

### NOTICE

#### **Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product with a damp cloth only (fresh water).
- ▶ Only use fresh water with a temperature of no more than 65 °C for cleaning.
- ▶ If dirt cannot be removed, the product must be sent to an authorised Ottobock Service Center. Your contact is the O&P professional.

## 4.4 Information on the Power Supply/Battery Charging

### CAUTION

#### **Charging the product without taking it off**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ For safety reasons, take the product off before charging the battery.

### CAUTION

#### **Charging the product with damaged power supply unit/charger/charger cable**

Falling due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

### NOTICE

#### **Use of incorrect power supply unit/battery charger**

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

## 4.5 Battery charger information

### WARNING

#### **Storing/transporting the product near active implanted systems**

Interference with active implantable systems (e.g. pacemaker, defibrillator, etc.) due to the product's magnetic field.

- ▶ When storing/transporting the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.

- ▶ Make sure to observe any operating conditions and safety notices stipulated by the manufacturer of the implant.

#### NOTICE

##### **Penetration of dirt and humidity into the product**

Lack of proper charging functionality due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

#### NOTICE

##### **Mechanical stress on the power supply/battery charger**

Lack of proper charging functionality due to malfunction.

- ▶ Do not subject the power supply/battery charger to mechanical vibrations or impacts.
- ▶ Check the power supply/battery charger for visible damage before each use.

#### NOTICE

##### **Operating the power supply unit/charger outside of the permissible temperature range**

Lack of proper charging functionality due to malfunction.

- ▶ Only use the power supply unit/charger for charging within the allowable temperature range. The section "Technical data" contains information on the allowable temperature range (see page 66).

#### NOTICE

##### **Independent changes or modifications carried out to the battery charger**

Lack of proper charging functionality due to malfunction.

- ▶ Have any changes or modifications carried out only by Ottobock authorised, qualified personnel.

#### NOTICE

##### **Contact of the battery charger with magnetic data storage devices**

Wiping of the data storage device.

- ▶ Do not place the battery charger on credit cards, diskettes, audio or video cassettes.

#### NOTICE

##### **Improper care of the housing**

Damage to the casing through the use of acetone, white spirit or similar solvents.

- ▶ Only clean the housing with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).

## **4.6 Information on Proximity to Certain Areas**

#### **CAUTION**

##### **Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm to HF communication devices is recommended.

**⚠ CAUTION**

**Operating the product in very close proximity to other electronic devices**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

**⚠ CAUTION**

**Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for people (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations...).
- ▶ If this cannot be avoided, make sure to at least have a safeguard when walking or standing (e.g. a handrail or the support of another person).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected changes in the damping behaviour of the product.
- ▶ In general, monitor the product for unexpected changes in the damping behaviour when electronic or magnetic devices are in the immediate vicinity.

**⚠ CAUTION**

**Entering a room or area with strong magnetic fields (e.g. magnetic resonance tomographs, MRT (MRI) equipment...)**

- > Falling due to unexpected restriction of the product's range of motion caused by metallic objects adhering to the magnetised components.
- > Irreparable damage to the product due to the effect of strong magnetic fields.
- ▶ Take off the product before entering a room or area with strong magnetic fields and store the product outside this room or area.
- ▶ Damage to the product caused by exposure to strong magnetic fields cannot be repaired.

**⚠ CAUTION**

**Remaining in areas outside the allowable temperature range**

Falling due to malfunction or the breakage of load-bearing product components.

- ▶ Avoid remaining in areas with temperatures outside of the permissible range (see page 66).

## 4.7 Information on Use

**⚠ CAUTION**

**Walking up stairs**

Falling due to foot placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking up stairs and place most of the area of the sole of your foot on the stair surface.
- ▶ Particular caution is required when carrying children up stairs.

**⚠ CAUTION**

**Walking down stairs**

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking down stairs and roll over the edge of the step with the middle of the shoe.
- ▶ Observe the warning/error signals (see page 70).
- ▶ Be aware that resistance in the flexion and extension direction can change in case of warning and error signals.
- ▶ Particular caution is required when carrying children down the stairs.

**⚠ CAUTION**

**Overheating of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill)**

- > Falling due to unexpected behaviour of the product because of switching into overheating mode.
- > Burns due to touching overheated components.
- ▶ Be sure to pay attention when pulsating vibration signals start. They indicate the risk of overheating.
- ▶ As soon as these pulsating vibration signals begin, you must reduce your level of activity so the hydraulic unit can cool down.
- ▶ Full activity may be resumed after the pulsating vibration signals stop.
- ▶ If the activity level is not reduced in spite of the pulsating vibration signals, this could lead to the hydraulic element overheating and, in extreme cases, cause damage to the product. In this case, the product should be inspected for damage by an O&P professional. If necessary, they will forward the product to an authorised Ottobock Service Center.

**⚠ CAUTION**

**Overloading due to unusual activities**

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by a prosthetist. If necessary, the prosthetist will forward the product to an authorised Ottobock Service Centre.

**⚠ CAUTION**

**Improper mode switching**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.
- ▶ Take the weight off the product and correct the switching, if required.

**⚠ CAUTION**

**Improper use of the stance function**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Make sure that you are standing safely when using the stance function and check the lock of the knee joint before placing your full weight on the prosthesis.
- ▶ Make sure you have been instructed in the correct use of the stance function by the O&P professional and/or therapist. Information on the stance function see page 56.

**⚠ CAUTION**

**Quickly pushing the hip forward with the prosthesis extended (e.g. serve while playing tennis)**

- > Falling due to unexpected activation of the swing phase.
- ▶ Note that the knee joint may flex unexpectedly when the hip is pushed forward quickly while the prosthesis is extended.
- ▶ Therefore, familiarise yourself with swing phase activation in such situations under secure conditions (e.g. while holding on to parallel bars) and with the instruction of trained, qualified personnel.
- ▶ For sports where this movement pattern may occur, use a corresponding pre-configured MyMode. For further information about the MyModes, see the section 'MyModes' (see page 61).

**⚠ CAUTION**

**Overloading due to changed body weight when carrying heavy objects, backpacks or children**

- > Falling due to unexpected behaviour of the product.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Note that the behaviour of the product can change due to increased weight. The swing phase may not be triggered, or triggered at the wrong time.
- ▶ Make sure that the maximum permissible body weight is not exceeded due to the additional weight.

#### 4.8 Notes on the safety modes

**⚠ CAUTION**

**Using the product in safety mode**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see page 70) have to be observed.
- ▶ Particular caution is necessary when using a bicycle without a freewheel (with a fixed gear).

**⚠ CAUTION**

**Safety mode cannot be activated due to malfunction caused by water penetration or mechanical damage**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Do not continue using the defective product.
- ▶ Consult the O&P professional promptly.

### **⚠ CAUTION**

#### **Safety mode cannot be deactivated**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ If safety mode cannot be deactivated by recharging the battery, a permanent error has occurred.
- ▶ Do not continue using the defective product.
- ▶ The product must be inspected by an authorised Ottobock Service Center. The O&P professional is your contact.

### **⚠ CAUTION**

#### **Safety signal occurs (ongoing vibration)**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Observe the warnings/error signals (see page 70).
- ▶ Do not continue using the product after the safety signal has been emitted.
- ▶ The product must be inspected by an authorised Ottobock Service Center. The O&P professional is your contact.

## **4.9 Instructions for use with an osseointegrated implant system**

### **⚠ WARNING**

#### **High mechanical loads due to normal or unusual situations, such as falling**

- > Overloading of the bone, which can lead to pain, loosening of the implant, death of bone tissue or bone fracture, among other things.
- > Damage or breakage of the implant system or its components (safety components...).
- ▶ Verify compliance with the fields of application, conditions of use and indications according to the information of the manufacturers, both for the knee joint and for the implant system.
- ▶ Note the instructions of the clinical personnel that indicated the use of the osseointegrated implant system.
- ▶ Note changes in your state of health that result in restrictions or doubt regarding the use of the osseointegrated connection.

## **5 Scope of Delivery and Accessories**

### **5.1 Scope of delivery**

- 1 pc. 3B5-3=ST Genium X3 (with threaded connector) or
- 1 pc. 3B5-3 Genium X3 (with pyramid) each with preinstalled 4X900 or 4X193-1 Genium X3 Protective Cover
- 1 pc. prosthesis passport
- 1 pc. instructions for use (user)
- 1 pc. 2R19 AXON tube adapter
- 1 pc. 757L16-4 power supply
- 1 pc. 4E60\* inductive charger
- 1 pc. 646C107 Bluetooth PIN card

### **5.2 Accessories**

The following components are not included in the scope of delivery and may be ordered separately:

- Genium X3 Protective Cover: 4X900
  - Genium X3 Protective Cover: 4X193-1
  - USB charging adapter: 757L43
- To connect the 757L43 USB charging adapter to the respective charger, observe the information in the USB charging adapter instructions for use.

• **“Cockpit” user app: 4X441-V\*=\***

for download from the app stores (Apple App Store, Google Play, ...). To do so, either enter the search terms “Ottobock”, “Cockpit” or scan the QR code.

Further information about the app and how it works can be found either in the link in the description of the app stores or in the installed app.

## 6 Charging the battery

The following points must be observed when charging the battery:

- Use the 757L16-4 power supply unit and 4E60\* battery charger to charge the battery.
- With average use, the capacity of the fully charged battery is sufficient for about 5 days.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.
- Note the permissible temperature range for charging the battery (see page 66).
- The distance between the battery charger and the receiver on the product must not exceed 2 mm.

### 6.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place (see fig. 1).
- 2) Connect the round, **three-pin** plug of the power supply to the receptacle on the battery charger (see fig. 2) so that the plug locks into place.

**INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.**

- 3) Plug the power supply unit into the outlet (see fig. 3).
  - The green LED on the back of the power supply lights up.
  - The LED ring (status indicator) on the rear of the charger lights up green to indicate the correct connection to the power supply.

→ If the green LED on the power supply and the LED ring on the battery charger do not light up, there is an error (see page 70).

## 6.2 Charging the prosthesis battery

### INFORMATION

When the Protective Cover is installed, the battery charger cable has to point to the upper closure. A correct knee joint charging process is only ensured with this alignment.

- 1) Connect the inductive charger to the receiver of the charging unit on the back of the product. The charger is held in place by a magnet.
  - The LED ring on the rear of the charger pulsates purple (4-second cycle).
  - If the LED ring lights up in a different colour, this indicates an error (see page 70).
- 2) The charging process starts.
  - Once the product battery is fully charged, all LEDs on the side of the battery charger light up.
- 3) After the charging process is complete, hold the prosthesis still and remove the inductive charger from the receiver.
  - A self-test is performed as well as an inspection of the product for impermissible changes to the product parameters (cybersecurity). The product is ready for operation only after corresponding feedback (see page 72).

## 6.3 Display of the current charge level

### 6.3.1 Display of battery charge level without additional devices

### INFORMATION

The charge level cannot be queried during the charging process or when MyMode is activated, e.g. by turning the prosthesis over. The product is in charging mode.

- 1) Turn the prosthesis by 180° (the sole of the foot must face up).
- 2) Hold still for 2 seconds and wait for beep signals.

Beep signal	Battery charge level
5x short	more than 80%
4x short	60% to 80%
3x short	40% to 60%
2x short	20% to 40%
1x short	less than 20%

## INFORMATION

If the parameter **Volume** is set to '0' in the user app or if mute mode is activated (silent mode), the beep signals are not emitted.

### 6.3.2 Display of the current charge level during the charging process

During the charging process, the current battery charge level is indicated by the number of LEDs lit on the side of the charger.

Quantity	Charge level
0	0 %-10 %
1	10 %-30 %
2	30 %-50 %
3	50 %-70 %
4	70 %-100 %
5	100 %

## 7 Use

### INFORMATION

#### Knee joint movement noise

When using exoprosthetic knee joints, servomotor, hydraulic, pneumatic or brake load dependent control functions can cause movement noises. This kind of noise is normal and unavoidable. It generally does not indicate any problems. If movement noise increases noticeably during the lifecycle of the knee joint, the knee joint should be inspected by the O&P professional immediately.

### 7.1 Standing

Knee control through high hydraulic resistance and static alignment.

A stance function can be enabled by the O&P professional. Please see the following section for further information on the stance function.

#### 7.1.1 Stance function

### INFORMATION

To use this function, it needs to be enabled by the O&P professional. It also has to be activated using the user app.

The stance function (standing mode) is a functional supplement to the basic mode (mode 1). This function makes it easier, for example, to stand on an inclined surface for a longer time. The joint is fixed in the flexion direction.

The stance function must be enabled by the O&P professional. In addition, the type of joint locking (deliberate/intuitive) has to be established by the O&P professional. The locking type cannot be changed via the user app.

#### Intuitive locking of the joint

The intuitive stance function recognises any situation that puts strain on the prosthesis in the flexion direction but where flexion is not permitted. Examples of this include standing on uneven or sloping surfaces. The knee joint is always locked in the flexion direction when the prosthetic leg is

not fully extended, under some amount of load and at rest. Upon forward or backward rollover or extension, the level of resistance is immediately reduced to stance phase resistance again. The knee joint is not locked when the above conditions are met and a sitting position is assumed (for example while driving).

**Deliberate locking of the joint**

- 1) Assume the desired knee angle.
  - 2) Do not remove the entire load from the prosthesis.
  - 3) Do not change the knee angle for a brief period (1/8 second). This time period prevents unintentional activation of the stance function while walking.
- The blocked joint can now be loaded in the flexion direction.

**Deliberate unlocking of the joint**

- By deliberately extending or unloading the knee joint, is it unlocked again.

**INFORMATION**

**Stance function with hip disarticulation amputation level**

Due to personal abilities and experiences with prostheses, these users may encounter difficulties with activating/deactivating the stance function. If these users wish to stand with the knee joint bent and locked for longer periods of time, the O&P professional can configure a MyMode that can be switched on/off using the user app.

**7.2 Walking**

Initial attempts at walking with the prosthesis always require the instruction of trained, qualified personnel. The hydraulics stabilise the knee joint in the stance phase and release the knee joint in the swing phase so that the leg can swing forward freely. Switching to the swing phase requires that the prosthesis roll over to the front out of the stride position.

**7.3 Running short distances ("walk-to-run" function)**

For covering short distances quickly, the knee joint detects a transition from walking to running in basic mode and automatically changes the following settings:

- The swing phase angle is increased
- Preflexion of 4° at heel strike (PreFlex) is reduced to 0°

The requirements to automatically switch to the running motion are fast forward movement of the prosthetic leg and high dynamic load on the knee joint. When stopping from the running motion, the changed settings are set back to the standard values.

**INFORMATION**

For covering longer distances, the O&P professional can configure a "Running" MyMode (see page 61).

## 7.4 Sitting down

The resistance in the prosthetic knee joint while sitting down ensures even bending into the sitting position.

The O&P professional can configure whether the sitting process is to be supported or not.

- 1) Place both feet side by side at the same level.
- 2) While sitting down, weight should be distributed evenly between both legs and the arm supports used if available.
- 3) Move the buttocks in the direction of the back support and lean the upper body forward.

## 7.5 Sitting

### INFORMATION

While sitting, the knee joint also switches to energy saving mode. This energy saving mode is activated regardless of whether the sitting function is activated or not.

If the patient is in a sitting position for more than two seconds (i.e. the thigh is close to horizontal and there is no load on the leg), the knee joint switches the resistance to a minimum in the extension direction.

A sitting function can be enabled by the O&P professional. For more information about the sitting function, see the following section.

### 7.5.1 Sitting function

### INFORMATION

To use this function, it needs to be enabled by the O&P professional. It also has to be activated using the user app.

In the sitting position, the resistance in the flexion direction is reduced in addition to the reduction of resistance in the extension direction. This makes it possible to swing the prosthetic leg freely.

## 7.6 Standing up

The flexion resistance is continuously increased while standing up.


- 1) Place the feet at the same level.
- 2) Lean the upper body forward.
- 3) Place the hands on arm supports, if available.
- 4) Stand up with support from the hands while keeping weight evenly distributed on the feet.

## 7.7 Walking up stairs step-over-step

### INFORMATION

To use this function, it needs to be enabled by the O&P professional. It also has to be activated using the user app.

Although the knee joint is passive, which means it cannot actively initiate movements, negotiating stairs step-over-step is possible. This function must be practised and executed consciously.

- 1) Lift the extended prosthesis off the floor.
- 2) Immediately after lifting the extended leg off the floor, extend the hip briefly and then abruptly flex it. This requires a sufficiently secure hold in the socket and a certain level of residual limb strength.  
→ This whip motion flexes the knee, because the knee joint automatically recognises the movement and sets the flexion resistance to minimum.

**INFORMATION: Take note of people behind you before executing the whip motion.**

- 3) When sufficient knee flexion has been achieved, the knee joint increases extension resistance so that there is enough time to position the foot on the next step before the knee joint is extended again.
- 4) Set the foot onto the next step.  
The support area for the foot on the step must be large enough that the heel does not extend back too far over the edge. With too little support area, the lower leg would extend too early and position the leg too far backwards. In this phase, the knee joint has already set the flexion resistance to maximum (blocked). The knee joint can no longer be flexed but only extended. This ensures that the leg does not buckle if the hip strength is not sufficient for the extending motion.
- 5) Support yourself with your hand on the contralateral side. A smooth wall will also work. This lateral support is intended to prevent the residual limb from twisting in the socket. Twisting can lead to unpleasant surface tension between the skin and the socket. Lateral support also improves balance.
- 6) Bring the knee into extended position. When the knee joint is fully extended, the initial position has been reached.
- 7) You can climb the next step or continue walking normally.

## 7.8 Overcoming obstacles

### INFORMATION

To use this function, it needs to be enabled by the O&P professional. It also has to be activated using the user app.

The stair function can also be used to cross obstacles:

- 1) Lift the extended prosthesis off the floor.
- 2) Briefly extend the hip.
- 3) Quickly flex the hip. This causes the knee to flex.
- 4) With the knee flexed, step over the obstacle.  
With sufficient knee flexion, the extension resistance is increased to allow enough time for crossing the obstacle.

## 7.9 Walking down stairs

This function must be consciously practised and executed. Only when the sole is properly positioned can the knee joint react correctly and permit controlled flexion.

- 1) Hold the handrail with one hand.
- 2) Position the leg with the prosthesis on the step so that the foot projects halfway over the edge of the step.  
→ This is the only way to ensure a secure rollover.
- 3) Roll the foot over the edge of the step.  
→ This flexes the prosthesis slowly and evenly in the knee joint.
- 4) Place the foot of the other leg onto the next step.
- 5) Place the leg with the prosthesis on the next step after that.

## 7.10 Walking down a ramp

Under increased flexion resistance, permit controlled flexion of the knee joint which lowers the body's centre of gravity.

## 7.11 Bluetooth

### 7.11.1 Establishing the Bluetooth connection

The Bluetooth function provides a wireless connection between the component and various devices. Bluetooth must be switched on at the component in order to establish a connection.

Bluetooth can be switched on in the following ways:

- Hold the prosthesis with the prosthetic foot down and then rotate it by 180° so that the prosthetic foot faces up. An acoustic signal and a vibration signal are emitted.
- Connect the battery charger to the component and disconnect it again after about 5 seconds.

### 7.12 Mute mode (silent mode)

Activating mute mode (silent mode) turns off the audible feedback signals and the vibration signals. However, warnings in case of component errors are still generated (see page 70).

Mute mode can be activated/deactivated via the user app.

#### INFORMATION

Connecting the battery charger automatically deactivates mute mode again.

## 7.13 Deep sleep mode

#### INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.





#### INFORMATION

If the **Volume** parameter is set to '0' in the user app, there are no acoustic signals.

When this mode is activated, the knee joint can be placed into what is known as a deep sleep, in which power consumption is reduced to a minimum. The knee joint offers no functionality in this mode. The safety mode resistance values are activated.

Deep sleep mode can be ended by using the user app or by connecting the battery charger. Waking from deep sleep mode using the user app can take up to 30 seconds. After ending deep sleep mode, the knee joint is in basic mode again.

## 8 MyModes

In addition to basic mode, the O&P professional can enable and configure up to 5 MyModes. These can be selected via the user app. Only the first 3 MyModes can be selected using movement patterns. Switching via movement patterns must be enabled by the O&P professional.

### 8.1 Running function as configured MyMode

For a prolonged walking movement, a **"Running"** MyMode can be configured by the O&P professional, which can be switched on via the user app or via a movement pattern.

In this mode, every step will be performed as a running step with a larger swing phase angle and no pre-flexion at heel strike (Preflex).

#### INFORMATION

The running function requires specialised running feet, the 1E95 Challenger, or prosthetic feet with axial compression such as the 1C61 Triton Vertical Shock. Feet without axial compression are generally not suited for running. Contact your O&P professional for further information.

### 8.2 Switching MyModes using motion patterns

#### INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

#### INFORMATION

If the 'Volume' parameter is set to '0' in the user app, there are no acoustic signals. Observe the vibration signal in this case.

#### Information on switching

- Switching and the number of movement patterns must be enabled by the O&P professional.
- Before the first step, always check whether the selected mode corresponds to the required motion type.


#### Requirements for successful switching using motion patterns

The following points must be observed to carry out switching successfully:

- Switching using movement patterns has to be enabled by the O&P professional.
- Position the prosthetic leg back slightly and bounce on the forefoot with the leg extended while maintaining constant contact with the floor.

- Weight must be placed on the forefoot while bouncing.
- Do not take the weight off fully during unloading while bouncing.

### Switching process

- 1) Position the prosthetic leg slightly to the rear (step position).
- 2) While maintaining constant contact with the floor, bounce on the forefoot with the leg extended a number of times in one second depending on the desired MyMode (MyMode 1 = 3 times, MyMode 2 = 4 times).
- 3) Fully unload and keep the prosthetic leg still in this position (step position).  
→ A beep and vibration signal will sound to confirm that the movement pattern has been recognised.

**INFORMATION: If this beep and vibration signal does not sound, the requirements were not observed when bouncing the foot or mute mode (silent mode) is activated. For more information about mute mode, see the section "Mute mode" (see page 60).**

- 4) After the beep and vibration signal sounds, keep the prosthetic leg extended and still for 1 second.

→ A confirmation signal will sound to indicate that the prosthesis has successfully switched to the corresponding MyMode (2 times = MyMode 1, 3 times = MyMode 2).

**INFORMATION: If this confirmation signal does not sound, the leg with the prosthesis was not correctly kept still or mute mode (silent mode) is activated. Repeat the process for correct switching. For more information about mute mode, see the section "Mute mode" (see page 60).**

### 8.3 Switching from a MyMode back to basic mode

#### Information on switching


- Regardless of the configuration of the MyModes, it is always possible to switch back to basic mode (mode 1) using a movement pattern.
- It is always possible to switch back to basic mode (mode 1) by connecting/disconnecting the battery charger.
- Before the first step, always check whether the selected mode corresponds to the required motion type.

#### Requirements for successful switching using motion patterns

The following points must be observed to carry out switching successfully:

- Position the prosthetic leg back slightly and bounce on the forefoot with the leg extended while maintaining constant contact with the floor.
- Weight must be placed on the forefoot while bouncing.
- Do not take the weight off fully during unloading while bouncing.

### Switching process

- 1) Position the prosthetic leg slightly to the rear (step position).
- 2) While maintaining constant contact with the floor, and with the leg extended, bounce on the forefoot three or more times.
- 3) Fully unload and keep the prosthetic leg still in this position (step position).  
→ A beep and vibration signal will sound to confirm that the movement pattern has been recognised.

**INFORMATION: If this beep and vibration signal does not sound, the requirements were not observed when bouncing the foot or mute mode (silent mode) is activated. For more information about mute mode, see the section "Mute mode" (see page 60).**

- 4) After the beep and vibration signal sounds, keep the prosthetic leg extended and still for 1 second.

→ A confirmation signal will sound to indicate that the prosthesis has successfully switched over to basic mode.

**INFORMATION: If this confirmation signal does not sound, the leg with the prosthesis was not correctly kept still or mute mode (silent mode) is activated. Repeat the process for correct switching. For more information about mute mode, see the section "Mute mode" (see page 60).**

## 9 Additional operating states (modes)

### 9.1 Empty battery mode

Beeps and vibration signals are emitted if the available battery charge level is 5 % (see page 70). During this time, the flexion resistance is set to the safety mode values. This may be low or high depending on the setting configured by the O&P professional. The prosthesis is then switched off. You can switch back to basic mode (mode 1) from empty battery mode by charging the product.

### 9.2 Mode for charging the prosthesis

The product is non-functional during charging.

The product is set to the flexion resistance of safety mode. This may be low or high depending on the setting configured by the O&P professional.

### 9.3 Safety mode

The product automatically switches to safety mode if a critical fault occurs (e.g. failure of a sensor signal). Safety mode remains in effect until the error has been rectified.

The switch to safety mode is indicated by beeps and vibration signals immediately prior to switching (see page 70).

Safety mode can be disabled by connecting and disconnecting the battery charger. If the product switches into safety mode again, this means a permanent error exists. The product must be inspected by an authorised Ottobock Service Centre.

Depending on the type and severity of the error, different remaining functionality is offered in safety mode. The swing phase control and stance phase extension resistance may be available or not depending on the type of error. This makes limited walking possible for the user depending on the type of error.

#### **The following remaining functionality is available:**

- **Moderate error**, e.g. tube adapter not connected: Continuous stance phase flexion resistance is set, with the ability for stance release.
- **Severe error**: Safety mode flexion resistance is configured. This flexion resistance may be low or high depending on the setting configured by the O&P professional. Depending on the type of error, the product may also be locked completely in the flexion direction.

#### **The following functions are deactivated in safety mode:**

- Stairs and obstacles function

- Stance function
- Sitting function

## 9.4 Overheating mode

### INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

When the hydraulic unit overheats due to uninterrupted, increased activity (e.g. extended walking downhill), the flexion resistance is increased along with the rising temperature in order to counteract the overheating. When the hydraulic unit cools down, the product switches back to the settings that existed prior to overheating mode.

Overheating mode is not activated in the MyModes.

Overheating mode is indicated by a long vibration every 5 seconds.

#### The following functions are deactivated in overheating mode:

- Sitting function
- Display of the battery charge level without additional equipment
- Switching to a MyMode

## 10 Storage and bleeding

Air may accumulate in the hydraulic unit if the product is stored for longer periods and not in an upright position. This is noticeable through sounds and irregular damping behaviour.

The automatic bleeding mechanism ensures that all functions of the product are again intact after approximately 10 - 20 steps.

### Storage

- Before storing the knee joint, the knee head has to be extended. The knee head must not be flexed!
- Avoid extended disuse of the product (use the product regularly).

## 11 Cleaning

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

### INFORMATION

Please note that the weight of dirt adhering to the prosthesis can affect the gait pattern.

## 12 Maintenance

Regular maintenance (service inspections) at 12-monthly intervals is mandatory in the interest of your own safety and in order to maintain operating reliability and protect the warranty, to maintain basic safety and the essential performance characteristics, and to ensure safety in regards to EMC.

When maintenance is due, this is indicated by feedback signals after disconnecting the battery charger (see "Operating states / error signals section" see page 69).

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

Components required for maintenance or repair:

Prosthesis, battery charger and power supply unit.

## 13 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 13.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 13.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 13.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The full text of the EU Declaration of Conformity is available at the following Internet address: <http://www.ottobock.com/conformity>

### 13.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

#### **Caution: Exposure to Radio Frequency Radiation.**

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s)..

Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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**Caution: Exposure to Radio Frequency Radiation.**

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population.

**Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.**

## 14 Technical data

<b>Ambient conditions</b>	
Transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage in the original packaging (≤3 months)	-20 °C/-4 °F to +40 °C/+104 °F Max. 93% relative humidity, non-condensing
Long-term storage in the original packaging (>3 months)	-20 °C/-4 °F to +20 °C/+68 °F Max. 93% relative humidity, non-condensing
Transport and storage between applications (without packaging)	-25 °C/-13 °F to +70 °C/158 °F Max. 93% relative humidity, non-condensing
Operation	-10 °C/+14 °F to +60 °C/+140 °F Max. 93% relative humidity, non-condensing
Time for warming to the operating temperature after storage between applications, from -25 °C/-13 °F at an ambient temperature of +20 °C/+68 °F	30 minutes
Time for cooling to the operating temperature after storage between applications, from +70 °C/+158 °F at an ambient temperature of +20 °C/+68 °F	30 minutes
Charging the battery	+10 °C/+50 °F to +45 °C/+113 °F

<b>Product</b>	
Reference number	3B5-3*/3B5-3=ST*
Mobility grade according to MOBIS	3 and 4
Maximum body weight including additional weight	150 kg
Protection rating	IP66/IP68 Maximum water depth: 3 m Maximum time: 1 hour
Waterproofness	Waterproof, corrosion-resistant, protected against penetration from jets of water

<b>Product</b>	
Range of Bluetooth connection to mobile device	Max. 10 m
Weight of the prosthesis without tube adapter, with Protective Cover	Approx. 1700 g
Information on the product's ruleset and firm-ware version	Accessible using the user app
Expected lifetime given compliance with pre-scribed maintenance intervals	6 years
Test procedure	ISO 10328-P6-150 kg/3 million load cycles

<b>Data transfer</b>	
Wireless technology	Bluetooth Smart Ready
Range	approx. 10 m / 32.8 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK, $\pi/4$ DQPSK, 8DPSK
Data rate (over the air)	2178 kbps (asymmetrical)
Maximum output power (EIRP):	+8.5 dBm

<b>Tube adapter</b>	
Reference number	2R19
Weight	190 g–300 g
Material	Aluminium
Max. body weight	150 kg
Protection rating	IP66/IP68 Maximum water depth: 3 m Maximum time: 1 hour
Water resistance	Waterproof, corrosion-resistant, resistant to penetration from jets of water
Lifetime	6 years

<b>Prosthesis battery</b>	
Battery type	Li-ion
Charging cycles (charging and discharging cycles) after which at least 80 % of the original battery capacity remains available	500
Charge level after 1 hour charging time	30 %
Charge level after 2 hours charging time	50 %
Charge level after 4 hours charging time	80 %
Charge level after 8 hours charging time	Fully charged
Product behaviour during the charging process	The product is without function.
Operating time of the prosthesis with new, fully charged battery at room temperature	Approx. 5 days with average use

<b>Power supply unit</b>	
Reference number	757L16-4
Type	FW8001M/12

<b>Power supply unit</b>	
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing
Operation	0 °C/+32 °F to +50 °C/+122 °F Max. 95% relative humidity Air pressure: 70–106 kPa (up to 3,000 m without pressure equalisation)
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V ===

<b>Battery charger</b>	
Reference number	4E60*
Storage and transport in original packaging	-25 °C to 70 °C/-13 °F to 158 °F
Storage and transport without packaging	-25 °C to 70 °C/-13 °F to 158 °F Max. 93% relative humidity, non-condensing
Operation	5 °C to 40 °C/41 °F to 104 °F Max. 93% relative humidity, non-condensing
Protection rating	IP40
Input voltage	12 V ===
Wireless technology	Proprietary protocol
Frequency range	270 kHz to 450 kHz
Modulation	ASK, load modulation
Maximum output power (EIRP)	-12.7 dBμA/m @ 10 m

## 15 Appendices

### 15.1 Symbols Used

	In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with regulations in your country can be harmful to health and the environment. Please observe the instructions of your national authority pertaining to return and collection procedures.		Non-ionising radiation
	Manufacturer		Compliance with the requirements of "FCC Part 15" (USA)
	Type BF applied part	IP40	Protection against penetration of solid foreign objects with a diameter greater than 1 mm, no protection against water
	Complies with the requirements of the Radiocommunications Act (AUS)	IP66	Dust-tight, protected against strong jets of water
		IP68	Dust-tight, protection against continuous submersion. Maximum depth: 3 m Maximum time: 1 hour
			Protect from wetness

<b>CE</b>	Declaration of conformity according to the applicable European directives
<b>SN</b>	Serial number (21)YYYYWWNNNN YYYY – year of manufacture WW – week of manufacture NNN – sequential number
<b>MD</b>	Medical device
<b>LOT</b>	Lot number (10)PPPPYYYYWW PPPP – plant YYYY – year of manufacture WW – week of manufacture
<b>UDI</b>	UDI number (Unique Device Identifier)
<b>REF</b>	Article number

	Data matrix code
<b>GTIN</b>	Global Trade Item Number
	Caution, hot surface
	Please note the instructions for use
	Limits for temperature
	Limits for atmospheric pressure
	Limits for relative humidity

## 15.2 Operating states/error signals

The prosthesis indicates operating states and error messages through beeps and vibration signals.

### 15.2.1 Signals for operating states

#### Battery charger connected/disconnected

Beep signal	Vibration signal	Event
–	3x long	Charging mode started (3 sec. after connecting the battery charger)
1x short	1x short	Self-test completed successfully, product is operational

#### Mode switching

##### INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

##### INFORMATION

If the **Volume** parameter is set to '0' in the user app, there are no acoustic signals.

Beep signal	Vibration signal	Additional action performed	Event
1x short	1x short	Mode switching using the user app	Mode switching is performed using the user app.
1x long	1x long	Bouncing on the forefoot followed by weight taken off the prosthetic leg	Bouncing pattern recognised.
1x short	1x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to basic mode (mode 1) carried out.

Beep signal	Vibration signal	Additional action performed	Event
2x short	2x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to MyMode 1 (mode 2) carried out.
3x short	3x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to MyMode 2 (mode 3) carried out.



## 15.2.2 Warnings/error signals

### Error during use



Beep signal	Vibration signal	Event	Required action
-	1x long at intervals of approx. 5 seconds (if mute mode (silent mode) is activated, this signal is not generated)	Overheated hydraulics	Reduce activity.
-	3x long	Battery charge level under 25%	Charge battery soon. Remaining operating time approx. 24 hours
-	5x long	Battery charge level under 10%	Charge battery soon. Remaining operating time approx. 6 hours
5x long	5x long repeated every 60 seconds	<b>Error of moderate severity (see page 63)</b> e.g. a sensor is not operational	Walking possible with restrictions. Note the change in flexion resistance. The product must be inspected by an O&P professional immediately.
10x long	10x long	Charge level 5% After the beep and vibration signals, the product switches to empty battery mode and then switches off.	Charge the battery.
30x long	1x long, 1x short repeated every 3 seconds	<b>Severe error/indication of safety mode activation (see page 63)</b> e.g. one or more sensors are not operational	Attempt to reset this error by connecting/disconnecting the battery charger. If the error persists, use of the product is prohibited. The product must be inspected by an O&P professional immediately.

Beep signal	Vibration signal	Event	Required action
-	Continuous	<b>Total failure</b> Electronic control no longer possible. Safety mode active or undetermined valve state. Undetermined product behaviour.	Attempt to reset this error by connecting/disconnecting the battery charger. If the error persists, use of the product is prohibited. The product must be inspected by an O&P professional immediately.

### Error while charging the product

LED on power supply	Status LED on battery charger	Error	Resolution
○		Country-specific plug adapter not fully engaged on power supply	Check whether the country-specific plug adapter is fully engaged on the power supply.
		Non-functional socket	Check socket with another electric appliance.
		Defective power supply	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●		No connection between battery charger and power supply	Check whether the charging cable plug is fully engaged on the battery charger.
		Defective battery charger	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.

	Status LED	Charging status indicator (5 LEDs)	Error	Resolution
	The LED ring is lit in weak purple	No LED is lit	Distance between battery charger and receiver of the charging unit on the prosthesis too great. If the distance is more than 2 mm, the prosthesis cannot be charged.	Reduce distance between battery charger and receiver of the charging unit.

	Status LED	Charging status indicator (5 LEDs)	Error	Resolution
	The LED ring is lit up yellow	LED 2 and 4 are lit up	Charger excess temperature	Check whether the specified ambient conditions for charging the battery are met (see page 66).
		LED 1, 3 and 5 are lit up	Excessively high or low temperature of the prosthesis	
		LED 3 is lit up	The prosthesis is not being charged Distance between battery charger and receiver of the charging unit too great.	Connection may be improved by reducing the distance between the charger and receiver of the charging unit.
	The LED ring is lit up green		The battery charger is functional but not connected to the receiver, or the distance between the battery charger and receiver of the charging unit is too great.	Connect the battery charger or reduce the distance between the charger and the receiver of the charging unit on the prosthesis.
	The LED ring is flashing red		The prosthesis is not being charged Defective battery charger.	Reset the error by disconnecting and connecting the power supply. If the error persists, the battery charger and power supply must be inspected by an authorised Ottobock Service Center.

### 15.2.3 Status signals

#### Battery charger connected

LED on power supply	Status LED on battery charger	Event
		Power supply and battery charger operational. Battery charger not yet connected to receiver.
		Battery charger connected to receiver, good connection. This display turns off automatically after approximately one minute to avoid bothersome light at night. The charging process is not affected.

## Battery charger disconnected

Beep signal	Vibration signal	Event	Resolution
1x short	1x short	Self-test and check for impermissible changes to the product parameters (cybersecurity) completed successfully. Product is ready for operation.	
3x short	3x short	Maintenance note: e.g.: Maintenance interval has been exceeded, temporary disruption of a sensor signal	<ul style="list-style-type: none"> <li>• Use the user app to check the next maintenance date for the prosthesis. If the date is within the next month, make a maintenance appointment with the O&amp;P professional. For this appointment, the O&amp;P professional has to be provided with the battery charger and power supply in addition to the prosthesis and tube adapter.</li> <li>• Conduct the self-test again by connecting/disconnecting the battery charger.</li> <li>• If the beep is emitted again and the maintenance date has not been reached or exceeded, you should visit the O&amp;P professional soon. If necessary, the O&amp;P professional will forward the prosthesis to an authorised Ottobock Service Center.</li> <li>• The product can be used without restrictions. However, vibration signals may not be generated.</li> </ul>

### Battery charge level

During the charging process, the current battery charge level is indicated by the number of LEDs lit on the side of the charger.

LEDs	0	1	2	3	4	5
Charge level	0 %-10 %	10 %-30 %	30 %-50 %	50 %-70 %	70 %-100 %	100 %

## 15.3 Directives and manufacturer's declaration

### 15.3.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Observe the safety notes in the section "Information on proximity to certain areas" (see page 45).

## Electromagnetic emissions

Interference measurements	Compliance	Electromagnetic environment directive
HF emissions according to CISPR 11	Group 1/class B	The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighbouring electronic devices is unlikely.
Harmonics according to IEC 61000-3-2	Not applicable – power below 75 W	–
Voltage fluctuations/flicker according to IEC 61000-3-3	Product meets the requirements of the standard.	–

## Electromagnetic interference immunity

Phenomenon	EMC basic standard or test procedure	Interference immunity test level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air,
High-frequency electromagnetic fields	IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz
Magnetic fields with rated power frequencies	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Electrical fast transients/bursts	IEC 61000-4-4	± 2 kV 100 kHz repetition rate
Surges Line against line	IEC 61000-4-5	± 0.5 kV, ± 1 kV
Conducted interference induced by high-frequency fields	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage drops	IEC 61000-4-11	0% $U_T$ ; 1/2 period At 0, 45, 90, 135, 180, 225, 270 and 315 degrees
		0% $U_T$ ; 1 period and 70% $U_T$ ; 25/30 periods Single phase: at 0 degrees
Voltage interruptions	IEC 61000-4-11	0% $U_T$ ; 250/300 periods

## Interference resistance against wireless communication devices

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	1.8	0.3	28
710	704 to 787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 to 960	GSM 800/90- 0, TETRA 800, iDEN 820, CDMA 850, GSM 800/90- 0, LTE band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1,720	1,700 to 1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1,845						
1,970						
2,450	2,400 to 2,570	Bluetooth WLAN 802.1- 1 b/g/n, RFID 2450 LTE band 7	Pulse modulation 217 Hz	2	0.3	28
5,240	5,100 to 5,800	WLAN 802.1- 1 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5,500						
5,785						

## Immunity to magnetic fields in close range

Test frequency	Modulation	Interference immunity test level [A/m]
30 kHz	CW	8
134.2 kHz	Pulse modulation 2.1 kHz	65
13.56 MHz	Pulse modulation 50 kHz	7.5

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