
1C70 Evanto

DE	Gebruiksaanweisung	7
EN	Instructions for use	13
FR	Instructions d'utilisation	20
IT	Istruzioni per l'uso	27
ES	Instrucciones de uso	34
PT	Instruções de utilização	41
NL	Gebruiksaanwijzing	48
SV	Bruksanvisning	55
DA	Brugsanvisning	61
NO	Bruksanvisning	67
FI	Käyttöohje	74
PL	Instrukcja użytkowania	80
HU	Használati útmutató	87
CS	Návod k použití	94
RO	Instrucţiuni de utilizare	100
HR	Upute za uporabu	107
SL	Navodila za uporabo	114
SK	Návod na použitie	120
BG	Инструкция за употреба	127
TR	Kullanım kılavuzu	134
EL	Οδηγίες χρήσης	140
RU	Руководство по применению	148
UK	Інструкція щодо використання	155
JA	取扱説明書	162
ZH	使用说明书	168
KO	사용 설명서	174

1

Bench Alignment, transtibial

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INFORMATION

Date of last update: 2025-11-20

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.

- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The 1C70 Evanto prosthetic foot has two exchangeable heel wedges. The heel wedges influence the behaviour of the prosthetic foot in different phases of the gait cycle. This allows the prosthetic foot to be individually adapted.

2 Safety-related information

CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, treat the product with care.
- ▶ Observe the information on product combinations in the instructions for use for the products.
- ▶ Observe the maximum lifetime of the product (see page 15).
- ▶ Observe the product's operating conditions (see page 14).
- ▶ Do not expose the product to prohibited environmental conditions (see page 14).
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not use the product if it is damaged or if you are not sure whether it is working properly. Have the product inspected, repaired or replaced by the manufacturer or a specialist workshop.
- ▶ The product is intended for use by one person only and may not be reused by other persons.

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover behaviour) or delamination of the spring are indications of loss of functionality. Unusual noises can indicate a loss of functionality.

3 Intended use

3.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

3.2 Operating conditions

- **The prosthetic foot is not suitable for sports prostheses.**
Certain movement patterns in sports (e.g. soccer, handball or tennis) can damage the prosthetic foot: pronounced stopping when running and spontaneous changes in direction

The product is recommended for:

- Mobility grade 2: restricted outdoor walker
- Mobility grade 3: unrestricted outdoor walker
- Mobility grade 4: unrestricted outdoor walker with particularly high demands

Body weight [kg]	Spring stiffness
Up to 58	1
59 to 72	2
73 to 86	3
87 to 100	4
101 to 125	5

TIP: Also try the next higher stiffness so the user can experience the difference.

3.3 Environmental conditions

Storage and transport
Storage temperature: -20 °C to +60 °C, relative humidity: 20 % to 90 %, no vibrations or impacts

Allowable environmental conditions
Temperature range: -10 °C to +45 °C
Chemicals/liquids: Fresh water, soap solution, chlorine water, salt water ≤ 3.5 % salt content
Moisture: Submersion: max. 1 h in 2 m depth, relative humidity: no restrictions
Solids: dust, occasional contact with sand
Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and tear and damage (see page 19).

Prohibited environmental conditions
Solids: highly liquid-binding particles (e.g. talcum), dust in high concentrations (e.g. construction site), intensive contact with sand
Chemicals/liquids: perspiration, urine, acids, continuous use in liquid media

3.4 Product combinations

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Limited combination options for Ottobock components

Combination with prosthetic knee joints	
From stiffness 4 and foot size 28	Do not combine with 3R85 Dynion
From stiffness 5 and foot size 26	Do not combine with 3R55

The prosthetic foot generates high moments in the ankle area. Use structural components with higher weight limits:

Body weight [kg]	73 to 86	87 to 100	101 to 125
Foot size [cm]	From 28	From 26	From 28
Structural component weight limit [kg]	≥ 125	≥ 125	≥ 150

3.5 Lifetime

Depending on the user's activity level, the maximum lifetime of the product is 3 years.

Footshell, protective sock

The products are wear and tear parts that are subject to normal degeneration.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	647G1431=*
1	Quick start guide	647G1809=*
1	Prosthetic foot	1C70=*
1	Protective sock	Size 22 to 25: 2U3=1-7 Size 26 to 30: 2U3=2-7
Additional heel wedges		
1	Upper heel wedge	2F70=*
1	Lower heel wedge	2F71=*

5 Accessories

Designation	Reference number
Cover	Size: 22 to 23: 2Z685=2223-1 Size: 24 to 25: 2Z685=2425-1 Size: 26 to 27: 2Z685=2627-1 Size: 28 to 30: 2Z685=2830-1
Sealing	Size: 22 to 23: 2Z840=2223 Size: 24 to 25: 2Z840=2425 Size: 26 to 27: 2Z840=2627 Size: 28 to 30: 2Z840=2830
Screw for strap tensioner	2Z684=1 2Z684=2 2Z684=3 2Z684=4 2Z684=5 2Z684=6

6 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- Observe the alignment, assembly and adjustment instructions.

NOTICE

Altering the prosthetic foot or footshell

Premature wear resulting from damage to the product

- Do not alter the prosthetic foot or footshell.

6.1 Applying/removing the footshell

CAUTION! Always use the prosthetic foot with a protective sock and footshell.

- 1) Pull the protective sock over the prosthetic foot to protect the prosthetic foot against abrasive particles and prevent noises in the footshell.
- 2) Apply the footshell as described in the footshell instructions for use.

Removing the footshell

- Remove the footshell as described in the footshell instructions for use.

6.2 Bench Alignment

Bench alignment, transtibial prosthesis

> **Required materials and tools:** 662M4 goniometer, 743S12 heel height measuring device, 743A80 50:50 gauge, alignment apparatus (e.g. 743A220 PROS.A. Assembly)

- 1) Position the prosthetic foot in the alignment apparatus.
 - **Heel height: effective heel height of everyday shoe + 5 mm**
- 2) **Sagittal plane:** Align the prosthetic foot with the alignment reference line.
 - Position the alignment reference line on the markings of the footshell.
- 3) **Frontal plane:** Align the prosthetic foot with the alignment reference line.
 - Position the alignment reference line on the markings of the footshell: **external rotation approx. 5°.**
- 4) Install the structural components between the prosthetic foot and prosthetic socket.

- 5) **Sagittal plane:** Align the prosthetic socket with the alignment reference line.
 - Determine the centre of the prosthetic socket proximally and distally with the 50:50 gauge and mark the centre line.
 - Mark the socket reference point on the centre line: **For transtibial prostheses, at the height of the centre of the patella.** The alignment reference line runs vertically through the socket reference point.
 - Adjust the socket flexion by rotating around the socket reference point: **individual residual limb flexion + 5°**
- 6) **Frontal plane:** Align the prosthetic socket with the alignment reference line.
 - Take the abduction or adduction of the residual limb into account.

Bench alignment, thigh/knee disarticulation prosthesis

- Observe the instructions for use for the prosthetic knee joint.

6.3 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the 3D L.A.S.A.R. Posture and adapting it as needed.
- **Check the length of the prosthesis:** Compression of the prosthetic foot and other components (e.g. a vacuum pump) may require correction of the length.

Optionally, the alignment recommendations (transfemoral modular leg prostheses: **646F219***, transtibial modular leg prostheses: **646F336***) can be requested from Ottobock as needed.

6.4 Dynamic Trial Fitting

The optimal gait pattern is developed during dynamic fitting. The prosthesis is aligned optimally in the frontal plane and the sagittal plane for this purpose.

Walking and standing with Evanto looks and feels different. The user typically immediately feels increased flexibility when walking and standing.

Recommended action

- Let the user walk for a few minutes to get used to the new feel and flexibility of the Evanto. Observe the process of heel strike, rollover and toe-off.
- Have the user walk up ramps. Compared to the previously used prosthetic foot, this should be perceived as equally difficult or easier.
- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal planes is achieved after the leg begins to bear weight following the heel strike. Avoid medial movement of the knee joint.
 - If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction.
 - If the knee joint moves in the medial direction in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.

6.5 Adjusting stiffness

The stiffness of the prosthetic foot can be adjusted by replacing the two heel wedges:

- **Upper heel wedge:** Determines the heel stiffness.
- **Lower heel wedge:** Determines the forefoot stiffness and heel stiffness.

Heel stiffness affects the knee flexion effect:

- **Soft heel wedge:** less knee flexion effect.
- **Hard heel wedge:** more knee flexion effect).

Forefoot stiffness affects the rollover of the prosthetic foot.

Examples

Situation	Measures
The user feels a stepped rollover (second half of the stance phase).	First insert a harder lower heel wedge. To make the heel softer again, insert a softer upper heel wedge.

Situation	Measures
The user feels that there is too much resistance when walking up a slope.	Move the prosthetic foot gradually in the posterior direction.
The user feels a double rocking during loading response.	Insert a harder upper heel wedge.
The user experiences a dead point during the rollover (first half of the stance phase).	Insert a harder upper heel wedge.
Some users (mobility grade 3 and particularly mobility grade 2) prefer a slightly more anterior foot position. It reduces the knee flexing effect and provides a secure heel strike and stable rollover.	

Replacing heel wedges

Heel wedge degrees of hardness: The heel wedges are marked with numbers. The smaller the number, the softer the heel wedge.

CAUTION

Using the prosthetic foot without seal and cover

Risk of injury due to penetration of abrasive particles into the prosthetic foot

- Reinstall the seal and cover before using the foot.

The strap tensioner is used to tighten the strap looped around the foot. It consists of a screw and a thrust piece, which is pressed against the strap by the screw. The strap is relaxed to replace the heel wedges. First, the cover and seal must be removed, which protect against dust and other particles.

- 1) Lift off the cover of the prosthetic foot towards the heel (see fig. 3).
- 2) **If present:** Take off the seal above the strap tensioner (see fig. 4).
- 3) Fully unscrew the strap tensioner screw (see fig. 5).
- 4) **NOTICE! Make sure the strap stays in position. This prevents damage to the base spring of the prosthetic foot.**

Remove the lower heel wedge (see fig. 6).

- 5) **Optional:** Replace the upper heel wedge, aligning the recess in the heel wedge upwards (see fig. 7, see fig. 8).

INFORMATION: As soon as the upper heel wedge has been removed, the seat of the strap tensioner can be easily checked from the bottom.

- 6) Insert the lower heel wedge (see fig. 9).
- 7) Check whether the strap lies correctly in the guide:
 - The strap is positioned all around the guide.
 - **If the strap is reinforced:** The sewn-on strap layers are aligned inwards and sit over the surfaces to the left and right of the strap tensioner (see fig. 10, see fig. 11).
 - The strap tensioner is not at an angle (see fig. 12).
- 8) Tighten the strap tensioner screw **4 Nm** (see fig. 13).
- 9) **If present:** Reposition the seal above the strap tensioner (see fig. 14). The recess on the underside of the seal must sit exactly on the belt tensioner screw.
- 10) Put the cover of the prosthetic foot back on (see fig. 15 **1**, **2**).

TIP: The cover sits tightly above the seal. Pull the cover apart slightly to the side when applying it so that you do not move the flexible seal.
- 11) Investigate whether the seal is correctly positioned under the cover:
 - An even semicircle is visible between the cover and the transition to the pyramid (see fig. 16).
 - The seal does not protrude anywhere under the cover.

7 Cleaning

> **Allowable cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)

- 1) Pull the footshell and protective sock off the prosthetic foot to clean all products individually.
- 2) **NOTICE! To avoid product damage, only use the permissible cleaning agents.**
Clean all products with clear fresh water and a pH-neutral soap. Some dirt can also be removed with compressed air.
- 3) **If present:** Use a toothpick to remove dirt from the contours that drain off the adapters or footshell and rinse.
- 4) Rinse the soap away with clear fresh water. Rinse the footshell several times until all dirt has been removed.
- 5) Dry the products with a soft cloth or blow off the water with compressed air.
- 6) Allow to air dry in order to remove residual moisture.

If noises occur or dirt is difficult to remove:


Visit the Evanto YouTube playlist. You will find video instructions for intensive cleaning of the prosthetic foot there.

8 Maintenance

- ▶ **After the first 30 days of use:** Perform a visual inspection and functional test of the product and the prosthesis.
- ▶ **Every 12 months:** Check the function of the product and the prosthesis and inspect for signs of wear.

8.1 Strap tensioner – replacing the screw

The strap tensioner consists of a thrust piece and a screw. To compensate for fabrication tolerances of the strap, the screw is available in 6 different lengths. During the fabrication of the prosthetic foot, the appropriate screw is selected for the strap length.

Point indentations on the screw head indicate the length of the screw. The more indentations (1 to 6), the longer the screw.

- ▶ Replace a defective screw with a new screw of the same length.
- ▶ If you are not sure of the correct length of the screw: Contact Ottobock.

9 Disposal

Do not dispose of the product in household waste. Improper disposal is harmful to health and the environment. Observe the instructions of the local authority regarding return, collection and disposal.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for

damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

A series of horizontal lines for writing, spanning the width of the page.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Ottobock SE & Co. KGaA
Max-Näder-Straße 15 · 37115 Duderstadt · Germany
T +49 5527 848-0 · F +49 5527 848-3360
healthcare@ottobock.de · www.ottobock.com