8E38=6, 8E39=6

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	14
FR Instructions d'utilisation (personnel spécialisé)	25
IT Istruzioni per l'uso (personale tecnico specializzato)	37
ES Instrucciones de uso (personal técnico especializado)	48
PT Instruções de utilização (pessoal técnico)	59
NL Gebruiksaanwijzing (deskundig personeel)	70
SV Bruksanvisning (fackpersonal)	81
DA Brugsanvisning (faguddannet personale)	92
HU Használati útmutató (szakszemélyzet)	103
TR Kullanım kılavuzu (uzman personel)	114
EL Οδηγίες χρήσης (Τεχνικό προσωπικό)	125
JA 取扱説明書(有資格担当者)	137

1


2

3

4

5

INFORMATION

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INFORMATION

Date of last update: 2026-02-24

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The products 8E38=6*, 8E39=6* are referred to as the product/prosthesis/terminal device in the following.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Function

The System Electric Hand DMC plus is equipped with the proportional DMC plus control device.

Two independent measurement and control systems proportionally control both the gripping speed and gripping force. The gripping speed and gripping force are determined by the strength of the muscle signal. This makes the gripping process significantly more physiological and fragile objects can be gently grasped with ease.

The DMC plus control device is installed in a black housing with a black adjuster. The coaxial bushing and hand cable are marked with a brown ring.

The System Electric Hands with DMC plus control device are operated with the 13E125, 13E200 or 13E202 electrodes.

Two different operating modes can be selected with a function plug. The System Electric Hand DMC plus is delivered with the function plug connected in DMC mode.

2.1.1 Wrist joint versions

The products 8E38=* and 8E39=* are distinguished by different variants of the prosthetic wrist joint:

8E38=* (quick-disconnect wrist)

Makes it easy to separate the terminal device from the prosthetic socket. The terminal device can be removed quickly when needed via a 360° rotating movement and replaced by a different terminal device with the same quick-disconnect unit.

8E39=* (lamination ring)

Low-profile connection for users with a long forearm or transcarpal amputation. The hand can be rotated against a constant friction, which can be adjusted during the fitting process.

The 9S110=* lamination ring has to be ordered separately.

The 13E190 or 13E190=150 switch block is required when this terminal device is used.

2.2 Product combinations

This product can be combined with the following Ottobock components:

Power supply (rechargeable battery)

- MyoEnergy Integral 757B35=1, 757B35=3, 757B35=5
- 757B20/757B21 EnergyPack

Battery chargers

The following battery chargers can be used depending on the rechargeable battery:

- 757L20 battery charger (including 757L16=* power supply unit)
- 757L35 battery charger (including 757L16=* power supply unit)

Lamination rings

- 10S1=* lamination ring (for 8E38=*)
- 9S110=* lamination ring (for 8E39=*)

Active rotation

- 10S17 Electric wrist rotator
- 13E205 MyoRotronic

Passive rotation

- 9E169 coaxial plug
- 10S4 coupling piece

Elbow joints

- 12K44=* ErgoArm Hybrid Plus
- 12K50=* ErgoArm Electronic plus
- DynamicArm, DynamicArm Plus 12K100N=*/12K110N=*

3 Intended use

3.1 Indications for use

The product is intended **exclusively** for upper limb exoprosthetic fittings.

3.2 Conditions of use

The product is intended **exclusively** for use on **one** user. Use of the product by other persons is **prohibited** by the manufacturer.

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports which involve excessive strain and/or shocks to the wrist (pushups, downhill racing, mountain biking, etc.) or extreme sports (free climbing, paragliding, etc.). Furthermore, the product should not be used to operate motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

See the technical data for the permissible environmental conditions (see page 23).

3.3 Indications

- Amputation level below-elbow, above-elbow and shoulder disarticulation
- For unilateral or bilateral amputation
- Dysmelia of the forearm or upper arm
- The user must be able to understand usage instructions and safety notices and put them into practice.
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the sections on "Safety" (see page 16) and "Intended use" (see page 15).

3.5 Qualification

Treatment may only be carried out by trained, qualified personnel.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 WARNING
The heading describes the source and/or the type of hazard The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible: > E.g.: Consequence 1 in the event of failure to observe the hazard > E.g.: Consequence 2 in the event of failure to observe the hazard ► This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

 WARNING
Using the prosthesis while operating a vehicle Accident due to unexpected prosthesis behaviour. ► The prosthesis should not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).
 WARNING
Using the prosthesis while operating machines Injury due to unexpected prosthesis behaviour. ► The prosthesis should not be used to operate industrial machines or motor-driven equipment.

⚠ WARNING

Operating the prosthesis near active implanted systems

Interference with active implantable systems (e.g. pacemaker, defibrillator, etc.) due to electromagnetic emissions of the prosthesis.

- ▶ When operating the prosthesis in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety notes stipulated by the manufacturer of the implant.

⚠ WARNING

Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

⚠ WARNING

Skin contact with leaking lubricants due to mechanical defects

Injury due to skin irritation.

- ▶ Avoid contact between leaking lubricants and the mouth, nose and eyes.
- ▶ The product must be inspected by an authorised Ottobock Service Centre.

⚠ CAUTION

Signs of wear on the product

Injury due to faulty control or malfunction of the product

- ▶ In the interest of user safety and in order to maintain operating reliability, the terminal device has to be inspected by an authorised Ottobock Service Centre in case of noticeable restrictions of functionality.
- ▶ Note that the functionality of the terminal device may be limited if the battery charge level gets too low.

⚠ CAUTION

Independent manipulation of the product

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Use of a damaged product

Injury due to loss of product functionality.

- ▶ Prior to use, conduct an external visual inspection to verify that all parts of the product are undamaged.
- ▶ In case of damage, have the product repaired promptly.

⚠ CAUTION

Penetration of product with dirt and humidity

Injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

4.4 Information on Alignment/Adjustment

⚠ CAUTION

Operator errors when adjusting settings using the adjustment software

Injury due to unexpected product behaviour.

- ▶ Participation in an Ottobock product training course is required prior to initial use. Additional product training courses may become necessary to qualify for software updates.
- ▶ Transfer the setting changes to the terminal device before you check the settings on the user.
- ▶ Use the online help function integrated into the software.

⚠ CAUTION

Incorrect electrode settings/electrode assignment

Injury due to unexpected product behaviour.

- ▶ The electrodes are to be placed on intact skin only and with as much skin contact as possible. In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary. If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact the Ottobock branch responsible for your country.
- ▶ Set the electrode sensitivity as low as possible in order to reduce interference from powerful electromagnetic radiation (e.g. visible or concealed theft prevention systems at the entrance/exit of stores), metal detectors/body scanners for persons (e.g. in airports) or other sources of strong electromagnetic interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ Make sure that the electrode connection positions correspond to physiological opening and closing for the corresponding muscle groups.
- ▶ Make sure that selecting a control program that is unsuitable for the user may lead to injuries.

⚠ CAUTION

Incorrect electrode settings due to muscle fatigue

Injury due to faulty control or malfunction of the product.

- ▶ Allow the user to rest during the adjustment of the electrodes.

⚠ CAUTION

Failure to observe the instructions for use of all prosthesis components used

Injury due to unexpected product behaviour.

- ▶ Observe the instructions for use of the prosthesis components used.

⚠ CAUTION

Use of unsuitable prosthetic components

Injury due to unexpected product behaviour.

- ▶ Use the product only in combination with components listed in the section "Combination possibilities" (see page 15).

⚠ CAUTION

Battery not connected correctly

Injury due to faulty control or malfunction of the product.

- ▶ Make sure the battery connections are not switched.
- ▶ Check the battery connections by pulling slightly on the cable.

⚠ CAUTION

Use of silicone spray when putting on the prosthetic glove

Injury due to a loss of grip resulting from insufficient adhesion of the prosthetic glove on the hand.

- ▶ Only use 640F18 donning spray to put on the prosthetic glove. Note the instructions for use for the prosthetic glove.

4.5 Information on Proximity to Certain Areas

⚠ CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm from HF communication devices is recommended.

⚠ CAUTION

Operating the product in very close proximity to other electronic devices

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify intended use in the existing setup.

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

4.6 Information on Use

⚠ CAUTION

Improper handling

Injury due to faulty operation or malfunction of the product.

- ▶ Instruct the user in the proper use of the product.

CAUTION

Improper product care

- > Injuries due to faulty control/malfunction of the product or damage to the mechanical components
- > Damage or breakage due to brittleness of plastics caused by the use of acetone or petrol, for example.
- ▶ Clean the product only as described in the "Cleaning" section (see page 22).
- ▶ Do not clean the product under running water.
- ▶ When using a prosthetic glove, also note the instructions for use of the glove.

CAUTION

Grasping objects with incorrect gripping forces

Injury due to unexpected product behaviour.

- ▶ Note that the gripping has to be controlled manually depending on the consistency (soft/hard) of the object being grasped.

CAUTION

Risk of pinching between the fingertips

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the fingertips when using the product.
- ▶ Ensure that no body parts are between the fingertips when closing the hand.
- ▶ When closing the hand, ensure that fingers and other body parts are not in the area of the finger joints.
- ▶ Make sure the product is switched off for cleaning.

CAUTION

Unintentional unlocking of the terminal device

Risk of injury due to the terminal device disconnecting from the forearm (e.g. while carrying objects).

- ▶ Ensure that the connection is carried out correctly when connecting the terminal device to the prosthetic socket or components.
- ▶ When using a quick-disconnect wrist, position the terminal device so that slight twisting cannot disconnect the terminal device from the forearm.

CAUTION

Use of unapproved accessories

- > Injury due to product malfunction as a result of reduced resistance to interference.
- > Interference of other electronic devices due to increased emissions.
- ▶ Only use the product in combination with the accessories, signal converters and cables listed in the sections "Combination possibilities" (see page 15), "Scope of delivery" (see page 20) and "Accessories" (see page 21).

5 Scope of Delivery and Accessories

5.1 Scope of delivery

- 1 pc. 8E38=6* System Electric Hand DMC plus
or

- 1 pc. 8E39=6* System Electric Hand DMC plus
- 1 pc. Instructions for use (qualified personnel)
- 1 pc. instructions for use (users)

5.2 Accessories

- USB charging adapter: 757L43 (only in combination with MyoEnergy Integral)
To connect the 757L43 USB charging adapter to the respective charger, observe the information in the USB charging adapter instructions for use.

6 Preparing the product for use

6.1 Charging the battery

The following information can be found in the instructions for use of the respective batteries or the elbow components:

- Handling the battery
- Querying the charge level
- Feedback (beeps and vibration signals)

6.2 Adjusting the electrodes

The ideal position of the electrodes is determined with the 757M11=X-CHANGE MyoBoy, see 647G265=ALL_INT instructions for use.

Adjust the electrodes so the user can hold the respective electrode signal above the HIGH value for about 2 seconds (see fig. 4).

6.3 Control programs

The DMC plus control device encompasses the conventional DMC control device and a DMC control device with a higher activation threshold (DMC plus). The desired control variant is selected with the 13E185 function plug. The DMC plus control device is installed in a black housing with a black adjusting ring. The coaxial bushing and hand cable are marked with a brown ring.

Two different operating modes can be selected with the function plug of the product with the DMC plus control device.

6.3.1 Operating mode of digital control device: Function plug connected

Function plug connected: see fig. 1.

The gripping speed and gripping force are determined by the strength of the muscle signal. When the strength of the muscle signal changes, the gripping speed and gripping force immediately adjust proportionally to the changed muscle signal.

6.3.2 Operating mode of dual-channel control device: Function plug not connected

Function plug not connected: see fig. 2.

The gripping speed and gripping force are determined by the strength of the muscle signal. When the strength of the muscle signal changes, the gripping speed and gripping force immediately adjust proportionally to the changed muscle signal.

After gripping with maximum force, the activation threshold in the OPEN direction is raised to a higher value. This increase reduces the risk of opening the hand with unintended muscle signals.

The reliability of the grip is improved, for example, for eating.

To change the operating mode, remove the lock ring of the system inner hand. Partially pull down the system inner hand and remove the function plug. (see fig. 3).

Pull the system inner hand and the lock ring back over. In order for the system to recognise the new operating mode, the rechargeable battery must be briefly removed and reinserted (= reset).

7 Use

7.1 Terminal device on/off

The ON/OFF button is located directly under the palm of the hand.

To avoid unintentional activation, the ON/OFF button has to be pressed for about 1 second to activate and deactivate the prosthesis. (see fig. 5)

Press the ON/OFF button only when you want to maintain a grip for a longer period of time, while holding cutlery, writing or during longer passive pauses, for example, while travelling by rail or air. This prevents unintentional opening of the hand caused by inadvertent tensing of the muscles or extreme electrical interference. The operating time of the prosthesis battery is also extended. The hand switch is easy to operate with your other hand or by pressing against the leg or the edge of a table.

The ON/OFF button can be operated by pressing on the prosthetic glove.

Area	Function
Back of the hand	ON/OFF (--- FEHLENDER LINK ---)

7.2 Opening the terminal device in an emergency

The integrated slip clutch allows the terminal device to be opened regardless of the active control signals.

8 Storage

If the terminal device will not be used, store the terminal device in the open position to protect the thumb sensor and mechanism.

9 Cleaning

- 1) Switch off the product before cleaning.
- 2) Clean the product with a damp cloth and mild soap when needed.
Make sure that no liquids get into the product and product components.
- 3) Dry the product with a lint-free cloth and allow it to air dry fully.

INFORMATION

When using a prosthetic glove, note the cleaning information in the instructions for use for the prosthetic glove.

10 Maintenance

To prevent injuries and maintain the quality of the product, it is recommended that regular maintenance (service inspection) is performed every 24 months by an authorised Ottobock service centre.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

11 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

11.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

11.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The full text of the EU Declaration of Conformity is available at the following Internet address: <http://www.ottobock.com/conformity>

12 Technical data

Environmental conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85 % relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90 % relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95 % relative humidity, non-condensing

Reference number	8E38=6*, 8E39=6* (size 7)
Opening width	79 mm
Proportional speed	15–300 mm/s
Proportional gripping force	0–90 N
Weight (including system inner hand)	355 g
Operating voltage	Approx. 7.4 V
Static current	1 mA
Lifetime	5 years

Rechargeable battery reference number	757B20	757B21	757B3-5=1	757B3-5=3	757B3-5=5
Temperature range for charging [°C]	+5 to +40				
Capacity [mAh]	900 ¹ /950 ₁	680 ¹ /800 ₁	600	1150	3450
Rechargeable battery lifetime [years]	2				
Product behaviour during the charging process	The product is non-functional				
Operating time of the product with fully charged battery [grip cycles]	Approx. 2500–3000	Approx. 2000–2500 (depending on the battery capacity)	Approx. 1000–2000	Approx. 4000	Approx. 10,000
Charging times (rechargeable battery fully discharged) [hours]	Approx. 3.5	Approx. 3	Approx. 2.5		Approx. 3
Nominal voltage [V]	Approx. 7.4				

Rechargeable battery reference number	757B20	757B21	757B3-5=1	757B3-5=3	757B3-5=5
Battery technology	Li-ion		Li-poly	Li-ion	
Approved battery chargers	757L20		757L35		

¹ See labelling on rechargeable battery

13 Appendices

13.1 Symbols Used

	In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with regulations in your country can be harmful to health and the environment. Please observe the instructions of your national authority pertaining to return and collection procedures.
	Serial number (21)YYYYWWNNNN YYYY – year of manufacture

	WW – week of manufacture NNN – sequential number
	Medical device
	Manufacturer
	Declaration of conformity according to the applicable European directives
	UDI number (Unique Device Identifier)

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