1C50 Taleo, 1C53 Taleo LP

DE	Gebruuchsanweisung	6
EN	Instructions for use	13
FR	Instructions d'utilisation	21
IT	Istruzioni per l'uso	29
ES	Instrucciones de uso	37
PT	Instruções de utilização	45
NL	Gebruiksaanwijzing	53
SV	Bruksanvisning	61
DA	Brugsanvisning	68
NO	Bruksanvisning	75
FI	Käyttöohje	82
PL	Instrukcja użytkowania	90
HU	Használati útmutató	98
CS	Návod k použití	105
RO	Instrucţiuni de utilizare	113
HR	Upute za uporabu	121
SL	Navodila za uporabo	128
SK	Návod na použitie	135
BG	Инструкция за употреба	143
TR	Kullanım kılavuzu	151
EL	Οδηγίες χρήσης	159
RU	Руководство по применению	167
UK	Інструкція щодо використання	176
JA	取扱説明書	185
ZH	使用说明书	192
KO	사용 설명서	198

ISO 10328 - P n - m kg*)


*) Body mass limit not to be exceeded!

Stiffness	Weight [kg]	Label text
1	51	ISO 10328-P3-51-kg
2	58	ISO 10328-P3-58-kg
3	67	ISO 10328-P3-67-kg
4	77	ISO 10328-P4-77-kg
5	88	ISO 10328-P4-88-kg
6	100	ISO 10328-P5-100-kg
7	115	ISO 10328-P5-115-kg
8	130	ISO 10328-P6-130-kg
9	150	ISO 10328-P7-150-kg

*Actual test load reflect maximum body mass.

1

Bench alignment, transtibial

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INFORMATION

Date of last update: 2026-02-11

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- Please keep this document for your records.

The 1C50 Taleo and 1C53 Taleo Low Profile prosthetic feet are suitable for walking on various surfaces and at a wide cadence range.

Carbon and polymer spring elements enable perceptible plantar flexion, a natural rollover and high energy return. An exchangeable heel wedge provides absorption at heel strike.

2 Safety-related information

CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, treat the product with care.
- ▶ Observe the information on product combinations in the instructions for use for the products.
- ▶ Observe the maximum lifetime of the product (see page 15).
- ▶ Observe the product's operating conditions (see page 14).
- ▶ Do not expose the product to prohibited environmental conditions (see page 15).
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not use the product if it is damaged or if you are not sure whether it is working properly. Have the product inspected, repaired or replaced by the manufacturer or a specialist workshop.
- ▶ The product is intended for use by one person only and may not be reused by other persons.

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover) or delamination (detachment of layers) of the spring are indications of loss of functionality. Unusual noises can indicate a loss of functionality.

3 Intended use

3.1 Indications for use

The product is to be used **solely** for lower limb exoprosthetic fittings.

3.2 Operating conditions

The product is recommended for:

- Mobility grade 2: restricted outdoor walker
- Mobility grade 3: unrestricted outdoor walker
- Mobility grade 4: unrestricted outdoor walker with particularly high demands

The table below shows the suitable spring stiffness of the prosthetic foot for the user's body weight and activity level.

Spring stiffness relative to body weight and load level		
Body weight [kg]	Normal load	High load
up to 51	1	2
52 to 58	2	3
59 to 67	3	4
68 to 77	4	5
78 to 88	5	6
89 to 100	6	7
101 to 115	7	8
116 to 130	8	9
131 to 150	9	–

3.3 Environmental conditions

Storage and transport
Storage temperature: -20 °C to +60 °C, relative humidity: 20 % to 90 %, no vibrations or impacts
Allowable environmental conditions
Temperature during use: -10 °C to +45 °C
Chemicals/liquids: Fresh water, salt water, perspiration, urine, soapsuds, chlorine water
Moisture: Submersion: max. 1 h in 3 m depth, relative humidity: no restrictions
Solids: Dust, occasional contact with sand
Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and tear and damage (see page 19).
Prohibited environmental conditions
Solids: Highly liquid-binding particles (e.g. talcum), continuous contact with sand
Chemicals/liquids: Acids, continuous use in liquid media

3.4 Product combinations

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

3.5 Lifetime

Prosthetic foot

Depending on the user's activity level, the maximum lifetime of the product is **3 years**.

Footshell, protective sock

The product is a wear part, which means it is subject to normal wear and tear.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Prosthetic foot	–
1	Heel wedge set	2F50=*
1	Protective sock	Up to size 24: 2U3=1-7 From size 25: 2U3=2-7 From size 28: SL=Spectra-Sock-7

5 Accessories

Spare parts/accessories	
Designation	Reference number
Footshell	2C8=*, 2C15=*
Connection cover	For 2C15=* (normal): 2C19=*
	For 2C15=* (slim): 2C20=*
	For 2C8=* (normal): 2C21=*
Screw cover for 1C50	2F51=*
Screw cover for 1C53	2F52=*

6 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- Observe the alignment, assembly and adjustment instructions.

6.1 Applying/removing the footshell

CAUTION! Always use the prosthetic foot with a protective sock and footshell.

- 1) Pull the protective sock over the prosthetic foot to protect the prosthetic foot against abrasive particles and prevent noises in the footshell.
- 2) Apply the footshell as described in the footshell instructions for use.

Removing the footshell

- Remove the footshell as described in the footshell instructions for use.

6.2 Bench Alignment

INFORMATION: Use the prosthetic foot only with an installed heel wedge.

INFORMATION: If the user wants more stability: Move the prosthetic foot approx. 1 cm further in the anterior direction. This lengthens the forefoot and increases the knee-securing effect.

- **Illustration of bench alignment:** see fig. 1

- > **Required materials and tools:** 662M4 goniometer, 743S12 heel height measuring device, 743A80 50:50 gauge, alignment apparatus (e.g. 743A220 PROS.A. Assembly)
- > **Transfemoral prosthesis or knee disarticulation prosthesis:** Observe the instructions for use for the prosthetic knee joint.
- 1) Position the prosthetic foot in the alignment apparatus.
 - **Heel height: effective heel height of everyday shoe + 5 mm**
- 2) **Sagittal plane:** Align the prosthetic foot with the alignment reference line.
 - Position the alignment reference line on the markings of the footshell.
- 3) **Frontal plane:** Align the prosthetic foot with the alignment reference line.
 - Position the alignment reference line on the markings of the footshell: **external rotation approx. 5°**
- 4) Install the structural components between the prosthetic foot and prosthetic socket.
- 5) **Sagittal plane:** Align the prosthetic socket with the alignment reference line.
 - Determine the centre of the prosthetic socket proximally and distally with the 50:50 gauge and mark the centre line.
 - Mark the socket reference point on the centre line: **For transtibial prostheses, at the height of the centre of the patella.** The alignment reference line runs vertically through the socket reference point.
 - Adjust the socket flexion by rotating around the socket reference point: **individual residual limb flexion + 5°**
- 6) **Frontal plane:** Align the prosthetic socket with the alignment reference line.
 - Take the abduction or adduction of the residual limb into account.

6.3 Static Alignment

► Illustration of static alignment: see fig. 2

- Ottobock recommends checking the alignment of the prosthesis using the 3D L.A.S.A.R. Posture and adapting it as needed.

Optionally, the alignment recommendations (transfemoral modular leg prostheses: **646F219=***, transtibial modular leg prostheses: **646F336=***) can be requested from Ottobock as needed.

6.4 Dynamic Trial Fitting

The optimal gait pattern is developed during dynamic fitting. The prosthesis is aligned optimally in the frontal plane and the sagittal plane for this purpose.

► TT fittings: Make sure that physiological knee movement in the sagittal and frontal planes is achieved after the leg begins to bear weight following the heel strike. Avoid medial movement of the knee joint.

- If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction.
- If the knee joint moves in the medial direction in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.

► After completing the dynamic fitting and the walking exercises: Remove the plastic part that protects the area under the pyramid from scratching during the fitting.

6.4.1 Optimising the heel characteristics

NOTICE

Altering the prosthetic foot or footshell

Premature wear resulting from damage to the product

- Do not alter the prosthetic foot or footshell.

INFORMATION: Use the prosthetic foot only with an installed heel wedge.

The behaviour of the prosthetic foot at heel strike and during heel contact in the mid-stance phase can be adapted by replacing the heel wedge. Heel wedges of various hardness are included in the scope of delivery.

Heel wedge hardness: The colour of the heel wedge indicates the hardness (see fig. 5). Ottobock recommends starting with the preinstalled heel wedge.

- 1) Slightly bend open the spring at the heel of the prosthetic foot and remove the installed heel wedge.

- 2) Align the other heel wedge so the Ottobock lettering is right side up and the tip points in the anterior direction.
- 3) Insert the heel wedge into the prosthetic foot (see fig. 6). First insert the tip of the heel wedge into the holder ❶ , then push down at the back ❷ .

7 Cleaning

- > **Permissible cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)
- 1) **NOTICE! To avoid product damage, only use the permissible cleaning agents.**
Clean the product with clear fresh water and a pH neutral soap.
- 2) Remove dirt from water drainage contours using a toothpick and rinse.
- 3) Rinse the soap away with clear fresh water. Rinse the footshell until all dirt has been removed.
- 4) Dry the product with a soft cloth.
- 5) Allow to air dry in order to remove residual moisture.

8 Maintenance

- **After the first 30 days of use:** Perform a visual inspection and functional test of the product and the prosthesis.
- **Every 12 months:** Check the function of the product and the prosthesis and inspect for signs of wear.

9 Disposal

Do not dispose of the product in household waste. Improper disposal is harmful to health and the environment. Observe the instructions of the local authority regarding return, collection and disposal.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

10.3 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

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