

7. TECHNICAL DATA



- Size: L 150 mm x ø 61 mm
- Bacteria Retention: 7 log steps *Brev. dim.*
- Pore size: 0.2 µm (microfiltration)
- Operating time*: up to 3 months
- Chlorine resistance**: ≤ 10 ppm
- Flow rate: 22 l/min (initially, at 5 bar, 20 °C)
- Max. operating pressure***: 5 bar
- Operating temperature****: ≤ 40 °C
- Connection: Quick-release coupling

* The service life is dependent on the user-specific water quality, pressure conditions, frequency of use and requisite flow rate.

** Continuous dosage of ≤ 10 ppm over the useful life, short-term (1h) high dosage (400,000 ppm h) for chemical disinfection

*** The filter must not be connected to low-pressure systems (e.g. boiler).

**** Application at higher temperatures only after consultation with Aqua free.

8. STORAGE

The Baclyser® filter must be stored correctly, i.e. dry and hygienically, as well as protected from frost and heat.

9. PRODUCT MARKINGS



Manufacturer



Do not use if package is damaged



Medical device



Expiry date



Batch code



Store in a dry place



Consult instructions for use



Non-sterile

Aqua free is EN ISO 13485 certified.



Aqua free Solutions Ltd.

The Media Centre, 7 Northumberland Street
Huddersfield, HD1 1RL
Great Britain
E-Mail : info@aqua-free.co.uk
Phone: +44 (0)1484 483045



Headquarters:

Aqua free GmbH

Winsbergring 31
22525 Hamburg
Germany
E-mail: info@aqua-free.com
Phone: +49 (0)40 46899990
Fax: +49 (0)40 46899999

Aqua free

Instructions for use

Baclyser® WG



- Inline hollow-fibre membrane filter for birthing pool water hygiene



Instructions-for-use-Baclyser WG-03/2018-V002



1. Product

The Baclyser® WG Filter is a CE-marked medical device that has been developed specifically for inline installation in the water supply pipes of birthing pools. The filter offers effective protection from waterborne bacteria utilising a hollow fibre membrane with a bacterial retention of 7 Log (*Brev. Dim.*).

The couplings with integrated flow stop valves offer safe, secure connection and quick filter changes.

2. Preparing for installation

- Turn off the water supply of the birthing pool prior to installation of the filter.
- Open the pool supply connection / filter shaft.
- Disconnect the water pipe at the point where the filter is to be inserted. *Caution!* Residual water can leak out of the pipes.
- Fit the appropriate couplers (see accessories) (seal using Teflon tape or similar to avoid leaks). Avoid contamination.
- Turn on the water supply and check the connectors for leaks.




3. Installing the filter

Attention: Before installing the filter, ensure that the birthing pool pipes and outlets are spotlessly clean. This can be ensured by overhauling the unit and replacing any pipes and hoses etc.


Note: Disinfect your hands prior to fitting the filter and wear disposable gloves. Turn off the water supply. Disinfect the outer surfaces of the connecting pipework using an alcohol based surface disinfectant. Only use filters which are in original packaging.

- Take the filter out of the packaging. Ensure that the filter outlet does not get contaminated.
- Note the flow direction (arrow) of the filter during installation.

- Insert the filter outlet side into the prepared coupling (A).



- Insert the filter inlet side into the prepared coupling (B).



- Turn on the water supply once the filter is installed. Check the flow rate to the unit is adequate and that all connections are free from leaks.

4. Operation

- The filter must be changed max. 3 months after installation. If the water flow is so low as to be insufficient for operation the filter should be changed immediately.
- Note the next date for changing the filter on an adhesive label and stick it to the filter where it is clearly visible.
- The filter housing can be wiped clean on the outside using a standard surface disinfectant.

- We recommend turning off the water supply to the filter if the unit isn't being used for a long period of time and to check the system tightness at regular intervals.
- Please refer to "Technical data" for information on the max. operating pressure and water temperature.

In the event of a deviation from the specification, the product should be returned to Aqua free using the "Complaints Report" form (which can be downloaded from www.aqua-free.com).

5. Dismantling and changing the filter

Note: Disinfect your hands prior to fitting the filter and wear disposable gloves.

- Turn off the water supply.
- Release the quick coupler at the inlet side (B) (the flow stop valve in the filter connection and coupling is activated though some residual water may leak out).
- Then release the quick coupler on the outlet side (A). Attention: Some filtered residual water may drip out of the filter (collect in a bag or container).
- Dispose of the used filter in the appropriate way.
- Install the new filter as described in 3.

6. Accessories

The standard scope of supply does not include any accessories. Couplings with the required connecting threads for standard birthing pools offered by Apal, Vivipar, Febromed etc. can be ordered/supplied separately.

Details of other accessories are available on request.