MGH00460

Valid from / en cours de validité depuis le / en espera de validación desde el / in corso di validità dal / gültig seit dem

25/04/2019

**Product Sheet (EN) 2**

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**Eigenschaften (DE) 10**

PRODUCT SHEET			valid from 25/04/2019	
<b>IDENTIFICATION OF THE MEDICAL DEVICE</b>	Type	Rechargeable battery		
	Commercial designation	Batterie medical SAP AGILIA de Fresen 6V 2500mAh PHR4		
	Reference	MGH00460		
	EAN	3660766505286		
	Brand	NX		
	Compatible / Original battery	Compatible		
Packaging	Unitaire			
<b>RECOMMENDED USAGE</b>				
<i>Follow the instructions and recommendations specific to each model, using the technical instructions and document resources from the devices in which the battery is used</i>				
<b>Brands</b>	<b>Equipment</b>	<b>Models</b>		<b>PN</b>
Fresenius - Vial	Batterie médicale pompe	Injectomat Agilia		
<b>Identification</b>	<b>GENERAL TECHNICAL CHARACTERISTICS</b>	Chemistry		NiMH
		Type	AA 5S1P	
		IEC designation	5 HRH 15/49	
		Rated voltage	6V	
		Nominal capacity	2.5Ah	
		Internal resistance $\Omega$	125m $\Omega$	
<i>The voltage and the actual capacity in use can be affected by several factors, especially the temperature, the discharge current, the pack's history (ex:use, storage), etc</i>				
<b>ELECTRICAL CHARACTERISTICS</b>	<b>CHARGE</b>	Maximum charging voltage	8.5V	
		Standard charging current (15h)	250mA	
		Fast charging current (2,5h)	1000mA	
	<b>DISCHARGE</b>	Range of operating voltage	5.5V at 7V	
		Min tension in discharge	5V	
		Max discharge current	1000mA	
		Lifespan 80% DOD (0,5 C)	500 Cycles	
	<b>MAINTENANCE</b>	Frequency of maintenance charges at 20°C	3 Months	
	<b>CONTROL ELECTRONICS</b>	Electrical protection	Yes	
		Low voltage power cut	No	
High voltage power cut		Yes		
Max power cut voltage		Yes		
<i>These devices not only designed to protect the pack in case of an equipment failure. They must not be used to control the discharge. The protection circuits have a response time of a few milliseconds.</i>				
<b>MECHANICAL CHARACTERISTICS</b>	Dimensions (+/- 2mm)	Length	72mm	
		Width	14.5mm	
		Depth	52.5mm	
	Weight (+/- 5g)	143Gr.		
	Mechanical protection	Shrink sleeve		
	Wire length (+/- 10mm)	60mm		
Terminal	JST PHR4			
<b>CONDITIONS OF USE, STORAGE, AND TRANSPORT</b>	<b>CONDITIONS OF USE</b>	Charging temperature	0 at 40°C	
		Discharge temperature	0 at 50°C	
	<b>CONDITIONS OF STORAGE</b>	Storage temperature	-20 at 30°C	
		Level of humidity	65.00 %	
		Max storage time	2 Years	
	<b>TRANSPORT</b>	UN code	Class 9	
		ADR/RID classification	Class 9	
		IMDG classification	Class 9	
IATA classification		Class 9		

<b>INSTRUCTIONS</b>	<b>COMMISSIONING</b>	<ul style="list-style-type: none"> <li>• Check the batteries and the connectors: wires not damaged, battery not swollen, burnt smell, oxidation of the connectors, leak...</li> <li>• Respect the polarity</li> <li>• Do a full charge with the adequate charger before the first use</li> </ul>
	<b>CHARGE</b>	<ul style="list-style-type: none"> <li>• Use an adequate charger</li> <li>• The battery is warmer during the charge: during the first charge, check that the battery's temperature stays in the temperature operating ranges.</li> <li>• In case of an abnormal heating, stop the charge by unplugging the charger within the realms of possibility, remove the battery from the equipment de l'équipement, have the equipment, the battery and the charger checked by a technician.</li> </ul>
	<b>CASE OF NON-WATERPROOF BATTERIES</b>	<ul style="list-style-type: none"> <li>• It is normal to observe a release of gas during the charge and use. Do not smoke. Place in suitable premises.</li> <li>• Open batteries need regular maintenance carried out by a qualified technician.</li> </ul>
	<b>CASE OF LITHIUM ION BATTERIES</b>	There is a fire hazard with lithium ion batteries in the following cases: overload, short circuit, charge and use outside the voltage and temperature ranges.
	<b>WARNINGS</b>	<ul style="list-style-type: none"> <li>• Read the instructions of your device.</li> <li>• Only use in compatible devices.</li> <li>• Respect the load and storage conditions.</li> <li>• Do not use if the battery is damaged, do not burn, do not pierce, do not dismantle or modify. The protection circuits protect the battery and the equipment: do not deactivate them.</li> </ul>

<b>Marking</b>	<b>CE</b>
----------------	-----------

<b>EXPLANATION OF SYMBOLS</b>	<b>LABELLING</b>	
		Catalogue reference
		Lot number
		manufacturer's address
		Conform to the directive CE 93/42 DM class I
		To recycle in a suitable salvage and recycling structure
	Read the product sheet and the instruction manual	

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