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ENGLISH

Thank you for purchasing this **LRV7 laboratory food pH electrode**.

For best results with your electrode:

- **read** these operating instructions carefully,
- **comply** with the precautions for use.

The CE marking indicates conformity with European directives, in particular LVD and EMC.

The rubbish bin with a line through it indicates that, in the European Union, the product must undergo selective disposal in compliance with Directive WEEE 2012/19/UE. This equipment must not be treated as household waste.

DELIVERY CONDITION

The food pH electrode is delivered in a cardboard box with:

- a storage bottle,
- an S7-BNC cable approximately 1 metre long or a S7-tight 8-point DIN cable approximately 1 metre long, depending on the order,
- a shim,
- a bottle of storage solution.

For the accessories and spares, consult our web site:

www.chauvin-arnoux.com

USE


This laboratory electrode is used to measure the pH of liquid, semi-solid, and solid foods.

- Unscrew the protective cap.
- Connect the electrode to the pH-meter using the cable supplied with it, S7-BNC or S7- DIN 8. Screw the S7 connector of the cable to the S7 connector of the electrode.
- Remove the storage bottle.
- Calibrate the electrode with buffer solutions, rinsing with de-ionized water between measurements.
- To make a measurement, dip the end of the electrode into the solution to be measured; make sure than the junction is immersed. With a semi-solid or solid food, insert the tip of the electrode into the food as far as the junction. To avoid damaging the electrode, drill a hole in a solid food before inserting the electrode.
- Wait for the measurement to stabilize.
- Rinse the electrode with de-ionized water or a suitable cleaning solution between measurements.
- When the measurements are over, put the electrode back in the storage bottle.
- Unscrew the connecting cable from the pH-meter and screw the protective cap back on.

To avoid damaging the electrode, do not make measurements in the any of the following solutions:

- hydrofluoric acid solution,
- acid and alkaline solutions at high concentrations (> 2 mol/L),
- clogging solution (resin, paint, etc.).

PRESENTATION

CHARACTERISTICS

Reference: polymer.
Junction: ceramic and single pore.
pH measurement range: 0 to 14.
Temperature range: 0 to 60°C.

Body of the electrode: PVC
Dimensions: 132 x Ø 16 mm
Weight: approximately 30 g
Connector: S7

CLEANING

The presence of crystallized white deposits on the electrode is a normal consequence of the evaporation of the storage solution. To eliminate a deposit, simply rinse the electrode with de-ionized water, then wipe dry with a clean, lint-free cloth.

Do not use aggressive and/or abrasive products (hydrofluoric acid, acetone, etc.). The cleaning operation must not scratch the glass membrane.

Contamination of the junctions is the commonest cause of measurement problems.

- Remove any impurities deposited on the surface of the glass membrane and on the junctions using de-ionized water, or a soft cloth if the deposit is large.
- The commonest cause of clogging is proteins (dairy products, meat). If rinsing is not enough, use a solution of pepsin + 0.1 mol/L HCl.

MAINTENANCE

Keep your electrode in its bottle in a clean 3 mol/L KCl electrolytic solution. Change it as soon as it becomes cloudy.

Do not keep an electrode dry or in distilled or de-ionized water. The membrane may be damaged and the life of the electrode shortened.

Keep the connector clean and dry. If necessary, clean it using a soft cloth wet with soapy water. Rinse using a damp cloth and dry quickly with a dry cloth or forced air.

Regularly check the condition of your electrode and the cleanliness of the junction and of the glass membrane.

WARRANTY

Electrodes are wear parts. The life of your electrode depends on how you use and maintain it.

Except as otherwise stated, our warranty is valid for **12 months** starting from the date on which the equipment was sold. Extract from our General Conditions of Sale, provided on request.

The warranty does not apply in the following cases:

- traces of corrosion on the inside of the connector;
- traces of clogging agents (glue, paint, resin, etc.) on the body, the glass ball or the junctions;
- breakage of the glass;
- damage to the body of the electrode (deformation, discoloration).

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