FR - Notice de fonctionnement
EN - User's manual
DE - Bedienungsanleitung
IT - Manuale d'uso
ES - Manual de instrucciones
SV - Bruksanvisningar
FI - Käyttöohje

 CHAUVIN
ARNOUX

XRPTST1

Électrode ORP combinée / T°
Combined ORP / T° electrode
Kombinierte Redox-Elektrode / T°
Elettrodo ORP combinato / T°
Electrodo de ORP combinado / T°
Kombinerad ORP-elektrod / T°
Yhdistelmäelektrodi / T°

Measure up

English	6
Deutsch	10
Italiano	14
Español	18
Svenska	22
Suomalainen	26

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ENGLISH

Thank you for purchasing this **XRPTST1 combined ORP electrode with integrated temperature measurement**.

- For best results with your electrode:
- **read** these operating instructions carefully,
- **comply** with the precautions for use.

This electrode is for use only with the C.A 10101 pH meter.

The CE marking indicates conformity with European directives, in particular LVD and EMC.

The rubbish bin with a line through it indicates that, in the European Union, the product must undergo selective disposal in compliance with Directive WEEE 2012/19/UE. This equipment must not be treated as household waste.

DELIVERY CONDITION

The ORP electrode is delivered in a cardboard box with:

- a shim,
- a storage bottle.

For the accessories and spares, consult our web site:

www.chauvin-arnoux.com

USE

The ORP electrode serves to measure the oxidation reduction potential of a solution.


ORP = Oxidation Reduction Potential

- Remove the storage bottle.
- Connect the electrode to the measuring instrument.
- Calibrate the electrode with a buffer solution, then rinse it with de-ionized water.
- Plunge the end of the electrode in the solution to be measured; make sure that the platinum and the junction are submerged.
- Wait for the measurement to stabilize.
- Rinse the electrode with de-ionized water between measurements.
- When the measurements are over, put the electrode back in the storage bottle.

So as not to damage the electrode, do not make measurements in any of the following solutions:

- hydrofluoric acid solution,
- acid and alkaline solutions at high concentrations (> 2 mol/L)
- clogging solution (resin, paint, etc.)

PRESENTATION

CHARACTERISTICS

Reference: Ag/AgCl gel
Junction: ceramic.
mV measurement range: ± 1999 mV.
Temperature measurement range: 0 to 60°C.
Temperature measurement: PT1000 probe.

Dimensions: 190 x \varnothing 18 mm
Body of the electrode: 122 x (dia.) 12 mm, polycarbonate
Length of the lead: 1 m
Weight: approximately 120 g
Connector: tight 8-point

CLEANING

The presence of crystallized white deposits on the electrode is a normal consequence of the evaporation of the storage solution. To eliminate a deposit, simply rinse the electrode with de-ionized water, then wipe it dry with a clean, lint-free cloth.

Do not use aggressive and/or abrasive products (hydrofluoric acid, acetone, etc.).

Contamination of the diaphragm is the commonest cause of measurement problems.

- Remove any impurities deposited on the metal and on the outside of the diaphragm with de-ionized water, or a soft cloth if the deposit is large.
- If rinsing is not enough, it is possible to use various chemical cleaning solutions, depending on the nature of the contamination:

Type of clogging	Cleaning solution
Silver sulphide	Solution of thiourea in 0.1 mol/L hydrochloric acid (HCl)
Silver chloride	Concentrated ammonia solution
Proteins	Solution of pepsin + 0.1 mol/L HCl
Other clogging	Ultrasound bath with water or 0.1 mol/L HCl

MAINTENANCE

Keep your electrode in its bottle in an electrolytic 3 mol/L KCl solution.

Do not keep an electrode dry or in distilled or de-ionized water. The membrane may be damaged and the life of the electrode shortened.

Regularly check the condition of your electrode and the cleanliness of the diaphragm and of the glass membrane.

WARRANTY

Electrodes are wear parts. The life of your electrode depends on how you use and maintain it.

Except as otherwise stated, our warranty is valid for **12 months** starting from the date on which the equipment was sold. Extract from our General Conditions of Sale, provided on request.

The warranty does not apply in the following cases:

- traces of corrosion on the inside of the connector;
- traces of clogging agents (glue, paint, resin, etc.) on the body, the glass ball or the ceramic junction;
- breakage of the glass;
- damage to the body of the electrode (deformation, discoloration).

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