Intended use

The ArthroFit brace is portable, non-sterile, non-invasive medical device intended for long-term immobilization / support of the knee joint while injured. The device enables the controlled mobilization of the joint by setting limits for flexion and extension. It supports moderate to severe levels of osteoarthritis.

Indications

- ✓ Conservative treatment of advanced unilateral osteoarthritis in combination with complex band instability;
- ✓ Post-operative pressure relief & compartmental protection;
- ✓ Functional pre- or post-operative care of knee instabilities.

Contraindications

If undue pain, unpleasant sensation or allergy develops while you are wearing the brace, please take it off immediately and consult your healthcare professional.

Fitting instructions

- Open all the belts and pull the numbered end from the belt eyelets
Tip: to prevent the belt from getting stuck, it is advisable to close the Velcro tab in the attachment.
- Sit on the edge of a chair and spread the ArthroFit over the aected leg. Slightly bend your leg. Place the ArthroFit in such a position that the wedges visible from the front are at the level of the upper edge of the kneecap. Please take care that the ArthroFit is not the wrong way round on your knee.
- Tighten the lowest belt on the front side (belt 1) and the top belt at the back (belt 2). Press the joints gently backwards.
- Tighten the lower belt at the back of the thigh (belt 3). Check that the belt is fitting properly. The lateral joints should be positioned slightly backwards (more exactly: behind the lateral middle line of the knee).
- Tighten the front upper lower leg belt (belt 4), followed by the back lower leg belt (belt 5) and nally the front thigh belt (belt 5) and then the front thigh belt (belt 6).
- Take a few steps or bend your knee a few times to check the fit. If necessary, the belt can be readjusted in the same sequence to guarantee a good fit.

Setting the flexion/extension limits

(only to be done by an orthopaedic technician)

The orthosis is equipped with a 0° extension wedge. The orthosis includes additional wedges.

- ⚠ The extension and flexion limits should only be limited in accordance with the instructions from the responsible doctor.

To avoid hyperextension, extension wedges must be used. The wedges exactly specify the radius of motion of the orthosis. The situation in soft tissue may cause compression, making it necessary to use a larger wedge to restrict the movement of the knee to within the desired radius. Select the required wedge.

Extension wedges are marked with an E and flexion wedges with an F. Insert the wedge into the joint and lock it by clamping it with help of the frame ends. The wedge then slides automatically into the correct position and does not rotate. Tighten the wedge using the screw. It is sufficient to do this by hand. No great force it necessary.

- ⚠ The medial (inside) and lateral (outside) limit settings must always be the same.

Setting the varus-valgus decompression

Get the responsible technician to give you a thorough explanation about how to adjust the ArthroFit. This will help you to master the handling rapidly and to benefit from optimal support from the orthosis. Sit down and adjust the varus-valgus adjustable joints to the neutral position, using an Allen wrench. You have reached the neutral position when the padding of the ArthroFit joint is covering the surface but producing little or no pressure. Stand up and set the tension to a quarter cycle (on both adjusting screws). This is a good initial setting for the first time you wear the ArthroFit. If you want, you can adjust the ArthroFit to the breadth of you knee with the various joint pads enclosed. Take care that the ArthroFit is sitting properly. For conservative therapy (without an operation), the ArthroFit should not be worn for more than two hours per day. You can then wear the ArthroFit for longer periods, depending on how well you tolerate it. It is most important to wear the ArthroFit for activities that normally cause pain. By using the Allen wrench, the decompression can gradually be adjusted over time, or during a course of a day. For postoperative therapy, the ArthroFit must always be word when there is stress. Please ask your doctor for instructions.

- ⚠ Extreme settings have no therapeutic value, they only cause additional discomfort

Care instructions

The inner pads and the condyle pads can be washed by hand with a mild soap/detergent.

- ⚠ Do not wring out. Do not iron.
- ✖ ✖ ✖ ✖

When used normally the brace is maintenancefree. However, after using the brace in water, the hinges can be lubricated using Teflon™ spray. Check regularly that the joint screws are secure, tighten them if necessary.

Storage instructions

Store in dry place and do not expose to direct sunlight and moisture.

Warranty

Currently valid legal conditions applicable.

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⚠ **Important notes**

- The brace is recommended for one user only.
- The brace is not intended to be applied on breached skin.
- Respect the cleaning instructions.
- Remove device before showering.

If allergic skin rashes occur, please consult your physician. Any serious incidents that occurred in relation to this device should be reported to the manufacturer and the competent authority of the member state in which the user is established.

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