KINETEC 6080™

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1.	General information.....	2
1.1.	Definition	2
1.2.	Clinical Benefits	2
1.3.	Indications	2
1.4.	Contraindications	2
2.	Warning and safety instructions	3
3.	Regulatory compliance	4
4.	Unpacking and packing.....	4
5.	Installing the device	4
6.	Description.....	5
7.	Assembly	5
8.	Electrical connection: safety first	5
9.	Starting the unit	6
10.	START / STOP / REVERSE function	6
11.	Procedure to stop the machine.....	6
12.	Adjusting the SPEED	6
13.	Adjusting the limits.....	6
14.	Use of the KINETEC® hygienic pads.....	7
15.	Setting up the patient	7
16.	Possible movements	8
17.	Product information.....	9
17.1.	Maintenance	9
17.2.	Breakdown	9
17.3.	Cleaning	9
17.4.	Disposal and recycling.....	9
17.5.	Technical specifications	10
17.6.	Symbols used	10
17.7.	Warranty	10
17.8.	Information about electromagnetic compatibility (EMC)	11
17.9.	List of cables.....	11
17.10.	Recommended separation distances.....	11
17.12.	Electromagnetic emissions.....	12
17.13.	Magnetic and electromagnetic immunity	12
17.14.	Electromagnetic immunity, radio-frequency portable equipment.....	13

1. General information

1.1. Definition

The Kinetec 6080™ is used for passive mobilization of the Elbow, providing the following motion:

- Flexion-extension of the Elbow from 0° to 135°.
- Simultaneous 90° pronation and 90° supination with flexion-extension

1.2. Clinical Benefits

- Effectively breaks the vicious circle of: trauma → immobility → effusion → atrophy.
- Improves the nutrition of joint surfaces.
- Speeds the recovery of post-operative range of motion.
- Speedy recovery of the motor pattern for the limb operated on.
- Promotes joint cartilage healing.
- Provides immediate post-operative continuous passive motion.
- Reduces hospitalization time.
- Reduces the need for pain medication.
- Visual feedback for users.
- Realization of postures (resting, stretching muscles).).

1.3. Indications

- Intra-articular fractures of the elbow
- Metaphyseal fractures, osteosynthesised in the elbow area
- Arthrolysis for post-traumatic stiffness with limitation of elbow movement
- Elbow prostheses
- Synovectomies

1.4. Contraindications

- Rheumatoid arthritis in the inflammatory phase,
- Algodystrophy in the inflammatory phase (hyperpainful),
- Para-osteo-arthroplasty,
- Unhealed infected wounds,
- Bone cancer,
- Infectious arthritis,
- Deformed joint surfaces,
- Paralysed limbs (atonic or spastic),
- Non-stabilised fractures.
- The machine is not suitable for patients over 1.95 m (6'7") or under 1.45m (4'7") tall.

2. Warning and safety instructions

- WARNING: The machine must be installed and commissioned according to the information provided in this manual.
- WARNING: If you need any assistance in the assembly, use or maintenance of the device, please contact your KINETEC® distributor.
- WARNING: The practitioner determines the protocol and ensures its proper implementation (settings, session duration and frequency of use).
- WARNING: Run a cycle with the device unloaded before installing the patient on the machine.
- WARNING: For optimum safety, always give the hand control to the patient before starting the system. The patient must know the start/stop/reverse function on the hand control (see page 6).
- WARNING: Danger, risk of explosion: Do not use the machine with anaesthetic gas or in an environment that is rich in oxygen.
- WARNING: For Type B Class I devices, and to avoid all risks of electric shock, the machine should only be connected to a power supply that has protective earthing, see page 5.
- WARNING: Before using this machine, always check that the electrical socket is in good condition and is suitable for the splint power supply cord. Only use the original cable supplied with the machine. Check that the cables remain free around the device so that they do not get damaged.
- WARNING: Before using this machine, always check that the machine is not damaged, in particular the protective housings.
- WARNING: Please do not touch the fixed or moving parts while the unit is running: pinching or crushing risk. Keep children and pets away from the machine.
- WARNING: Modifying the machine in any way is strictly forbidden.
- WARNING: Only the accessories, spare parts and supplies described in this manual should be used with this machine.
- WARNING: Do not connect the device to other devices not described in this manual.
- WARNING: If unforeseen events or malfunctions occur, please contact your KINETEC® distributor.
- WARNING: In the event of reciprocal, electromagnetic or other interference with other devices, move the device.
- WARNING: The device must not be used next to other devices or stored with the latter because this can cause incorrect operation. If such use is necessary, the device and other machines should be monitored to ensure normal operation.
- WARNING: The use of accessories and cables other than those specified or provided by KINETEC SAS can increase electromagnetic emissions or reduce immunity of the device and cause inappropriate operation.
- WARNING: Portable radio-frequency devices should not be used (including peripherals such as coil cables and external coils) within 30 cm (12 inches) of any part of the device including cables specified by Kinetec. Should the opposite occur, the performance of these devices could be altered.

3. Regulatory compliance

The Kinetec 6080™ device was designed, manufactured and distributed in accordance with the requirements of Directive 93/42/CEE applicable to medical devices. In this respect, the Kinetec 6080™ device is CE marked.

The Kinetec 6080™ device complies with the following standards:

- IEC 60601-1 relating to requirements and tests in terms of electrical safety,
- IEC 60601-1-2 relating to electromagnetic compatibility,
- IEC 60601-1-11 which defines requirements for electro-medical devices used in the home care environment.

The Kinetec 6080™ device also complies with:

- the essential requirements of Directive Machine n°2006/42/CE,
- directive 2011/65/UE relating to the limitation of the use of certain hazardous substances in electrical and electronic equipment,
- directive 2002/96/CEE relating to electrical and electronic equipment waste (DEEE).

4. Unpacking and packing

Unpacking

When you unpack the machine, don't forget that you may need to pack it up again. We recommend that you keep the packaging materials, boxes and plastic bags.

Recommendations for plastic bags: do not put them over the head as there is a risk of suffocation, and keep them out of the reach of children.

Be careful with small-sized pieces: they could be swallowed by a child.

Be careful with cables and wires: risk of strangulation.

Packing

To prevent any problems when the machine is transported, only pack it using its original packaging.

- Stop the unit at 5° of flexion.
- Consult the package diagram of the device.

5. Installing the device

The device Kinetec 6080™ machine is designed to be used in hospitals, clinics or in doctors' offices.


The machine is delivered with its own chair. It must not be used with another chair.

The machine must be installed on a flat surface that is wide enough to accommodate the entire device; you should leave a free space of 1 m (39 in.) on each side.

6. Description

The device Kinetec 6080™ consists of the following components:

- A. Upper arm support
- B. Forearm support
- C. Telescopic column
- D. Motor
- E. Swiveling handle
- F. Control box
- G. Base

7. Assembly

- Attach the telescopic column (C) to the red base (G) with the lock nut and the washer.
- Place the wire basket (1) over the telescopic column (C).
- Attach motor (D) to the telescopic column with the toothed link (2).
- Place the control box (F) on the wire basket (1).

IMPORTANT: To avoid the risk of instability of the Kinetec 6080™ when moving it, we recommend to carry the control box separately (F) of the rest of the device.


8. Electrical connection: safety first

The Kinetec 6080™ is a Type B Class I machine.

Before plugging in the machine, make certain that the available voltage corresponds to the voltage indicated on the back of the control box (100-110V~ or 240V~ 50/60Hz).

Connect the cords to the back of the control box as follows:

- Attach the hand control cord (3) to terminal (4).
- Attach the cord (5) from the control box (6) to the motor (D).
- Attach the power cord (7) to the control box (8).

IMPORTANT: Cords (3) and (5) must be securely fastened. Check to be sure that the plug fits the socket being used

IMPORTANT

For Type B Class I devices, and to avoid all risks of electric shock, the machine should only be connected to a power supply that has protective earthing.

To connect the power supply, only use the original cable supplied with the machine.

Check that the cables remain free around the device so that they do not get damaged.

Check that the machine is not damaged, in particular the protective housings.

9. Starting the unit

Turn on the power switch (10). The green light will turn on.


10. START / STOP / REVERSE function

As with all KINETEC® systems, the device Kinetec 6080™ is equipped with a START/STOP/REVERSE function.

- Press the button once, the machine starts.
- Press the button twice, the machine stops
- Press the button a third time, the machine reverses direction.

IMPORTANT

for maximum safety, the hand control should always be given to the patient.

11. Procedure to stop the machine

To stop the machine's movement: press the hand control switch (11).

To switch the unit off: press the ON / OFF switch (10).

12. Adjusting the SPEED

The SPEED is set by turning the button (12). Minimum speed is all the way to the left;


Maximum speed is all the way to the right.

The speed can be modified while the machine is running or stationary.

13. Adjusting the limits

The physician positions the cursors (13) and (14) on the graduated scale and locks them with the safety nuts (15) and (16).

The graduated scale placed on the screw drive casing gives only an approximate measurement.

14. Use of the KINETEC® hygienic pads

The KINETEC® patient pad kit is designed for rapid fitting, optimal hygiene and maximum patient comfort.

FOR OPTIMAL HYGIENE, A NEW SET OF PADS SHOULD BE USED FOR EACH PATIENT.


Place the pads (17) and (18) on the adhesive parts of the forearm and arm support, with the sponge side up.

The straps (19) and (20) are used for holding patient's arm on the machine. Do not use strap (19) when machine is used for pronation-supination.

For glove use, see the drawing.

Part number to order pads: 4650000646

Part number to order glove: 4650000654


CLEANING:

- Disinfecting the straps: Wash at 30°C, using a disinfectant solution during the rinse cycle. Example of product that can be used: Solution "Bac linge" at 0.125 % or "Souplanios" at 0.125% from ANIOS Laboratories. A complete list of distributors in your country is available on request.

15. Setting up the patient

- Make sure that clean pads are on the machine.
- Place the patient comfortably in a sitting or reclining position.
- Place arm into the machine, in the maximum tolerable flexion.
- Adjust the height of the machine on the telescopic column (C) and the incline of the toothed link (2).
- Adjust the joint axis of elbow and the joint axle of the machine.
- Lock casters of base (G).
- Place the pronation-supination hand grip (E) in the palm of hand with hand in supine position.
- Fix the top of arm with the strap (20).

Flexion-Extension of the elbow (H).

The Kinetec 6080™ allows flexion-extension of the elbow from 0° to 135° with fixed PRONO-SUPINATION (H). Choose the hand position between full PRONATION to 90° and full SUPINATION to 90°.

- Position the cleat (25) on the LEFT of its slide (Fig.9•1).
- Lock the cleat (25) with the locking nuts (21 and 22; Fig.9•1).
- Wind the rope according to the indicated direction on the primary pulley (23), and put the rope into the flexion pulley (24) and into the cleat (Fig. 9•2 or 9•3).
- Tighten the rope to bring the hand to the chosen position.
- Cinch the rope into the cleat (25)..

Fig.9•2: representation for a left arm; Fig.9•3: representation for a right arm.

AT THE TIME OF FLEXION-EXTENSION, THE HAND-GRIP WILL STAY AT THE INITIAL ANGLE.

PRONATION-SUPINATION in conjunction with Flexion-Extension of the elbow (H+I).

PRONATION-SUPINATION in conjunction with Flexion-Extension of the elbow (H+I)

The Kinetec 6080™ allows PRONO-SUPINATION 90° to 90° in conjunction with flexion-extension of the elbow from 0° to 135° (H+I).

- Position the cleat (25) on the RIGHT of its slide (Fig.9•4).
- To obtain pronation-supination associated less than 180°, move the cleat to the left on the slide.
- The more the cleat is moved to the left, the less pronation-supination occurs.

Fig.9•5: representation for a left arm; Fig.9•6: representation for a right arm.

17. Product information

17.1. Maintenance

Direct a jet of compressed air into the protection screw casing.

Do not clean the screw with cotton. Move plug that holds up the nut, and grease the screw.

The motor does not need maintenance. However, compressed air can be used to remove carbon buildups. After 4,000 hours of operation, returning the machine to Kinetec SAS for maintenance is recommended.

When the machine is no longer in working condition, please return it to us, together with its accessories, for destruction.

17.2. Breakdown

If you need a spare parts catalog and/or service manual, please request it from your KINETEC® distributor.

- After connecting the power cord on the Kinetec 6080™ machine and switching ON the machine (10), if:
 - o The green light is off,
 - o Try another electrical device to make sure that the power supply is uninterrupted.
 - o Check the fuse of the socket (9) by replacing it with a fuse of the same type and value.
 - o If the green light is still off, contact your KINETEC® specialist.
- The Kinetec 6080™ fails to operate but the green light remains on,
 - o Press the START/STOP switch two more times.
 - o If your Kinetec 6080™ still fails to operate, contact your KINETEC® specialist.
- The red light (26) is on,
 - o Verify that the hand control is firmly connected. If so, contact your KINETEC® specialist.

17.3. Cleaning

Before carrying out any cleaning operation, SWITCH OFF the unit and disconnect the power supply. In order to ensure optimal hygiene, you are advised to clean the machine for each new patient.

Cleaning should be carried out in the environmental conditions specified in the "Technical Specifications" section below.

Use a DISINFECTANT product (alcohol-free or <5% alcohol solution) in spray.

In order to ensure optimal hygiene, you are advised to clean the covers for each new patient. All the consumables enable hazard-free disposal.

17.4. Disposal and recycling

- a. **Packaging:** The packaging must be separated into plastic and paper / cardboard components and taken to special recycling sites.
- b. **KINETEC® patient pad kit:** Clean with a disinfectant product then take it to special recycling sites.
- c. **Unit:** It contains electronic components, cables, aluminium, steel and plastic parts. When the machine is no longer operational, disassemble it, separate it into different types of material and take these to authorised recycling centres or return the machine to Kinetec SAS for destruction. Or contact the local authorities to determine the appropriate method of disposal for parts and accessories that are potentially hazardous to the environment.

17.5. Technical specifications

Product:

Lifespan of the machine:	12 years
Weight:	21 kg
Splint dimensions:	72cm x 65cm x 130cm
Angular limits:	from 0° to 135° of flexion with 90° pronation to 90° supination
Speeds:	from 35 to 135° per minute
Patient sizes:	from 1,40m to 1,90m
Maximum weight of the user:	135 kg
Acoustic pressure:	<70dB

Electricity:

Power supply:	240V~ or 100-110V~ (US - Japan)
Frequency:	50/60Hz
Power consumption:	50VA
Class:	Device of Type B Class I
Protection class (device):	IP 20 (protected against solid objects greater than 12.5mm, but not protected against liquids)
Protection class (carrying case):	IP 01 (non-protected against solid foreign objects, protected against vertically falling water drops)
Fuse:	T 315mA 250V 5x20mm (Ref. KINETEC®: 4610002351) T 630mA 250V 5x20mm (Ref. KINETEC®: 4610002351) (US - Japan)

Environment:

Storage/transport conditions:	Temperature: -25 to 70°C / -13 to 158°F. Relative humidity: up to 93% without condensation.
Operating conditions:	Temperature: 5 to 40°C / 41 to 104°F. Relative humidity: 15% to 93% without condensation. Atmospheric pressure: 700 hPa to 1060 hPa.

17.6. Symbols used

	STOP (power off)
	Maximum SPEED
	Follow the instructions for use
	Contains electric and electronic components; do not throw away with household refuse.
	Keep dry during storage or transport
IP20 IP01	See "Protection class" in section "Technical specifications"

	ON (power on)
	Minimum SPEED
	Warning or CAUTION (consult the accompanying documentation)
	Alternating current
	Right way up when box is stored
	Humidity limit during storage or transport

	"ON" / "OFF" (2 stable positions)
	Do not push
	TYPE B device (protection against electric shocks)
	Fuse
	Fragile
	Temperature Limit during storage or transport

17.7. Warranty

The KINETEC® warranty is strictly limited to the replacement, free of charge, or to factory repairs of part(s) recognised as defective.

Kinetec SAS guarantees its continuous passive motion systems for 1 year against all defects of manufacture from the date of purchase by the consumer.

Kinetec SAS is the only organization able to assess the application of the warranty to its systems.

The warranty will be considered null and void if the device has been used abnormally or under conditions of use other than those indicated in the user's manual.

The warranty will also be considered null and void in the event of deterioration or an accident due to negligence, inappropriate surveillance or inappropriate maintenance, or due to transformation of the equipment or an attempt to repair the equipment.

17.8. Information about electromagnetic compatibility (EMC)

All the information shown below is from requirements which the manufacturers of electro-medical devices are subject to, in the context of standard IEC60601-1-2.

Radio-frequency emissions from the Kinetec 6080™ device are very low and is therefore not likely to cause interference with electronic equipment installed close by (radios, computers, telephones, etc.). Nevertheless, the user shall ensure that any electromagnetic interference does not pose an additional risk, such as radio-frequency emitters or other electronic devices.

The Kinetec 6080™ device is designed to tolerate foreseeable disturbances resulting from electrostatic discharges, magnetic fields of the mains power supply or radio-frequency emitters. Nevertheless, certain types of mobile telecommunication devices such as mobile phones are likely to interfere with the medical device. The separation distances recommended in this chapter must therefore be scrupulously followed.

In this chapter you will find information about installing and using your medical device under the best conditions in terms of electromagnetic compatibility.

WARNING: In the event of reciprocal, electromagnetic or other interference with other devices, move the device away.

WARNING: The Kinetec 6080™ device must not be used next to other devices or stored with the latter because this can cause incorrect operation. If such use is necessary, the device and other machines should be monitored to ensure normal operation.

17.9. List of cables

WARNING: The use of accessories and cables other than those specified or provided by KINETEC SAS can increase electromagnetic emissions or reduce immunity of the device and cause inappropriate operation.

List of cables likely to alter the compliance of the Kinetec 6080™ device in relation to emission requirements and electromagnetic immunity:

Designation	Kinetec reference	Maximum length
Europe mains lead	4610006626	3.5m
US/Japan mains lead	4610006634	4m
Brazil mains lead	4610008953	3.5m
UK mains lead	4610007400	3.5m
Australia mains lead	4610007161	3.5m
Swiss mains lead	4610009422	3.5m
Control box/Motor cable	4610003747	1.7m
Hand control cable	4610000975	2.2m

17.10. Recommended separation distances

The Kinetec 6080™ device is anticipated to be used in an electromagnetic environment in which radiated radio-frequency disturbances are controlled.

WARNING: Portable radio-frequency devices should not be used (including peripherals such as coil cables and external coils) within 30 cm (12 inches) of any part of the Kinetec 6080™ device including cables specified by Kinetec. Should the opposite occur, the performance of these devices could be altered.

17.12. Electromagnetic emissions

The Kinetec 6080™ device is anticipated to be used in the electromagnetic environment described in the table below.

The user and the installer should ensure that the Kinetec 6080™ device is used in the environment described below.

Emission tests	Conformity	Electromagnetic environment - Notes
Disturbance from electromagnetic radiation (radiated emissions) CISPR 11	Group 1	The Kinetec 6080™ device uses radiofrequency energy only for its internal operation. Consequently, its radiofrequency emissions are very low and are not likely to create any interference with neighbouring equipment.
Disruptive voltage at power supply terminals (conducted emissions) CISPR 11	Class B	The Kinetec 6080™ device is suitable for use in the health care at home environment and in the environment of a professional health care institution.
Emissions from harmonic currents IEC 61000-3-2	Class A	
Variations of voltage, fluctuations of voltage and flickering IEC 61000-3-3	Compliant	

17.13. Magnetic and electromagnetic immunity

The Kinetec 6080™ device is anticipated to be used in the magnetic and electromagnetic environment described in the table below.

The user and the installer should ensure the conformity of the electromagnetic environment.

Immunity tests	Test level according to IEC60601	Level of conformity	Electromagnetic environment - Notes
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV in contact ±15kV to air	±8 kV in contact ±15kV to air	Suitable for use in the health at home care environment and in the environment of a professional health care institution.
Rapid / in salvo electrical transitionals IEC 61000-4-4	±2 kV for electrical power supply lines ±1 kV for signal ports	±2 kV for electrical power supply lines ±1 kV for signal ports	The quality of the electrical power supply network must be equivalent to that of a health at home care environment and of an environment of a professional health care institution.
Shock waves IEC 61000-4-5	±1 kV between phases ±2 kV between phase and earth	±1 kV between phases ±2 kV between phase and earth	The quality of the electrical power supply network must be equivalent to that of a health at home care environment and of an environment of a professional health care institution.
Magnetic field at the assigned industrial frequency IEC 61000-4-8	30A/m	30A/m	The intensity of the magnetic field must be of the level found in a health at home care environment and in an environment of a professional health care institution.
Voltage hollow IEC 61000-4-11	0% U_T for 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T for 1 cycle 70% U_T for 25 cycles at 50Hz for 30 cycles at 60Hz single phased at 0°	0% U_T for 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T for 1 cycle 70% U_T for 25 cycles at 50Hz for 30 cycles at 60Hz single phased at 0°	The quality of the electrical power supply network must be equivalent to that of a health at home care environment and of an environment of a professional health care institution. If the use of the device requires the continuation of operation during power cuts, it is recommended powering the medical device using a separate power source (UPS, etc.).
Voltage interruptions IEC 61000-4-11	0% U_T for 250 cycles at 50Hz for 300 cycles at 60Hz	0% U_T for 250 cycles at 50Hz for 300 cycles at 60Hz	

NOTE: U_T corresponds to the alternative network voltage before the application of the test level.

17.14. Electromagnetic immunity, radio-frequency portable equipment

The Kinetec 6080™ device is anticipated to be used in the magnetic and electromagnetic environment described in the table below.

The user and the installer should ensure the conformity of the electromagnetic environment.

WARNING: Portable radio-frequency devices should not be used (including peripherals such as coil cables and external coils) within 30 cm (12 inches) of any part of the Kinetec 6080™ device including cables specified by Kinetec. Should the opposite occur, the performance of these devices could be altered.

Immunity tests	Test level	Level of conformity	Electromagnetic environment - Notes
Radiated radio-frequency electromagnetic fields IEC 61000-4-3	3 V/m from 80 MHz to 2.5 GHz 10 V/m from 80 MHz to 2.7 GHz 80% MA to 1kHz	3 V/m from 80 MHz to 2.5 GHz 10 V/m from 80 MHz to 2.7 GHz 80% MA to 1kHz	The Kinetec 6080™ device is suitable for use in the health care at home environment and in the environment of a professional health care institution.
Proximity fields emitted by radio-frequency wireless communication devices IEC 61000-4-3	9 V/m 710 MHz, 745 MHz, 780 MHz, 5240 MHz, 5550 MHz, 5785 MHz 27 V/m 385 MHz 28 V/m 450 MHz, 810 MHz, 870 MHz, 930 MHz, 1720 MHz, 1845 MHz, 1970 MHz, 2450 MHz	9 V/m 710 MHz, 745 MHz, 780 MHz, 5240 MHz, 5550 MHz, 5785 MHz 27 V/m 385 MHz 28 V/m 450 MHz, 810 MHz, 870 MHz, 930 MHz, 1720 MHz, 1845 MHz, 1970 MHz, 2450 MHz	The Kinetec 6080™ device is suitable for use in the health care at home environment and in the environment of a professional health care institution.
Disturbances conducted, induced by radio-frequency fields IEC 61000-4-6	3 V from 150 kHz to 80 MHz 6 V in ISM band and band between 0.15 MHz and 80 MHz, amateur radio band included 80% MA to 1 KHz	3 V from 150 kHz to 80 MHz 6 V in ISM band and band between 0.15 MHz and 80 MHz, amateur radio band included 80% MA to 1 KHz	The Kinetec 6080™ device is suitable for use in the health care at home environment and in the environment of a professional health care institution.
The intensities of electromagnetic fields for fixed radio-frequency emitters, as determined by an electromagnetic environment measurement (a) must be less than the conformity level for each range of frequency.  Interference can occur close to equipment identified by the following symbol:			
NOTE: These specifications may not apply in all situations. Electromagnetic propagation is affected by the absorption and by the reflection of structures, objects and people. (a) The intensities of electromagnetic fields for fixed radio-frequency emitters such as base stations for mobile phones (cell/wireless), mobile radios, amateur radios, AM/FM radio emissions and TV emissions cannot be determined accurately by theory. To evaluate the electromagnetic environment due to fixed radio-frequency emitters, an electromagnetic environment measurement must be made. If the intensity measured of the radio-frequency field in the immediate environment of use of the product exceeds the level of radio-frequency conformity specified above, it is necessary to test the performance of the product to check that it is compliant with the specifications. If abnormal performance is noticed, additional measurements may be necessary, such as re-directing or moving the product.			

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