REDUCE PAIN & SWELLING

WITH THE LATEST COLD COMPRESSION TECHNOLOGY

IMPORTANT NOTES

This product is recommended for one patient use only. If it is used for treating more than one patient, the manufacturer's product liability according to the Medical Devices Act will become invalid.

The final decision concerning the usage and selection of the appliance shall be made by the doctor or medical specialist.

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Manual CCT Shoulder 05/03/2018

dr. **aktive**
LIVE LIFE

Shoulder Support
Cold Compression Therapy

Handleiding
Instructions for use
Mode d'emploi
Gebrauchsanweisung

Easy Use

Durable

Lightweight

INSTRUCTIONS FOR USE ENGLISH

Your doctor has prescribed the **CC Shoulder Pad** for you. Please read through these instructions carefully to ensure that you fit your brace correctly.

IMPORTANT

This product is recommended for one patient use only. If it is used for treating more than one patient, the manufacturer's product liability according to the Medical Devices Act will become invalid.

The final decision concerning the usage and selection of the appliance shall be made by the doctor or medical specialist.

INDICATIONS

- Sprains
- Sport injuries
- Rehabilitation
- Arthritis
- Muscle tiredness

CONTRA-INDICATIONS

No contra-indications are known until today.

WASHING INSTRUCTIONS




Store the product under well ventilated and dried condition and avoid direct sunlight.

Only use mild detergents and warm water to wash the product. Do not use strong detergents!

Using the washing machine and the spin-drier to wash and dry the product is not recommended.

FITTING INSTRUCTIONS STEP BY STEP

- 1 Check the parts first. (1) compressor (2) gel pack (3) one strap (4) pump with tube.
- 2 Refrigerate (2) at least 2 hours before usage.
- 3 Put (2) back on the (1) (Picture 1) after refrigeration.
- 4 Connect the tube to the ON/OFF switch and turn "ON" the switch.
- 5 Pull the strap with black Velcro through the square hook (from 'IN' to 'OUT') then fasten it. (Picture 1).
- 6 Aim at the center of shoulder shin (A) and attach the Velcro first Red then Blue on the proper position. (Picture 2). Picture 2 applies for the left shoulder. If for the right, use the opposite direction.
- 7 Pass the strap (3) through your left or right armpit and pull it through the square hook ('IN' to 'OUT') . Then adjust with Velcro until comfortable (Picture 3).
- 8 If you want to change to the right shoulder, only loose the Velcro of strap (3) in front of chest. Take it off and wear it directly to the right shoulder, then adjust the strap in front of chest until comfortable.
- 9 When turning ON/OFF the switch, please fix it with one hand and turn it with the other to ensure the switch is on the appropriate position and then pump.
- 10 To get the best effect: ON for 5-10 minutes, OFF for 5-10 minutes to let the skin rest and then repeat.
- 11 Empty it before taking off and then refrigerate (2) for next usage.

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