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
**Andratec**  
Advanced Interventional Technologies

AndraTec GmbH  
Simmerner Str. 70  
56075 Koblenz  
Germany  
ph: +49 (0) 261 29675172  
fax: +49(0) 261 29675173  
www.andratec.com

#### Instructions for Use:

#### **Exeter-Snare and Exeter-Snare micro**

The Exeter-Snare and the Exeter-Snare micro is only available as set with the suited Introducer Sheath.

#### **Device description:**

The Exeter-Snare is a flexible snare for interventional retrieval of foreign bodies out of blood vessels or body cavities. The Exeter-Snare consists of a super elastic Nitinol wire with a special head. Due to the special material properties and the variable head of the snare, it is reconfiguring into its original shape. According to the position how much the snare is advanced out of the introducer sheath the retrieval snare has a different shape and a varying catching area. The Exeter-Snare comes preloaded with the AndraTec Introducer Sheath. The interventional retrieval should only be performed by physicians with experience in interventional techniques.

#### **Indications:**

The Exeter-Snare is indicated for use in patients for interventional retrieval of displaced foreign bodies. The snare can be used as well for catching guide wires and for special interventional guide wire loop procedures. The Exeter-Snare is available in different lengths and variable diameters.

#### **Contra-indications**

- Foreign bodies, which are not retrievable with a snare.
- Foreign bodies, which are hooked with the surrounding tissue.
- Patients, where is a risk of vessel perforation or damage of organs during the retrieval procedure.
- Patients with a contra-indication for antiplateletted/ anticoagulant therapy.

#### **Adverse effects**

Possible adverse effects include, but are not limited to the following:

- Arrhythmias, including ventricular fibrillation.
- Bleeding complications (particular haemorrhage or haematoma at puncture site).
- Cerebral vascular accidents
- Embolism
- Infection
- Injury or rupture of surrounding vessel structures
- Loss of foreign body
- Vessel occlusion

#### **Method:**

Prior to use the AndraTec introducer sheath of adequate size must be placed for vessel access. Select an appropriate snare size for foreign body removal. The introducer sheath must be placed near the site of the foreign body. With the help of a small plastic tubing, the snare head is folded together and can be advanced easily into the introducer sheath. According how far the snare is advanced out of the sheath it is reconfiguring with variable catching areas and shapes for being pushed over the foreign body. The movements of the snare can be accurately controlled with the help of a torquer. When the snare is around the foreign body the introducer sheath can be advanced with fixed torquer for closing the snare loop. The foreign body can be fixed to the tip of the introducer sheath with slight tension and can then be removed. For removal of bigger foreign bodies, larger sheaths may be required or a surgical cut down.

#### **Preoperative treatment:**

The physician himself is responsible for the drug regimen.

#### **Warning:**


Prior to use the package, the introducer sheath and the snare have to be inspected to exclude any damage. The interventional removal of foreign bodies has to be performed under X-ray visualization or endoscopic view and the (angiographic) result should be documented. In emergency situations the interventional removal has to be performed with surgical stand-by or if there is any

increased risk of vessel perforation. If there is any increased resistance during foreign body retrieval the cause has to be detected immediately to undertake adequate reactions. If necessary the removal has to be stopped or has to be performed surgically. Application of extended force may cause damage to surrounding vascular structures. Prior to use the introducer sheath has to be flushed with a sterile isotonic (heparin zed) (or NaCl 0.9%) solution.

#### **Disclaimer of warranty:**

Although this medical product has been carefully designed, manufactured, and tested prior to sale, it may fail to perform its intended function satisfactorily for a variety of reasons. Handling, storage, cleaning and sterilization of this instrument as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond the manufacturer's control directly affect the instrument and the results obtained from its use. The manufacturer, therefore, disclaims all warranties, both expressed and implied, with respect to the product. The manufacturer shall not be liable for any incidental or consequential damages by any use, defect, or failure of the product, whether the claim is based on warranty, contract or otherwise. The manufacturer is not liable for all incidental, special, or consequential loss, damage, or expense direct or indirect, from the use of its product. The manufacturer assumes no liability with respect to instruments reused, reprocessed or resterilized. Descriptions or specifications in printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

**Sterile and nonpyrogenic unless package is opened or damaged.  
Check for damage before use.**

**Andratec**  
Advanced Interventional Technologies

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