# VascuFlex®

- ⒹⒺ Selbstexpandierender, peripherer Gefäßstent
- ⒼⒷ Self Expandable Peripheral Vascular Stent
- ⒻⓇ Stent vasculaire périphérique, auto-expandible
- ⒺⓈ Stent vascular periférico autoexpandible
- ⒾⓉ Stent vascolare periferico autoespandibile
- ⓈⒺ Självexpanderande, perifer kärlstent
- ⒻⒾ Itselaajeneva perifeerinen stentti
- ⓃⒻ Zelfexpanderende, perifere bloedvatstent
- ⓅⓉ Endoprótese vascular periférica autoexpansível
- ⒸⓏ Samoexpandující, periferní cévní stent
- ⓈⓀ Samoexpandovateľný, periférny cievny stent
- ⓅⒻ Stent obwodowy samorozprężalny
- ⒷⒼ Саморазгъващ се периферен съдов стент
- ⒽⓁ Öntáguló, perifériás ér-stent
- ⓉⓇ Kendinden genişleyen periferik damar stenti
- ⓇⓄ Stent vascular periferic, cu auto-expansiune

## VascuFlex®

**Produktbeschreibung**

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## VascuFlex®

### Description

VascuFlex® is a self-expanding stent system consisting of two main components: (a) a self-expanding Nitinol stent with Tantal markers, mounted on (b) a co-axial delivery system. The delivery system is introduced into the body through an introducer sheath over a guide wire (see table 1) and advanced up to the stenosis or the occlusion. After the retraction of the exterior catheter of the delivery system, the self-expanding mechanism of the stent is triggered by the influence of the body temperature. Thus the stent reaches its preshaped diameter and squeezes due to its radial force excess material into the vessel wall, restoring the patency.

Table 1 – Overview over the VascuFlex® variants

System	Recommended introducer sheath	Recommended guide wire	Usable length of the delivery system
VascuFlex® 5F	5F (1.7 mm)	.018"/.021"	80 cm or 120 cm
VascuFlex® 5F transradial/transbrachial	5F (1.7 mm)	.018"	180 cm
VascuFlex® 6F	6F (2 mm)	.035"	80 cm or 130 cm
VascuFlex® 6F LRF	6F (2 mm)	.035"	80 cm or 130 cm

### Delivery system

The delivery system (10) consists of an interior (9) and an exterior catheter (4). The lumen of the interior catheter is the guide wire lumen. The stent is fitted at the distal end of the delivery system between the interior and the exterior catheter. In order to facilitate the positioning of the stent two radiopaque ringmarkers indicate the proximal (3) and the distal (3) end of the stent. An atraumatic tip is molded (2) to the distal end of the interior catheter.

The interior catheter (9) is proximally connected to a stainless steel tube (8). This runs through a Y-shaped haemostatic valve (5) and ends proximally in a Luer-Lock-

### Stent

The stent is cut from a Nitinol tube by a laser into the required size. During multiple manufacturing steps the stent is given a so called "memory" shape and in a further step the Tantal markers are applied to the ends of the stent. Exposed to body temperature, the stent expands to its preshaped diameter. The stent features a high level of flexibility and its high radial forces enable the patency of the vessel.

connector (9) through which the guide wire exits. The exterior catheter is proximally connected to the haemostatic valve. The haemostatic valve is dispatched in a not completely locked status. It has to be completely locked before the application and during the insertion of the system in order to prevent a premature release of the stent. The haemostatic valve is locked by a clockwise turn of the end (7) and it is opened by an anti-clockwise turn. During the release of the stent the metal shaft (8) is to be fixated whilst the haemostatic valve is to be retracted simultaneously with the exterior catheter in proximal direction over the metal shaft.


Figure 1 – VascuFlex® Stentsystem

### Indication

VascuFlex® is indicated for the following applications:

- atherosclerotic diseases of the peripheral arteries, e.g. A. iliaca, A. femoralis and A. femoralis superficialis and in other vessels which correspond in their anatomic dimensions to the available sizes of stents (see table 2 "Selection of appropriate stent" further below).
  - treatment of bile duct strictures.
- VascuFlex® 5F transradial/transbrachial is not indicated for treating bile duct strictures.

### Contraindications

In general the same contraindications are valid for the application of stents in the blood vessel system as for PTA (percutaneous transluminal angioplasty). The implantation of a stent in the vascular system as well as in the bile duct is contraindicated under the following circumstances, however not exclusively:

- for patients with severely calcified stenosis which cannot be treated by PTA
- for patients whose target lesion borders on a larger conglomeration of acute or sub-acute thrombus
- for patients with a proven vessel perforation
- if there is an aneurysm immediately before or behind the target lesion
- for patients with known and untreated blood-coagulation disorders
- for patients who cannot receive an anticoagulation or antiplatelet therapy
- for functionally relevant obstruction of the inflow path, poor distal outflow or defect distal lumen connection
- for existing fresh, non-organized thrombus or embolic material
- for lesions of bifurcations in which a stent implantation would block the access to the main branches which in future could profit from a PTA
- in intracranial arteries as well as in the A. carotis
- in case of ascites
- in the case of a perforated bile duct in which because of the stent the leakage from the bile duct would increase through the stent meshes
- if the stenosis cannot be passed with a guide wire
- for extensively metastasis-affected liver or liver carcinoma
- in the case of main bifurcations which might be used for diagnostic and therapeutic purposes at a later stage

### Relative Contraindications

- for patients with a risk for an allergic reaction to contrast media
- for pregnant patients

### Warnings and Precautions

Please remember that failure to observe the following may result in serious patient complications.

### In general

- The stent system should only be used by physicians who are experienced in interventional techniques as PTA or stent procedures or under guidance of physicians experienced in this field. The physicians should be accustomed to the complications, the side effects and the dangers that are commonly attributed to stenting.
- The stent system is only suited for single use. Reprocessing and resterilisation of the stent system may lead to major faults such as stent system malfunctions, mechanical properties can no longer be guaranteed, material properties can no longer be guaranteed, risk of infection

and risk of contamination.

- The stent system may not be used, if the sterile package has been opened or damaged.
- The stent system may not be used if it is visibly or tangibly damaged.
- The stent system must be used before the expiry date.
- The non-compliance of sterile precautions can lead to infectious complications.
- The VascuFlex® 5F stent with a 180 cm application instrument is recommended for the transradial / transbrachial approach.
- The VascuFlex® stent is not intended for the employment in the heart, in the central circulation system or central nervous system (according to EC-directive 93/42/EEC, annex IX, definitions, point 1.7. and 1.8.).
- The VascuFlex® stent should not be implanted if other stents from different materials affect the VascuFlex® stent or are in the vicinity. In this case, electrolytic corrosion and potentially thrombosis formation can be caused.
- The application system must not be attached to high pressure injector.
- The application system must not come into contact with organic solvents (e.g. alcohol).

### The delivery system

- The delivery system is **not** designed for high pressure injections.
- Do not use with Lipiodol or Ethiodol contrast media.
- Do not expose the delivery system to organic solvents (e.g. alcohol).

### During the Implantation

- The stent system must always be inserted, moved or retracted over a suitable guide wire (see table 1).
- As soon as the release respectively the partial release of the stent has been started it **may not be retracted into the delivery system or repositioned.**
- As soon as the system has been inserted into the vascular or biliary system it may only be manipulated under fluoroscopic monitoring. For the X-ray visualization only high-resolution equipment should be used.
- The stent system may not be advanced against strong resistances. Before the procedure can be continued the cause of the resistance must be determined under fluoroscopic monitoring and measures for the retraction must be considered.
- The blockage of vital lateral branches by the stent is to be avoided.
- A faulty positioning technique can lead to a failure in the release of the stent and thereby to an incorrect implantation.

### The Stent

- The stent should be easy to release, i.e. the exterior catheter should be retractable with ease. Under no circumstance try to release the stent by exercising increased force. Remove the complete system and replace it by a new system.
- As soon as the stent is partially released it may not be retracted into the delivery system. Also it may not be repositioned.
- If several lesions are to be treated during one procedure, first the distal lesion and then the proximal lesion is to be stented. Stenting in this order avoids having to pass already positioned stents again, thereby possibly dislocating them.

- The renewed passage of a stent should only be exercised if absolutely necessary and then with greatest care.
- In order to prevent an overirritation of the vessel wall and to maintain the flexibility of the vessel it is best advice to position the stents end to end and not overlapping.
- If the size of the stent is not matched exactly to the native vessel the stent might migrate if the diameter has been chosen too small due to the insufficient radial force. In the case of a diameter chosen too large the intima can be damaged due to too high radial force. The stent length is to be chosen so as to completely cover the lesion. For the correct choice of stent please confer the selection table.

#### **Adverse reactions**

Procedures with percutaneous catheters may only be exercised by physicians who are acquainted with the potential complications. Complications can occur at any time during or after the procedure. If due to a complication the surgical removal of the stent becomes necessary the clinically accepted surgical standard method is to be adopted. Among others the following potential complications – known today – are listed:

#### **related to vascular application**

- Difficulties accessing or dilating the stenosis
- Spasm of the vessel wall
- Different extents of injuries of the vessel wall including perforation, rupture, or dissection
- Haemorrhage which necessitates a transfusion
- Thrombus formation, distal embolism or total occlusion
- Pseudo-aneurysm or haematoma in the access area
- Wrong placement of the stent
- Migration of the stent with subsequent embolism
- Infection following contamination
- Allergic reaction to one or several used materials, to the contrast medium, or to the administered drugs
- Emergency surgery of an arterial bypass
- Renal failure
- Amputation
- Pneumothorax
- CVA (cerebral vascular accident) – stroke
- Coronary ischemia
- Arrhythmia

#### **related to bile duct applications**

- A stent protruding too far into the duodenum can lead to a damage or an occlusion of the bowel tract
- Perforation of the bile duct
- External biliary fistula
- Parenchymatous or intraperitoneal haemorrhage
- Cholangitis
- Cholecystitis
- Biliary peritonitis
- Haemobilia
- Liver abscess
- Tissue necrosis
- Subcapsular liver haematoma
- Stent obstruction as a consequence of tumor ingrowth
- Wrong placement of the stent
- Stent migration
- Proliferation of the stent ends
- Occlusion by sludge

#### **Implantation Procedure**

Artificial stents must always be implanted by surgeons

that are trained and familiar with the specific techniques of such operations. Please also observe the following points when implanting the VascuFlex® stent:

#### **Preparation**

- The percutaneous placement of the stent in a stenotic or obstructed vessel is to be exercised under fluoroscopic monitoring in an angiographic cath-lab.
- The preparation of the patient as well as the sterile measures to be taken are to be the same as for any angioplastic or transhepatic catheter treatment.
- The medication is set by the treating physician, as well as the measures to be adopted during and after the treatment to reduce coagulation considering the latest scientific guidelines and suited to the individual patient.

#### **Angiography/Cholangiography**

- The extent of the lesion(s) and the collateral circulation must be established employing angiography respectively cholangiography in the case of the bile duct.
- If there is a thrombus formation or a suspected thrombus formation a thrombolysis should be exercised using a clinically customary and proven standard method.
- The access vessels must be sufficiently patent or re-channelled, before the intervention can be continued.
- The lesion should be identified and marked under fluoroscopy and under consideration of the distal part of the stenosis or stricture.

#### **Pre-dilation of the Lesion**

- It is advised to pre-dilate the lesion by standard PTA technique in order to be able to determine the correct size of the stent.
- The balloon catheter is removed and the access to the lesion is maintained by the remaining guide wire.
- **Please note:** With malignant biliary strictures a pre-dilation is not practiced as a standard procedure. However, if the practicing physician considers it necessary a pre-dilation can be exercised using a standard PTA technique.

#### **Selection and preparation of the stent system and compatibility with the accessories**

##### **1. Selection of stent size**

- In order to determine the necessary stent length, the length of the target lesion has to be measured.
- The lesion is to be completely covered by the stent.
- **Please note:** Should more than one stent be necessary to cover the lesion the distal stent must be implanted first, followed by the stent(s) to be placed proximally.
- In order to avoid an irritation of the vessel wall and to maintain the flexibility of the vessel it is best advice to position the stents end to end and not overlapping.
- In order to determine the correct diameter of the stent, the diameter of the reference vessel (proximally and distally to the lesion) must be measured. It is necessary to select a stent which expanded nominal diameter is at least one millimeter bigger than the largest diameter of the reference vessel so that the stent fixes itself securely into the vessel wall.
- For the correct selection please take note of the following table 2:

**Table 2 – Stent-Selection Table**

Diameter of the vessel or bile duct lumen	Nominal diameter of the stent in expanded status
3.0 – 4.0 mm	5 mm
4.0 – 5.0 mm	6 mm
5.0 – 6.0 mm	7 mm
6.0 – 7.0 mm	8 mm
7.0 – 8.0 mm	9 mm
8.0 – 9.0 mm	10 mm

##### **2. Preparation of the delivery system**

- First of all please make sure that the sterile package has not been damaged or opened.
- Open the sterile pouch and remove the blister under sterile conditions.
- Take out the delivery system from the blister. Check it for visual or palpable damages. Check whether the stent is mounted on the delivery system.
- **Please note:** If there is reasonable doubt that the sterility or the functionality is impaired the system may not be used.
- The system may not be used, if the stent is partially deployed or if the space between exterior catheter and catheter tip is more than 2 mm. By opening the haemostatic valve (7) and slightly pulling the interior catheter in the proximal direction the gap should be closed. If this is not the case discard the system.
- The haemostatic valve (7) is locked on delivery. Before further manipulations of the delivery system are exercised, it should be ensured by turning the valve clockwise that it is completely locked. Thereby a premature release of the stent is avoided.
- Before the use of the delivery system it must be completely rinsed with heparinized saline solution through both proximal accesses.
- For this a syringe with heparinized saline solution is fitted to the lateral exit of the Y-shaped haemostatic valve (6) and pressure applied. Ensure that the haemostatic valve (7) is locked and rinsing fluid continuously drips from the catheter tip between interior and exterior catheter.
- Then the central guide wire lumen is rinsed by connecting the syringe to the proximal Luer-Lock connector (9).

#### **Stentimplantation**

##### **1. Insertion of the delivery system**

- Using common techniques, such as the Seldinger technique, a percutaneous access is created. Then an introducer sheath compatible with the delivery system (see table 1) is inserted.
- Through the introducer sheath a guide wire of suitable gauge and length (see table 1) is advanced up to and beyond the stenosis/stricture.
- First of all – if necessary – a balloon catheter is inserted over the guide wire and the lesion pre-dilated. Then the balloon catheter is removed. The guide wire remains in place to maintain access to the lesion.
- **Please note:** All manipulations with the delivery system or other devices must always be exercised over the guide wire.
- Always use an introducer sheath with a haemostatic

valve. For contra-lateral accesses it is advised to use a long introducer sheath specially designed for contra-lateral accesses which simplifies the crossing of the bifurcation.

- **Manufacturer recommendation:** In order to manipulate the pelvis or leg vessels on the opposite side using the cross-over technique, a long bend-resistant sheath should be used of which the tip is pushed into the arteria iliaca communis or externa, so as to ensure stability, and a Nitinol guide wire.
- In order to remove air bubbles from the guide wire lumen (9) and the interior catheter (4) these are once again rinsed with heparinized saline solution.
- By turning the haemostatic valve (7) – which connects the interior and the exterior catheter – clockwise, it is ensured that an accidental premature release of the stent is prevented.
- The delivery system (10) is advanced to the lesion over the guide wire (1) under permanent fluoroscopic control. The exact positioning of the stent in the stenosis/stricture is achieved by radiopaque markers (3) between which the stent is placed.
- **Warning:** The stent system must not be advanced against strong resistance. Before the procedure can be continued the reason for the resistance must be determined under fluoroscopic monitoring, and measures for the retraction must be considered.

##### **2. Release of the stent**

- The delivery system is advanced so far that both the radiopaque markers on the interior catheter are positioned proximally and distally to the target lesion. Only then the correct positioning of the stent is ensured. The haemostatic valve (7) which connects the interior and the exterior catheter can then be opened by an anti-clockwise turning.
- By fixing the stainless steel tube (8), which is connected to the interior catheter, in a stable position with one hand and slowly retracting the exterior catheter (4) at the haemostatic valve (7) over the stainless steel tube in direction of the proximal end with the other hand, the stent is positioned and can expand under the influence of the body temperature.
- **Please note:** During the complete release the stainless steel tube (8) is to be held in a fixed position.
- In order to achieve the best possible stent positioning attention must be given to the fact that the proximal marker (3) remains in the same position during the release.
- **Caution:** Do not advance the stainless steel tube (8) in distal direction. This can lead to a faulty positioning of the stent and damage the vessel/bile duct.
- **Warning:** The stent should be easy to release. Interrupt the procedure if an unusual effort is necessary. A change of system is necessary.
- The release is completed when the haemostatic valve (7) touches the proximal Luer Lock.

##### **After stent implantation**

- After the complete release of the stent the exterior catheter is moved again forward over the interior catheter. The complete delivery system is removed slowly and carefully under fluoroscopic control over the guide wire out of the body. Thereby the guide wire remains still in position.
- **Please note:** Should the tip of the delivery system get caught in a distal stent segment it is to be freed from it

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
Figure 1 – VascuFlex® Stentsystem

by careful movement.

- By conventional execution of a post angiography or post cholangiography it is checked whether the stent has completely expanded.
- Has the stent not completely expanded at any position within the lesion, an additional balloon dilation (PTA technique) can be exercised. For this a balloon catheter of suitable size is chosen and the lesion dilated in common technique. The inflated diameter of the balloon has to correspond with the reference vessel.
- **Caution:** Never expand the stent with a balloon which is larger in diameter than the expanded diameter of the stent.
- Has the procedure been successful, the guide wire and the introducer sheath can also be removed and the haemostasis of the access site can be induced.
- Should the result of the stenting not yet be satisfactory, further measures to be taken must be decided.

### MR – Compatibility

The stent is MR-compatible. It will not be influenced by MR and causes only low till none image artifacts. Nevertheless, it is recommended to carry out a MR examination only about six weeks after the implantation to avoid a dislocation. The current standard of knowledge regarding magnetizing force and adjustment to the magnetic field should be considered to receive optimal results.

### Storage and packaging

VascuFlex® is a single-use product that has been sterilized with ethylene oxide. It is supplied in sterile packaging and will remain sterile as long as the package is not opened or damaged. The product is free of pyrogen.

The stent system must be kept in dark, dry and cool place. Direct sun exposure must be avoided – heat could have a considerable, negative influence on the functionality of the product.

Do not expose the stent system to ionizing radiation or ultra-violet light.

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