# Antibiotic Disc

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### DESCRIPTION

Antibiotic Disc are paper discs with special features, that are impregnated with antibiotic and used for the susceptibility test according to the Kirby-Bauer antibiotic testing (KB testing or disk diffusion antibiotic sensitivity testing).

Antibiotic Disc are available in a large variety of configurations. Each configuration is available in packages of 50 and 250 tests.

### CONTENTS OF THE PACKAGES

The 50-test box contains 1 cartridge with 50 discs packed in a desiccant envelope.

The 250-test box contains 5 cartridges of 50 discs, each cartridge individually packed in a desiccant envelope.

Each package also contains a transparent resealable bag.

### METHOD PRINCIPLE

The discs are applied to the surface of a culture medium inoculated with the broth culture containing pure colonies of the microorganism under examination. After incubation, the plates are examined, the inhibition halos around each disc are examined and compared with the standard inhibition haloes: in this way the microorganisms are defined as being susceptible, intermediate or resistant to the tested antimicrobial agents.

### COMPOSITION

Liofilchem antimicrobial susceptibility test discs are made of high-quality paper in compliance with WHO and FDA specifications.

The discs are manufactured under the quality systems UNI EN ISO 9001:2008 and EN ISO 13485:2012, and to DIN specification for potency, i.e. the concentration of each antibiotic is within 90-125% of the concentration stated on the disc.

### GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility test are taken up by culture media that have been previously swabbed with the sample under examination. In the case of mixed colonies the bacterial strains must be purified before they are swabbed on the plates for the susceptibility test.

### TEST PROCEDURE

For the details of the test procedure, refer to the bibliography indicated below or to microbiology manuals.

Here the Kirby-Bauer method is shortly summarized:

1. Allow unopened envelop to come to room temperature before opening it, for minimizing condensation on the discs as such condensation could affect long-term stability;
2. Swab 4 to 5 well isolated and morphologically similar colonies with a culture medium and suspend them in 5 ml of a suitable broth culture such as Tryptic Soy broth;
3. Incubate the broth culture at  $36\pm 1^\circ\text{C}$  until 0.5 McFarland turbidity is reached;
4. Immerse a sterile swab in the broth culture or in a diluted form thereof and squeeze it on the wall of the test tube to eliminate excess liquid. Drag it along the surface of the susceptibility test medium contained on the plate so as to produce even growth;
5. Position the antibiotic discs within 15 minutes from inoculation of the plates, pressing them with a sterile pliers on the surface of the agar and within 15 minutes of depositing of the discs leave the plate to incubate at  $36\pm 1^\circ\text{C}$  for 18-24 hours.
6. Reinsert the cartridge containing unused discs into the desiccant envelope and then into the resealable bag.

### EVALUATING THE RESULTS / QUALITY CONTROL

At the end of the incubation period, measure the inhibition halos and interpret according to the current reference standards:

[www.liofilchem.net/pdf/disc/disc\\_interpretative\\_table.pdf](http://www.liofilchem.net/pdf/disc/disc_interpretative_table.pdf)

### CLINICAL INTERPRETATION

The susceptibility test carried out *in vitro* cannot exactly reproduce *in vivo* conditions. Nevertheless, it shows the effect of the concentration of the antibiotic, which varies in the culture medium in relation to the growth of the microbial population.

The final choice of antibiotic to administer to the patient is the responsibility of the clinician who possesses all the information on the patient.

### LIMITS

Diffusion susceptibility tests use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion susceptibility test. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

### PRECAUTIONS

The Antibiotic Disc cannot be classified as being hazardous according to current legislation but fall within the specific field of application where a safety data sheet must be supplied because they can cause phenomena of sensitisation in sensitive subjects if they come into contact with the skin.

Antibiotic Disc are disposable products. Antibiotic Disc are only for diagnostic *in vitro* use and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

### STORAGE

The unopened package of Antibiotic Disc can be stored in most cases at  $-20^\circ\text{C}$  to  $+8^\circ\text{C}$  till the expiry date. Some products have to be stored at  $-20^\circ\text{C}$  as maximum storage temperature. The recommended temperature limits can be found both on the product envelop and on the label box. Leftover discs from an opened cartridge should be stored for no more than 7 days. Return unused discs to the refrigerator/freezer as soon as the application of the discs has been completed. Dispose of expire discs.

### ELIMINATING USED MATERIAL

After use, Antibiotic Disc and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

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RUO: For Research Use Only, not CE marked.

\* **Packaging** Each pack contains 5 cartridges of 50 discs (5x50 discs = 250 discs).

Single cartridges of 50 discs are available: add /1 to the catalogue ref. no. to indicate the relevant item.

Example: ref. 9045/1 indicates Vancomycin 30 µg in one single cartridge of 50 discs.

**CE Marking** Unless otherwise indicated, all products intended to be used for clinical applications are CE marked according to the European Directive 98/79/EC for *In Vitro* Diagnostic Medical Devices.

**Note that testing and reporting antimicrobial agents for which there are no interpretive criteria are the responsibility of the chief microbiologist and such decisions should be made with input from the infectious disease clinicians.**

For full details on specific organism/agent combinations refer to current CLSI and EUCAST recommendations.

Also see the EUCAST guidance at [http://www.eucast.org/clinical\\_breakpoints/when\\_there\\_are\\_no\\_breakpoints/](http://www.eucast.org/clinical_breakpoints/when_there_are_no_breakpoints/)

<b>Antibiotic discs in cartridge</b>			
<b>Veterinary</b>			
<b>Description</b>		<b>µg</b>	<b>Ref.*</b>
Ceftiofur	FUR	30	9251
Enrofloxacin	ENR	5	9233
Florfenicol	FFC	30	9234
Neomycin	N	30	9032
Spiramycin	SP	100	9088
Tiamulin	T	30	9094
Tylosin	TY	30	9082

- \* **Packaging** Each pack contains 5 cartridges of 50 discs (5x50 discs = 250 discs).  
 Single cartridges of 50 discs are available: add /1 to the catalogue ref. no. to indicate the relevant item.  
 Example: ref. 9082/1 indicates Tylosin 30 µg in one single cartridge of 50 discs.

#### BIBLIOGRAFIA / BIBLIOGRAPHY / BIBLIOGRAFIA / BIBLIOGRAPHIE / ΒΙΒΛΙΟΓΡΑΦΙΑ / BIBLIOGRAFIA

1. CLSI. Performance Standards for Antimicrobial Susceptibility Testing; 27th ed. CLSI Supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.
2. CLSI. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard - Twelfth Edition. CLSI document M02-A12. Wayne, PA: Clinical and Laboratory Standards Institute; 2015.
3. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 7.1, 2017. <http://www.eucast.org>
4. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST, Version 7.0, 2017. <http://www.eucast.org>
5. EUCAST Disk Diffusion Method for Antimicrobial Susceptibility Testing - Version 6.0 (January 2017).
6. EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 2.0 (July 2017).
7. CLSI. Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated From Animals, 3rd Edition. CLSI supplement VET01S. Wayne, PA: Clinical and Laboratory Standards Institute; 2015.
8. CLSI. Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated From Animals; Approved Standard - 4th Edition. CLSI document VET01. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
9. FDA (1978) Codes of Fed.Rebs. 21.Part 460.
10. WHO (1977) Tech rep.Ser.n°610.
11. DIN 58940-2 Medical microbiology - Susceptibility testing of microbial pathogens to antimicrobial agents - Part 2: Active substance carriers for the agar diffusion test; 2007-10.

For more information on Antimicrobial Susceptibility Discs,  
 please go to [www.liofilchem.net/antibioticdisc/](http://www.liofilchem.net/antibioticdisc/)

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

**LIOFILCHEM® s.r.l.**

Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
 Tel. +39 0858930745 Fax +39 0858930330 [www.liofilchem.net](http://www.liofilchem.net) [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

F01511  
 Rev.13 / 20.11.2017

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© Copyright LIOFILCHEM 2017