# ANTIFUNGAL DISC

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# ANTIFUNGAL DISC

Discs for antifungal susceptibility test of yeasts

## DESCRIPTION

**ANTIFUNGAL DISC** are paper discs with special features that are impregnate with antimicrotics, used for the susceptibility test of yeasts according to the CLSI method. Each antimicrotic is available in packages of 50 and 250 discs.

## CONTENTS OF THE PACKAGES

The 50-test package contains 1 cartridge with 50 discs packaged in blister with silica gel and an instruction sheet.

The 250-test package contains 5 cartridges with 50 discs, packaged in blister with silica gel and an instruction sheet.

## METHOD PRINCIPLE

The discs are applied to the surface of a culture medium inoculated with the suspension containing pure colonies of the microorganism under examination. After incubation, the plates are examined, the inhibition halos around each disc are examined and compared with the standard inhibition haloes: in this way the microorganisms are defined as being susceptible, susceptible-dose dependent, resistant or non susceptible to the tested antimicrobial agents.

## COMPOSITION

The discs are made of high-quality paper in compliance with WHO and FDA specifications.

Each disc is impregnated with a standardised amount of antimicrotic.

## GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility test are taken up by culture media that have been previously swabbed with the sample under examination. In the case of mixed colonies the strains must be purified in blood medium or Sabouraud Dextrose Agar before they are swabbed on the plates for the susceptibility test.

## TEST PROCEDURE

For procedural details of the test, see the attached bibliography or the microbiology literature. A brief description of the CLSI method follows.

1. Take the blister of cartridges from the refrigerator and leave it for about an hour at ambient temperature in order to prevent condensation forming on the discs when it is opened as such condensation could affect long-term stability.
2. Swab 4 to 5 well isolated and morphologically similar colonies with a culture medium and suspend them in 5 ml of a sterile physiological solution until 0.5 McFarland turbidity is reached for *Candida* spp. and 1 Mcfarland for *C. neoformans*.
3. Immerse a sterile buffer in the suspension and squeeze it on the wall of the test tube to eliminate excess liquid. Drag it along the surface of the susceptibility test medium of yeasts (M.HINTON+GLUCOSE+METHYLEN BLUE AGAR, Ref. 11510).
4. Position the antifungal discs within 15 minutes from inoculation of the plates, pressing them with a sterile pliers on the surface of the agar and within 15 minutes of depositing of the discs leave the plate to incubate at  $36\pm 1^{\circ}\text{C}$  for 20-24 hours; read at 48 hours only when insufficient growth is observed after 24 hours of incubation.

## EVALUATING THE RESULTS

At the end of the incubation period, measure the inhibition halos and compare them with the diameters of the standard CLSI inhibition halos.

Pinpoint microcolonies at the zone edge or large colonies within a zone are encountered frequently and should be ignored.

## CLINICAL INTERPRETATION

The susceptibility test carried out *in vitro* cannot exactly reproduce *in vivo* conditions. Nevertheless, it shows the effect of the concentration of the antimicrotic, which varies in the culture medium in relation to the growth of the microbial population.

The final choice of antimicrotic to administer to the patient is the responsibility of the clinician who possesses all the information on the patient.

## QUALITY CONTROL

Each batch of **ANTIFUNGAL DISC** is subjected to precise and thorough checks in compliance with CLSI standards using the following strains: *Candida albicans* ATCC 90028; *Candida parapsilosis* ATCC 22019; *Candida tropicalis* ATCC 750; *Candida krusei* ATCC 6258.

## LIMITS

Diffusion susceptibility tests use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion susceptibility test. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

## PRECAUTIONS

The **ANTIFUNGAL DISC** cannot be classified as being hazardous according to current legislation but fall within the specific field of application where a safety datasheet must be supplied because they can cause phenomena of sensitisation in sensitive subjects if they come into contact with the skin. **ANTIFUNGAL DISC** are disposable products. **ANTIFUNGAL DISC** are for *in vitro* diagnostic use only and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

## STORAGE

Store the blister at  $-20$  to  $8^{\circ}\text{C}$  till the expiry date. Leftover discs from an opened package should be stored at  $2-8^{\circ}\text{C}$  avoiding humidity until the expiry date. To further reduce humidity, store the discs inside the disc dispenser which contains desiccant. In such conditions, **ANTIFUNGAL DISC** can be used until the expiry date shown on the label. Do not use after this date. Dispose of if they show signs of deterioration.

## ELIMINATING USED MATERIAL

After use, **ANTIFUNGAL DISC** and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

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**BIBLIOGRAFIA/REFERENCES/BIBLIOGRAPHIE/REFERENCIAS**

- CLSI document M44-A2—Method for Antifungal Disk Diffusion Susceptibility Testing of Yeasts; Approved Guideline—Second Edition (August 2009).
- CLSI document M44-S3—Zone Diameter Interpretive Standards, Corresponding Minimal Inhibitory Concentration (MIC) Interpretive Breakpoints, and Quality Control Limits for Antifungal Disk Diffusion Susceptibility Testing of Yeasts; Third Informational Supplement.

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**LIOFILCHEM® S.r.l.**