

Warm up

Our warm up program prepares your body for the workout. This ensures a good start to the training session.

Training	Description	Frequency	Pulse length	Duty Cycle	Duration (min.)
∑∑∑ Warm Up	Ready, set, go! Warm up and get your muscles ready for the next training session. The frequency of the stimuli slowly increases during the process: from 7 to 40 stimulations per second. You can choose the period of time during which this increase will take place.	7-40Hz	350µs	-	3, 6

Cool down

With the cool down program, you help your muscles to "power down" slowly after training.

Training	Description	Frequency	Pulse length	Duty Cycle	Dura- tion (min.)
}}} Cool Down	Training completed? Great, so relax – and start with the first regeneration for your muscles. The frequency of stimulation slowly decreases during the process: from 40 to 7 stimulations per second. You can select a period of 3 or 6 minutes.	7-40Hz	350µs	-	3, 6