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10031106 10031107

KLAR  **FIT**

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Dear Customer,

Congratulations on purchasing this equipment. Please read this manual carefully and take care of the following hints to avoid damages. Any failure caused by ignoring the mentioned items and cautions mentioned in the instruction manual are not covered by our warranty and any liability.

Technical Data

Item number	10031106, 10031107
Power supply	12 V
Power consumption	36 W

Safety Instructions

General Hints

- In order to ensure safe & proper use of the product, please refer to the manual before using.
- This product is not suitable to the people, who has no healthy body, mental disorder, or who lack of massage experience or knowledge (including children), to use, unless there is one responsible person for their safety supervision and guide them to use it.
- Watch the children, to ensure them not play with it.
- The foot massager with heater, surfaces will be very hot, the one who is not sensitive to heat must pay attention when using it.
- If the power cord damage, in order to avoid danger, it must be the manufacturer or its designated maintenance organization, or professional maintenance personnel to change it.
- Far away from sharp spines, if the cloth cover are broken, stop to use it to avoid toes roll into the foot massager, and must send it to dealerships or designated maintenance organization for replacement the cover or repairing.
- Before cleaning or after massaging, please pull out the power plug from the socket.
- Please do not put more cover by yourself when massaging, in case of leading the motor temperature too high and cause failure or damage.
- During massage. once feel uncomfortable, please immediately stop to use it.
- Once machine can't work well, please contact the service center for inspection and maintenance.
- Please use the manufacturers recommend spare parts only.
- Please don't put anything into the opening mouth of the foot massagers.
- Please do not use your wet hand to operate the hand controller.
- Don't let water touch the hand controller.
- If you want to close the power . must close all the working procedure firstly, then pull out the plug from the socket.
- This product is not for the self diagnosis or treatment, such as the need to diagnosis or treatment, please go for the professional hospital.

Hints on Location

- Please do not use the foot massage in humid and too much dust environment.
- Don't use the foot massage in room . where the temperature is higher than 40°C.
- Don't put the foot massages outdoor.
- Do not use foot massager nearby the electric heaters, stoves or under the direct shining sunshine.

The below person, please consult the doctor before using:

- Someone who are under the treatment by the doctor or whom feel unwell.
- Patients with high blood pressure,heart disease, cerebral arteriosclerosis, cerebral thrombosis.
- Patients with malignant tumor.
- Pregnant women and women in menstrual period.
- Patients with bone osteoporosis.
- Patients who are using a pacemaker or implanted medical apparatus.
- Patients who with acute disease.

Hints on the Duration of Use

- Use 1-2 times a day, 5-10 minutes per time, effectively eliminate sore muscles & improve blood circulation.Massage 5 to 15 minutes before going to bed.can effectively improve sleeping quality.

Hints on Use

- Don't sleep when using the foot massage.
- Don't use the foot massage when get drunk or don't feel well.
- Don't massage the foot more than 30 minutes , warm foot less than 2 hours per time.
- Once the fabric cover destroyed , stop to use it.
- Don't using it , one hour before eating.
- If the product used by low temperature storage, move to the warm environment, it may cause the parts inside water vapor and influence the operation , Even cause malfunction. So when moving, the product first brought to normal room temperature open in an hour, then move to the warm room accordingly.

Product Description

- 1 Time
- 2 Mode
- 3 Intensity
- 4 Heating
- 5 Power
- 6 Massage Bag

Time: to control the massage time period.

- **Blue:** massage time period is 10 minutes.
- **Red:** massage time period is 15 minutes.
- **Blue and Red:** massage time period is 20 minutes.

Use and Operation

Mode Key

- Mode 1 is the defaults mode, press the key, enter into mode 2, press one more key, enter into mode 3, in this cycle.
- Mode 1: blue lights, mode 2: red light, mode 3: alternating red and blue lights flashing.
- Model 1 and 2 are foot turn and gasbag massage, mode 3 is for airbag massage.

Intensity Key

- Defaults to low intensity, press once, enter into the strength of a mid-range, press one more time, enter into the strength of high-grade, in this cycle.
- Low Intensity: blue lights; press once to middle grade Intensity: green lights; press once more to the high grade intensity: red lights.

Heating Key

The default is in the opening status, the sole heating function automatically open, light blue; press the heating key, can close the heating function, once the light close, the heating function stop.

Power Switch

- The foot massage connect to the power supply, is on the ready mode, red light flashing.
- Press the power button, the foot massage start working (heating + in the air waving + mode 1), All lights light up, the default time last 10 minutes.
- Working time finish, or the setting time reaches, foot massage returns to the ready mode. red light flashing.
- Finish using the foot massage, please close the power by getting out the plug from the socket

Cleaning and Care

Storage

- When no using, put the massage in the box and keep it in the safe, dry, airy and cool place.
- Don't tweak the wire.
- Do not use wire hanging products.
- Avoid sharp points scraping the surface broken.
- Please do not place the product in direct sunlight or high temperature place.

Cleaning

- Before cleaning, please make sure the power switch is closed, and pull out the power plug.
- Please use dry cloth to clean the massage surface. Do not use corrosive cleanser, such as Benzene or thinner cleanser.

Care

- Keep the product clean and far away from the wet place.
- If long time not use this product, please use dust cover trapping the massagers, prevent the dust fall into.

Troubleshooting

Problem	Possible Cause	Solution
Can not work.	No power.	Access to the power supply and check wheter the power supply is good.
	The function keys can't open.	Press the ON/OFF (Power) button.
Sudden stop while operating.	The power off.	Check the power interface is connected well or not.
	Massaging time is more than 15 Minutes.	Re-start the ON/OFF button.
	Continuous working too long, overheating.	Re-start the switch button (Power) after cooling down.
	Force is too large, massage ball blocked.	Relax your feet soles when massaging.

Hints on Disposal

According to the European waste regulation 2012/19/EU this symbol on the product or on its packaging indicates that this product may not be treated as household waste. Instead it should be taken to the appropriate collection point for the recycling of electrical and electronic equipment. By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product. For more detailed information about recycling of this product, please contact your local council or your household waste disposal service.

Declaration of Conformity

Producer: Chal-Tec GmbH, Wallstraße 16, 10179 Berlin, Germany.

This product is conform to the following European Directives:
2014/30/EU (EMC)
2011/65/EU (RoHS)

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Produttore: Chal-Tec GmbH, Wallstraße 16, 10179 Berlino, Germania.