SUPERLETIC

Atilla

Akupressurmatte

Acupressure Mat

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Tapis d'acupression

Tappetino per agopressione

10037713

10039844

10039845

10039846

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Dear customer,

Congratulations on the purchase of your appliance. Please read the following instructions carefully and follow them to prevent potential damage. We accept no liability for damage caused by disregarding the instructions or improper use. Please scan the QR code to access the latest operating instructions and for further information about the product.

CONTENTS

Safety Instructions 16

Care instructions 16

Basics of acupressure and acupuncture 17

The acupressure set in detail 18

Possible applications 19

MANUFACTURER & IMPORTER (UK)

Manufacturer:

Chal-Tec GmbH, Wallstraße 16, 10179 Berlin, Germany.

Importer for Great Britain:

Berlin Brands Group UK Limited

PO Box 42

272 Kensington High Street

London, W8 6ND

United Kingdom

SAFETY INSTRUCTIONS

The use of the acupressure set is not recommended for:

- **Inflammation:** In case of acute inflammation, you must avoid the affected area. These include wounds, febrile illnesses or vascular diseases.
 - **Skin diseases:** If you suffer from skin diseases, do not massage the irritated skin area. This precaution also applies to skin burns, such as sunburn.
 - **Poor general health:** Never use the acupressure set for serious illnesses such as tumour diseases, infectious diseases or heart diseases.
 - **Pregnancy:** During pregnancy, you should first talk to your attending doctor about the possible use of the acupressure set.
-

CARE INSTRUCTIONS

Note: Never wash the acupressure mat and cushion in the washing machine or tumble dry as this will permanently damage the tips of the acupressure pads.

- If necessary, you can remove the covers with the acupressure pads and wash them by hand in lukewarm water with a mild liquid detergent.
- Use a soft brush to thoroughly clean the acupressure pads as well.
- Allow the wet garments to drip dry and dry slowly while lying on a clothes horse. Do not expose the covers to direct sunlight, otherwise the pads may be damaged and the colour of the covers may fade.
- The soothing scent of the lavender pillow may fade over time. This is a typical characteristic of authentic lavender flowers. To release the essential oils in the flowers, take out the small pillow and press and rub it several times. Afterwards, the pleasant feel-good scent can unfold again.

BASICS OF ACUPRESSURE AND ACUPUNCTURE

Both acupuncture and acupressure developed from the theory of Traditional Chinese Medicine (TCM) and are based on a disturbed energy flow (also called Qi flow), which is the trigger for any pain. The individual pathways of energy are called meridians and their access points are the acupuncture points that are stimulated by therapy, thereby releasing the blockages and pain. In acupuncture, needles are used for stimulation, which trigger reactions in the autonomic nervous system through these energy points and alter the muscle tension.

How acupressure works

In contrast to acupuncture, acupressure stimulates the access points of the meridians via pressure points. The extensive treatment loosens hardening and tension and stimulates the circulation. In addition, the body releases happiness hormones (endorphin, oxytocin) through the treatment. These relieve pain, reduce stress, support self-healing and lead to an overall better sense of well-being. In most cases, the first results are visible after just a few applications. However, regular acupressure treatment is a condition for successful self-therapy.

When acupressure can help

Acupressure can be seen as invigorating relaxation and increases the body's oxygen intake. It not only helps with the most common symptoms such as headaches, backaches and neck pain, but can also help regulate blood pressure, alleviate circulatory problems, stimulate digestion, prevent or reduce disc problems, cure sciatic pain, improve cellulite and can also provide relief from stress and anxiety.

By stimulating the reflex zones of the feet, there is a connection to the internal organs and therefore these can also be positively influenced by acupressure. Therefore, the treatment with the acupressure mat can have a holistic positive effect on health.

THE ACUPRESSURE SET IN DETAIL

The acupressure set consists of a mat made of mite-repellent coconut fibres and a cushion filled with buckwheat husks. The covers are made of a skin-friendly cotton-linen blend. The total of 278 acupressure pads made of High Impact Polystyrene (HIPS), each with 30 tips, are very resistant due to the particularly hard plastic and are also completely safe in terms of health. A small lavender cushion is included as an extra. The soothing lavender scent supports self-therapy.

POSSIBLE APPLICATIONS

Neck and shoulder

Neck pain is one of the most common side effects of migraines and headaches. You can relieve tension in your neck and shoulders while lying on your back on the mat and pillow by moving your head slowly and in a controlled way from top to bottom or from side to side. To treat the neck more intensively, you can also roll up the mat and then lie on it with your neck and move your head.

Effect after **10** minutes:

- Immediate relief of pain in the neck and shoulder area.
- Increased release of endorphins, helps with migraines and headaches.
- Promotes blood circulation.

Back

The total of 8,340 tips stimulate the pressure points on the back. Position yourself on the mat in a supine position and rest your head on the floor. The legs are slightly bent. If required, you can place your head on the pillow. As soon as your back is used to the treatment, you can use your arms to change and intensify the treatment by pulling your knees closer to your body with your arms. If the pressure feels too intense for you, you can also place your legs stretched out on the floor first.

Effect after **20** minutes:

- Immediate relief of pain in the entire back area.
- Muscular relaxation.
- Improved blood circulation.
- Relief of the spine.

Lumbar vertebrae

Back pain often originates in the lumbar spine and there is a close connection to the intervertebral discs. Regular treatment of the area is already beneficial for the prevention of disc damage. Carefully lay down with your back on the mat and put your feet up. After a few minutes you can increase the pressure on the lumbar region by pulling your feet to your abdomen with your arms.

Alternatively, you can place the cushion on the mat so that your lower back is arched over the mat.

Effect after **15** minutes:

- Prevention of herniated discs.
- Relief from back pain.
- Holistic relaxation.

Belly

A treatment on the abdomen has an activating effect and vitalises. In addition, digestion is stimulated and stomach/intestinal disorders can be positively influenced. To do this, lie on your stomach facing the mat. Take care to avoid the breast, as the tissue around the breast is very sensitive. You can put your head on your arms, but also on the pillow and still treat your jaw at the same time.

Effect after **10 minutes**:

- Stimulates digestion.
- Prevention of panic attacks.
- Relaxation in the chest area.

Waist

To treat the waist, position yourself sideways on the mat. In order to be able to apply the pressure more specifically, it is recommended to place a rolled-up towel underneath. You can rest on your elbow, but you can also rest on your entire side. You will feel a pleasant warmth which has a relaxing effect on your muscles.

Effect after **10** minutes:

- Stimulates digestion.
- Prevention of panic attacks.
- Relaxation in the chest area.

Arms and hands

After a long day at the desk, the forearms and hands are often limp and feel heavy. With the help of the acupressure mat you will be able to re-energise them. Use your hands to familiarise yourself with the treatment and place them on the mat. Once you get used to the feeling, you can additionally place the forearms on the mat and treat both the front and the back.

Effect after **15** minutes:

- Prevention of numbness in the hands and forearms.
- Relaxation after a hard day's work.
- Lighter fingers, hands and forearms.

Legs

To treat the waist, position yourself sideways on the mat. In order to be able to apply the pressure more specifically, it is recommended to place a rolled-up towel underneath. You can rest on your elbow, but you can also rest on your entire side. You will feel a pleasant warmth which has a relaxing effect on your muscles.

Effect after **10** minutes:

- Stimulates digestion.
- Prevention of panic attacks.
- Relaxation in the chest area.

Feet

All organs and body parts are found in the foot reflex zones. With the acupressure mat, these can be stimulated and positively influenced. To do this, you can place your feet on the mat either in a sitting or standing position. In the standing position, the pressure is more intense and can be further intensified by standing on one leg. You can do this exercise, for example, while working at the computer or watching TV in the evening with no additional effort.

Effect after **5** minutes:

- Reflexology treatment for the whole body.
- Promotes blood circulation.
- More strength and energy for the day.

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