


KYOTO[®]

The World's Most
Advanced Massage Chairs[®]



User Manual

**Kokoro M888
4D Massage Chair**

Staying Safe

Important Safety Instructions

Read all instructions before using this massage chair.

When using this massage chair, please take the following precautions:

DANGER – To reduce the risk of electric shock:

1) Always unplug the massage chair from the electrical outlet immediately after using and before cleaning.

WARNING – To reduce the risk of burns, fire, electric shock, or injury to persons:

- 1) Massage chair should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
- 2) Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
- 3) Close supervision is necessary when this massage chair is used by, on, or near children, invalids, or disabled persons.
- 4) Use this massage chair only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- 5) Never operate this massage chair if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the massage chair to a service center for examination and repair.
- 6) Do not carry this massage chair by supply cord or use cord as a handle.
- 7) Keep the cord away from heated surfaces.
- 8) Never operate the massage chair with the air openings blocked. Keep the air openings free of lint, hair, and the like.
- 9) Never drop or insert any object into any opening.
- 10) Do not use outdoors.
- 11) Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- 12) To disconnect, turn all controls to the off position, then remove plug from outlet.
- 13) Use heated surfaces carefully. May cause serious burns. Do not use on insensitive skin areas or in the presence of poor circulation. The unattended use of heat by children or incapacitated persons may be dangerous.
- 14) Do not stand on or in the massage chair. Use only while seated.
- 15) Keep children away from extended foot support (or other similar parts).
- 16) Connect the massage chair to a properly grounded outlet only. See Grounding Instructions.
- 17) Do not use the massage chair in close proximity to loose clothing or jewelry.
- 18) Keep long hair away from massage chair while in use.

SAVE THESE INSTRUCTIONS.

Not Suitable For Below Persons

We recommend those with any of the following conditions consult your physician before using.

- Those who are using electronic medical instruments embedded in the body, such as a cardiac pacemaker.
- Those undergoing treatment, especially those who are not feeling well.
- Those with malignant tumors or malignant abscesses. (Need to be supervised)
- Women who are pregnant.
- Those with osteoporosis or a fractured spine.
- Those with skin diseases or those with skin injuries.
- Those with a high fever.
- This massage chair can be used by children aged 8 years and above and persons with reduced physical, sensory, or mental capabilities, or lack of experience and knowledge, if they have been given supervision or instruction concerning the use of the massage chair in a safe way and understand the hazards involved.
- Children shall not play with the massage chair.
- Cleaning and user maintenance shall not be made by children without supervision.

Best Circumstances

- Do not use this massage chair when ambient temperature is over 104 F.
- Women who are pregnant should not use this chair.
- Do not use this massage chair in humid or dusty environment, or when any part of the body is in contact with plumbing or any similar ground.
- Please place the product on a flat level floor.
- Do not store near heat or open flame.

Control Buttons Precautions



- Do not operate the control buttons with wet hands.
- Do not spill water or other liquids onto the control buttons to avoid fault.
- Do not lay heavy objects on the controller.

Cautions Regarding the Adapters and Cord



- Always unplug the massage chair from the electrical outlet immediately after using and before cleaning.
- Massage chair should never be left unattended when plugged in. Unplug from outlet when not in use.
- Unplug from outlet during thunder or lightening storms, and before putting on or taking off parts.
- Check that the power voltage stated on the rating label corresponds to your local electricity supply.
- Do not use a damp hand to unplug the cord to avoid electric shock.
- Do not drag the cord, reach to the front plastic end to unplug it.
- Do not squeeze, over bend, pull hard, or twist the cord.
- To disconnect, turn all controls to the off position, then remove plug from outlet.
- The massage chair is only to be used with the power supply unit provided with the massage chair.
- It must only be supplied at safe extra low voltage corresponding to the marking on the massage chair.
- If the supply cord is damaged, it must be replaced by the manufacturer or its service agent or a similarly qualified person in order to avoid a hazard.

Best Practices for Safe Use




- Only for users who weigh below 300 lbs.
- Massage chair should only be used by one person. Do not allow two or more persons to sit in the chair at the same time to avoid damage to the product.
- Do not use the product right after a meal, it may cause discomfort.
- Check the upholstery under the cushions and other spots, if any damages are present, do not use and refer to a qualified service professional for servicing.
- If you feel any abnormal discomfort during the operation of this product, turn off the power immediately and do not use the product.
- Use of the massage chair for 20 minutes each time is highly recommended.
- Make sure there are no obstacles behind the chair and there is enough clearance before adjusting the backrest or legrest angle.
- Make sure there is no other excessive load pressed on the backrest or legrest besides the user.
- Make sure there are no children, pets or other objects under the backrest or legrest.
- Please be seated properly while using the product.
- Do not drag or push the chair, or roll it on its caster on wood or carpet flooring to avoid damage.
- Do not fall asleep while the product is operating.
- Avoid using the product under the influence of alcohol.
- If there is any smoke or unusual smell while using, stop using immediately, and unplug the power cord.
- Do not disassemble this product yourself to avoid failure.
- Do not dismantle, fix or modify the product yourself. Please contact authorized dealer or service agent.

Grounding Instructions

This product must be grounded. If it should malfunction or breakdown, grounding provides a path of least resistance for electric current to reduce the risk of electric shock. This product is equipped with a cord having an equipment-grounding conductor and a grounding plug. The plug must be plugged into an appropriate outlet that is properly installed and grounded in accordance with all local codes and ordinances.


DANGER - Improper connection of the equipment-grounding conductor can result in a risk of electric shock. Check with a qualified electrician or serviceman if you are in doubt as to whether the product is properly grounded. Do not modify the plug provided with the product - if it will not fit the outlet, have a proper outlet installed by a qualified electrician.

This product is for use on a nominal 120-V circuit and has a grounding plug that looks like the plug illustrated. Make sure that the product is connected to an outlet having the same configuration as the plug. No adapter should be used with this product.




Your Chair

External Structure



Remote Control



Welcome to

KYOTA®

Thank you, and congratulations on purchasing your very own Kyota Massage Chair product, you've made a great buying decision. Welcome to the Kyota family.

Inside this manual, you will find everything you need to setup and safely enjoy your chair. Please review the entire manual before assembling or using your chair. If you are assembling the chair yourself, please see the Assembly Guide next and reference the online assembly video.

For convenience, we've also created a Quick Start Guide. We suggest you keep it in the remote pocket as a tool to get to know all the features of your new chair.

Sit back, relax, and enjoy your new Kyota massage chair!

Assembly Video

Before assembling your chair, watch an assembly video at:
www.kyotamassagechairs.com/assembly-tutorials

Alternatively, **scan the QR code** to go directly to the URL above.



SCAN CODE

Warranty Registration

Register your product at **www.kyotamassagechairs.com/warranty/registration** to ensure it is in our system. This will simplify the warranty process should you ever need it.

Alternatively, **scan the QR code** to go directly to the URL above.



SCAN CODE

Chair Features

- **Tri-zone Deep Penetrating Heat**
- **Full-Body Scanning Technology**
- **Zero Wall Technology**
- **10 Auto Programs**
- **Customizable Massage**
- **Air Ionizer**
- **Oscillating Calf Massage**
- **Arm and Hip Air Compression**
- **Chromotherapy mood LED lighting**
- **Shoulder Air Compression**
- **Three Zero Gravity Positions**
- **Wireless Controller**
- **Quick Access Control Panel**
- **Total Sole Reflexology™**

Table of Contents


Staying Safe	2	Getting Started	
Preparation		Airbag Massage	18
Installation Site	8	Airbag Region	18
Clearance Space for the Chair	8	Airbag Intensity	18
Floor Protection	8	All Airbags Off	18
Move Methods	9	Foot Rollers	19
Connect to Power	9	Heating	19
		Air Ionizer	19
Chair Setup		Settings	
Main Power Switch	10	Setting Tool	20
Adjust Pillow	10	Language	20
Increase the Intensity of your Massage - Manually	10	LED Lights	20
Connecting the Wireless Remote Control to Bluetooth	11	Timing	20
Quick Access Control Panel	12	Stop Massage	20
Body Scanning	13	Bluetooth Audio	20
Getting out of the Chair	13	Intelligent Voice Control	21
Auto Programs		Fixing Issues	
Auto Programs (AUTO)	14	Troubleshooting	22
Remote Control Screen Idle	14	FCC Information	22
		Frequently Asked Question	23
		Specifications	23
Create a Massage			
Manual Programs	16		
Massage Mode	16		
Massage Position	16		
Massage Width	17		
Massage Speed	17		
4D Massage Intensity	17		
Footrest	17		

Preparation

Installation Site

Clearance Space for the Chair


- 4" clearance for the recline of the backrest; 20" clearance for the extension of footrest and sliding forward of the chair.
- Keep it 3 feet away from the TV, radio or other audio & video appliance to avoid signal interference



- Keep the distance between the power socket and power switch within 60" to avoid damage of power cord while massage chair sliding forward.
- Straighten the power cord while in use to avoid tangle or damage.
- When the massage chair sliding forwards or back, please check to make sure there is no children, pets or other obstacles around to avoid accident.

Floor Protection


We recommend to place a piece of thin carpet or soft mat underneath it to avoid damaging the floor.



- Do not use the massage chair in high moisture environment such as near swimming pool or bathroom to prevent electric shock accident.
- Please keep the massage chair away from heated surfaces to prevent fire hazard or damage to the product.

Move Methods

- Apply pressure to the back of the massage chair from behind to pivot on the Castors. Have a second person up front to guide the Footrest until in the correct position. Maintain clearances as mentioned above. Gently lower into place. Do not pull on the side arms.




Castors



- Do not move the chair when there is a user inside.
- To carry this product, do not lift by the footrest or armrest, hold the chassis and steady the backrest (have 4 or more people to carry the product).

Connect to Power

Plug into a 110V grounded outlet.



This product must be grounded. If it should malfunction or breakdown, grounding provides a path of least resistance for electric current to reduce the risk of electric shock. This product is equipped with a cord having an equipment-grounding conductor and a grounding plug. The plug must be plugged into an appropriate outlet that is properly installed and grounded in accordance with all local codes and ordinances.




- Ensure there is no damage to the supply cord before switching on the massage chair.
- Ensure that the main power switch is turned off before connecting to power

Chair Setup

Getting Comfortable

Main Power Switch

To power up the chair, flip up the main power switch located on the back of the base of the chair to **ON**.



- Children should be supervised to ensure that they do not play with the chair or touch the power box.
- Always unplug the massage chair from the electrical outlet immediately after use.
- Unplug from outlet during thunder or lightening storms, and before putting on or taking off parts.
- Turn off the power switch and unplug the massage chair from the electrical outlet while not in use to avoid children power on the massage chair accidentally.
- When not in use for a long period of time, we recommend that the chair be covered, the power cord be coiled and that the unit be placed in an environment free of dust and moisture.




Adjust Pillow

Take a seat in the chair and take note of where the neck pillow rests on your neck. Adjust the pillow so that it rests comfortably in the arch of your neck. Do this by removing it from the velcro, aligning it, and then reattaching to the velcro once it is in the desired position.

Increase the Intensity of your Massage - Manually

To make the massage more intense, remove the **neck pillow** and **shoulder pad** cushioning between you and the massage mechanism in the seat back. To make the foot rollers more intense, remove the foot pads that came with your chair from the foot bed of the footrest.



Connecting the Wireless Remote Control to Bluetooth


Before you pair your wireless remote control with the massage chair, make sure the chair is powered on. To power on the chair without your wireless remote, press the **ON/OFF** button on the quick access control panel located on the top of the right arm panel.

To pair your wireless remote control, simultaneously hold down the **MENU** and **OK** buttons for 4 seconds. The screen will indicate that Bluetooth pairing has begun. If the connection is successful, the Bluetooth icon will light up in the top right corner of your remote control screen.

After pairing successfully, turn the remote control off and then on again to complete pairing procedure. Your wireless remote control is now ready to use to operate your massage chair.

Please note: Your wireless remote control should be charged upon receipt of your chair. If it won't turn on or has low battery, simply charge it by plugging in the provided USB cable into the USB port on the bottom of your remote control.


Controller Reset: If the remote control disconnects or experiences pairing issues, use a pin or similar object to press the reset button located on the bottom of your remote control, and then repeat the pairing process above.



Chair Setup

Quick Access Control Panel

Another feature of your chair is a built in quick access control panel. The quick access control panel is located on the top of the right arm panel. This control panel allows you to turn your chair ON/OFF, emergency shutoff the chair, cycle through 4 Auto programs, adjust the legrest and backrest, or toggle through zero gravity positions.



On/Off: Short press; turn chair ON/OFF

Emergency Shutoff: Short press; immediately stops all chair functions.

AUTO: Short press; toggles through auto programs shown. Program icon will light up when active.

Adjust legrest: Press and hold; to raise or lower legrest.


Adjust backrest: Press and hold; to incline or recline backrest.

Zero Gravity: Short press; press once to enter ZG position 1. Press again for ZG position 2. Press a 3rd time for ZG position 3. Press a 4th time to return chair to upright position.

Body Scanning

Navigate to any auto program and press the **OK** button. The chair will automatically begin scanning your body, mapping your back, shoulders, and neck. During this process it is important to remain seated back into the chair without lifting your head.

The scan can be viewed on the remote control screen. If fine tuning is needed just after your body scan, wait for the **shoulder height adjust** prompt and then use the **up** and **down** buttons to move the massage rollers up or down for the perfect shoulder kneading position.



Getting out of the Chair


To easily exit the chair, press the **power** button on the remote control or the quick access control panel to restore the chair to the upright position.

Auto Programs

Auto Programs (AUTO)

When the massage chair is operating, press the **AUTO** button on the remote control to enter the Auto Program menu display. Choose from **Health, Relief, Stretch, Music Sync, Renew, and Relaxation** programs. Additionally, select the Special option from the list to access 6 more auto programs: **Neck & Shoulder, Waist & Spine, Back Rehab, Active, Relax, and Healing.**

Press **OK** to select the program you desire. Selected program will show on the in-use screen.



Programs

Health

Pain Relief

Stretch

Music Sync

Experience

Relaxation

Special

Neck & Shoulder

Waist & Spine

Back Rehab


Active

Relax

Healing

Remote Control Screen Idle

Whilst running an Auto Program or in Manual mode, the remote control screen will return to the in-use screen (as shown) **after 10 seconds of being idle**. To exit the in-use display, press the AUTO or MENU button on your remote control



What for?

Description

This program is the best choice to relax tired muscles and increase physical vigor.

Comprehensive massage session using deep tissue kneading, and tapping.

This is a great program for people who work at a desk or travel a lot.

The most intense of the massage programs. This strong kneading and shiatsu massage effectively decompresses the entire body.

Best for winding down and re-energizing yourself after a long day.

Yoga like stretching motions grab, hold and gently extend the body, while deep shiatsu massage helps to boost the vitality of muscles.

Designed for morning or midday recharging.

Targets the lower back, buttocks and feet with kneading and tapping. Massage follows the rhythm of music.

A quick way to feel all the massage techniques and chair movements in a brief 8 minute session.

Rapid demonstration of all massage techniques and chair motions.

Reduce built up tension carried in the body from everyday stress.

Total body relaxation with shiatsu and rolling massage techniques, and air cell massage.

Soothing pain and tension in the neck and shoulders.

A sequence of neck, shoulder, and upper back massage techniques.

This program is recommended to relieve lower back pain.

Focuses on the lumbar vertebra with targeted kneading and tapping massage to relieve tension.

Relieving lower and upper back pain.

Assists in relieving lower and upper back pain through a series of massage techniques.

Recharge and re-vitalize your tired muscles after physical exertion.

Uses moderate pressure to decrease recovery time after physical exertion, relieve the body of stress and improve circulation.

Designed for relaxing. Try this program as you prepare your body to rest for the night.

The lowest intensity massage. Shiatsu and gentle rolling massage techniques sooth and relax the body.


Helps to relax and recover tired and sore muscles. Whole body massage with a focus on the main back muscles along your spine.

Featuring neck, back, and waist kneading this massage relaxes muscles, improves circulation and aids in recovery.

Create a Message


Manual Programs

Press **MENU** button on remote control. Select **Manual** from Menu display, then use up or down buttons to navigate among: **Mode**, **Position**, **Width**, **Speed**, **4D**, and **Footrest**.



Message Mode

Navigate to **Mode**, press right button, then use the up or down button to select your desired mode: **Kneading**, **Tapping**, **Knocking**, **Knead & Tap**, or **Shiatsu**. Press **OK** button to select mode.



Message Position


Navigate to **Position**, press right button, then use the up or down button to select your desired message position: **Whole**, **Partial**, or **Point**. Press **OK** button to select position.

Whole travels the entire length of the message track.

Partial travels a few inches above and below the current position.


Point stops the massage mechanism in current position.

Partial and Point location can be fine tuned using the arrow up or arrow down buttons.




Massage Width

Navigate to **Width**, press right button, then use the up or down button to select your desired roller width: **Narrow**, **Middle**, or **Wide**. Press **OK** button to select width. (Massage Width only works in Tapping, Knocking, and Shiatsu techniques).



Massage Speed

Navigate to **Speed**, press right button, then use the up or down button to select your desired speed: **1 = slowest**, **2, 3, 4, 5 = fastest**. Press **OK** button to select the speed setting. (Massage speed can only be adjusted when in manual mode).




4D Massage Intensity

To adjust the intensity of your massage, navigate to **4D**, press right button, then use the up or down button to select **4D+**, or **4D-**. Press **OK** button to select the 4D intensity setting. Current 4D intensity will display in the top left of the in-use screen.

4D+ will increase the massage intensity.

4D- will decrease the massage intensity.




Footrest

To adjust the extension of the footrest, navigate to **Footrest**, press right button, then use the up or down button to select **Out**, or **In**. Press **OK** button to select footrest setting.

Out will extend the footrest outwards away from the chair.


In will retract the footrest inwards towards the chair.



Getting Started


Airbag Massage

Press **MENU** button on remote control. Select **Airbag** from Menu display, press up or down button to navigate among: **Position**, **Intensity**, and **All Off**.




Airbag Region

Navigate to **Position**, press right button, then use the up or down button to select your desired airbag region: **Full Body**, **Upper Body**, **Waist & Seat**, and **Legs & Feet**. Press **OK** button to select




Airbag Intensity

Navigate to **Intensity**, press right button, then use the up or down button to select your desired intensity: **1 = light**, **2, 3, 4, 5 = firm**. Press **OK** button to select the airbag intensity.



All Airbags Off

Navigate to **All Off**. Press **OK** button to turn all airbag massage OFF.




Foot Rollers

Press **MENU** button on remote control. Select **Roller** from Menu display, then use the up or down button to navigate among: **Quick**, **Soft**, and **Stop**. Press **OK** button to select the speed level you desire.

Quick - Faster speed

Soft - Slower speed


Stop - deactivates foot rollers.



Heating

Press **MENU** button on remote control. Select **Heating** from Menu display, then press up or down button to navigate among desired heat settings: **Back**, **Waist**, **Foot**, **TEMP+**, and **TEMP-**.

Select **Back**, **Waist**, or **Foot** to toggle ON/OFF heating in that particular zone. Press **OK** button to select. Once selected, activated heating zones will display in the in-use screen as an icon next to the respective region.



To adjust the temperature of the heating, select **TEMP+**, or **TEMP-**. Current temperature is displayed in the top right of the in-use screen.

TEMP+ will increase the temperature.

TEMP- will decrease the temperature.

Air Ionizer

While the **In-Use screen** is displayed on your remote control, press the **OK** button to toggle ON/OFF the Air Ionizer. Once turned on, the air ionizer icon will appear on the screen.



Settings

Setting Tool

Press **MENU** button on remote control. Select **Settings** from Menu display, then use the up or down button to navigate among: **Language**, or **LED Lights**.

Language


Navigate to **Language**, press right button, then use the up or down button to select your desired language: **Chinese, English, Spanish, or Korean**. Press **OK** button to select language.

LED Lights

Navigate to **LED Lights**, press right button then press up or down button to select your desired LED light setting: **On** or **Off**. Press **OK** to select LED light setting.


Timer

Press **MENU** button on remote control. Select **Timer** from Menu display, then use up or down button to select your desired timer: **10 minutes, 20 minutes, or 30 minutes**. Press **OK** to select your desired massage timer length. 20 minutes is the default setting. The massage chair will turn off automatically and reset after the timer is up.



Stop Massage

Press **MENU** button on remote control. Select **Stop** from Menu display to stop the current massage and return to the main interface.



Bluetooth Audio

Your massage chair includes Bluetooth compatibility so you can play and enjoy your favorite music. To connect your Bluetooth enabled audio device, go into your audio device settings, navigate to Bluetooth settings and pair with the device named **Kyoto M888#####**.

Intelligent Voice Control

Your chair comes with built in intelligent voice control. This feature allows you to control your chair and massage experience completely hands-free by using one of the many voice commands.

To activate the intelligent voice control function, simply say **“Turn on speech recognition mode”**. If your voice command is recognized, the speaker will respond after 1 second. You may now give your chair any of the commands below **within 3 minutes**. Make sure to speak clearly and slowly. After 3 minutes the Intelligent Voice Command function will turn off automatically.

Voice Recognition Mode

At any time, Simply say **‘Voice off’** to turn voice recognition off. Otherwise, this mode stays active to receive commands for 3 minutes. After 3 minutes you will need to reactivate by saying the activation phrase again.

Activation Phrase: **“Turn on speech recognition mode”**.

Commands

“Health massage”	“Experience”
“Pain relief”	“Relaxation”
“Body stretch”	“Massage chair shutdown”
“Zero gravity start”	“Voice off”

Fixing Issues

Trouble Shooting

Symptoms	Possible Causes	Possible Solution
Noise is heard while in use	These noises are normal and arise from friction between the massage rollers and the leather cover	No need to take any measures
Massage Chair does not start	The power cord is not secured into the socket Not switched ON The fuse is blown	Insert the power plug properly into the power socket Turn on the switch found on chair side. Then push power button on remote Replace only with a fuse of same type and rating
No music comes from the speaker	The volume on your device is turned to the lowest setting	Ensure the volume on your device is turned up and the Bluetooth pairing was successful
Backrest or footrest can not raise or recline	Excessive load has been pressed on the backrest or legrest	Reduce the load and try again
Airbag malfunction	The air hose is blocked	Straight the air hose and clear out the blockage

FCC Information

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this device.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Frequently Asked Questions

Q: I find the chair is too intense to be comfortable. What can I do?

A: Your body needs to become acclimated to therapeutic massage. To be effective, the massage techniques must compress the tissue in your body. We recommend the following steps to help get started. Then over time as your body acclimates you can increase the intensity by reducing the layers between you and the chair.

1. Use the neck pillow and shoulder pad to reduce upper body intensity.
2. Use the foot bed pads that came with your chair.
3. Press **MENU** button on remote control. Select **Airbag** from Menu display, press up or down button to navigate to **Intensity**, then select **level 1**.
- If a specific airbag is too intense, turn off full body and turn on only the airbags you want.
4. Adjust the recline on the remote. More upright will be less intense. Zero Gravity will be more intense

Q: I'm not a foot massage person, what should I do?

A: Turn off foot rollers, Press **MENU** button on remote control. Select **Roller** from Menu display, then use the up or down button to **Stop**. Press **OK** button to select. You may also want to turn off the **Legs & Feet** airbags or reduce the air intensity to level 1.

Q: My chair is not performing as expected. What should I do?

A: Your Kyota massage chair is an amazing combination of programming and robotics. On occasion the chair just may need to be rebooted. To do this power down the chair and unplug from the wall for 5 minutes then plug back in and power back up to see if it self-healed. If you are still having an issue, give Kyota a call at 800-242-0761.

Specifications

Model	Kokoro M888 4D
Dimension of Main Body (in)	58 x 34 x 51 (L x W X H)
Necessary Clearance to Wall (in)	4"
Net Weight (lbs)	270
Shipping Weight (lbs)	435
Max. User Weight (lbs)	300
Voltage	120V
Fuse	3.15 AMP
Power Consumption	230W
Safety Standards	Type 1

KYOTA[®] | **The World's Most
Advanced Massage Chairs[®]**

72 Stard Road, Seabrook, NH 03874

Toll Free #: (800) 242-0761

www.kyotamassagechairs.com

Version 2.3