



The quality of Leonisa products far exceeds that of other brands...,

Dr. Elie Levine, MD
Plastic Surgery & Dermatology of NYC, New York, NY

REVIEWS

"It is recommended that post-surgical garments are kept on as much as possible during the prescribed timeframe by doctors. For this reason, EMCAS doctors recommend the use of Leonisa's high quality Post-Surgical Garments for their comfort to ensure our clients have the best care in their recovery."

- Dr. Phạm Xuân Khiêm, EMCAS Ho Chi Minh City, Vietnam

"I have used Leonisa garments for my patients for over nine years. Leonisa's customer service is top notch and their garments are durable, longlasting and aid in optimal results. In fact, many of my patients continue to wear them long after their recovery."

Craig Staebel, MD Georgetown Plastic Surgery Georgetown, TX

"After trying Leonisa... I have left all other surgical garment brands in the past. The quality of Leonisa products far exceeds that of other brands...Not only are my patients receiving what I believe to be the best surgical garments on the market today, but they are doing so at a fraction of what it once cost me."

Dr. Elie Levine, MD Plastic Surgery & **Dermatology of NYC** New York, NY

"The Leonisa post-surgical garments...are both "We are very happy we switched over to Leonisa" for our post-procedure compression garments. They are so comfortable and securing for the patients, a lot of them come back for more! They offer a variety of options and colors to fulfill each of our patients' needs as well."

Office of Dr. Tim Zelko, D.O. Zelko Aesthetics Greenfield, WI

WE'RE PROUD TO EXHIBIT IN:

- Atlanta Breast Surgery Symposium
- The Aesthetic Meeting
 - Plastic Surgery: The Meeting

What Makes Leonisa So Unique?

Decades' experience

in compression garment manufacturing



Made in-house, in

the Americas, with exclusive European machinery



Breathable. moisture-wicking, fast-drying

warp knit fabric



Zippers and hooks backed by a barrier to

protect incisions



Beautiful designs combine aesthetics and function



Comfortable enough

to wear even after recovery

SIZE CHART				
SIZE	WAIST (IN)	HIPS (IN)	DRESS/PANTS SIZE	
XS	25-27"	36-37.5"	0-2	
S	28-30"	38-40.5"	4-6	
М	31-33"	41-43"	6-8	
L	34-36"	43.5-46"	8-10	
XL	37-39"	46-48"	10-12	
XXL	40-42"	48.5-50.5"	14-16	
3XL	43-45"	51-53"	18-20	
4XL	46-48"	53.5-55.5"	22-24	





SIZES:

32B-42B 32C-42C

34D-42D

34DD-42DD





802 700 NUDE BLACK

SUGGESTED USES:

BREAST AUGMENTATION, REDUCTION, LIFT, RECONSTRUCTION





STYLE: 054004

SIZES: S/M-L/XL



852 SOFT NATURAL









Powerslim®





2 side zippers with inner hooks for easy wear

STYLE: 018489

SIZES:

S-M-L-XL-2XL



880 700 NATURAL TAN BLACK



SUGGESTED USES:

LIPOPLASTY, HIPS LIPOSUCTION, THIGHPLASTY



SIZES:

S-M-L-XL-2XL



880 NATURAL TAN



Powerslim®



SUGGESTED USES:

BBL, **BUTTOCK AUGMENTATION**

STYLE: 018478

SIZES:

S-M-L-XL-2XL





Powerslim®









SIZES:

S-M-L-XL-2XL








700 BLACK

SUGGESTED USES:

LIPOSUCTION, THIGHPLASTY





SIZES:

S/M - L/XL





852 700 SOFT NATURAL BLACK

SUGGESTED USES:

BRACHIOPLASTY



• **SkinFuse®** compression bands comfortably support

• 3 pieces with **easy-use** Velcro to attach and detach

- Soft, stretchy material is comfortable for all-day everyday wear
- Easy to fit and adjust according to size

SIZES:

UNIVERSAL FIT




852 SOFT NATURAL

SUGGESTED USES:

CHIN LIPOSUCTION, FACELIFT, NECKLIFT, JAW LIPOSUCTION, CHEEKBONE SURGERY, CHIN AUGMENTATION



Hook closure Adjust each belt for a perfect fit



Can be used with catheter to drain fluids by inserting tube between belts

STYLE: E15792

SIZES:

XS/S-M/L-XL/2XL





857 700 NATURAL BLACK

Skinfuse®

SUGGESTED USES:

SMARTLIPO, LASER LIPOSUCTION



Adjustable velcro closure

> Adjustable hook-and-eye gusset

STYLE: 012885

SIZES:

S-M-L-XL





802 NUDE

BLACK

Powerslim®

STYLE: 012807M

SIZES:


XS/S-M-L-XL/2XL-3XL





NATURAL

Skinfuse®



SUGGESTED USES:

SMARTLIPO, LASER LIPOSUCTION

STYLE: 012727

SIZES:

S/M-L/XL-2XL-3XL





852 BEIGE

Skinfuse®



Targeted compression smooths out torso and legs

> Sheer mesh in the back for no flattening

Mesh cutouts for comfort and flexibility







Hook-and-eye front closure for easy wear

Firmly compresses torso and upper thighs to promote healing and improved results

STYLE: 018491

SIZES:

S-M-L-XL-2XL-3XL





880 NATURAL TAN BLACK

Powerslim®

SUGGESTED USES:

ABDOMINOPLASTY, LIPOPLASTY, TUMMY TUCK, FLANK AND STOMACH, LIPOSCULPTURE, SMARTLIPO



Thick, adjustable straps

Covers bust without flattening

> Easy-use hook-and-eye front closure

Firmly compresses torso and thighs to promote healing and improved results



Microfiber on back for no flattening

STYLE: 018687

SIZES:

S-M-L-XL-2XL





880 NATURAL TAN

Powerslim®



4-position adjustable straps

> Front zipper with inner hooks for easy wear

> > Firmly compresses abs, chest, back, and thighs to help recover

Open gusset



STYLE: 038000

SIZES:

S-M-L-XL-2XL



880 NATURAL TAN



Powerslim®

SUGGESTED USES:

HIP LIPOSUCTION, ABDOMINOPLASTY, BARIATRIC, TUMMY TUCK, GYNECOMASTIA, LIPOSTRUCTURE, SMARTLIPO, HERNIA SURGERY

SIZE CHART LEO MEN'S SHAPEWEAR

SIZE	CHEST (INCHES)	HIPS
S	33.5 - 36.5	35.5 - 37.5
M	37 - 40	37.5 - 39.5
L	40.5 - 43.5	39.5 - 41.5
XL	44 - 47	41.5 - 43.5
2XL	47.5 - 50.5	43.5 - 45.5
3XL	51 - 54.5	45.5 - 47.5
4XL	55 - 58	47.5 - 49.5

SIZES:

S-M-L-XL-2XL-3XL-4XL





864 700 NATURAL BLACK

Powerslim®



SUGGESTED USES:

SMARTLIPO, LASER LIPOSUCTION

STYLE: 035020

SIZES:

S-M-L-XL-2XL



700 BLACK

Powerslim[®]



Leonisa POST-SURGICAL

Compression Garments

LEONISAWHOLESALE.COM

