


Congratulations on the purchase of your new Violife Slim Sonic toothbrush!
Slim Sonic operates at 22,000 brush strokes per minute for cleaner, healthier teeth and gums.

Parts

- A. Vented-Cap
- B. Battery Cover/Base
- C. Brush Head
- D. Power Button
- E. Battery Compartment
- F. Battery



How to Use

1. Push the power button firmly until vibration begins.
2. Brush teeth as recommended by your dental professional.
3. To turn Slim Sonic off, press power button firmly until vibration stops.

Changing the Replaceable Brush Head


Caution: Do not pull the brush head straight off.

1. Hold the base of the toothbrush firmly in one hand, and with the other hand, turn brush head counter clockwise and pull up.
2. Align new brush head and turn clockwise until brush head clicks into place.

Note: Notches on neck should align with inner notches of toothbrush head to lock in place.

Ordering Extra Brush Heads

Please check with your place of purchase or visit www.violife.com for replacement brush heads.




Care Instructions

1. After each use, rinse the brush head thoroughly with water.
2. Shake excess water from bristles. Wipe the handle dry after each use. Dry your brush head before covering and storing.
3. Remove batteries if not used for an extended period of time.
4. To clean your Slim Sonic, wipe surface with a damp cloth. Make sure the battery cover is closed securely to prevent water from getting into the battery compartment.

Caution

- Do not immerse your Slim Sonic in water.
- Keep out of reach of children and pets as this product may be dangerous for them.
- **Warning: Choking hazard** - this product contains small parts. Seek immediate medical attention if swallowed or inhaled.
- Do not use this product for other purposes.
- Do not expose this product to any chemical solvent, direct sunlight or high temperature.
- If an unexpected event happens when operating, please report it to www.violife.com/customerservice.



line on back of base aligns with dot on body to lock in place

Replacing the Battery

Alkaline batteries are recommended.

1. Hold the battery cover with one hand and grip firmly around both sides of the power button section with your other hand.
2. Turn the battery cover clockwise and pull down, exposing the battery compartment.
3. Remove the old battery. Insert a new AAA battery.
4. Replace the battery cover and turn counter clockwise until it is securely in place.

Note: When battery is low, vibration will weaken and product may stop working. This is not a malfunction. You will need to replace your battery.

Battery Disposal

Follow applicable laws and regulations for transport, shipping and disposal of batteries. Contact your local recycling service.

Dispose of at public collection point in the EU countries – 2002/96/EC WEEE Directive.





Warranty & Return/Exchange/Refund Policy

EN

Questions?

Email customerservice@violife.com
or go to www.violife.com/customerservice.

ACTIVATE YOUR WARRANTY!

Visit www.violife.com/register to register your Slim Sonic Toothbrush and receive discounts and special offers on Violife products!

Limited 1 Year Warranty With Proof of Purchase

Violife LLC warrants to the original purchaser of this product that it is free from defects in materials and workmanship for a period of 1 year from date of purchase.

If the unit becomes defective within the 1-year period, please email us at customerservice@violife.com or go to www.violife.com/customerservice for replacement or troubleshooting assistance. If your request is approved, you may be required to submit the following:

- Original proof of purchase with a valid date even if the item was a gift
- Defective product

At Violife's discretion we will:

- Issue a credit for the full purchase value
- Replace the unit with a comparable product, free of charge

Note: Warranty does not cover scratches, abrasions, discolorations, fading, damage from misuse, abuse, alterations, neglect, mold, other incidentals resulting from normal wear and tear, brush heads and replacement batteries.

Warranty Exchanges for Third Party Purchases

(Items not purchased directly through Violife)

Any warranty exchanges for third party purchases will require approval from Violife.

Please follow the above instructions on how to contact Violife.

This warranty gives you specific legal rights and you may also have other rights, which vary from state to state.

Other Inquiries:

Please email us at customerservice@violife.com with as much information as possible about your order, starting with your first and last name, order number, email address, phone number and nature of inquiry. Customer Service will reply within 2-3 business days.

Return/Exchange/Refund Policy

Email customerservice@violife.com or go to www.violife.com/customerservice for further assistance.

All return/refund/exchange requests need to be approved.

Item(s) without a Return Authorization Number (RA#) cannot be processed.

Returns/refund/exchange within 30 days for United States and International:

You have 30 days to decide if the item you purchased is right for you. Violife will exchange, refund or credit the purchase value of any item(s) within 30 days of purchase. Exchanges are not possible at this time for international orders. Proof of purchase is required even if it is a gift. The original shipping and handling fees are not refundable.

Return/Exchange Policy – Third Party Purchases

(Items not purchased directly from Violife):

Violife reserves the right to exempt third party purchases from return/exchange policy. If item is not directly purchased from Violife no refund will be given. Customers will be required to submit proof of purchase, even if it is a gift, in order to determine eligibility. The original shipping and handling fees are not refundable.

If it is determined that your return/exchange request is ineligible, customers will be instructed to contact their original place of purchase.

Redeeming Violife Store Credit:

If you have received a Violife store credit, you can redeem it online at www.violife.com.

Return/Exchange/Refund Status:

Please allow 2-3 business days from the time you emailed us for your inquiry to be processed.